
# FLAEM DUO

AEROSOL + DOCCIA NASALE

**I** MANUALE ISTRUZIONI D'USO **GB** INSTRUCTIONS FOR USE MANUAL

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# Duo

Mod. P1207EM F400

## Complete system for the treatment of the entire respiratory tract

Congratulations on your purchase and thank you for choosing our unit. Our goal is to fully satisfy our customers by offering them cutting-edge products for treating ailments of the respiratory system. **Please read these instructions carefully and keep them handy for future reference. Use this device only as described in this manual. This is a home medical device for the nebulization and administration of medication prescribed or suggested by your doctor.**

Please visit [www.flaem.it](http://www.flaem.it) to see our entire range of Flaem products

### THIS DEVICE COMES EQUIPPED WITH:

#### Ⓐ Aerosol Therapy Device

- A1 - Switch
- A2 - Air outlet port
- A3 - Air filter
- A4 - Nebulizer port
- A5 - Carrying handle
- A6 - Power cord
- A7 - Cable compartment

#### Ⓑ Connection tube (main unit / nebulizer)

#### Ⓒ Accessories

- C1 - RF7 Dual Speed Nebulizer
  - C1.1 - Bottom piece
  - C1.2 - Nozzle
  - C1.3 - Top piece
  - C1.4 - Speed selector
- C2 - Child SoftTouch mask
- C3 - Adult SoftTouch mask
- C4 - Mouthpiece
- C5 - Adult nasal prong

#### C6 - Nasal shower RHINO CLEAR®

- C6.1 - Nasal adaptor with a soft rubber tip, suitable for children up to 3 years of age
- C6.2 - Nasal adaptor with a soft rubber tip, suitable for children from 3 to 10 years of age
- C6.3 - Nasal adaptor with soft rubber tip, suitable for 10 years of age and over
- C6.4 - Separator
- C6.5 - Atomizer
- C6.6 - Shower body with ergonomic handle
  - C6.6.1 - 10 ml nebulization chamber
  - C6.6.2 - Collection chambre
- C6.7 - Dispensing button
  - C6.7.1 - Button
  - C6.7.2 - Rubber cap
- C7 - Nasal shower-tube connect fitting

#### Ⓓ Carrying bag

### ⚠ IMPORTANT SAFEGUARDS

- Before using the device for the first time and periodically during the product's life, check the integrity of its structure and of the power cable, making sure that there is no damage; if it is damaged, do not insert the plug and take the product immediately to an authorized technical assistance centre or to your local dealer.
- The average lifespan of accessories is 1 year; however, we recommend you replace the nebulizer every 6 months in intensive use (or earlier, if the nebulizer is blocked) to always ensure maximum therapeutic effect.
- Children and disabled persons should always use the unit under strict supervision of an adult who has read this manual.
- Several parts of the device are small enough to be swallowed by children; therefore keep the device out of reach of children.
- Do not use the tubes and the cords of the device for purposes other than those specified as they may cause strangulation hazard; utmost attention must be paid to children

and persons with special needs; usually, such persons are not capable of assessing the hazard correctly.

- The unit should not be used in the presence of an anaesthetic mixture inflammable with air, oxygen or nitrous oxide.
  - Keep the cord away from heated surfaces.
  - Keep the power supply cord away from animals (for example rodents) as they may damage the insulation of the power supply cord.
  - Do not handle the unit with wet hands. Do not use the device in humid environments (e.g. while taking a bath or shower). Do not immerse the device in water; should this accidentally occur, unplug the device immediately. Do not remove or touch the device immersed in water before pulling out the plug. Take the device immediately to an authorized service center or to your local dealer.
  - Only use the device in dust-free environments; otherwise, therapy may be compromised.
  - The casing of this unit is not protected against the penetration of liquids. Do not wash the device under running water or submerge in water. Keep away from splashing water or other liquids.
  - Do not expose the device to high temperatures.
  - Do not put the device near heat sources, direct sunlight, or in hot environments.
  - Do not obstruct or insert any objects in the filter or in its housing.
  - Do not obstruct the air slots placed on both sides of the unit.
  - During use, always place the unit on a firm surface clear of obstacles.
  - Check that there is no material obstructing the air slots before each use.
  - Do not place any objects inside the air slots.
  - The average duration provided for compressor ranges is: F400: 400 hours, F700: 700 hours, F1000: 1000 hours, F2000: 2000 hours.
  - Repairs must be performed by authorized personnel only. Any unauthorized repairs will void the warranty and may pose a safety hazard for the user.
- WARNING:** Do not modify this equipment without authorization of the manufacturer
- The Manufacturer, the Retailer and the Importer are responsible for the safety, reliability and performance of the device only if: a) the device is used in accordance with the operating instructions b) the electrical system of the premises where the device is used is in compliance with current laws.
  - Interactions: The materials used that come into contact with medications have been tested with a wide range of medications. However, due to the variety and continuous evolution of medication, the possibility of a chemical interaction cannot be excluded. We recommend that once the medication has been opened you should use it as soon as possible and avoid prolonged contact of the medication with the nebulizer.
  - The manufacturer should be contacted for reporting problems and/or unexpected events related to device operation.

## **NASAL SHOWER RHINO CLEAR® INSTRUCTIONS FOR USE**

**Before each use, wash your hands thoroughly and sanitize RHINO CLEAR® as described in section “CLEANING, SANITIZATION, DISINFECTION AND STERILIZATION”. During the application it is recommended to be adequately protected against dripping.**

1. Remove the adapter (C6.1, C6.2, C6.3) and the separator (C6.4).
2. Place the solution (max 10 ml) in the nebulizing chamber (C6.6.1) after having inserted the atomizer (C6.5).
3. Apply the separator (C6.4) verifying that the arrow on the shower body is aligned with

the separator (Assembly diagram). Then, insert the correct adapter for your age, or the most appropriate to fit the size of your nostrils (C6.1, C6.2, C6.3).

4. Connect the RHINO CLEAR® device through the tube (B) and the fitting (C7).
5. Turn on the device using the switch (A1)
6. Place the device close to the nose and gently insert the adapter nosepiece into one nostril.
7. During treatment it is preferable to hold the other nostril closed, where the shower is not applied, by applying light finger pressure.
8. Run the shower holding down the manual dispensing button (C6.7) during the inhalation phase and release it before exhaling; at this stage, to help expel the mucus, blow gently through the nostril treated directly into the RHINO CLEAR®.
9. Repeat from Step 6 to Step 8 for the other nostril, trying to evenly distribute the solution between the two nostrils.
10. When used on children or other dependent people, or in any case when there are problems in coordinating the various phases of paragraph 8., we recommend to run the application alternating the nostrils about 20 seconds each, with pauses required to expel the mucus.
11. The operations described in the preceding paragraphs is useful for nasal cleaning using physiological solutions. The treatment with drugs is advisable after a first cleaning with physiological solution in order to achieve a more effective distribution of the drug on the nasal mucosa.
12. After application, turn the device off and sanitize it following the instructions given in **"CLEANING, SANITIZATION, DISINFECTION AND STERILIZATION"**.

## **RF7 DUAL SPEED NEBULIZER INSTRUCTIONS FOR USE**

**Before each use, clean hands thoroughly and clean the device as described in the section on "CLEANING, SANITIZATION, DISINFECTION AND STERILIZATION". During use, it is advisable to protect yourself from any dripping. It is recommended that each person use their own nebulizer cup and accessories to prevent risk of infection due to contamination.**

1. Plug the power cord (A6) into a power socket corresponding to the voltage of the unit. This must be positioned so that it is not difficult to make the disconnection from the mains.

2. Open the nebulizer by turning the top piece (C1.3) anticlockwise.

3. Pour the medicine prescribed by your physician in the bottom piece (C1.1). Close the nebulizer by turning the top piece (C1.3) clockwise.

4. Connect the accessories as shown in the "Assembly diagram".

5. Sit in a comfortable position holding the nebulizer in your hand. Place the mouthpiece on your mouth or use the nose piece or mask. Should you opt for the mask accessory, place it on your face as shown in the picture (with or without using the elasticated strap).

6. Turn on the device using the switch (A1) and then breathe in and exhale deeply.

7. When the treatment is finished, shut off the device and unplug the power cord.

**WARNING:** If after the therapy session moisture accumulates in the tube (B), remove the tube from the nebulizer and dry it using the device's compressor fan; this operation will prevent mold from forming inside the tube.


## How to use the “RF7 Dual Speed” nebulizer with the speed selector

Professional and fast, this device is suitable for administering all types of medications, including more expensive ones, even in patients with chronic diseases. The geometry of the internal lines of the RF7 Dual Speed nebulizer ensures the ideal granulometry for active treatment all the way down to the lower respiratory tract.


To speed up the inhalation therapy, move the speed selector button (C1.4) by pressing on MAX with your finger.


To increase the effectiveness of the inhalation therapy, press with your finger on the wording MIN of the speed selector (C1.4). In this position, the speed selector acts as a valve and allows nebulizing the optimal amount of drug to the lower respiratory tract, reducing its dispersion in the environment.

## SoftTouch Mask

Soft biocompatible material  
Dispersion Limiting Device

**SoftTouch** masks have an outer edge made of **soft biocompatible material** that ensures excellent adherence to the face, and is also equipped with an **innovative Dispersion Limiting Device**. These distinctive elements that distinguish it allow greater sedimentation of medication in the patient and also **limit dispersion**.

During the inspiratory phase, the **tab**, acting as a **Dispersion Limiting Device**, bends inwards towards the mask

During the expiratory phase, the **tab**, bends outwards from the mask.

## CLEANING, SANITIZATION, DISINFECTION AND STERILIZATION

*Turn off the unit before any cleaning and disconnect the power cord from the mains.*

### UNIT AND OUTER SURFACE OF THE TUBE

#### CLEANING

Only use a cloth dampened with antibacterial detergent (non-abrasive and free of solvents of any kind).

#### ACCESSORIES SANITIZATION, DISINFECTION AND STERILIZATION

We recommend a personal use of the accessories to avoid any risk of infection or contagion.

Disassemble the accessories as described below:

**C1 - RF7 Dual Speed Nebulizer** Open the nebulizer by turning the top (C1.3) counter-clockwise, then remove the nozzle from the top (C1.3), and disassemble it as shown in the “Assembly diagram” section C1.

**C7 - Disassemble the NASAL SHOWER RHINO CLEAR®** as shown in the “Assembly diagram” section C6.

#### SANITIZATION

Before and after each use, sanitize the nebulizer and accessories according to one of the methods described below.

**(Method A): The accessories that can be sanitized with the method A are (C1.1-C1.2-C1.3-C1.4-C2-C3-C4-C5-C6.1-C6.2-C6.3-C6.4-C6.5-C6.6-C6.7.1-C6.7.2-C7)**

Sanitize the accessories in warm (approximately 40° C) drinking water with mild dishwashing liquid (non abrasive) or in the dishwasher with hot cycle

**(Method B): The accessories that can be sanitized with the method B are (C1.1-C1.2-C1.3-C1.4-C2-C3-C4-C5-C6.1-C6.2-C6.3-C6.4-C6.5-C6.6-C6.7.1-C6.7.2-C7)**

Sanitize the accessories by immersion in a solution consisting of 50% water and 50% white vinegar.

Finally rinse thoroughly with warm (approximately 40° C) drinking water.

**(Method C): The accessories that can be sanitized with the method C are (C4-C5-C6.1-C6.2-C6.3-C6.4-C6.5-C6.6-C6.7.1-C6.7.2-C7)**

Sanitize the accessories by boiling in water for 20 minutes; it is preferably use use of demineralized or distilled water to avoid scale deposits.

After sanitization shake the accessories vigorously and place them on a paper towel, or alternatively, dry them with a jet of hot air (e.g. using a hair dryer).

## **DISINFECTION**

**The accessories that can be disinfected are (C1.1-C1.2-C1.3-C1.4-C2-C3-C4-C5-C6.1-C6.2-C6.3-C6.4-C6.5-C6.6-C6.7.1-C6.7.2-C7)**

The disinfection procedure described in this paragraph is to be performed before using the accessories and is effective only if all the prescriptions herein are met and all the components involved are previously sanitized.

The disinfectant to be used must be of electrolytic oxidizing chlorine type, (active ingredient: sodium hypochlorite), specific for disinfection and is available in all pharmacies.

### **Procedure:**

- Fill a container, suitable to contain all the individual components to be disinfected, with a solution of drinking water and disinfectant, while respecting the proportions indicated on the packaging of the disinfectant itself.
- Completely immerse each individual component in the solution, taking care to avoid the formation of air bubbles in contact with the components. Leave the components immersed for the period of time indicated on the packaging of the disinfectant, and associated with the concentration chosen to prepare the solution.
- Remove the components now disinfected and rinse thoroughly with lukewarm drinking water.
- After disinfection shake the accessories vigorously and place them on a paper towel, or alternatively, dry them with a jet of hot air (e.g. using a hair dryer).
- Dispose of the solution as recommended by the disinfectant manufacturer.

## **STERILIZATION (ONLY FOR NASAL SHOWER)**

**The components that can be sterilized are (C6.1-C6.2-C6.3-C6.4-C6.5-C6.6-C6.7.1-C6.7.2-C7)**

The sterilization procedure described in this paragraph is effective only if all the prescriptions herein are met, all the components involved are previously sanitized and the procedure itself is validated according to ISO 17665-1 standard.

**Unit:** Steam sterilizer with fractionated vacuum and overpressure according to EN 13060 standard.


**Procedure:** Package each component to be treated in sterile barrier system or packaging complying with EN 11607 standard. Place the packed components in the steam sterilizer. Run the sterilization cycle according to the instructions manual of the device, selecting a temperature of 134° C and a duration of 10 minutes.

**Storage:** Store the sterilized components as per the instructions of the sterile barrier package or system chosen.

## REPLACING THE FILTER

The device has a suction filter (A3) that should be replaced when it is dirty or changes color. Do not wash or reuse the same filter. Regular replacement of the filter is required to ensure correct performance of the compressor. Check the filter on a regular basis. Contact your local dealer or authorized service center for spare filters.

## SYMBOLS

Class II device

Important: check the operating instructions

Switch off

Switch on

Complies with: European standard EN 10993-1 “Biological Evaluation of Medical Devices” and European Directive 93/42/EEC “Medical Devices” Phthalates free. In conformity with Reg. (EC) no. 1907/2006

SN

Serial number of device

Manufacturer

Type BF applied part

Alternating current

Hazard: electrocution. Consequence: Death.  
Do not use device while taking a bath or shower

0051

CE Medical Marking ref. Dir. 93/42 EEC and subsequent amendments

TÜV Certification ref. EN 60601-1 3° Ed. - EN 60601-1-11

**IP21**

Enclosure protection rating: IP21.

(Protected against solid bodies over 12 mm. Protected against access with a finger; Protected against vertically falling water drops.)

## AVAILABLE SPARE PARTS

### Description

RF7 Dual Speed Nebulizer + Mouthpiece  
 Child SoftTouch mask + Adult SoftTouch mask  
 Connection tube  
 Nasal shower RHINO CLEAR®

### Code

ACO557P  
 ACO462P  
 ACO35  
 ACO558P

## TECHNICAL DATA

### Mod. P1207EM F400

Voltage:	230V~ 50Hz 130VA	115V ~ 60Hz	220V ~ 60Hz
Safety certifications:	 <div style="display: inline-block; vertical-align: middle; margin-left: 10px;"> <p>Type Approved</p> <p>www.tuv.com ID 200000000</p> </div>		

Max Pressure:	1.8 ± 0.3 bar
Compressor air output:	9 l/min approx
Noise level (at 1 m):	54 dB (A) approx
Operation:	Continuous use
Operating conditions:	
Temperature:	min 10°C; max 40°C
Air humidity:	min 10%; max 95%
Atmospheric pressure:	min 69KPa; max 106KPa
Storage conditions:	
Temperature:	min -25°C; max 70°C
Air humidity:	min 10%; max 95%
Atmospheric pressure:	min 69KPa; max 106KPa
Dimensions (L)x(P)x(H):	17.5x19.5x10.5 cm
Weight:	1.4 Kg

### APPLIED PARTS

Type BF applied parts include: (C2,C3,C4,C5,C6.1,C6.2,C6.3)

#### RF7 Dual Speed nebulizer

Medication minimum capacity: 2 ml  
 Medication maximum capacity: 8 ml  
 Operating pressure (with neb.): 0.65 bar approx.

	<b>Speed selector C1.4 in pos Max</b>	<b>Speed selector C1.4 in pos Min</b>
<sup>(1)</sup> <b>Delivery:</b>	0.50 ml/min approx.	0.22 ml/min approx.
<sup>(2)</sup> <b>MMAD:</b>	4.58 $\mu\text{m}$	3.78 $\mu\text{m}$
<sup>(2)</sup> <b>Breathable fractions &lt; 5 <math>\mu\text{m}</math> (FPF):</b>	54.4%	63%

<sup>(1)</sup> data detected according to Flaem I29-P07.5 internal procedure  
<sup>(2)</sup> In vitro characterization certified by TÜV Rheinland Product Safety GmbH - Germany in compliance with the new European Standard for aerosol therapy units, Standard EN 13544-1, ANNEX CC. Further details are available on request.

### Nasal shower Rhino Clear®

Medication maximum capacity: 10 ml

<sup>(1)</sup> <b>10 ml dispensing time</b>	2' 10"
<sup>(2)</sup> <b>MMD</b>	48 $\mu\text{m}$
<sup>(2)</sup> <b>% &lt;10 <math>\mu\text{m}</math></b>	3,2

<sup>(1)</sup> Tests carried out with 0,9% NaCl physiological solution according to Flaem internal procedure  
<sup>(2)</sup> Data gathered with Laser Malvern Mastersizer system (in a certified laboratory)

#### DEVICE DISPOSAL

In conformity with Directive 2012/19/CE, the symbol shown on the device to be disposed of indicates that it is considered as waste and is therefore subject to "sorted waste collection". The user must therefore deliver (or have delivered) the above waste to an appropriate waste collection centre provided by the local authorities, or deliver it to the dealer when purchasing a new appliance of the same type. Pre-sorted waste collection and the subsequent treatment, recovery and disposal operations favor the production of appliances made of recycled materials and limit the negative effects to the environment and public health due to incorrect waste management. The unlawful disposal of the product by the user could result in administrative fines as provided by the laws transposing Directive 2012/19/CE of the European member state or of the country in which the product is disposed of.

#### ELECTROMAGNETIC COMPATIBILITY

This device has been designed to satisfy requirements currently required for electromagnetic compatibility (EN 60 601-1-2:2007). Electrical medical devices require special care. During installation and use with respect to EMC requirements, it therefore required that they be installed and/or used according to the manufacturer's specification. Potential risk of electromagnetic interference with other devices, in particular with other devices for analysis and treatment. Radio and mobile telecommunications devices or portable RF (mobile phones or wireless connections) may interfere with the operation of electrical medical devices. For further information visit [www.flaemnuova.it](http://www.flaemnuova.it). Flaem reserves the right to make technical and functional changes to the product without notice.

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Rivenditore (timbro e firma)

The warranty terms here specified are valid only in Italy for Italian residents. In all other countries, the warranty will be provided by the local dealer that sold you the unit, in accordance with the applicable laws.

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© 2013 FLAEM NUOVA®  
All rights reserved  
cod. 15508A0-0 Rev. 01/2015  
(TÜV 4)

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FLAEM NUOVA S.p.A.