SINCE 1966
MADE IN ITALY

Universal

Plus

by **FLAEM**®

ITALIANO
Pag. 3

ENGLISH
pg. 10

FRANCAIS
Pag. 17

NEREDLANDS
Pag. 24

DEUTSCH
Pag. 31

ESPAÑOL
PÁG. 38

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MOD. UNIVERSAL PLUS

Ultrasound aerosoltherapy unit

Congratulations on your purchase and thank you for choosing our unit.

*Our goal is to fully satisfy customers by offering them cutting-edge systems for the treatment of respiratory tract ailments. Carefully read these instructions and keep them for future reference. Use the unit only as described in this instruction manual. **This is a medical device for home use and may only be used with medication prescribed or recommended by a doctor.** Visit our Internet site www.flaem.it to view the whole range of Flaem products.*

THE UNIT STANDARD FITTINGS INCLUDE:

- | | |
|--|---|
| <p>(A) Nebulizer
A1) ON/OFF button
A2) ON/OFF operation LED
A3) Upper part
A4) Flow adjustment knob</p> | <p>(E) Child mask
(F) Extension tube
(G) Power supply unit
(H) Carry bag
(I) Table support
(L) Spare parts kit
(nebulizer upper part O ring and drug anti-dispersion valve)</p> |
| <p>(B) Mouthpiece
(C) Non-invasive nasal prong
(D) Adult mask</p> | |

IMPORTANT CAUTIONARY NOTES

- Before using the unit for the first time, and periodically throughout its life, check the integrity of its structure and the power cable, making sure that there is no damage; if it is damaged, do not insert the plug and take the product immediately to an authorized FLAEM technical assistance centre or to your local dealer. Take the appliance immediately to an authorised service centre or your dealer.
- This unit is especially suitable for administering drugs in solution form and can also work with drugs in suspension form. However, in the case of some glycol containing drugs the instructions for use advise against using ultrasound devices, since these drugs cannot be nebulised with them.
- Interactions: the materials used that come into contact with drugs have been tested with a wide range of medications. However, due to the variety and continuous evolution of medication, the possibility of a chemical interaction cannot be excluded. Once opened, the medication should be used as soon as possible, and prolonged contact between drug and nebuliser chamber should be avoided.
- The appliance is not suitable for essential oils or balsamic substances.
- Do not run the appliance if it does not contain any solution.
- If the upper part is not correctly fitted in, the appliance does not work; If the upper part is removed during the operation, the appliance switches off automatically.
To start it again place the upper part in its housing and press button (A1).
- Use the unit with a dose of at least 4 ml to be nebulised; if necessary, dilute the drug with physiological solution.
- It is normal for there to be a non-nebulised drug residual of around 0.5 – 0.8 ml at the end of treatment; therefore, when nebulisation becomes discontinuous, stop the treatment as it is pointless to insist on nebulising the remainder, thus causing the unit to overheat.
- Always keep the unit in an upright position to prevent the medication from spilling out.
- Do not use the tubes and the cords of the device for purposes other than those specified as they may cause strangulation hazard; utmost attention must be paid to children and persons with special needs; usually, such persons are not capable of assessing the hazard correctly.
- The average service life of the device is about 600 hours.
- Children and disabled persons should always use the unit under strict supervision of an adult who has read this manual.
- Several parts of the device are small enough to be swallowed by children; therefore keep the device out of reach of children.
- The unit should not be used in the presence of anaesthetic mixture inflammable with air, oxygen or nitrous

oxide.

- Keep the power unit away from any liquid substance.
- Do not expose the device to high temperatures.
- Do not put the device and power adaptor near heat sources, direct sunlight, or in hot environments.
- The power unit provided has been specifically designed for use with the unit. Never use the power unit with other appliances or for uses other than those indicated in this manual and never use with other power units.
- Never use adapters for mains voltages different from those shown on the power supply unit (plate details).
- Always keep the power supply and power cord away from heated surfaces.
- Do not handle the unit with wet hands. Do not use the device in humid environments (e.g. while taking a bath or shower). Do not immerse the device in water; should this accidentally occur, unplug the device immediately. Do not remove or touch the device immersed in water before pulling out the plug. Take the device immediately to an authorized service center or to your local dealer.
- Repairs must be performed by authorized personnel only. Any unauthorized repairs will void the warranty and may pose a safety hazard for the user.

WARNING: Do not modify this equipment without authorization of the manufacturer.

- The Manufacturer, the Retailer and the Importer are responsible for the safety, reliability and performance of the device only if: a) the device is used in accordance with the operating instructions b) the electrical system of the premises where the device is used is in compliance with current laws.
- The manufacturer should be contacted for reporting problems and/or unexpected events related to device operation.

OPERATING INSTRUCTIONS

Before each use, the nebuliser and accessories should be cleaned according to the instructions in the section “SANITIZATION AND DISINFECTION”. A personal use of the nebuliser and its accessories is recommended in order to avoid any risks of infection.

1. Remove the upper part (A3) from the unit body (fig. 1)


fig. 1

2. Insert the medication dose (max 6 ml) in the relevant tray (if a syringe is used, remove the needle) (fig. 2); The medication is best diluted to at least 4 ml using saline solution to optimise nebulisation.




fig. 2

3. Reinsert the unit upper part and connect the nebulizer to the power mains using the special power unit as indicated in figure 3.


4. Sit comfortably, holding the nebulizer in your hands, place the mouth-piece (B) over your mouth or use the nasal prong (C) or mask (D) or (E) (fig. 4).


While using the nasal prong (C) we recommend breathing in strongly to prevent the formation of condensation which would prevent correct nebulisation.

5. To give the dose of medicine to the patient while he/she is lying down (fig. 5), use the extension tube (F), one end of which must be connected to the mask (D) or (E) while the other must be connected to the appliance suitably housed in the special table stand (I).

6. Press the ON/OFF button for a moment (fig. 6); the nebulizer led (A2) will show a steady green light and the nebulised medication will start to come out.

7. Regulate the supply flow by means of the special knob (A4); turn it clockwise (position +) to treat the upper airways, or anticlockwise (position -) for the lower airways (fig. 7).
To increase the effect of the therapy, after inspiration, hold your breath for a second so the drops of inhaled aerosol are able to deposit along the airways. Then expire slowly.


fig. 7

8. To stop nebulisation press button (A1).

9. When the medication is about to run out, nebulisation becomes intermittent. Switch off the application by pressing the button (A1). Despite the unit having a sensor indicating the presence of medication, it is always best not to use the appliance dry (without medication).

AUTOMATIC SWITCH-OFF: The appliance switches off automatically after about 10 minutes; nonetheless, if you wish to continue with the therapy, simply press the button again and the appliance will start operating.

SANITIZATION AND DISINFECTION

- Switch off the device before all SANITIZATION AND DISINFECTION procedures and unplug the power cable from the socket.
- Separate the unit (A) from the power, disconnecting the power cable (G)
- Remove the top part (A3) from the body (Fig. 1) and the connected accessory. .

SANITIZATION

Before and after each use, sanitise the components (A-B-C-D-E-F) by choosing one of the methods described below.

(Method A): Wash only the unit and its accessories under warm drinking water and using a sponge soaked in a mild washing-up liquid (non-abrasive).

(Method B): Position only the unit and its accessories inside a plastic, porcelain or glass container (not metal); pour a 50% water and 50% white vinegar solution directly on top of these parts. Finally, rinse thoroughly with drinking water;

DISINFECTION

The disinfectable components are (A-B-C-D-E-F)

The disinfection procedure described in this paragraph is to be performed before using the accessories and is effective only if all the prescriptions herein are met and all the components involved are previously sanitized.

The disinfectant to be used must be of electrolytic oxidizing chlorine type, (active ingredient: sodium hypochlorite), specific for disinfection and is available in all pharmacies.

Procedure:

- Fill a container, suitable to contain all the individual components to be disinfected, with a solution of drinking water and disinfectant, while respecting the proportions indicated on the packaging of the disinfectant itself.
- Completely immerse each individual component in the solution, taking care to avoid the formation of air bubbles in contact with the components. Leave the components immersed for the period of time indicated on the packaging of the disinfectant, and associated with the concentration chosen to prepare the solution.
- Remove the components now disinfected and rinse thoroughly with lukewarm drinking water.

- After disinfection shake the accessories vigorously and place them on a paper towel, or alternatively, dry them with a jet of hot air (e.g. using a hair dryer).
- Dispose of the solution as recommended by the disinfectant manufacturer.

After the CLEANING AND / OR DISINFECTION operations dry the area of the transducer, without exerting unnecessary pressure on the transducer itself, dry the accessories and the outside of the appliance with a paper towel or even better with a jet of hot air (e.g. hair dryer). Finally, check that the “power plug” area of the device is completely dry.

SPARE PARTS

Description	Code
- Extension tube	ACO104T
- Mouthpiece and Nasal prong	ACO252P
- PVC adult mask	ACO90
- PVC child mask	ACO242

SYMBOLS

	Class II unit		Power supply unit for home use
	Important: check the operating instructions		Safety transformer with cut-out
	Type BF applied part		Thermofuse
	Direct current		Approvals
	Hazard: electrocution. Consequence: Death. Do not use the unit while bathing or showering.		Plug for very low safety voltage
IP21	Enclosure protection rating: IP21. (Protected against solid bodies over 12 mm. Protected against access with a finger; Protected against vertically falling water drops.)		Compliant with Directives
	CE Medical Marking ref. Dir. 93/42 EEC and subsequent amendments Complies with: EN 60601-1-11		Important: check the operating instructions
	Manufacturer		Alternate current
	In compliance with European standard EN 10993-1 Biological evaluation of medical devices and 93/42/EEC Medical Devices Directive”		
SN	Serial number of the device		

TROUBLESHOOTING

PROBLEM	REMEDY
1. While pressing the handset button the green led is not lit.	- Make sure that the handset and the power supply unit are properly plugged in.
2. In the case of poor nebulisation	- Make sure that the drug is not too cold – drug temperature should be close to room temperature.
3. If the unit does not nebulise	<ul style="list-style-type: none">- Make sure that the upper part is correctly fitted in.- Check that the nebulizer contains medication.- Make sure that quantity of medication is between 4 and 6 ml max.- Check that the medication is diluted sufficiently with saline solution 0.9%.- Make sure the adjustment knob and the disc underneath are perfectly clean.- Make sure the medication is not glycol based.
4. If the appliance interrupts nebulisation and the green LED flashes, this means there is no more medication	- Put some more medication into the tray or proceed to clean and disinfect.
5. Failure to clean the nebulizer from drug deposits may endanger proper operation of the unit.	- Carefully follow the instructions for cleaning and disinfection.

Should the unit again fail to operate after following the above instructions, please contact your local dealer or authorised after-sales centre.

NOTICE (ref. A Directive 2007/47/EC)

Flaem Nuova SpA has been manufacturing top-quality electro-medical appliances for the past 40 years or more, using carefully selected and tested materials in compliance with the requirements of reference directives and standards, with patient/user safety as number one priority. Because of this, we are able to state that the use of the Flaem device and original accessories purchased by you does not involve any risk whatsoever for the user's health.

The specific tests conducted at accredited laboratories (as per UNI EN ISO 10993-17: 2004 rule on Biological Assessment Part 17 Definition of allowable limits for leachable substances) actually demonstrate that the BPA used is released in amounts much lower than the maximum allowed for the safety threshold.

TECHNICAL SPECIFICATIONS

Mod. UNIVERSAL PLUS

NEBULIZER

Power:	15V 
Absorption:	850 mA
Frequency:	2.4 MHz
Operation:	30 min ON/ 30 min OFF
Automatic switch-off :	after 10 minutes
Max tray capacity :	6 ml.
Max. operating speed ⁽¹⁾ :	0.7 ml/min approx
MMAD ⁽²⁾ :	1.9 µm

⁽¹⁾Tests performed with phys. sol. 0.9% NaCl according to Flaem internal procedure.

⁽²⁾Certified with TSI Aerosizer 3220 computerized laser system according to Flaem internal procedure.

Noise (at 10 cm):	< 50 dB(A)
Dimensions:	9x5x12h cm
Weight:	200 g
In conformity with:	 0051
Dimensions of bag:	20x13x13h cm
Weight of bag:	1 Kg
(complete with appliance and accessories)	Temperature: min 10°C; max 40°C
Operating conditions:	Air humidity: min 10%; max 95% RH
	Temperature: -25°C; max 70°C
	Air humidity: min 10%; max 95% RH
	min. 690 hPa; max. 1060 hPa P
Storage conditions:	
Operating/storage atmospheric pressure:	

POWER SUPPLY UNIT

Use only power supply unit code 10529

Primary:	230V~ 50Hz
Secondary:	15V  850mA
Dimensions:	10x7x9 cm
Weight:	510 g

DEVICE DISPOSAL

In conformity with Directive 2012/19/EC, the symbol shown on the appliance to be disposed of indicates that this is considered as waste and must therefore undergo "sorted waste collection". The user must therefore take (or have taken) the above waste to a pre-sorted waste collection centre set up by the local authorities, or else give it back to the dealer when purchasing a new appliance of the same type. Pre-sorted waste collection and the subsequent treatment, recovery and disposal operations favour the production of appliances made of recycled materials and limit the negative effects of any incorrect waste management on the environment and public health. The unlawful disposal of the product by the user could result in administrative fines as provided by the laws transposing Directive 2012/19/EC of the European member state in which the product is disposed of.

ELECTROMAGNETIC COMPATIBILITY

This device has been designed to satisfy requirements currently required for electromagnetic compatibility (EN 60 601-1-2:2007). Electrical medical devices require special care. During installation and use with respect to EMC requirements, it therefore required that they be installed and/or used according to the manufacturer's specification. Potential risk of electromagnetic interference with other devices, in particular with other devices for analysis and treatment. Radio and mobile telecommunications devices or portable RF (mobile phones or wireless connections) may interfere with the operation of electrical medical devices. For further information visit www.flaemnuova.it. Flaem reserves the right to make technical and functional changes to the product without notice.

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