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④ Measuring sensor

Important Safety Instructions

- The instrument may be used only for measuring body temperature!
- The minimum measurement time until the beep is heard must be maintained without exception!
- Ensure that children and disabled persons do not use the instrument unsupervised; some parts are small enough to be swallowed.
- Protect the instrument from impact and dropping!
- Avoid ambient temperatures above 60 °C. NEVER boil the instrument!
- Use only the commercial disinfectants listed in the section «Cleaning and Disinfecting» to clean the instrument.

Batteries and electronic devices must be disposed of in accordance with the locally applicable regulations, not with domestic waste. Illegal dumping of the product by the user entails the administrative sanctions stated in article 50 et seq of Legislative Decree n. 22/1997.

Important: check the operating instructions

Type BF applied part

Turning on the Thermometer

To turn on the thermometer, press the ON/OFF button ① a short beep signals «thermometer ON». A display test is performed. All segments should be displayed. Then at an ambient temperature of less than 32 °C, an «L» and a flashing «°C» appear at the display field ②. The thermometer is now ready for use.

Function Test

Correct functioning of the thermometer is tested automatically each time it is turned on. If a malfunction is detected (measurement inaccuracy), this is indicated by «ERR» on the display, and a measurement becomes impossible. In this case, the thermometer must be replaced.

Using the Thermometer

Choose the preferred measuring method. When taking a measurement, the current temperature is continuously displayed and the «°C» symbol flashes. If the beep is heard 10 times and the «°C» is no longer flashing, this means that the measured increase in temperature is less than 0.1 °C in 16 seconds and that the thermometer can now be read. To prolong the battery life, turn off the thermometer by briefly pressing the ON/OFF button ①. Otherwise the thermometer will automatically turn off after about 10 minutes.

Storage of Measured Values

If the ON/OFF button ① is pressed for more than 3 seconds when turning on the thermometer, the automatically stored maximum temperature during the last measurement will be displayed. At the same time, a «M» for memory will appear on the display. About 2 seconds after the button is released, the temperature value disappears and the thermometer is ready for measurement.

Technical Specifications

Type:

Maximum thermometer

Measurement range:

32.0 °C to 43.9 °C

Temp. < 32.0 °C: display «L» for low (too low)

Temp. > 43.9 °C: display «H» for high (too high)

± 0.1 °C between 34 °C and 42 °C

Measurement accuracy:

10 - 40 °C; 15-95 % relative maximum humidity

Operating temperature:

-25 - +60 °C; 15-95 % relative maximum humidity

Storage temperature:

1.5V; LR 41

Battery:

EN 12470-3, clinical thermometers;

ASTM E1112; IEC 60601-1; IEC 60601-1-2 (EMC)

This device complies with the requirements of the Medical Device Directive 93/42/EEC.

Technical alterations reserved.

WARRANTY: The warranty terms here specified are valid only in Italy for Italian residents. In all other countries, the warranty will be provided by GB the local dealer that sold you the unit, in accordance with the applicable laws.

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www.flaem.it

Microlife AG, Esplanstrasse 139
9443 Widnau / Switzerland

Onbo Electronic (Shenzhen) Co., Ltd
No. 497, Ta Laneg Nan Road,
Ta Laneg Street, Baoan District,
Shenzhen China

CE 0044

Cod.16328 rev. 06/2014

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**FLAEM NUOVA S.p.A. - Via Colli Storici, 221