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# Lella the Ladybug

## COMPLETE SYSTEM FOR AEROSOL THERAPY

### Mod. P0406EM F700

Congratulations on your purchase and thank you for choosing our unit. Our goal is to fully satisfy consumers by offering them cutting-edge systems for treating ailments of the upper airways. Carefully read these instructions and keep them in a safe place for future reference.

**Use the unit only as described in this instruction manual.**

**This is a home medical device designed to administer nebulizer medications prescribed or recommended by your physician.** Please visit [www.flaem.it](http://www.flaem.it) to see the entire range of Flaem products.

### THIS DEVICE COMES EQUIPPED WITH:

#### A - Aerosol Therapy Device

- A1 - Switch
- A2 - Air outlet port
- A3 - Air filter
- A4 - Nebulizer port
- A5 - Socket for power cord
- A6 - Power cord

#### B - Connection tube (main unit / nebulizer)

#### C - Accessories

- C1 - RF6 Basic<sup>2</sup> nebuliser
  - C1.1 - Lower part
  - C1.2 - Complete nozzle
  - C1.3 - Upper part
  - C1.4 - Cap
- C2 - Mouthpiece
- C3 - Adult mask (PVC)
- C4 - Child mask (PVC)

#### C5 - RHINO CLEAR® Nasal shower

- C5.1 - Nasal transparent adapter "A", suitable for children up to 3 years of age
  - Nasal white adapter "B", suitable for children 3 to 10 years old
  - Nasal blue adapter "C", suitable for children 10 years old and up

#### C5.2 - Separator

#### C5.3 - Atomizer

#### C5.4 - Shower body with ergonomic handle

#### C5.4.1 - 10 ml nebulization chamber

#### C5.4.2 - Collection chambre

#### C5.5 - Dispensing button

#### C5.5.1 - Button

#### C5.5.2 - Rubber cap

#### C6 - Nasal shower-tube connect fitting

#### E - Carrying bag

#### F - "Lella the Ladybug" fairy tale

### IMPORTANT SAFEGUARDS

- Before using the device for the first time and periodically during the product's life, check the integrity of its structure and of the power cable, making sure that there is no damage; if it is damaged, do not insert the plug and take the product immediately to an authorized FLAEM technical assistance centre or to your local dealer.
- The average lifespan of accessories is 1 year; however, we recommend you replace the nebulizer every 6 months in intensive use (or earlier, if the nebulizer is blocked) to always ensure maximum therapeutic effect.
- Do not allow to be used as a toy. Close attention is necessary when used by or near children.
- Children and disabled persons should always use the unit under strict supervision of an adult who has read this manual.
- Several parts of the device are small enough to be swallowed by children; therefore keep the device out of reach of children.
- Do not use the tubes and the cords of the device for purposes other than those specified as they may cause strangulation hazard; utmost attention must be paid to children and persons with special needs; usually, such persons are not capable of assessing the hazard correctly.
- The unit should not be used in the presence of anaesthetic mixture inflammable with air, oxygen or nitrous oxide.
- Keep the cord away from heated surfaces.
- Keep the power supply cord away from animals (for example rodents) as they may damage the insulation of the power supply cord.
- Do not handle the unit with wet hands. Do not use the device in humid environments (e.g. while taking a bath or shower). Do not immerse the device in water; should this accidentally occur, unplug the device immediately. Do not remove or touch the device immersed in water before pulling out the plug. Take the device immediately to an authorized service center or to your local dealer.
- Only use the device in dust-free environments; otherwise, therapy may be compromised.

- The casing of this unit is not protected against the penetration of liquids. Do not wash the device under running water or submerge in water. Keep away from splashing water or other liquids.
- Do not expose the device to high temperatures.
- Do not put the device near heat sources, direct sunlight, or in hot environments.
- Do not obstruct or insert any objects in the filter or in its housing.
- Do not obstruct the air vents placed on both sides of the unit.
- During use, always place the unit on a firm surface clear of obstacles.
- Check that there is no material obstructing the air vents before each use.
- Do not place any objects inside the air vents.
- The average duration provided for compressor ranges is: F400: 400 hours, F700: 700 hours, F1000: 1000 hours, F2000: 2000 hours.
- Repairs must be performed by authorized personnel only. Any unauthorized repairs will void the warranty and may pose a safety hazard for the user.  
WARNING: Do not modify this equipment without authorization of the manufacturer
- The Manufacturer, the Retailer and the Importer are responsible for the safety, reliability and performance of the device only if: a) the device is used in accordance with the operating instructions b) the electrical system of the premises where the device is used is in compliance with current laws.
- Interactions: The materials used that come into contact with medications have been tested with a wide range of medications. However, due to the variety and continuous evolution of medication, the possibility of a chemical interaction cannot be excluded. We recommend that once the medication has been opened you should use it as soon as possible and avoid prolonged contact of the medication with the nebulizer.
- The manufacturer should be contacted for reporting problems and/or unexpected events related to device operation.

## INSTRUCTIONS FOR USE

### RF6 Basic<sup>2</sup> NEBULISER

**Before each use, the nebuliser and accessories should be cleaned according to the instructions in the section "CLEANING, SANITIZATION, DISINFECTION". A personal use of the nebuliser and its accessories is recommended in order to avoid any risks of infection.**

During application, it is advisable to suitably protect yourself from any dripping.

1. Plug the power cord (A6) into the socket (A5) and plug the power cord into a power socket corresponding to the voltage of the unit. This must be positioned so that it is not difficult to make the disconnection from the mains.
2. Assemble the nozzle as indicated in the "Connection diagram" in point C1.2 Insert the complete nozzle in the bottom part (C1.1). Insert the Upper ventilation cover (C1.4) in the upper part (C1.3) as shown in the "Connection diagram" in point C1.
3. Pour the medication prescribed by the physician into the lower part (C1.1). Close the nebuliser by turning the upper part (C1.3) clockwise.
4. Connect accessories as indicated in the "Connection diagram" on the cover.
5. Sit comfortably holding the nebuliser in your hand, place the mouthpiece onto your mouth or alternatively use the nose piece or mask. Should you opt for the mask accessory, place it on your face as shown in the picture (with or without using the elasticated strap).

The cap (C1.4) provided with the nebuliser RF6 Basic<sup>2</sup> makes it possible to deliver the medication in 2 modes:

**"standard mode"**  
for optimum delivery and minimum waste of medication

**"high speed mode"** for a quicker treatment

6. Turn on the device using the switch (A1) and then breathe in and exhale deeply.
7. When the treatment is finished, shut off the device and unplug the power cord.


**WARNING:** If after the therapy session moisture accumulates in the tube (B), remove the tube from the nebulizer and dry it using the device's compressor fan; this operation will prevent mold from forming inside the tube.

## NASAL WASH RHINO CLEAR®

**Before each use, wash your hands thoroughly and sanitize RHINO CLEAR® as described in section “CLEANING SANITIZATION DISINFECTION”. The NASAL WASH SYSTEM and accessories are for personal use only to avoid potential risks of transmitting infection.**

During the application it is recommended to be adequately protected against dripping.

- 1) Insert the atomiser (C5.3) into the nebulisation chamber (C5.4.1).
- 2) Place the solution (max 10 ml) in the nebulizing chamber (C5.4.1).
- 3) Apply the separator (C5.2) verifying that the arrow on the shower body is aligned with the separator (Pic. B). Then, insert the correct adapter for your age, or the most appropriate to fit the size of your nostrils (C5.1A/C5.1B/C5.1C).
- 4) Connect the RHINO CLEAR® device through the tube (B) and the fitting (C6).
- 5) Turn on the device using the switch (A1).
- 6) Place the device close to the nose and gently insert the adapter nosepiece into one nostril.
- 7) During treatment it is preferable to hold the other nostril closed, where the shower is not applied, by applying light finger pressure.
- 8) Run the shower holding down the manual dispensing button (C5.5) during the inhalation phase and release it before exhaling; at this stage, to help expel the mucus, blow gently through the nostril treated directly into the RHINO CLEAR®.
- 9) Repeat from Step 6 to Step 8 for the other nostril, trying to evenly distribute the solution between the two nostrils.
- 10) When used on children or other dependent people, or in any case when there are problems in coordinating the various phases of paragraph 8), we recommend to run the application alternating the nostrils about 20 seconds each, with pauses required to expel the mucus.
- 11) The operations described in the preceding paragraphs is useful for nasal cleaning using physiological solutions. The treatment with drugs is advisable after a first cleaning with physiological solution in order to achieve a more effective distribution of the drug on the nasal mucosa.
- 12) Once the application is complete, switch off the device and clean it, following the instructions provided in the “CLEANING SANITISATION DISINFECTION” paragraph.

## CLEANING SANITISATION DISINFECTION

*Switch the device off and disconnect it from the power supply before proceeding with any cleaning operation.*

### CLEANING THE DEVICE AND THE TUBE'S OUTER SURFACE

Only use a cloth with antibacterial detergent (non-abrasive and with no solvents of any kind)

### ACCESSORIES

We recommend personal use of accessories to avoid any risk of spreading the infection.

**Disassemble the accessories as described below:**

**C1** - Open the nebulizer by turning the upper part (C1.3) anticlockwise, remove the nozzle from the upper part (C1.3) and then disassemble it as shown in the “Assembly diagram” in section C1.2.

**C5** - Disassemble it as shown in the “Assembly diagram” in section C5.

## SANITISATION AND DISINFECTION

### SANITISATION

Before and after each use, sanitise the nebuliser cup and the accessories, choosing one of the methods described below.

**method A:** Sanitise accessories **C1.1-C1.2-C1.3-C1.4-C2-C3-C4-C5.1A-C5.1B-C5.1C-C5.2-C5.3-C5.4-C5.5.1-C5.5.2-C6** under potable hot water (approximately 40°C) with a gentle, non abrasive dish detergent.

**method B:** Sanitise accessories **C1.1-C1.2-C1.3-C1.4-C2-C3-C4-C5.1A-C5.1B-C5.1C-C5.2-C5.3-C5.4-C5.5.1-C5.5.2-C6** in the dishwasher with a hot cycle.

**method C:** Sanitise accessories **C1.1-C1.2-C1.3-C1.4-C2-C3-C4-C5.1A-C5.1B-C5.1C-C5.2-C5.3-C5.4-C5.5.1-C5.5.2-C6** by immersing them in a solution of 50% water and 50% white vinegar, then rinse thoroughly under potable hot water (approximately 40 °C).

If you want to also perform the cleaning for DISINFECTION, jump to the DISINFECTION paragraph.

After having sanitised the accessories, shake them vigorously and place them on a paper towel. Alternatively, dry them with a jet of hot air (for example, a hair dryer).

### DISINFECTION

After sanitising the nebuliser cup and the accessories, disinfect them choosing one of the methods described below.

**method A:** Accessories **C1.1-C1.2-C1.3-C1.4-C2-C3-C4-C5.1A-C5.1B-C5.1C-C5.2-C5.3-C5.4-C5.5.1-C5.5.2-C6** can be disinfected.

The disinfectant must be an electrolytic chloroxidizer (active principle: sodium hypochlorite) specific for disinfecting, which is available in any pharmacy.

**Implementation:**

- Fill a container big enough to hold all of the parts to disinfect with a solution of potable water and disinfectant, according to the proportions indicated on the packaging of the disinfectant.
- Completely immerse each part in the solution, taking care to avoid the formation of air bubbles on the parts. Leave the parts immersed for the amount of time indicated on the packaging of the disinfectant associated with the concentration chosen for the solution.
- Remove the disinfected parts and rinse abundantly with warm potable water.
- Dispose of the solution following the instructions provided by the disinfectant manufacturer.

**method B:** Disinfect the accessories **C1.1-C1.2-C1.3-C1.4-C2-C5.1A-C5.1B-C5.1C-C5.2-C5.3-C5.4-C5.5.1-C5.5.2-C6** by boiling them in water for 10 minutes; use demineralised or distilled water to prevent calcium deposits.

**method C:** Disinfect the accessories **1C1.1-C1.2-C1.3-C1.4-C2-C5.1A-C5.1B-C5.1C-C5.2-C5.3-C5.4-C5.5.1-C5.5.2-C6** with a hot steam steriliser for baby-bottle (not the microwave type). Perform the process faithfully following the instructions of the steriliser. To ensure that the disinfection is effective, choose a steriliser with an operating cycle of at least 6 minutes.

After having disinfected the accessories, shake them vigorously and place them on a paper towel. Alternatively, dry them with a jet of hot air (for example, from a hair dryer).


## AIR FILTERING

The unit is equipped with an air filter (A3) which should be replaced when dirty or when its colour changes. Do not wash or re-use the used filter. The filter must be regularly replaced to help ensure efficient compressor performance. The filter should be regularly checked. For spare filters contact your dealer or an authorised service centre.

### To replace the filter:

Lift the lid and remove the filter using a screwdriver (see fig.). The filter is designed to remain fixed in its housing.

**Use only Flaem original accessories.**

## SYMBOLS

TÜV approval ref. EN 60601-1 3rd Ed. - EN 60601-1-11

EC Marking medical ref. Dir 93/42 EEC and subsequent updates

Class II device

Attention: check the instructions for use

"OFF" for part of Equipment

"ON" for part of equipment

In compliance with: European Standard EN 10993-1 "Biological Evaluation of medical devices" and European Directive 93/42/EEC "Medical Devices."

Phthalate-free. In compliance with: Reg. (EC) no. 1907/2006

Device serial number

Manufacturer

Type BF applied part

Alternating current

Risk: electrocution.  
Consequence: Death.  
Using the device while taking a bath or shower is prohibited

Keep dry

## TECHNICAL DATA

### Mod. P0406EM F700

|          |                  |             |             |
|----------|------------------|-------------|-------------|
| Voltage: | 230V~ 50Hz 140VA | 115V ~ 50Hz | 220V ~ 50Hz |
|----------|------------------|-------------|-------------|

Safety certifications:

|                        |                  |
|------------------------|------------------|
| Max Pressure:          | 2.6 ± 0.4 bar    |
| Compressor air output: | 10 l/min approx  |
| Dimensions (L)x(P)x(H) | 21x16x13 cm      |
| Weight                 | 1.600 Kg         |
| Noise level (at 1 m):  | 55 dB (A) approx |
| Operation:             | Continuous use   |

### APPLIED PARTS

Type BF applied parts include: patient accessories (C2,C3,C4,C5.1A-C5.1B-C5.1C)

### RF6 Basic<sup>2</sup> Nebulizer

|                                |         |
|--------------------------------|---------|
| Medication minimum capacity    | 2 ml    |
| Medication maximum capacity:   | 8 ml    |
| Operating pressure (with neb.) | 0.8 bar |

|   | Mod. High speed     | Mod. standard       |
|---|---------------------|---------------------|
| <sup>(1)</sup> <b>Delivery:</b>                             | 0.42 ml/min approx. | 0.32 ml/min approx. |
| <sup>(2)</sup> <b>MMAD:</b>                                 | 3.6 µm              | 3.6 µm              |
| <sup>(2)</sup> <b>Breathable fractions &lt; 5 µm (FPF):</b> | 67%                 | 66%                 |

<sup>(1)</sup> data detected according to Flaem I29-P07.5 internal procedure  
<sup>(2)</sup> In vitro characterization certified by TÜV Rheinland Product GmbH - Germany in compliance with the new European Standard for aerosol therapy units, Standard EN 13544-1. Further details are available on request.

### NASAL SHOWER SPECIFICATIONS

#### Mod.: Rhino Clear®

|                                      |                    |
|--------------------------------------|--------------------|
| Capacity:                            | drug 10 ml maximum |
| 10 ml dispensing time <sup>(1)</sup> | 2'                 |
| MMD <sup>(2)</sup>                   | 48 µm              |
| % <10 µm <sup>(2)</sup>              | 3.2%               |

<sup>(1)</sup> Tests carried out with 0,9% NaCl physiological solution according to Flaem internal procedure

<sup>(2)</sup> Data gathered with Laser Malvern Mastersizer system (in a certified laboratory)

#### Operating conditions:

Temperature: min 10°C; max 40°C  
Air humidity: min 10%; max 95%  
Atmospheric pressure: min 69KPa; max 106KPa

#### Storage conditions:

Temperature: min -25°C; max 70°C  
Air humidity: min 10%; max 95%  
Atmospheric pressure: min 69KPa; max 106KPa

### DEVICE DISPOSAL

 In conformity with Directive 2012/19/EC, the symbol shown on the device to be disposed of indicates that it is considered as waste and is therefore subject to "sorted waste collection". The user must therefore deliver (or have delivered) the above waste to an appropriate waste collection centre provided by the local authorities, or deliver it to the dealer when purchasing a new appliance of the same type. Pre-sorted waste collection and the subsequent treatment, recovery and disposal operations favor the production of appliances made of recycled materials and limit the negative effects to the environment and public health due to incorrect waste management. The unlawful disposal of the product by the user could result in administrative fines as provided by the laws transposing Directive 2012/19/EC of the European member state or of the country in which the product is disposed of.

### ELECTROMAGNETIC COMPATIBILITY

This device has been designed to satisfy requirements currently required for electromagnetic compatibility (EN 60601-1-2:2007). Electrical medical devices require special care. During installation and use with respect to EMC requirements, it therefore required that they be installed and/or used according to the manufacturer's specification. Potential risk of electromagnetic interference with other devices, in particular with other devices for analysis and treatment. Radio and mobile telecommunications devices or portable RF (mobile phones or wireless connections) may interfere with the operation of electrical medical devices. For further information visit [www.flaemnuova.it](http://www.flaemnuova.it). Flaem reserves the right to make technical and functional changes to the product without notice.

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© Copyright by FLAEM NUOVA S.p.A., 2017 cod. 10972F0-0 Rev. 08/2017 (TUV 4)