Rhino Clear Mobile

MADE
IN ITALY

**PHthalates
& BPA FREE**

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Nasal wash

These user instructions are provided for the device model P0907EM.

INTENDED USE

Medical device for the washing and treatment of the nasal cavities using saline solutions and/or medication, the treatment and medication of which must be prescribed by a doctor.

INDICATIONS FOR USE

Treatment, via the endonasal route, of upper airway pathologies (rhinitis, rhinosinusitis, sinusitis, and or as indicated by the treating physician).

Physiological washing of the nasal cavities (with saline, isotonic, hypertonic solutions, or thermal water). In children it helps with nasal hygiene, in adults it aids physiological nasal washing and can be used before upper airway therapy to assist drug administration.

CONTRAINDICATIONS

- The medical device should NOT be used for patients who are unable to breathe on their own or who are unconscious.
- Do not use the device in anaesthetic or assisted ventilation circuits.

INTENDED USERS

The devices are intended for use by legally authorised medical personnel/health workers (doctors, nurses, therapists, etc.). The device can be used directly by the patient.

TARGET GROUP OF PATIENTS

Adults, children of all ages, infants. Before using the device, the user manual must be read carefully and an adult responsible for safety must be present if the device is used by infants, children of any age or persons with limited abilities (e.g. physical, mental, or sensory). It is up to the medical personnel to assess the patient's condition and capabilities in order to determine, when prescribing therapy, whether the patient is able to operate the device safely on their own or whether the therapy should be administered by a responsible person.

Please refer to medical personnel for the evaluation of the use of the device on particular types of patients such as pregnant women, lactating women, infants, incapacitated persons or persons with limited physical capabilities.

OPERATING ENVIRONMENT

This device can be used in healthcare facilities, such as hospitals, outpatient clinics, etc., or even at home.

WARNINGS CONCERNING POSSIBLE MALFUNCTIONS

- Should your appliance fail to perform, please contact the authorised service centre for clarification.
- The manufacturer should be contacted to report problems and/or unexpected events relating to operation and if necessary for clarification of use and/or maintenance/hygienic preparation.
- Please also refer to the case history of faults and their resolution.

WARNINGS

- Use the device for therapeutic purposes only. This medical device is not intended as a life-saving device. Any other use is considered improper and can be dangerous. The manufacturer is not liable for any improper use.
- Always consult your general practitioner for identification of treatment.
- Follow your doctor's or therapist's instructions regarding the possible use of types of medication, dosage and indications for treatment.
- Pay attention to the indications provided with the drug and avoid using the devices with substances and dilutions other than those recommended.
- If you should experience allergic reactions or other problems while using the device, stop using it immediately and consult your doctor.
- Keep this manual safe for further reference.
- If the packaging is damaged or opened, contact the distributor or service centre.
- Do not expose the device to particularly extreme temperatures.
- The time required to switch from storage to operating conditions is about 2 hours.

- It is forbidden to access the device opening in any way. Repairs may only be carried out by personnel authorised by the manufacturer. Unauthorised repairs invalidate the warranty and may pose a danger to the user.
- Interactions: the materials used in the device are biocompatible materials and comply with statutory regulations, however possible allergic reactions cannot be completely excluded.
- Use the medicine as soon as possible once it has been opened and avoid leaving it in the ampoule body; once therapy has ended, do not leave the medicine inside the ampoule body and proceed with hygienic preparation.
- The nasal douche should not be used as a toy. The utmost attention should be paid when used by children.
- Do not use the device while driving, or in any other situation where you should not be distracted in order not to endanger yourself, people or pets nearby and the objects around you.
- **Suffocation risk:**
 - Some components of the device are small enough to be swallowed by children, therefore, keep the device out of the reach of children.
- **Fire risk:**
 - THIS device is not suitable for use in the presence of an anaesthetic mixture that is flammable with air, or with oxygen or nitrous oxide.
 - Do not place the device near sources of heat, light or in excessively hot and humid premises. Store as described in the chapter on environmental conditions.
- **Risk of electrocution:**
 - Before first use, and periodically during the product life, check the integrity of the device structure and the power cable of the charging base to ensure that there are no damages; if damaged, do not insert the plug and immediately take the product to an authorized service center or to your trusted dealer.
 - Never use adapters for power voltages other than those specified on the nasal shower charging base (nameplate data).
 - Keep the power cable away from animals (e.g., rodents), otherwise such animals could damage the power cable insulation.
 - Always keep the power cable away from hot surfaces.
 - Do not handle the plug and power cable of the charging base with wet hands.
 - Do not use the nasal shower in humid environments (e.g., while taking a bath or shower).
 - Do not immerse the charging base in water; if this happens, immediately unplug it. Do not remove or touch the immersed device, unplug it first. Take it immediately to an authorized FLAEM service center or to your trusted dealer.
 - Do not wash the device under running water or by immersion.
- **Risk of ineffectiveness of therapy:**
 - Do not obstruct the suction valve at the bottom of the douche compressor body (5); never put anything in its cover.
 - In the case of substances that are too dense, dilution with suitable saline solution may be necessary, as prescribed by a doctor.
 - Observe the hygienic preparation of the spray chamber as described under Hygienic Preparation.
 - Only use original Flaem spare parts, no liability is accepted if non-original parts are used
- **Risk of infection:**
 - Personal use of the nasal douche is recommended to avoid any risk of infection.
 - Follow hygienic preparation before each use. Ensure that the douche and its parts are not stored in the vicinity of other components or devices for different therapies (e.g. infusion devices).
 - Once the therapy is over, do not leave the liquid used (drug or saline solution) inside the ampoule body.
- **Risk of injury:**
 - Do not place the medical device on a soft supporting surface such as a sofa, a bed or a tablecloth.
 - Do not use the device outside of its intended use, the air generated by the compressor is available to the douche compressor body, this can cause a potential hazard.

WARNINGS ON INTERFERENCE RISKS DURING USE IN DIAGNOSTIC INVESTIGATIONS

This device is designed to meet the current requirements for electromagnetic compatibility. As far as EMC requirements are concerned, electromedical devices require special care when being installed and used, and are therefore required to be installed and/or used in accordance with the manufacturer's specifications. Risk of potential electromagnetic interference with other devices. Mobile or portable RF radio and telecommunication devices (mobile phones or wireless connections) may interfere with the operation of electromedical devices. The device may be susceptible to electromagnetic interference in the presence of other devices used for specific diagnosis or treatment. For more information visit www.flaemnuova.it.

TROUBLESHOOTING

Inconveniences	Solutions
If, after placing the shower compressor body on the charging base, the green light (5b) on the button does not turn on.	Check that the plug is plugged in securely. Try plugging another appliance into the outlet to see if it works. If the power cord or plug is damaged, don't use the charging base. Contact your trusted retailer or authorized service center.
If the nasal shower nebulizes little or for a short time.	The batteries may be low. Charge the nasal shower for at least 24 hours. Make sure the atomizer is properly inserted. Make sure the nasal shower is clean. If it's not, clean it following the instructions in the HYGIENIC PREPARATION section.
If the nasal shower does not nebulize.	The batteries may be low. Charge the nasal shower for at least 24 hours. Make sure the solution is poured correctly into the nebulization chamber of the ampoule body. Make sure the amount of liquid doesn't exceed 15 ml. If you're using a drug for aerosol therapy, make sure it's diluted sufficiently with saline solution (follow your doctor's advice). Make sure the nasal shower is clean. If it's not, clean it following the instructions in the HYGIENIC PREPARATION section.
If during operation the nasal shower interrupts nebulization (turns off).	The thermal protection built into the nasal shower has been activated due to a probable overload. Wait a few minutes for the unit to cool down. Then, press the button to turn the shower back on. It's best to wait at least 10 minutes before using it again. Make sure the nasal shower is clean. If it's not, clean it following the instructions in the HYGIENIC PREPARATION section.
If, after checking the above conditions, the device still doesn't work properly, we recommend that you contact your trusted retailer or the nearest FLAEM authorized service center. You can find a list of all service centers at the website http://www.flaemnuova.it/it/info/assistenza .	

DISPOSAL

Compressor body (5, 6)

 In accordance with Directive 2012/19/EC, the symbol on the equipment indicates that the equipment to be disposed of is considered as waste and must therefore be subject to 'separate collection'. Therefore, the user must deliver this waste (or have it delivered) to the separate collection centres set up by the local authorities, or hand it over to the retailer when purchasing a new appliance of an equivalent type. Separate waste collection and subsequent treatment, recovery and disposal operations promote the production of equipment from recycled materials and limit the negative environmental and health effects caused by improper waste management. Unauthorised disposal of the product by the user entails the application of the administrative sanctions provided for in the laws transposing Directive 2012/19/EC of the member state or country where the product is disposed of.

Components (1a, 1b, 1c, 2, 3, 4,)

to be disposed of as general waste after a sanitisation cycle.

Packaging

Product box

Heat-shrink film accessories.

Bag packaging pouch



NOTIFICATION OF SERIOUS EVENTS

Serious events occurring in connection with this product must be reported immediately to the manufacturer or the competent authority.

An event is considered serious if it causes or may cause, directly or indirectly, death or an unforeseen serious deterioration in a person's state of health.

COUNTRY	AUTHORITY
Ireland	Health Products Regulatory Authority Kevin O'Malley House, Earlsfort Centre, Earlsfort Terrace, IE - Dublin 2 E-mail: devicesafety@hpra.ie
Malta	Medical Devices Unit Medicines Authority Sir Temi Żammit Buildings, Malta Life Sciences Park San Ġwann SĠN 3000 E-mail : devices.medicinesauthority@gov.mt

SYMBOLS ON DEVICE OR PACKAGING

	Marcatura CE medicale rif. regolamento 2017/745 UE e successivi aggiornamenti		Serial number of the device
	Quality mark		Manufacturer
	Before use: Caution check instructions for use		Type BF applied part
	Phthalate- and bisphenol-free		Attention
	Model number		Protected against solid foreign bodies of \varnothing 1 mm or larger and against water jets
	Temperature limits		Medical device
	Atmospheric pressure limits		Unique device identifier
	Moisture limits		Production date
	Operating instructions		Alternating current
	Class II appliance		Direct current

INFORMATION ON RESTRICTIONS OR INCOMPATIBILITIES WITH CERTAIN SUBSTANCES.

- Interactions: The materials used in the device are biocompatible materials and comply with statutory regulations, however possible allergic reactions cannot be completely excluded.
- Use the medicine as soon as possible once it has been opened and avoid leaving it in the ampoule; once therapy has ended, do not leave the medicine inside the ampoule and proceed with hygienic preparation.

TECHNICAL SPECIFICATIONS

Model: P0907EM

Charging base voltage:	220-230V~ 50-60Hz / 18VA
Max. spray chamber capacity:	15 ml
Nebulisation time of 10 ml ⁽²⁾ :	about 2 minutes (1 ml in about 12 seconds)
Overmeasurement Cumulative Percentage in Mass: > 10 microns:	between 95 and 96 % ⁽¹⁾
(Non-breathable fraction useful for upper respiratory tract treatment)	

⁽¹⁾Values measured with a Malvern Spraytec computerised laser diffractometer. ⁽²⁾ Average values found with 0.9% NaCl saline solution.

Nasal douche dimensions:	7x6x15.5h cm
Nasal douche weight:	210 g
Operation:	Continuous
Charging base dimensions:	15x9x8h cm
Total weight (charging base complete with device and accessories):	480 g
Internal nasal douche batteries:	Ni-MH rechargeable 1.2V 620/650mAh
Type BF applied parts:	The 3 nasal adapters (1a, 1b, 1c)

ENVIRONMENTAL CONDITIONS

Operating conditions:

Ambient temperature	Between +5°C and +40°C
Relative air humidity	Between 15% and 90%
Atmospheric pressure	Between 69KPa and 106 KPa

Storage and transport conditions:

Ambient temperature	Between -10°C and +45°C
Relative air humidity	Between 10% and 75%
Atmospheric pressure	Between 69KPa and 106 KPa

SERVICE LIFE OF THE DEVICE

- (5, 6) Douche compressor body: approx. 1500 applications
- (1a, 1b, 1c) nasal adapters, (2) separator, (3) atomiser and (4) ampoule body: one year.

THE EQUIPMENT INCLUDES:

The device includes:	Information on materials
1a - Nasal adapter with soft rubber apex, suitable for children up to 3 years of age 1b - Nasal adapter with soft rubber apex, suitable for 3 to 10 years of age 1c - Nasal adapter with soft rubber apex, suitable after 10 years of age	Polypropylene + Thermoplastic Elastomers
2 - Separator	Polypropylene
3 - Atomiser	Polypropylene
4 - Removable and washable ampoule body consisting of 2 separate chambers: 4a - 15 ml spray chamber 4b - Recovery chamber	Polypropylene
5 - Douche compressor body with ergonomic handle	Acrylonitrile butadiene styrene
5a - Dispensing button	Silicone
5b Charging LED 5c Suction valve with cover	
6 - Recharger base with nasal adaptor compartment 6a - Adaptor compartment with lid 6b - Power supply cord compartment 6c - Power cord	Acrylonitrile butadiene styrene
7 - Fun child's colour sticker 8 - Handy carry bag	
IMPORTANT NOTE: There is an identification label on the packaging, remove it and affix it in the spaces provided on page 2. When replacing components, carry out the same procedure.	

1. CONNECTION AND CHARGING

- 1.1 - Always make sure that the mains voltage matches the information on the rating plate of the recharger base.
- 1.2 - Place the nasal irrigator complete with recharger base on a stable surface.
- 1.3 - Connect the power supply cord to the power socket (FIG. A).
- 1.4 - When the irrigator compressor body is placed in the recharger base, the green indicator light (5b) on the button will indicate the recharging status. If, after using the nasal irrigator, it is replaced on the recharger base, the indicator light will come on again. You cannot use the nasal wash while it is charging.
- 1.5 - Charging for first use takes about 24 hours and assures the utmost autonomy (6/7 applications of 10 ml). At the end of charging, disconnect the base from the power mains by unplugging the device from the socket, the position should ensure disconnection from the power mains can be carried out conveniently.
- 1.6 - After long periods of inactivity, recharge the unit by standing it on the recharger base for 24 hours, at least once every three months, to avoid damage to the battery inside. The internal batteries have an average life of 1 year and can be recharged about 200 times (approximately 1500 applications).

IMPORTANT INFORMATION TO KEEP THE BATTERY EFFICIENT

- We recommend keeping the unit on the charger base when it is not in use
- After long periods of inactivity, recharge the unit by standing it on the charger base for 24 hours, at least once every three months, to avoid damage to the battery.

2. INSTRUCTIONS FOR USE

- Before each use, wash your hands thoroughly and clean the nasal douche and its components as indicated in paragraph "HYGIENIC PREPARATION".
- Care should be taken during application because the solution may drip from the nose.

2.1 - Before use, insert the batteries as indicated in paragraph 1.

2.2 - For correct use of the nasal douche, the symbol of the person (♁) (4c), located on the ampoule body (4), must always be directed towards the patient. The ampoule can be attached with the position of the person on the same side as the dispensing button (5a) allowing operation of the nasal douche with the thumb, or with the person on the opposite side of the button for operation with the index finger.

In order to obtain one of the above conditions, the moment you start screwing the ampoule body (4) clockwise onto the douche compressor body (5) (fig. C), the person symbol must be on the opposite side to the position you want to obtain (fig. C).

2.3 - Pour the solution (max. 15 ml) into the atomisation chamber (4a) of the ampoule body, as shown in (fig. E). Also check the presence of the atomiser (3).

2.4 - Apply the separator (2) by aligning the dot (2a) with the person symbol (♁) (4c) on the outside of the ampoule body (4), as shown in (fig. F). Then insert the nasal adapter most suitable for the size of your nostrils (1a/1b/1c).

2.5 - Insert the hole of the nasal adapter into a nostril (Fig. G).

2.6 - Operate the nasal douche by holding down the dispensing button (5a) during the inhalation phase and release it before exhaling; during this phase, to help expel mucus, blow gently through the nostril treated with the nasal douche, keeping the other nostril closed.

2.7 - Even during atomisation, it is preferable to keep the other nostril where the spray is not applied, closed with a finger.

2.8 - Repeat steps 2.5 to 2.7 for the other nostril as well, trying to distribute the solution equally between the two nostrils.

2.9 - In the case of children or persons who are not self-sufficient, and in any case unable to co-ordinate the various steps in 2.6, it is recommended that the application be carried out alternating about every 15 seconds between the two nostrils, with pauses as necessary to expel the mucus.

2.10 - When administering drugs by aerosol, it is advisable to carry out an initial wash with saline solution to achieve more effective distribution of the drug on the nasal mucous membranes.

2.11 - When the application is finished, switch off the machine and carry out the hygienic preparation following the instructions given in paragraph "HYGIENIC PREPARATION".

HYGIENIC PREPARATION

Careful hygienic preparation is extremely important for the duration and proper functioning of the device. After use, disconnect the nasal adapter (1a-1b-1c), the separator (2), the atomizer (3) and unscrew the ampoule body (4) and empty the liquid present in the recovery chamber (4b). Then proceed to clean the parts by performing the following operations:

Cleaning the douche compressor body (5) and Charging Base (6)

⚠ Before carrying out cleaning operations on the Charging Base (6), unplug the power cord from the mains. Use only a damp cloth with an antibacterial detergent (non-abrasive and free of solvents of any kind).

ATTENTION: The douche compressor body is not protected against the penetration of liquids; therefore do not wash it under running water or by immersion.

Sanitisation.

Sanitise the components by choosing one of the methods provided in the table and described below.

Method A: under warm drinking water (approx. 40°C) with mild dishwashing detergent (non-abrasive).

Method B: by immersion in a solution of 50% water and 50% white vinegar. Finally, rinse thoroughly with warm drinking water (approx. 40°C);

After sanitising the components, shake them vigorously and place them on a paper towel, or alternatively dry them with a hot air jet (e.g. hair dryer).

Disinfection.

Disinfect the components by choosing one of the methods provided in the table and described below.

Method A: The disinfection procedure described in this paragraph is to be carried out before using

the components and is only effective on the components to be treated if every point is followed and only if the components to be treated are sanitised beforehand. The disinfectant must be of the electrolytic chloride type (active ingredient: sodium hypochlorite), specific for disinfection and available in all pharmacies.

Execution:

- Fill a container of a suitable size to hold all the individual components to be disinfected with a solution of drinking water and disinfectant, observing the proportions indicated on the packaging of the disinfectant.
- Completely immerse each individual component in the solution, taking care to avoid the formation of air bubbles in contact with the components. Leave the components immersed for the period of time indicated on the disinfectant packaging, and associated with the concentration chosen for the preparation of the solution.
- Recover the disinfected components and rinse them thoroughly with lukewarm drinking water.
- Dispose of the solution according to the disinfectant manufacturer's instructions.

Method B: by boiling in water for 10 minutes; demineralised or distilled water is preferable to avoid lime deposits.

After disinfecting the components, shake them vigorously and place them on a paper towel, or alternatively dry them with a hot air jet (e.g. hair dryer).

Sterilisation.

Sterilisation for clinical or hospital use is the responsibility of the facility performing the sterilisation and must be carried out according to the instructions in this manual.


The sterilisation procedure described in this paragraph is only effective on treated components if it is observed in all points and only if the components being treated are sanitised beforehand, and is validated in accordance with ISO 17665-1.

Equipment: Steam steriliser with fractionated vacuum and overpressure compliant with EN 13060.

Execution: Package each individual component to be treated in sterile barrier system or packaging in accordance with EN 11607. Place the packaged components in the steam steriliser. Perform the sterilisation cycle in accordance with the instructions for use of the equipment by selecting a temperature of 134°C and a time of 4 minutes first.


Storage: Store sterilised components according to the instructions for use of the system or sterile barrier packaging chosen.

SUCTION VALVE. The nasal douche is equipped with a suction valve (5c) to be cleaned when dirty. To access the valve, lift the battery compartment cover and remove the valve by prying it off with a pen; then wash the valve with water and detergent. Before reassembling, make sure the valve is completely dry.

ADDITIONAL INFORMATION FOR THE HYGIENIC PREPARATION of the spray chamber (except for the compressor body)

In order to prevent deposits from clogging the atomiser and spray chamber and thereby affecting the function of the nasal douche, it is recommended that hygienic preparation be carried out after each use, following the instructions in the instruction manual supplied with the nasal douche. After thoroughly rinsing the parts under lukewarm tap water, shake them to remove any remaining water droplets. Then place the pieces on a paper towel or dry them with a jet of hot air. **ATTENTION:** If your nasal douche does not atomise correctly after hygienic preparation, do not use sharp tools (e.g. screwdrivers, scissors, needles, hairpins, etc.) to try to clean or unclog the holes in the atomiser or spray chamber, as this may damage these parts and impair the proper functioning of the unit. Instead, perform additional hygienic preparation and sanitisation cycles as described in the instruction manual and recommended above, until the nasal douche begins to function properly again.

Table of planned methods / components						
Components Method						
	1a	1b	1c	2	3	4
HYGIENIC PREPARATION IN THE HOME						
Sanitisation						
method A	✓	✓	✓	✓	✓	✓
method B	✓	✓	✓	✓	✓	✓
Disinfection						
method A	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES
method B	✓ MAX 40 TIMES	✓ MAX 40 TIMES	✓ MAX 40 TIMES	✓ MAX 40 TIMES	✓ MAX 40 TIMES	✓ MAX 40 TIMES
HYGIENIC PREPARATION IN A CLINIC OR HOSPITAL						
Sterilisation	✓ MAX 40 TIMES	✓ MAX 40 TIMES	✓ MAX 40 TIMES	✓ MAX 40 TIMES	✓ MAX 40 TIMES	✓ MAX 40 TIMES
✓ : planned \: not planned						

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FLAEM NUOVA S.p.A. non risponde di ulteriori estensioni del periodo di garanzia assicurate da terzi.