User Manual  
Fingertip Pulse Oximeter

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Xuzhou Yongkang Electronic Science Technology Co., Ltd.  
#47 Building C8, 40 Jingshan Road, Economic and Technological  
Development Zone, 221000 Xuzhou,  
PEOPLE'S REPUBLIC OF CHINA

2. In menu 2. When the ""\*signal is shown on the "Sounds Setup", press the button (>3s) and enter into the menu 2 ,Press the button (<1s) can select item, then press the button (>3s) to setup data.choose ""\*or"" to plus or minus values.

Note: If the matching refers to clamp oximeter is basic, not equipped with parameter setting function

**Detailed descriptions of product functions:**

1. Display Type: OLED display
- 2.SpO2: Measurement range: 70%~99%
- Accuracy: 80%~99%±2% (Including 90%); 70%~79%: ±3%; Below 70% no requirement;
- 3.PR: Measurement range: 30BPM~240BPM
- Accuracy: ±1BPM or ±1% (the larger one)
- 4.Parameters of LED probe

	wavelength	radiation power
RED	660±2nm	1.8mW
IR	905±10nm	2.0mW

This parameter can be especially useful to clinicians.

5. PI index: measurement scope: 0~20 (optional)
6. Power: two AAA 1.5V alkaline batteries
7. Automatic standby: the product shuts off by itself when no finger is in the product about 8 seconds
8. Automatic startup: every 5 s instrument will automatically detect the signal, after the hole with my finger, timely automatically boot;(optional)
9. Accelerometer function: finger movement, the screen display will change with the accelerometer changes (optional)
10. Dimension: 58mm×36mm×3mm

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### Ver.1.0

#### Fingertip Pulse Oximeter

This Fingertip Pulse Oximeter is a kind of innovated medical device with non-invasive and continuous features for artery SpO2 and PR detection. Being portable, it is able to measure SpO2 and PR values quickly and precisely.

#### General Description

Haemoglobin Saturation is the percentage between the capacity of Oxyhemoglobin (HbO2) that compounded with oxygen and that of all combinative haemoglobin (Hb) obin (HbO2) in blood. In other words, it is the saturation of oxyhemoglobin in blood. It is a very important physiological parameter for Respiratory and Circulation Systems. Many respiratory diseases could reduce haemoglobin saturation in human blood. Moreover, factors such as Automatic Organic Regulation Malfunction caused by anaesthesia, trauma resulted from major operation and some medical examination can also cause problems in oxygen supply, which might reduce human haemoglobin saturation. As a result, such symptoms as megrim, vomiting and asthenia might appear to patients. Hence, it is very important to know haemoglobin saturation of patient timely in clinical medical aspects.

The fingertip pulse oximeter features in small volume, low power consumption, convenient operation and portability. It is only necessary for patient to put one finger into fingertip photoelectric sensor for diagnosis, and the display screen will directly show measured value of hemoglobin saturation. It has been proved in clinical experiments that it possesses rather high precision and repeatability.

#### Measurement principle

The principle of the oximeter is as follows: An experience formula of data process is established by exerting Lambert Beer Law according to Spectrum Absorption Characteristics of Reductive hemoglobin(R Hb) and oxyhemoglobin (O2 Hb) in glow and near-infrared zones.

Operation principle of the instrument is to combine Photoelectric Oxyhemoglobin Inspection Technology with Capacity Pulse Scanning and Recording Technology, so that two lights with different wavelength (660nm glow and 940nm near infrared light) can be focused onto human nail through perspective clamp finger-type sensor. Then measured signal can be obtained by a photosensitive element. This information acquired through which will be shown on two groups of LEDs through process in electronic circuits and microprocessor.

11. Operation environment:  
temperature: 5 °C ~ 40 °C  
Humidity: 15% ~ 80%  
Atmospheric pressure: 86kPa ~ 106kPa

Transport, storage environment:  
Temperature: -10 °C ~ 40 °C  
Humidity: <95%  
Atmospheric pressure: 50kPa ~ 106kPa, non-corrosive gas and well-ventilated atmosphere  
12. Declaration: EMC of this product comply with IEC60601-1-2 standard.

#### Classification

1. Management Class for Medical Devices (equipment)
2. Anti-electric Shock Type: normally powered equipment
3. Anti-electric Shock Degree: Type BF equipment

#### Maintenance and Preservation

1. Replace the batteries timely when low voltage lamp is on.
2. Clean the surface of fingertip pulse oximeter before it is used to diagnose patients.
3. Remove the batteries inside if you will not operate the Oximeter for a long time.
4. It would be better to preserve the product in -10~40 °C (14-104 °F) and humidity is 10%-80%.
5. It is recommended that the product should be kept dry anytime. A wet ambience might affect its lifetime and even damage the product.
6. Please follow the law of the local government to deal with used batteries.

#### Symbols and Definitions

	BF type application part		IP degree		Serial number
	Separate collection		Humidity range		Date of manufacture
	Reference manual		Keep dry		Manufacturer
	Cautions		Product certification		European union representative
	Standby		Temperature range		Avoid sunlight
	Up toward				

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### Ver.1.0

#### Fingertip Pulse Oximeter

This Fingertip Pulse Oximeter is a kind of innovated medical device with non-invasive and continuous features for artery SpO2 and PR detection. Being portable, it is able to measure SpO2 and PR values quickly and precisely.

#### General Description

Haemoglobin Saturation is the percentage between the capacity of Oxyhemoglobin (HbO2) that compounded with oxygen and that of all combinative haemoglobin (Hb) obin (HbO2) in blood. In other words, it is the saturation of oxyhemoglobin in blood. It is a very important physiological parameter for Respiratory and Circulation Systems. Many respiratory diseases could reduce haemoglobin saturation in human blood. Moreover, factors such as Automatic Organic Regulation Malfunction caused by anaesthesia, trauma resulted from major operation and some medical examination can also cause problems in oxygen supply, which might reduce human haemoglobin saturation. As a result, such symptoms as megrim, vomiting and asthenia might appear to patients. Hence, it is very important to know haemoglobin saturation of patient timely in clinical medical aspects.

The fingertip pulse oximeter features in small volume, low power consumption, convenient operation and portability. It is only necessary for patient to put one finger into fingertip photoelectric sensor for diagnosis, and the display screen will directly show measured value of hemoglobin saturation. It has been proved in clinical experiments that it possesses rather high precision and repeatability.

#### Product Accessories

1. One hang lace
2. One user manual

#### Guidance and manufacture's declaration-electromagnetic radiation-for other EQUIPMENTS and SYSTEMS

The Fingertip Pulse Oximeter is designed to be used in specified electromagnetic environment. Users of the Fingertip Pulse Oximeter must fit in the following environments.

Radiation Test	Compliance	Electromagnetic environment-guidance
RF interference CISPR 11	Group 1	RF signal of Fingertip Pulse Oximeter is simply carried by its internal function. Therefore, its RF interference is very low and is not likely to cause any interference to nearby electronic equipment.
RF interference CISPR 11	Class B	The Fingertip Pulse Oximeter applies to all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

#### Possible Problems and Resolutions

Problem	Possible reason	Solution
SpO2 or PR can not be shown normally	1. Finger is not plugged correctly 2. Patient's Oxyhemoglobin value is too low to be measured	1. Retry by plugging the finger 2. Try more times. If you can make sure there is no problem in the product, please go to hospital timely for exact diagnosis
SpO2 or PR is shown unsteady	1. The finger might not be plugged deep enough 2. Finger is trembling or the patient is moving	1. Retry by plugging the finger 2. Please remain at rest
The Oximeter can not be turned on	1. Inadequate power or power off 2. Batteries might be installed incorrectly 3. The Oximeter might be damaged	1. Please replace the batteries 2. Please reinstall the batteries 3. Please contact with local customer service centre
Indication lamps are suddenly off	1. The product automatically shuts off when no signal is detected in 8 seconds 2. Inadequate power	1. Normal 2. Replace the batteries

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

16. It may cause inaccurate reading when the positions of sensor and blood pressure cuff are on the same arterial catheter or intravascular line.
17. Hypotension, severe vasoconstriction, severe anemia, or hypothermia may cause inaccurate reading.
18. It may cause inaccurate reading by giving use of cardiotoxic to patient after his cardiac arrest or when he is in quiver.
19. Bright nail or painted nail may cause inaccurate SpO2 reading. Follow local ordinances and recycling instructions regarding to disposal or recycling of the device and device components, including batteries.

**Scope of application / Intended use**  
The fingertip pulse oximeter can be used to measure human haemoglobin saturation and pulse rate through finger , it can be used in hospitals, families, schools and medical centers.

#### contraindication: not found

#### note

1. The image in the instruction may have slight differences with the actual instruments.
2. The actual parameters and appearance change, without prior notice.

#### Product include: main machine and SpO2 sensor.

- ◆ OLED display
- ◆ Product adopts double color OLED display , can show the six different display mode
- ◆ If the hand movements, under the effect of accelerometer, the interface can save four different kinds of display mode (suitable for matching accelerometer function instrument)
- ◆ Low-power consumption, continuously work for more than six hours with two AAA batteries
- ◆ Low voltage indicator

- ◆ In the absence of signals, the product will be in after 8 seconds to enter standby state
- ◆ Small in volume, light in weight, and convenient to carry
- ◆ Instrument has 5s automatic signal detection function, when you insert finger, timely will automatically start;Automatic startup function instrument (applicable to Automatic startup function instrument)

#### Operation Instructions

1. Install two AAA batteries into battery cassette before covering its cover.
2. Plug one finger into rubber hole of the Oximeter (It is best to plug the finger thoroughly) before releasing the clamp with the nail upwards.

Reserves the right to technical change appearance, our products are subject to change without prior notice, please forgive me!

#### Statement:

1. Maintenance with data such as circuit diagram, components list, figure and the related rules of correction, injection, available only to the repair factory training qualified personnel and units.
2. The company can be in the form of email or other electronic files provide users with random files.
3. The instrument is not used for evaluation of blood oxygen probe pulse and pulse blood oxygen monitor accuracy.

#### After-sales service

#### Ensure that users

- ◆ Please read user manual before using the instrument;
- ◆ According to the requirement of the instruction manual for the operation and daily maintenance, and make sure the machine power supply, and environmental requirements

#### Maintenance time

Settings	Settings		
Alarm Setup	off	Sounds Setup	-
Beep	off	SpO2 Alm H	100
Demo	off	SpO2 Alm Lo	85
Restore	OK	PR Alm Hi	130
Brightness	OK	PR Alm Lo	50
Exit	3	Hi	Exit
Exit	+		

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