



SINCE 1966
MADE IN ITALY

Wi.Neb

by **FLAEM**®



**Powerful and compact
aerosol therapy unit**

Speed, Silence, Comfort


**New RF7 Dual Speed Plus
nebulizer**




Plus RF7

2 speed nebulizer with valve system developed based on the latest medical-scientific guidelines

Use with Mouthpiece




Maximum therapeutic effect is ensured by the valve system allowing optimisation of medication delivery and particle size to reach the lower respiratory tract.



Minimum dispersion of medication is ensured by the synchronised opening and closing of the valves. Exhaled air is directed towards the dedicated exit. The valve limits dispersion and the aerosol mixture is enriched for the next inhalation.

New SoftTouch Masks with dispersion limiting vent



During inspiration the vent bends inwards.





During expiration the vent bends outwards.

How to use

Delivery speed (2)

with valve max



⁽¹⁾ Data in compliance with EN 13544-1 and certified by TÜV Rheinland LGA Products GmbH - Germany

⁽²⁾ Data shown is in accordance with Flaem Nuova internal procedure I29-P07.5 and may vary depending on the patients respiratory capacity.



The new **Soft Touch Masks** feature an adaptive edge made of a soft biocompatible material ensuring an excellent adherence to the face.

The **Dispersion Limiting Vent** optimizes the treatment by reducing the amount of medication released into the environment. The aerosol mixture is enriched for the next inhalation.

Air filtration

Air filtration is ensured by the special filter. Easy to reach, it is located under the little lid on the side that together with the casing acts a **silencer**.



SINCE 1966
MADE IN ITALY
Wi.Neb
by **FLAEM**®



Why choose Wi.Neb



RF7 Dual Speed Plus nebulizer



high performance



extremely silent



adult/child SoftTouch masks



environmentally friendly
low energy consumption



usable in automobiles




convenient carrying bag

Wi.Neb with power supply unit or cigarette lighter cable

Compact size and powerful device. It can be powered by a 100/240V - 50-60Hz multi voltage power supply unit or by a cigarette lighter power cord.

Between 2 and 6 μm
for laryngotracheal
deposition

Between 0.5 and 3 μm
for alveolar deposition



MMAD (1): Mass Median Aerodynamic Diameter
FPF (1): Fine Particle Fraction (breathable fraction)
ml/min (2): delivery speed

(1) In vitro testing certified by
TÜV Rheinland LGA Products GmbH - Germany
in compliance with the **European Standard EN 13544-1**
for nebulizing systems

(2) Data shown is in accordance with **Flaem Nuova** internal
procedure **I29-P07.5**. Values shown for the aerosol delivery
may vary depending on the patients respiratory capacity

In compliance with:
European Standard **EN 10993-1**
"Biological evaluation of medical device"
and **93/42/EEC Medical Device Directive**



TECHNICAL FEATURES

12 Vdc powered rotary piston compressor (does not require lubrication and maintenance)	
Air filtering:	For better quality, the compressed air is sucked in from the external environment and filtered with a high efficiency filter system (density= 230 mg/m ³)
Compressor air output:	14 l/min' (approx.)
Max Pressure	2 ±0,3 bar
Noise level (at 1 m)	52 dB (A) (approx.)
Dimensions	13 (W) x 9,5 (D) x 5 (H) cm
Weight of device:	430 g
Power supply/Power:	Multivoltage PSU 100/240 V ~ 50/60 Hz/12V DC

APPROVALS

In compliance with 93/42/EEC Medical Device Directive	
TÜV approval in compliance with EN 60 601-1 / 3rd Ed.	 • BAUART GEPRÜFT • TYPE APPROVED
Licensed test mark: USA, CANADA	
In compliance with directive E/ECE/324-E/ECE/TRANS/505 Regulation No.10 - Add.9/ Rev.3	 10 R - 03 xxxx

NEBULIZER TECHNICAL FEATURES



RF7 Dual Speed Plus nebulizer	Speed selector	
	MIN with valve	MAX
Delivery ml/min' ⁽¹⁾	0.18	0.42
MMAD ⁽²⁾	2.95 µm	2.44 µm
Breathable fraction < 5 µm ⁽²⁾	74.7 %	80.7 %
Maximum fill volume ⁽¹⁾	8 ml	

(1) Data shown is in accordance with Flaem Nuova internal procedure I29-P07.5
Values shown for the aerosol delivery may vary depending on the patients respiratory capacity
(2) In vitro testing certified by TÜV Rheinland LGA Products GmbH - Germany in compliance with the European Standard EN 13544-1 for nebulizing systems

Wi.Neb^{go} WHITE
code WI00P00



INCLUDED ACCESSORIES

Code

1	RF7 Dual Speed Plus Nebulizer Kit composed of: Nebulizer, mouthpiece, nasal prong and Drug saving device	ACO436P
2	1 m connection tube	ACO35
3	SoftTouch Adult Mask	ACO437P
4	SoftTouch Pediatric mask	ACO438P
5	Multivoltage power supply unit	ACO458P
6	12VDC cigarette lighter power cord	ACO211P
7	Carrying bag	



Manufacturer:

FLAEM NUOVA S.p.A.
Via Colli Storici, 221
25015 S. Martino della Battaglia (BS) Italy
Tel. +39 030 9910168 - Fax +39 030 9910287
info@flaemnuova.it

FLAEM NUOVA S.p.A.
Company with certified quality system



Export Department of FLAEM NUOVA



3, Via del Passero - 20147 Milano - Italy
Phone (39) 02 4120769
Fax (39) 02 48376295
E-mail: info@medexpointernational.com