# Rhino Clear

**PHTHALATES  
& BPA FREE**

**IT** > MANUALE ISTRUZIONI D'USO  
**EN** > OPERATING INSTRUCTIONS MANUAL  
**FR** > MODE D'EMPLOI  
**DE** > BEDIENUNGSANLEITUNG  
**ES** > MANUAL DE INSTRUCCIONES DE USO  
**PL** > INSTRUKCJA OBSŁUGI I KONSERWACJI

**NL** > GEBRUIKSAANWIJZING  
**DE** > BETRIEBSANLEITUNG  
中文 > 使用说明书

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## Nasal douche

**These user instructions are provided for device model RC1-1. The nasal douche may be combined with a: Flaem compressor unit, non-Flaem compressor unit, centralised compressed air systems; in the latter two cases the working pressure must fall within the range of values defined in the paragraph TECHNICAL SPECIFICATIONS provided below in the manual.**

### INTENDED USE

Medical Device for washing and treatment of the nasal cavities using saline solutions and/or medicines, the treatment and medicines must be prescribed by a doctor.

### INDICATIONS FOR USE

Treatment of upper airway disorders, via the endonasal route (rhinitis, rhinosinusitis, sinusitis, and or as indicated by the treating physician).

Physiological washing of the nasal cavities (with saline, isotonic, hypertonic solutions, or hot spring water). In children it helps with nasal hygiene, in adults it aids physiological nasal washing and can be used before therapy of the upper airways to assist drug administration.

### **CONTRAINDICATIONS**

- The medical device should NOT be used for patients who are unable to breathe on their own or who are unconscious.
- Do not use the device in anaesthetic or assisted ventilation circuits.

### INTENDED USERS

The device is intended for use by legally authorised medical personnel/health workers (doctors, nurses, therapists, etc.). The device can be used directly by the patient.

### **TARGET GROUP OF PATIENTS**

Adults, children of all ages, infants. Before using the device, the user manual must be read carefully and an adult responsible for safety must be present if the device is used by infants, children of any age or persons with limited abilities (e.g. physical, mental, or sensory). It is the responsibility of the medical personnel to assess the patient's condition and capabilities in order to determine, when prescribing therapy, whether the patient is able to operate the device safely on their own or whether the therapy should be administered by a responsible person.

Medical personnel should be responsible for assessing the use of the device on specific types of patients such as pregnant women, lactating women, infants, incapacitated persons or persons with limited physical capabilities.

### OPERATING ENVIRONMENT

This device can be used in healthcare facilities, such as hospitals, outpatient clinics, etc.,

or even at home.

## **WARNINGS CONCERNING POSSIBLE MALFUNCTIONS**

- Should your device fail to perform, please contact the authorised service centre for clarifications.
- The manufacturer should be contacted to report problems and/or unexpected events relating to operation and if necessary for clarification on use and/or maintenance/cleaning.
- Please also refer to the troubleshooting section and problem resolution.

## **WARNINGS**

- Use the device for therapeutic purposes only. This medical device is not intended as a life-saving device. Any other use is considered misuse and may be dangerous. The manufacturer is not liable for any misuse.
- Always consult your general practitioner for identification of treatment.
- Follow your doctor's or therapist's instructions regarding the possible use of types of medicine, dosage and indications for treatment.
- Pay attention to the indications provided with the drug and avoid using the devices with substances and dilutions other than those recommended.
- If you should experience allergic reactions or other problems while using the device, stop using it immediately and consult your doctor.
- Keep this manual carefully for further reference.
- If the packaging is damaged or opened, contact the distributor or service centre.
- The time required to switch from storage to operating conditions is about 2 hours.
- Repairs may only be carried out by personnel authorised by the manufacturer. Unauthorised repairs invalidate the warranty and may pose a danger to the user.
- The nasal douche must not be used as a toy. The utmost attention should be paid when used by children.
- Suffocation risk:
  - Some components of the device are small enough to be swallowed by children, therefore, keep the device out of the reach of children.
- Risk of ineffectiveness of therapy:
  - In the case of substances that are too dense, dilution with suitable saline solution may be necessary, as prescribed by a doctor.
  - Comply with the cleaning instructions of the nasal douche as described in the hygienic preparation section.
  - Only use original Flaem spare parts, no liability is accepted if non-original parts are used
- Risk of infection:
  - Personal use of the nasal douche is recommended to avoid any risk of infection.
  - Carry out hygienic preparation before each use. Ensure that the douche and its parts are not stored in the vicinity of other accessories or devices for different therapies (e.g. infusion devices).
  - Once the therapy is over, do not leave the liquid used (drug or saline solution) inside the spray chamber and perform the hygienic preparation operations.

## TROUBLESHOOTING

Problems	Solutions	
If the nasal douche delivers insufficiently or for a short time.	Check that the atomiser (3) has been positioned correctly. Check that the nasal douche is clean and, if it is not, perform cleaning again by carefully following the instructions for hygienic preparation.	
If the nasal douche does not spray	Check that the solution has been correctly added into the spray chamber (4.1) of the ampoule body (4). Check that the amount of liquid does not exceed 10 ml. Check that the nasal douche is clean and, if it is not, perform cleaning again by carefully following the instructions for hygienic preparation.	
	The accessories are not properly connected to the compressor unit or flow of air from the centralised system	Check proper connection of the accessories with the compressor unit or centralised system air flow
	The compressor unit or centralised system does not deliver any air or the airflow is insufficient	Check that the compressor unit or centralised system delivers the airflow and that the air generates a working pressure within the range set out in the table provided in the section TECHNICAL SPECIFICATIONS below in the manual.
If, after checking the conditions described above, the device is still not working properly, we recommend contacting your trusted dealer or an authorised FLAEM service centre nearest to you. You can find a list of all Service Centres at <a href="http://www.flaemnuova.it/it/info/assistenza">http://www.flaemnuova.it/it/info/assistenza</a>		

## DISPOSAL

The nasal douche is to be disposed of as general waste after a sanitisation cycle.

### Packaging

Product box

Heat-shrink film

## NOTIFICATION OF SERIOUS EVENTS

Serious events occurring in connection with this product must be reported immediately to the manufacturer or the competent authority.

An event is considered serious if it causes or may cause, directly or indirectly, death or an unforeseen serious deterioration in a person's state of health.

COUNTRY	AUTHORITY
Ireland	<b>Health Products Regulatory Authority - Kevin O'Malley House, Earlsfort Centre, Earlsfort Terrace, IE - Dublin 2</b> <b>E-mail: <a href="mailto:devicesafety@hpra.ie">devicesafety@hpra.ie</a></b>
Malta	<b>Medical Devices Unit Medicines Authority - Sir Temi Żammit Buildings, Malta Life Sciences Park - San Ġwann SĠN 3000</b> <b>E-mail : <a href="mailto:devices.medicinesauthority@gov.mt">devices.medicinesauthority@gov.mt</a></b>

## SYMBOLS ON THE DEVICE OR ON THE PACKAGING

 Medical CE marking ref. regulation 2017/745 EU as amended

Phthalate- and bisphenol-free

Model number

Temperature limits

Atmospheric pressure limits

Moisture limits

Batch Number

Manufacturer

Attention

Type BF applied part

Medical device

Unique device identifier

Production date

Check the instructions for use

## INFORMATION ON RESTRICTIONS OR INCOMPATIBILITIES WITH CERTAIN SUBSTANCES.

- Interactions: The materials used in the device are biocompatible materials and comply with statutory regulations, however possible allergic reactions cannot be completely ruled out.
- Use the medicine as soon as possible once it has been opened and avoid leaving it in the nasal douche; once therapy has ended, do not leave the medicine inside the nasal douche and carry out hygienic preparation.

## TECHNICAL SPECIFICATIONS

**Model:** RC1-1

Maximum spray chamber capacity: 10 ml

The nasal douche does not work with ultrasonic or non-pneumatic aerosol therapy devices.

The nasal douche works with :  
1- Pneumatic aerosol therapy devices.  
2- Centralised air systems.

The nasal douche works with working pressure values within the range of 0.50 to 2.10 bar

Dispensing time of 10 ml: 2 minutes approx.

(Tests performed with saline solution 0.9% NaCl according to Flaem internal procedure)

MMD\*: > 10  $\mu\text{m}$

(\*) The MMD value shown was taken with the Malvern Spraytec instrument and refers to the use of saline (0.9% NaCl), so it may vary depending on the medicine used. The value indicated does not apply to medicines dispensed in suspension or with high viscosity. In that case, information must be obtained from the supplier of the drug.

Nasal douche dimensions: 7x6x15.5h cm

Nasal douche weight: 57 g

Type BF applied parts are the 3 nasal adapters (1A, 1B, 1C)

## ENVIRONMENTAL CONDITIONS

### Operating conditions:

Ambient temperature Between +10°C and +40°C

Relative air humidity Between 10% and 95%

Atmospheric pressure Between 69KPa and 106 KPa

### Storage and transport conditions:

Ambient temperature Between -25°C and +70°C

Relative air humidity Between 10% and 95%

Atmospheric pressure Between 69KPa and 106 KPa

## DURATION

**Model:** RC1-1

The service life of the device is one year

## EQUIPMENT AND INFORMATION ON MATERIALS

The supply includes:	Information on materials
1A - Transparent nasal adapter 'A', suitable for children up to 3 years of age 1B - White nasal adapter 'B', suitable for 3 to 12 years of age 1C - Blue nasal adapter 'C', indicated after 12 years of age	Polypropylene
2 - Separator	
3 - Atomiser	
4 - Douche body with ergonomic handle 4.1 - 10 ml spray chamber 4.2 - Recovery chamber	
5.1 - Dispensing button 5.2 - Rubber stopper	Acetal resin Elastomer
6 - Tube-Nose Douche connection	Polypropylene
IMPORTANT NOTE: Remove the identification label on the packaging and affix it in the appropriate spaces on page 2.	

## INSTRUCTIONS FOR USE

**Before each use, wash your hands thoroughly and clean the NASAL DOUCHE and accessories as described in the section Hygienic preparation.**

During application, it is recommended to adequately protect oneself from any dripping.

1. Insert the atomiser (3) into the spray chamber (4.1).
2. Place the solution (max. 10 ml) in the spray chamber (4.1)(Fig. A).
3. Attach the separator (2), ensuring that the arrow on the douche body is aligned with the arrow on the separator (fig. B). After that, insert the age-specific adapter, or the one most suitable for the size of your nostrils (1A/1B/1C).
4. Connect the Nasal Douche to the compressor unit or the central compressed air system using the tube (7 not supplied) and the connector (6).
5. Start up the nasal douche by switching on the compressor unit or supplying it with compressed air.
6. Place the adapter hole against one nostril.
7. During the treatment, it is preferable to close with one finger the other nostril where the NASAL DOUCHE is not applied.
8. Operate the NASAL DOUCHE by holding down the dispensing button (5) during inhalation and release it before exhaling; during this phase, to help expel mucus, blow gently through the treated nostril directly into the NASAL DOUCHE.
9. Repeat steps 6. to 8. for the other nostril as well, trying to distribute the solution evenly between the two nostrils.
10. In the case of children or impaired persons, and however unable to co-ordinate the various steps in 8., it is recommended to carry out the application alternating between the two nostrils every 20 seconds approximately, with the necessary pauses to expel the mucus.

11. What is described in the previous points is useful for nasal douching using saline solutions. It is advisable to carry out the drug treatment after initial washing with saline solution to achieve a more effective distribution of the drug on the nasal mucous membranes.
12. When application is over, switch off the machine and carry out cleaning following the instructions in paragraph Hygienic preparation.

## HYGIENIC PREPARATION

Disassemble the Nasal Douche as shown in the connection diagram.

Then proceed according to the instructions below.

### Sanitisation

Before and after each use, sanitise the components, choosing one of the methods described below.

**method A:** sanitise the components under warm drinking water (approx. 40°C) with mild dish-washing detergent (non-abrasive).

**method B:** sanitise the components in a dishwasher on a hot cycle.

**method C:** sanitise the components by soaking in a solution of 50% water and 50% white vinegar, then rinse thoroughly with warm drinking water (approx. 40°C).

If you also wish to carry out cleaning by DISINFECTION, skip to the DISINFECTION section.

After sanitising the components, shake them vigorously and place them on a paper towel, or alternatively dry them with a hot air jet (e.g. hair dryer).

### Disinfection

After sanitising the components, disinfect them using one of the methods described below.

**method A:** purchase an electrolytic chloride disinfectant (active ingredient: sodium hypochlorite), specifically for disinfection, available at all chemists'.

#### **Execution:**

- Fill a container of a suitable size to hold all the individual components to be disinfected with a solution of drinking water and disinfectant, complying with the proportions indicated on the packaging of the disinfectant.
- Completely immerse each individual component in the solution, taking care to prevent any air bubbles from forming in contact with the components. Leave the components immersed for the period of time indicated on the disinfectant packaging, and associated with the concentration chosen for the preparation of the solution.
- Retrieve the disinfected components and rinse them thoroughly with lukewarm drinking water.
- Dispose of the solution according to the disinfectant manufacturer's instructions.

**method B:** disinfect the components by boiling them in water for 20 minutes; use demineralised or distilled water to prevent hard water deposits.

**method C:** disinfect the components in a hot bottle steriliser (not microwave). Carry out the process strictly according to the hot bottle steriliser's instructions. For disinfection to be effective, choose a steriliser with an operating cycle of at least 6 minutes.

After disinfecting the components, shake them vigorously and place them on a paper towel, or alternatively dry them with a hot air jet (e.g. hair dryer).

### Sterilisation

Sterilisation for clinical or hospital use is the responsibility of the facility performing the sterilisation and must be carried out according to the instructions in this manual.

The sterilisation procedure described in this section is only effective on treated components if all the points are adhered to and only if the components to be treated are sanitised beforehand, and is validated in accordance with ISO 17665-1.

**Equipment:** Steam steriliser with fractionated vacuum and overpressure compliant with EN 13060.

**Execution:** Package each individual component to be treated in sterile barrier system or packaging in accordance with EN 11607. Place the packaged components in the steam steriliser. Perform the sterilisation cycle in accordance with the instructions for use of the equipment and select a temperature of 121°C and a time of 20 minutes.

**Storage:** Store the sterilised components according to the instructions for use of the system or sterile barrier packaging chosen.


#### **ADDITIONAL INFORMATION FOR CLEANING the spray chamber**

In order to prevent deposits from clogging the atomiser and spray chamber and thereby affecting the function of the nasal douche, it is recommended to carry out cleaning and sanitising after each use, following the instructions in the user manual supplied with the nasal douche. After thoroughly rinsing the parts under lukewarm tap water, shake them to remove any remaining drops of water. Then place the parts on a paper towel or dry them with a jet of hot air. **ATTENTION:** If your nasal douche does not atomise correctly after cleaning and sanitising, do not use sharp tools (e.g. screwdrivers, scissors, needles, hairpins, etc.) to try to clean or unclog the holes in the atomiser or spray chamber, as this may damage these parts and impair the proper functioning of the unit. Instead, perform additional cleaning and sanitisation cycles as described in the instruction manual and recommended above, until the nasal douche starts working properly again.

**Table of intended methods / components**

Components									
	Method	1A	1B	1C	2	3	4	5	6
<b>HYGIENIC PREPARATION IN THE HOME</b>									
<b>Sanitisation</b>									
method A	✓	✓	✓	✓	✓	✓	✓	✓	✓
method B	✓	✓	✓	✓	✓	✓	✓	✓	✓
<b>Disinfection</b>									
method A	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES
method B	✓ MAX 40 TIMES	✓ MAX 40 TIMES	✓ MAX 40 TIMES	✓ MAX 40 TIMES	✓ MAX 40 TIMES	✓ MAX 40 TIMES	✓ MAX 40 TIMES	✓ MAX 40 TIMES	✓ MAX 40 TIMES
method C	✓ MAX 40 TIMES	✓ MAX 40 TIMES	✓ MAX 40 TIMES	✓ MAX 40 TIMES	✓ MAX 40 TIMES	✓ MAX 40 TIMES	✓ MAX 40 TIMES	✓ MAX 40 TIMES	✓ MAX 40 TIMES
<b>HYGIENIC PREPARATION IN A CLINIC OR HOSPITAL</b>									
<b>Sterilisation</b>	✓ MAX 40 TIMES	✓ MAX 40 TIMES	✓ MAX 40 TIMES	✓ MAX 40 TIMES	✓ MAX 40 TIMES	✓ MAX 40 TIMES	✓ MAX 40 TIMES	✓ MAX 40 TIMES	✓ MAX 40 TIMES
✓: intended \: not intended									

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© 2023 FLAEM NUOVA®
**FLAEM NUOVA S.p.A. non risponde di ulteriori estensioni del periodo di garanzia assicurate da terzi.**
FLAEM NUOVA S.p.A.