Medical Devices

# AirFeel

IT > MANUALE ISTRUZIONI D'USO

EN > USER MANUAL

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## Medical device for respiratory gymnastics

**These operating instructions are provided for device model P0920EM.**

### **INTENDED USE**

Medical device for the prevention of increased accumulation of secretions and for therapy, in case of acute and chronic diseases of the lower airways. Prevention and/or treatment practices must be prescribed by professional medical personnel.

### **INDICATIONS FOR USE**

Diseases of the lower airways associated with increased accumulation of mucus and secretions. It allows you to perform PEP (Positive Expiratory Pressure) therapy and O-PEP (Oscillatory Positive Expiratory Pressure) therapy.

### **CONTRAINDICATIONS**

The device should not be used by persons with severe pneumothorax or those coughing up large quantities of blood (massive hemoptysis).

The following side effects may occur: Dizziness

If you notice any of these side effects, stop the treatment immediately. Contact your doctor or therapist to discuss how to proceed with the treatment.

Inform your doctor or therapist if you notice any side effects not listed in these instructions for use.

### **INTENDED USERS**

The devices are intended to be used by professional medical personnel. The device can be used directly by the patient.

### **GROUP OF RECIPIENT PATIENTS**

The device can be used:

- from children aged 4 to 10 years.
- by children aged 10 and over and by adults.

The frequency and duration of use are determined by professional medical personnel based on the individual needs of the patient. Before use, the device requires careful reading of the user manual and the presence of an adult responsible for safety if the use is intended for children or people with limited abilities (for example: physical, mental or sensory). It is up to the medical staff to evaluate the patient's conditions and abilities to establish whether the patient is able to operate independently or if the patient is not able to use the device safely independently and whether the prevention and/or therapy practices must be performed by a responsible person.

### **OPERATING ENVIRONMENT**

This device can be used in health care facilities, such as hospitals, clinics, etc., or even at home.

## WARNINGS

- Use the device only as intended. This medical device is not intended as a life saving device. Any different uses are to be considered improper and can be dangerous. The manufacturer is not responsible for any improper use.
- This product is to be used only as described in these instructions for use.
- The time required to pass from storage conditions to operating conditions is approximately 2 hours.
- Always consult professional medical personnel for identification of cure.
- Keep this manual carefully for any further consultation.
- When the device is used in combination with another medical device, the instructions for use of the same must also be followed.
- The materials used in the device are biocompatible materials and comply with mandatory regulations, however possible allergic reactions cannot be completely excluded.
- The device is made of an impact-resistant but not unbreakable material. Therefore, protect it from shocks and be careful not to drop it.
- The device does not administer drugs or substances.
- At the end of each use, after cleaning as per the HYGIENIC PREPARATION paragraph, store the device in a dry place away from dust.
- In the case of PEP therapy, no particular knowledge is required on the part of the user (who is normally the patient himself) except the ability to verify that the selected calibrated hole corresponds to the medical prescription, which occurs through a specific pressure gauge provided only to staff professional doctor.
- Risk of suffocation
- Some components of the device are so small that they can be swallowed by children; therefore keep the device out of the reach of children.
- Risk of ineffectiveness of therapy:
  - Use the device in the correct position (see instructions for use paragraph).
  - Use original Flaem spare parts, we decline all responsibility in the event non-original spare parts or accessories are used.
  - O-PEP treatment will be ineffective if the patient inhales/exhales too rapidly or has to force himself excessively to breathe.
- Risk of infection:
  - This medical device should only be used by a single patient.
  - Before and after each use, follow the hygienic preparation operations. Make sure the device is not stored near other accessories or devices for different therapies (e.g. infusions).

## TROUBLESHOOTING CASES AND RELATED SOLUTION

Problem	Cause	Solution
In case of PEP therapy: The respiratory effort is too heavy or too light	Incorrect setting of the hole on the ferrule	Check the correct position of the calibrated hole. If correct, consult your doctor / physiotherapist
In case of OPEP therapy: The sphere does not vibrate	The system is not assembled correctly	Check the correct assembly of the parts
	The ring nut is not positioned correctly	Position the ring nut in position "0"
	The cap is closed	Rotate the cap so that the indicator is aimed at the open padlock symbol

If, after having checked the conditions described above, the device still does not work properly, we recommend that you contact your trusted dealer or an authorized FLAEM service center closest to you. Find the list of all the Assistance Centers on the page <http://www.flaemnuova.it/it/info/assistance>

## DISPOSAL

All components of the product can be disposed of with domestic waste unless this is prohibited by the disposal regulations in force in the respective member countries.

### Packaging

Product box

Product packaging bag

Pouch packing bag

## COMMUNICATION OF SERIOUS EVENTS

Serious events that occur in relation to this product must be reported immediately to the manufacturer or the competent authority. An event is considered serious if it causes or may cause, directly or indirectly, death or an unexpected and serious deterioration in a person's state of health.

COUNTRY	AUTHORITY
Ireland	<b>Health Products Regulatory Authority</b> <b>Kevin O'Malley House, Earlsfort Centre, Earlsfort Terrace, IE - Dublin 2</b> <b>E-mail: devicesafety@hpra.ie</b>
Malta	<b>Medical Devices Unit Medicines Authority</b> <b>Sir Temi Żammit Buildings, Malta Life Sciences Park San Ġwann SĠN 3000</b> <b>E-mail : devices.medicinesauthority@gov.mt</b>

## SYMBOLS PRESENT ON THE DEVICE OR ON THE PACKAGE

Medical device

Manufacturer

Date of production

CE marking ref. regulation 2017/745 EU and subsequent updates

Free from phthalates and bisphenol

Temperature limits

Atmospheric pressure limits

Multiple use per single patient

Consult the instructions for use

Model number

Lot code

Unique identifier of the device

Humidity limits

Attention

Sphere lock

Sphere unlock

## TECHNICAL FEATURES

**Model:** P0920EM

Weight 0,110 Kg approx

## MEDICAL DEVICES USABLE IN CONJUNCTION WITH THE DEVICE

The device can be used in combination with the following medical device:  
Flaem nebulizer model RF7-1 or model RF7-2 (not supplied) for inhalation therapy and PEP / O-PEP therapy simultaneously.  
Combination therapy must be prescribed by medical personnel.

## ENVIRONMENTAL CONDITIONS

### Operating conditions:

Ambient temperature Between +10°C and +40°C

Relative humidity of the air Between 10% and 95%

Atmospheric pressure Between 69 KPa and 106 KPa

### Storage and transport conditions:

Ambient temperature Between -25°C and +70°C

Relative humidity of the air Between 10% and 95%

Atmospheric pressure Between 69 KPa and 106 KPa

## DURATION

**Model: P0920EM** The expected average life is 5 years

## EQUIPMENT COMPONENTS AND INFORMATION ON MATERIALS

Refer to the picture on page 1 of this instruction manual.

The device equipment includes:	Material information
1 - Main body	
2 - Mouthpiece cap	
3 - Mouthpiece	Polypropylene
4 - Valve cap	
5 - Valve	
6 - PEP module	
7 - Ring nut	
8 - Ball	
9 - PEP cap	
10 - Handle	
11 - Nose clip	

## INSTRUCTIONS FOR USE

**Before and after each application, the components of the product must be carefully cleaned following the instructions in the paragraph HYGIENIC PREPARATION.**

### PEP therapy

In PEP therapy, the expiratory resistance during expiration varies according to the different sizes of the holes in the ring nut (7) to which numbers ranging from 1 to 5 are associated, the number 1 (small hole) corresponds to a higher expiratory resistance, instead the lowest corresponds to number 5 (large hole). The increase in expiratory resistance has the effect of stabilizing the airways and stimulating the dissolution of the mucus.


The size of the hole must be determined individually for each patient by professional medical personnel.

### USING THE DEVICE IN PEP MODE

#### Preparation for therapy

Turn the PEP cap (9) to the symbol  (fig.1)

Make sure that the valve (5) is well assembled on the valve cap (4) and adheres perfectly to the internal surface of the cap (fig. 2).

#### Regulation of expiratory resistance

Lift the ring nut (7) and rotate it until the desired hole is aligned with the indicator  PEP module (6) (fig. 3).

Push the ring nut (7) downwards making sure that it fits correctly in the position key on the PEP module (6).

Turn the PEP Cap (9) clockwise until the position dot matches the symbol  (fig.4).

#### Carrying out the treatment

Proceed as follows to carry out the treatment:

- Before starting the treatment, make sure that all parts are firmly connected to each other.
- Make sure that the expiratory resistance is adjusted correctly (see: Regulation of expiratory resistance).
- If the breathing resistance seems too high or too low during a therapy session, stop the session and consult your doctor or therapist.
- Plug the nose with the supplied clip (11).

- Hold the mouthpiece (3) between your teeth, seal the mouthpiece tightly with your lips.
- Inhale slowly and deeply.
- Exhale through the mouthpiece (3).
- The exhaled air must pass through the hole in the adjustment ring.
- Cough to clear any secretions that may have moved during treatment.
- For hygienic reasons, avoid coughing into the AirFeel system.
- Continue the PEP therapy session for the period recommended by your doctor or physical therapist.

### **O-PEP therapy**

The patient repeatedly exhales into the device against a moving ball, causing a vibration that is transmitted to the lungs. This vibration tirs the mucus in the lower airways. During use, the airways are kept open and mucus clearance is enhanced.

If Airfeel is held horizontally, the frequency will fluctuate during use and equal to approximately 15 Hz. This frequency can be varied by tilting the device up (higher frequency) or down (lower frequency). This allows to select the ideal frequency for each individual patient.

### **USE OF THE DEVICE IN O-PEP MODE**

#### Preparation for therapy

Turn the PEP cap (9) to the symbol  (fig.1)


Make sure that the valve (5) is well assembled on the valve cap (4) and adheres perfectly to the internal surface of the cap (fig. 2).

#### **Regulation of expiratory resistance**

Lift the ring nut (7) and rotate it until the "0" symbol is aligned with the indicator  present on the PEP module (6). (fig.5).

Push the ring nut (7) downwards making sure that it fits correctly in the position key on the PEP module (6).

Make sure that the PEP Cap (9) is in position , if not, turn it anticlockwise until it is reached.

#### **Execution of the treatment**

Proceed as follows to carry out the treatment:

- Sit upright and relax.
- Hold the AirFeel handpiece (10) so that it is vertical to the ground.
- Hold the mouthpiece (3) between your teeth and seal it tightly with your lips.
- NOTE: Make sure the holes in the top are not blocked or covered by your hand.
- Inhale as slowly and deeply as possible through your nose.
- Exhale slowly and completely through your mouth into the device (as slowly and as long as possible). As you exhale, tense your cheek muscles and keep your cheeks tense.

- Repeat the inhalation and exhalation cycle following the instructions provided by the medical staff
- End the treatment session if the effort becomes excessive.

#### **Notes for optimal treatment:**

Place your free hand against your chest while inhaling and exhaling: when using the device correctly, you can feel an unmistakable regular vibration. The intensity of the vibration varies according to the alignment of the device on the ground: align the device in the position where the vibration is perceived to be strongest (tilt the device slightly upwards or downwards).


#### **Therapy in combination with a medical device model RF7-1 or RF7-2**

It is feasible both in the case of PEP therapy and O-PEP therapy (for all the instructions and restrictions on use of the device model RF7-1 or RF7-2, refer to the nebulizer's manual provided).

### **USE OF THE DEVICE IN COMBINED MODE WITH MEDICAL DEVICE MODEL RF7-1 or RF7-2**

Having defined the therapy to be applied PEP or O-PEP and carried out the operations described in the respective paragraphs "Preparation for the therapy" and "Adjustment of the expiratory resistance" and "Execution of the treatment" connect the output of the nebulizer model RF7-1 or RF7-2 in the valve cap (4) as shown in the figure.

In this mode, the mouthpiece to be used is the one supplied with the AirFeel device and not the one with the exhalation valve supplied with the Flaem nebuliser model RF7-1 or RF7-2. Do not use any masks supplied with the nebuliser model RF7-1 or RF7-2.

### **HYGIENIC PREPARATION**

Take apart the AirFeel system in its individual parts.

#### **Sanitization**

Sanitize the components by choosing one of the methods provided in the table and described below.

**method A:** Sanitize the accessories under hot drinking water (about 40°C) with a delicate dishwashing detergent (non-abrasive).

**method B:** Sanitize the accessories in the dishwasher with a hot cycle (70°C).

**method C:** Sanitize the accessories by immersing them in a solution of 50% water and 50% white vinegar, then rinse thoroughly with hot drinking water (about 40°C).

#### **Disinfection**

After having sanitized the components, disinfect them by choosing one of the methods provided in the table and described below. Each method can be used for a limited number

of times (see data in the table).

**method A:** Get a disinfectant of electrolytic chloroxidizing type (active ingredient: sodium hypochlorite), specific for disinfection, available in all pharmacies.

Execution:

- Fill a container of suitable size to contain all the individual components to be disinfected with a solution based on drinking water and disinfectant, respecting the proportions indicated on the disinfectant package.
- Completely immerse each single component in the solution, taking care to avoid the formation of air bubbles in contact with the components. Leave the components immersed for the period of time indicated on the disinfectant packaging, and associated with the concentration chosen for the preparation of the solution.
- Retrieve the disinfected components and rinse them thoroughly with warm drinking water.
- Dispose of the solution according to the directions provided by the disinfectant manufacturer.

**method B:** Disinfect the components by boiling in water for 10 minutes; use demineralized or distilled water to avoid limescale deposits.

**method C:** Disinfect the components with a hot steam type bottle sprayer (not microwave). Carry out the process by faithfully following the instructions of the vaporizer. For disinfection to be effective, choose a sterilizer with an operating cycle of at least 6 minutes

If you also want to carry out sterilisation, skip to the Sterilisation section.

After disinfecting the accessories, shake them vigorously and place them on a paper towel, or alternatively dry them with a hot air jet (e.g. hair dryer).

### **Sterilisation**

Device: Steam steriliser with fractionated vacuum and overpressure compliant with EN 13060.

**Execution:** Package each individual component to be treated in sterile barrier system or packaging in accordance with EN 11607. Place the packaged components in the steam steriliser. Perform the sterilisation cycle in accordance with the instructions for use of the equipment by selecting a temperature of 134°C and a time of 4 minutes first.

**Storage:** Store sterilised components according to the instructions for use of the system or sterile barrier packaging chosen.

The sterilisation procedure was validated in accordance with ISO 17665-1.

**At the end of each use, store the device complete with its accessories in a dry and dust-free place.**

Table of planned methods / patient accessories							
Accessory patient							
Method	1	2	3	4-5	6	7-8	9
<b>HYGIENIC PREPARATION IN THE HOME</b>							
<u>Sanitisation</u>							
method A	✓	✓	✓	✓	✓	✓	✓
method B	✓	✓	✓	✓	✓	✓	✓
method C	✓	✓	✓	✓	✓	✓	✓
<u>Disinfection</u>							
method A	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES
method B	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES
method C	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES
<b>HYGIENIC PREPARATION IN A CLINIC OR HOSPITAL</b>							
<u>Disinfection</u>							
method A	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES
method B	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES
<u>Sterilization</u>	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES
✓: expected    \: not expected							

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© 2023 FLAEM NUOVA® All rights reserved
**FLAEM NUOVA S.p.A. non risponde di ulteriori estensioni del periodo di garanzia assicurate da terzi.**
FLAEM NUOVA S.p.A.