
AirFeel IMT

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metodo A ✓ MAX 300 VOLTE

metodo B ✓ MAX 300 VOLTE

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Medical device for respiratory gymnastics

These operating instructions are provided for the device model P0124EM.

INTENDED USE

Medical device for training inspiratory muscles. Use must be prescribed by professional medical personnel.

DIRECTIONS FOR USE

This device is intended to relieve symptoms in major debilitating diseases such as COPD, asthma, heart disease, dyspnea or other medical indications.

Before use, the setting and any changes to the device settings must be performed by the Medical staff or by the patients as prescribed by their doctor/physiotherapist.

Notes for medical staff: The workload for respiratory gymnastics should be based on the patient's respiratory force, determined by measuring the maximum respiratory pressure (MIP) with a manometer. The MIP is the maximum negative pressure that the patient can generate during the inhalation phase.

CONTRAINDICATIONS

The device should not be used on:

- Patients who have recently undergone abdominal surgery and patients with abdominal hernia.
- Patients with asthma who have very low perception of symptoms and suffer from frequent and severe exacerbations or with abnormally low perception of dyspnea.
- Patients with a ruptured eardrum or any other ear disease.
- Patients with markedly elevated left ventricular end-diastolic pressure and volume.
- Patients with worsening signs and symptoms of heart failure after RMT/IMT.

The following conditions are also highlighted in order to request medical advice before using the device:

- People with severe pneumothorax or who cough up large amounts of blood (massive hemoptysis).
- Patients who have suffered from a ruptured eardrum or any other ear disease
- Patients with exceptionally high left ventricular volume and end-diastolic pressure
- Patients with worsening signs and symptoms of heart failure after the IMT session
- Pulmonary hypertension
- Large bullae on chest X-ray
- Marked osteoporosis with a history of rib fractures
- Desaturation during or after IMT (<94%)

The device is not recommended for use immediately after lung surgery and in subjects with a history of pneumothorax who require blood pressure monitoring.

If during therapy you experience severe tiredness, shortness of breath or a marked increase in heart rate, stop treatment and contact your doctor.

Inform your doctor or therapist if you notice any side effects not listed in these instructions for use.

INTENDED USERS

The devices are intended to be used by professional medical personnel. The device can be used directly by the patient.

GROUP OF TARGET PATIENTS

The device can be used by:

- Adults
- Persons under 16 years of age, only under the supervision of an adult

The device cannot be used with newborns.

The frequency and duration of use are determined by professional medical personnel based on the individual needs of the patient. Before using the device, the user manual must be carefully read and an adult responsible for safety must be present if the device is intended for use by children or persons with limited capabilities (for example: physical, mental or sensory). It is up to the medical personnel to evaluate the patient's conditions and capabilities to determine whether the patient is able to operate independently or is not able to safely use the device independently, and therefore whether the practices of use must be performed by a responsible person. The evaluation of the use of the device on particular types of patients such as pregnant women, breastfeeding women, subjects with limited physical and mental abilities is referred to professional medical personnel.

OPERATING ENVIRONMENT

This device can be used in healthcare facilities such as hospitals, clinics, etc., or even at home.

WARNINGS

- This product is not intended to diagnose, treat, cure, or prevent any disease.
- The following side effects may occur: Dizziness.
- Use the device only as per its intended use. This medical device is not intended to be a life-saving device. Any other use is considered improper and may be dangerous. The Manufacturer is not responsible for any improper use.
- This product must be used only as described in these instructions for use.
- The time required to pass from storage to operating conditions is approximately 2 hours.
- Please keep this manual carefully for future reference.

- When the device is used in combination with another medical device, the instructions for use of the latter must also be followed.
- The materials used in the device are biocompatible materials and comply with mandatory regulations, however it is not possible to completely exclude possible allergic reactions.
- The device does not administer drugs or substances.
- At the end of each use, after cleaning as per the HYGIENIC PREPARATION section, store the device in a dry place, protected from dust.
- It is recommended that you do not use the device if you suffer from cold, sinusitis or respiratory tract infection. In case of doubt, please consult your doctor.
- In case of cold, some users may experience mild ear discomfort during therapy, caused by inadequate equalization of pressure between the mouth and ears. If symptoms persist, consult your doctor.
- When using the device, you should feel a resistance to inhalation, but it should not be painful. If you feel any pain while using the device, stop training immediately and contact your doctor.

Choking Hazard:

- Some components of the device are small enough to be swallowed by children; therefore, keep the device out of the reach of children.

Risk of ineffectiveness of therapy:

- Use the device in the correct position (see instructions for use paragraph)
- Use original Flaem spare parts. Manufacturer declines any responsibility if non-original spare parts or accessories are used.

Risk of infection:

- Before and after each use, follow the hygienic preparation operations. Make sure that the device is not stored near other accessories or devices for different therapies (e.g. infusions).

FAILURE CASE STUDIES AND THEIR RESOLUTION

Problem	Cause	Remedy
The breathing effort is too heavy or too light	Incorrect respiratory load setting	Check the correct position of the indicator according to the medical prescription. If it is correct, consult your doctor/physiotherapist
The indicator is crooked	Incorrect assembly of the indicator on the internal screw	Disassemble the device and assemble it correctly.
When inhaling, no effort is felt.	Missing inhalation valve or sealing gasket	Check for their presence and, if necessary, reassemble the device with the missing parts.

If, after having checked the conditions described above, the appliance still does not function correctly, we recommend that you contact your trusted dealer or an authorised FLAEM service centre nearest you.

Find the list of all the Service Centers on the page <http://www.flaemnuova.it/it/info/associazione>

DISPOSAL

All components of the product can be disposed of with household waste unless this is prohibited by the disposal regulations in force in the respective member countries.

Packaging materials:

Product box

Product packaging bag

COMMUNICATION OF SERIOUS EVENTS

Serious events that occur in relation to this product must be reported immediately to the manufacturer or the competent authority.

An event is considered serious if it causes or may cause, directly or indirectly, death or an unexpected and serious deterioration in the state of health of a person.

SYMBOLS PRESENT ON THE DEVICE OR PACKAGING

Medical Device

Manufacturer

Production date

CE marking ref. regulation
2017/745 EU and subsequent
updates.

Free from phthalates and
bisphenol

Model number

Attention

See instructions for use

Lot code

Unique identifier of the
device

Atmospheric pressure limits

Temperature limits

Humidity limits

TECHNICAL FEATURES

Model: P0124EM

Weight: 0.05 Kg approx.

MEDICAL DEVICES THAT CAN BE USED IN COMBINATION WITH THE DEVICE

The device can be used in combination with the following medical device: AirFeel respiratory gymnastics device, model P0920EM, P0920EM-1 (not supplied) for IMT therapy and PEP / O-PEP therapy simultaneously.

Combination therapy must be prescribed by medical personnel.

The devices can be connected using the specific F.1 connection as illustrated in fig.F.

ENVIRONMENTAL CONDITIONS

Operating conditions:

Room temperature	Between +10°C and +40°C
Relative humidity of the air	Between 10% and 95%
Atmospheric pressure	Between 69 KPa and 106 KPa

Storage and transportation conditions:

Room temperature	Between -25°C and +70°C
Relative humidity of the air	Between 10% and 95%
Atmospheric pressure	Between 69 KPa and 106 KPa

DURATION

Model: P0124EM

The expected average life is 5 years.

APPLIANCE EQUIPMENT AND INFORMATION ON MATERIALS

Please refer to the image on page 1 of this instruction manual.

The equipment of the device includes: **Materials Information of applied parts**

1 - Main body	
2 – Central screw	
3 – Spring tensioner	
4 – Valve holder	
5 – Cap	
6 – Mouthpiece	Polypropylene
7 – Nose clip	
F.1 – Airfeel connector (optional)	
8 - Spring	

INSTRUCTIONS FOR USE

Before and after each application, the product components must be cleaned thoroughly following the instructions in the HYGIENIC PREPARATION section.

Inspiratory load adjustment:

To adjust the inspiratory load, act on the pin C.1 (fig.C): by turning the pin clockwise, the indicator ring moves closer to the mouthpiece, reducing the load; on the contrary, by turning the pin anti-clockwise, the indicator ring moves away from the mouthpiece, and

in this way the workload increases. The inspiratory load can also be adjusted after assembling the mouthpiece, by rotating it in its seat clockwise or anticlockwise as described for pin C.1.

Preparation for the session:

Insert the mouthpiece (6) inside the main body (1) as shown in fig.D.

It is recommended to wear the nose clip (7) on the nostrils, to maximize the therapy by inhaling and exhaling only through the mouth.

Make sure you are sitting or standing and feeling relaxed.

Holding the device in your hand, close the mouthpiece tightly between your lips, as shown in fig. E.

Inhale deeply so that the internal valve opens: when it intervenes, you can easily feel a flow of air passing through the device.

Continue therapy by inhaling and exhaling through the device, without removing it from your mouth.

Notes for therapy:

It is recommended to start with short training sessions, applying inspiratory pressure loading as directed by your doctor or physiotherapist.

Do not change this parameter without first consulting your doctor.

By performing the training consistently, it will be possible to increase the exercise time.

HYGIENIC PREPARATION

To open the device, rotate the cap (5) anticlockwise (fig.A) and remove it from the main body (1) as illustrated in fig.B. Dismantle the device into its individual parts as shown in the assembly diagram.

Warning: Do not remove the valve and the snap ring mounted on the valve holder (4). Do not remove the snap ring mounted on the main body (1).

These instructions also apply to the optional F.1 fitting.

Disinfection

Sanitize the components by choosing one of the methods described below:

method A: Sanitize the components under hot drinking water (about 40°C) with a mild dish detergent (non-abrasive).

method B: Sanitize the components in the dishwasher on a hot cycle (70°C).

method C: Sanitize the components by immersing them in a solution of 50% water and 50% white vinegar, then rinse thoroughly with hot drinking water (about 40°C).

Disinfection

After sanitizing the components, disinfect them by choosing one of the methods described below. Each method can be practiced for a limited number of times.

method A: Obtain an electrolytic chlor-oxidizing disinfectant (active ingredient: sodium hypochlorite), specifically for disinfection, available in all pharmacies.

Execution:

- Fill a container of a suitable size to contain all the individual components to be disinfected with a solution based on drinking water and disinfectant, respecting the proportions indicated on the disinfectant packaging.

- Completely immerse each individual component in the solution, taking care to avoid the formation of air bubbles in contact with the components. Leave the components immersed for the period of time indicated on the disinfectant packaging, and associated with the concentration chosen for the preparation of the solution.

- Recover the disinfected components and rinse them thoroughly with warm drinking water.

- Dispose of the solution according to the instructions provided by the disinfectant manufacturer.

method B: Disinfect the components by boiling in water for 10 minutes; use demineralized or distilled water to avoid limescale deposits.

method C: Disinfect the components with a hot steam bottle vaporizer (not microwave). Perform the process by faithfully following the instructions of the vaporizer. For disinfection to be effective, choose a sterilizer with an operating cycle of at least 6 minutes. If you also want to perform Sterilization, skip to the paragraph Sterilization.

After disinfecting the accessories, shake them vigorously and place them on a paper towel, or alternatively dry them with a jet of hot air (for example a hairdryer).

Sterilization

Equipment: Steam sterilizer with fractionated vacuum and overpressure compliant with EN 13060.

Execution: Package each individual component to be treated in a sterile barrier system or packaging compliant with EN 11607. Insert the packaged components into the steam sterilizer. Perform the sterilization cycle following the instructions for use of the equipment by selecting a temperature of 134°C and a time of 4 minutes.

Conservation: Store sterilized components as per the system instructions for use or in selected sterile barrier packaging.

The sterilization procedure has been validated in accordance with ISO 17665-1

At the end of each use, store the appliance, complete with accessories, in a dry place, protected from dust.

HYGIENIC PREPARATION IN THE HOME ENVIRONMENT

Disinfection

method A

✓

Method B

✓

method C

✓

Disinfection

method A

✓ MAX 300 TIMES

Method B

✓ MAX 300 TIMES

method C

✓ MAX 300 TIMES

HYGIENIC PREPARATION IN CLINICAL OR HOSPITAL ENVIRONMENT

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EN> The warranty terms here specified are valid only in Italy for Italian residents. In all other countries, the warranty will be provided by the local dealer that sold you the unit, in accordance with the applicable laws.

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All rights reserved
FLAEM NEW SpA