Medical Devices

Aspira

IT > MANUALE ISTRUZIONI
EN > OPERATING INSTRUCTIONS MANUAL
DE > BEDIENUNGSANLEITUNG
FR > NOTICE D'UTILISATION
NL > VERTALING VAN DE ORIGINELE GEBRUIKSAANWIJZING
ES > MANUAL DE INSTRUCCIONES
ΕΛ > ΕΓΧΕΙΡΙΔΙΟ ΟΔΗΓΙΩΝ ΧΡΗΣΗΣ
PL > INSTRUKCJA OBSŁUGI
CS > NÁVOD K POUŽITÍ
RU > РУКОВОДСТВО ПО ЭКСПЛУАТАЦИИ

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These operating instructions are provided for the FLAEM aspirator consisting of a depressor unit and the ancillary tubes and secretion collection vessel, in the available configurations as shown in the table below.

tubes		secretion collection vessel	depressor unit
#	#	#	#
A007	A009	ASP-4	P1211EM/11I, P1211EM/13.SI, P1211EM/20I
#	#	#	#
A005		ASP-2	P1211EM/11I, P1211EM/13.SI, P1211EM/20I

INTENDED USE. Medical Device, aspirator for medical and surgical use, for the aspiration of bodily fluids, pharyngeal aspiration and permanent tracheotomy. Its use must be prescribed by a doctor.

INDICATIONS FOR USE. Treatment of pathologies requiring the aspiration of body fluids (saliva, secretions, post tracheotomy fluids, etc.).

⚠ CONTRAINDICATIONS. The device must NOT be used in operating theatres, for chest drainage or on transport vehicles such as ambulances or emergency vehicles.

INTENDED USERS. The devices are intended for use by legally authorised medical personnel/health workers (doctors, nurses, therapists, etc.). The device can be used directly by the patient for home use, but under the guidance of medical personnel.

⚠ TARGET GROUP OF PATIENTS. Adults, children of all ages. Before using the device, the user manual must be carefully read and an adult responsible for safety must be present if the device is intended for use by children of all ages or people with limited capabilities (for example: physical, mental, or sensory). The evaluation on the use of the device on particular types of patients such as pregnant women, breastfeeding women, newborns, incapacitated subjects or those with limited physical capabilities is referred to medical personnel. All this must be done under the guidance of medical personnel. It is up to medical personnel to evaluate the patient's conditions and capabilities to establish, when prescribing the device, whether the patient is able to operate independently or whether the patient is unable to safely use the aspirator independently, and the therapy must be performed by a responsible person.

OPERATING ENVIRONMENT. All devices can be used in healthcare facilities, such as hospitals, outpatient clinics, etc., or even at home.

⚠ WARNINGS CONCERNING POSSIBLE MALFUNCTIONS. • Should your appliance fail to perform, please contact the authorised service centre for clarification. • The manufacturer should be contacted to report problems and/or unexpected events relating to operation and if necessary for clarifications regarding use and/or hygienic preparation. • Please also refer to the case history of faults and their resolution.

⚠ WARNINGS. • For correct operation and to prolong the life of the device, the operating instructions must be strictly followed.

- The operations for preparing the collection jar must be carried out on a stable support surface.
- If the packaging is damaged or opened, contact the distributor or service centre.
- The manufacturer, the seller and the importer are only responsible for safety, reliability and performance if: a) the appliance is used in accordance with the operating instructions b) the electrical installation in the room where the appliance is used is in accordance with the regulations and complies with the laws in force. • It is forbidden to access the device opening in any way. Repairs may only be carried out by authorised FLAEM personnel following the information provided by the manufacturer. Unauthorised repairs invalidate the warranty and may pose a danger to the user.

Suffocation risk: Some components of the appliance are small enough to be swallowed by children, so keep the appliance out of the reach of children.

Strangulation risk: Do not use the supplied pipes and cables outside of their intended use, they could cause a strangulation hazard, take special care for children and people with special difficulties, often these people are not able to correctly assess the dangers.

Fire risk: Do not expose the appliance to temperatures other than those indicated in the chapter "technical specifications"

- The aspirator is only intended for the collection of NON-flammable fluids, not suitable for use in the presence of flammable anaesthetic mixture with air, or with oxygen or nitrous oxide.

Risk of electrocution: - Before first use, and periodically during the life of the product, check the integrity of the structure of the appliance supplied and the power cord to ensure that there is no damage; if it is damaged, do not plug it in and immediately take the product to an authorised service centre or your dealer. - Do not carry out any hygienic preparation when the device is in use on a patient. - In the presence of children and care-dependent persons, the appliance must be used under the close supervision of an adult who has read these operating instructions. - Keep the power cable away from animals (e.g. rodents), otherwise such animals could damage the insulation of the power cable. - Always keep the power cable away from hot surfaces.

- Do not immerse the appliance in water; if this happens, pull the plug immediately. Do not pull out or touch the appliance while immersed in water, unplug it first. Take it immediately to an authorised FLAEM service centre or your dealer.

- Do not wash the appliance under running water or by immersion and keep it away from splashes of water or other liquids.

- The device is equipped with a safety fuse that can be easily inspected in the event of a fault. Before doing so, remove the power supply plug. - The supplied cable is fitted with a protective seal (17a) against the penetration of liquids into the device. Do not separate this seal from the cable.

Risk of ineffectiveness of therapy: - Only use original Flaem accessories or spare parts, no liability is accepted if non-original parts or accessories are used. - Ensure that the connections and closing of the collection vessel are carried out carefully to avoid aspiration leaks. - With the intervention of the protection device (10) aspiration ceases, switch off the appliance, empty the vessel (8/9) and carry out hygienic preparation.

• **Risk of infection:** - Personal use of the accessories, collection vessels and connecting tubes is recommended to avoid risks of infection.

- The manual aspirated flow control (15) is a sterile disposable product and must be replaced each time it is used.
- Check the expiry date on the original packaging of the manual aspirated flow control (15) or on the label on the packaging of the device and check the integrity of the sterile packaging, if the packaging is not intact please contact your Flaem authorised service centre or your dealer. FLAEM NUOVA will not be held liable for damage to the patient related to the deterioration of the above-mentioned sterile packaging due to manipulations carried out by third parties on the original packaging of the entire device.
- The hydrophobic/antiviral and antibacterial filter supplied with Flaem equipment is a disposable device and must be replaced after each use and if the filter becomes saturated. Check the expiry date on the original filter packaging or on the label on the aspirator packaging.
- Follow hygienic preparation before each use.
- Do not tip the vessel while it is connected to the aspirator in operation, as liquid can be sucked into the aspirator and thus damage the pump. If this happens, switch the aspirator off immediately and empty and clean the vessel (send it to an authorised FLAEM service centre).
- **Risk of injury:** - Do not place the medical device on a soft supporting surface such as a sofa, a bed or a tablecloth.
- Always operate it on a hard surface clear of any objects.

WARNINGS ON INTERFERENCE RISKS DURING USE IN DIAGNOSTIC INVESTIGATIONS. • This device is designed to meet the current requirements for electromagnetic compatibility. As far as EMC requirements are concerned, electromedical devices require special care when being installed and used, and are therefore required to be installed and/or used in accordance with the manufacturer's specifications. Risk of potential electromagnetic interference with other devices. Mobile or portable RF radio and telecommunication devices (mobile phones or wireless connections) may interfere with the operation of electromedical devices. The device may be susceptible to electromagnetic interference in the presence of other devices used for specific diagnosis or treatment. For more information visit www.flaemnuova.it.

- The Medical device is not suitable for use in MRI (Magnetic Resonance Imaging) environments.

TROUBLESHOOTING

Before performing any operation, switch off the appliance and unplug the mains cable from the socket.

PROBLEM	CAUSE	REMEDY
The device is not working	Power cable not correctly plugged into the appliance socket or in the mains power socket	- Plug the power cable correctly into the sockets
No aspiration	Collection vessel lid not properly inserted in the vessel	- Insert the lid correctly into the collection vessel
Aspiration failure caused by fluid spillage	Clogged filter	- Replace the filter
Blocked float	Fouling on the float	- Remove the lid from the vessel, dismantle the protective device and remove the float. Then carry out cleaning operations as per the paragraph "HYGIENIC PREPARATION"
Low and/or no vacuum power	Vacuum regulator all open	- Close the regulator completely and check vacuum power
	Protection filter blocked	- Replace the filter
	Obstructed, kinked or disconnected filter and device connection tubes	- Check the condition of the tubes, replace them if obstructed and connect them correctly as per the "CONNECTION DIAGRAM" in this manual
	Vessel lid overflow valve closed or blocked	- Unlock the overflow valve, hold the device upright
	Dirty or obstructed or damaged pump	- Take the device to your dealer or an authorised FLAEM service centre

If, after checking the conditions described above, the device still is not working properly, we recommend that you contact your trusted dealer or an authorised FLAEM service centre nearest to you. You can find a list of all Service Centres at <http://www.flaemnuova.it/it/info/assistenza>

DISPOSAL

Aspiration unit

 In accordance with Directive 2012/19/EC, the symbol on the equipment indicates that the equipment to be disposed of (excluding accessories) is considered as waste and must therefore be subject to 'separate collection'. Therefore, the user must deliver this waste (or have it delivered) to the separate collection centres set up by the local authorities, or hand it over to the retailer when purchasing a new appliance of an equivalent type. Separate waste collection and subsequent treatment, recovery and disposal operations promote the production of equipment from recycled materials and limit the negative environmental and health effects caused by improper waste management. Unauthorised disposal of the product by the user entails the application of the administrative sanctions provided for in the laws transposing Directive 2012/19/EC of the member state or country where the product is disposed of.

Accessories (Secretion collection vessel and piping)

Accessories to be disposed of as general waste after a sanitisation cycle.

Packaging

Product box

Product packaging bag

REPORTING OF SERIOUS INCIDENTS. Please report any serious incidents that occur in relation to this device to the manufacturer and the competent authority of the Member State in which you reside. An event is considered serious if it causes or may cause, directly or indirectly, death or an unexpected and serious deterioration in the state of health of a person.

⚠ INFORMATION ON RESTRICTIONS OR INCOMPATIBILITIES WITH CERTAIN SUBSTANCES.

Interactions: The materials used for contact with secretions are thermoplastic polymers with high stability and chemical resistance (PP, PC, SI). However, interactions cannot be excluded. Therefore, it is recommended to: a) Always avoid prolonged contact of the liquid with vessels or tubes, sanitise immediately after use. b) If abnormal situations occur, e.g. softening or cracking of accessories, quickly terminate operation and replace the materials used. Contact the authorised service centre specifying the method of use.

SYMBOLS ON DEVICE OR PACKAGING

Medical CE marking ref. regulation 2017/745 EU and subsequent updates	Model number	Serial number of the device	Moisture limits
Class II device	Temperature limits	Manufacturer	Unique device identifier
Before use: Caution check instructions for use	Atmospheric pressure limits	Type BF applied part	Greater vacuum
Switched 'ON' When switching off the unit, the switch stops compressor operation on only one of the two power phases.	Batch Code	Alternating current	Less vacuum
Switched 'OFF'	Single use	Production date	Medical device
Attention	Quality mark	See instructions for use	
IP22 Degree of protection of the casing: IP22. (Protected against the entry of solid objects larger than 12mm in diameter. Protected against dripping water at an angle within $\pm 15^\circ$).			

ASPIRATOR TECHNICAL SPECIFICATIONS

#	P1211EM/11 I	P1211EM/13.5 I	P1211EM/20 I
Power supply/power:	230V~ 50Hz 140VA	230V~ 50Hz 210VA	230V~ 50Hz 210VA
Fuse:	1 x T 2A 250V	1 x T 2A 250V	1 x T 2A 250V
Appliance dimensions:	32.5(L) x 12(D) x 23.5(H) cm	32.5 (L) x 12 (D) x 23,5 (H) cm	32.5 (L) x 12 (D) x 23.5 (H) cm
Weight:	2.3 kg	2.4 kg	2.4 kg
Usage:	Continued	Continued	Continued
Applied parts: Type BF applied parts are:	If necessary, a sterile cannula can be attached to the connection of the "MANUAL ASPIRATED FLOW CONTROL"		
Aspiration:	high vacuum / low flow	high vacuum / low flow	high vacuum / high flow
Adjustable vacuum level:	-0.05 to -0.80 bar approx. (accuracy class 2.5)	-0.05 to -0.86 bar approx. (accuracy class 2.5)	-0.05 to -0.82 bar approx. (accuracy class 2.5)
Max. air flow:	11 l/min approx	13.5 l/min approx	20 l/min approx
Noise level (at 1 m):	56 dB (A) (approx.)	56 dB (A) (approx.)	56 dB (A) (approx.)

ENVIRONMENTAL CONDITIONS

Operating conditions:

Ambient temperature	Between 5°C and +40°C
Relative air humidity	Between 10% and 95%
Atmospheric pressure	Between 69 KPa and 106 KPa

Storage and transport conditions:

Ambient temperature	Between -5°C and +45°C
Relative air humidity	Between 10% and 95%
Atmospheric pressure	Between 69 KPa and 106 KPa

DURATION

DEPRESSOR UNIT	SERVICE LIFE	ACCESSORIES	SERVICE LIFE
# 1211EM/11I	1000 hours	connecting tubes # A005	1 year
# P1211EM/13.5I	2000 hours	connecting tube # A007	
# P1211EM/20I	2000 hours	collection vessel # ASP-4	
		collection vessel # ASP-2	

APPLIANCE AND MATERIAL INFORMATION (DEPENDING ON THE MODELS SUPPLIED)

Nomenclature	Model	Information on materials
1) ASPIRATOR	P1211EM/11I P1211EM/13.5I P1211EM/20I	
2) SWITCH		
3) VACUUM GAUGE		
4) VACUUM ADJUSTMENT KNOB		
5) AIR INLET		
6a) CONNECTION PIPE - Ø 13 x 7.5 mm L 1300 mm	A005	Silicone
6b) CONNECTION PIPE - Ø 13 x 7.5 mm L 280 mm		
6c) CONNECTION PIPE - Ø 13 x 7.5 mm L 280 mm	A007	Silicone
6d) CONNECTION HOSE - Ø 13 x 7.5 mm L 2000 mm (Single use)	A009	PVC
7) HYDROPHOBIAL ANTIVIRAL/ ANTIBACTERIAL FILTER (SINGLE USE) (device not manufactured by FLAEM)		
8) COLLECTION VESSEL	ASP-4	Polycarbonate
10) PROTECTIVE DEVICE		Polypropylene+ Silicone
11) VESSEL CLOSURE CAP		Polypropylene + thermoplastic elastomers
12) VESSEL "VACUUM" SOCKET		
13) VESSEL "PATIENT" SOCKET		
14) LID PROTECTION CAPS		
9) COLLECTION VESSEL	ASP-2	Polycarbonate
10) PROTECTIVE DEVICE		Polypropylene+ Silicone
11) VESSEL CLOSURE CAP		Polypropylene + thermoplastic elastomers
12) VESSEL "VACUUM" SOCKET		
13) VESSEL "PATIENT" SOCKET		
14) LID PROTECTION CAPS		
15) MANUAL ASPIRATED FLOW CONTROL - (disposable) (device not manufactured by FLAEM)		
16) ASPIRATOR CANNULA. (NOT SUPPLIED WITH THE DEVICE). The cannula is not supplied and is not available as a spare part for the medical device. Any cannulae to be used with our aspirators must be of a shape adaptable to the 'MANUAL ASPIRATED FLOW CONTROL' connection supplied with the aspirator.		
17) POWER SUPPLY CABLE		
17a) PROTECTION SEAL		
IMPORTANT NOTE: There is an identification label on the packaging, remove it and affix it in the spaces provided on page 2. When replacing the Nebuliser and accessories, perform the same procedure.		


INSTRUCTIONS FOR USE. Before each use, the accessories must be carefully inspected, making sure that there is no dust, encrustation, lumps or liquid inside the connection pipe, the vessel and its cap. They must also be cleaned strictly in accordance with the instructions given in the section 'HYGIENIC PREPARATION'. Personal use of the accessories, collection vessels and connecting tubes is recommended to avoid the risk of infection.

Instructions for aspiration operation on the patient:

1. Connect the accessories, referring to the "CONNECTION DIAGRAM".

⚠ ATTENTION: the filter is hydrophobic and in addition to being antibacterial/antiviral, it also acts as a secondary protection and stops fluids that the primary protection device in the lid of the collection vessel accidentally fails to stop.

Follow the instructions below for the correct assembly of the filter:

2. Plug the power cord (17) into a mains socket corresponding to the voltage of the appliance, at the end of the use of the appliance pull the plug out of the socket, it must be positioned so that it is not difficult to disconnect from the mains.
3. Start up the appliance by actuating the switch (2), the appliance is intended for continuous use. Always use the appliance on a horizontal plane; this will allow the liquid ingress protection device in the aspiration pump to function properly.
4. For more comfortable suction, set the desired vacuum value (bar) via the vacuum regulator (4). Turning the knob clockwise gives a higher vacuum value and turning it anticlockwise a lower vacuum value; these values can be read on the "vacuum gauge" (3). In order to speed up the suction and simplify the HYGIENIC PREPARATION operations, it is recommended to introduce about 400 ml of water into the vessel (8/9).
5. Place your finger on the hole of the manual aspirated flow control (15) and, in pulsed mode, start aspiration on the patient through the cannula (NOT INCLUDED WITH THE DEVICE).
6. When the application is finished, switch off the appliance.
7. Disconnect the connecting pipes and if necessary close the "PATIENT" and "VACUUM" outlets with plugs 14a -14b. During aspiration, place caps 14a and 14b on fixing anchors 14d and 14c.
8. Disconnect the plug and empty and clean the vessel and connecting pipes.

HYGIENIC PREPARATION

Switch off the appliance before each cleaning operation and unplug the mains cable from the socket.

CLEANING THE APPLIANCE

Use only a cloth moistened with antibacterial detergent (non-abrasive and free of solvents of any kind).

COLLECTION VESSEL AND CONNECTION PIPES

• Detach the cannula (16) (NOT INCLUDED WITH THE DEVICE), the manual aspirated flow control (15) and the tube (6b/6d) from the vessel cap (11). Detach the filter (7) from the vessel cap (11), disconnect the tube (6b/6d) from both the filter (7) and the air inlet of the device (5), take the vessel from its stand in an upright position and empty it (at home in the toilet, in a hospital in the biological material container) and clean the vessel, disassembling it as shown below:

STERILISABLE VESSEL

Figure 3

1. Remove the protective device (10) from the vessel closure cap (11)

Figure 5

2. Pull the float (10a) out of the protective device (10)

Figure 6

SINGLE-PATIENT VESSEL

Figure 4

Sanitisation

Before and after each use, sanitise the collection vessel and connection pipes by choosing one of the methods as described below.

(method A): rinse individual components in hot drinking water (approx. 40°C) with mild dishwashing detergent (non-abrasive) or in a dishwasher on a hot cycle.

(method B): soak the individual components in a solution of 50% water and 50% white vinegar. Finally, rinse them thoroughly with warm (about 40°C) drinking water.

(method C): boil the individual components in water for 10 minutes; demineralised or distilled water is preferable to avoid lime deposits.

After sanitising the accessories, shake them vigorously and place them on a paper towel, or alternatively dry them with a hot air jet (e.g. hair dryer).

Disinfection

The disinfection procedure described in this section is to be carried out before using the accessories and is only effective on the components to be treated if every point is followed and only if the components to be treated are sanitised beforehand.

The disinfectant must be of the electrolytic chloride type (active ingredient: sodium hypochlorite), specific for disinfection and available in all pharmacies.

Procedure:

- Fill a container of a suitable size to hold all the individual components to be disinfected with a solution of drinking water and disinfectant, observing the proportions indicated on the packaging of the disinfectant.
- Completely immerse each individual component in the solution, taking care to avoid the formation of air bubbles in contact with the components. Leave the components immersed for the period of time indicated on the disinfectant packaging, and associated with the concentration chosen for the preparation of the solution.
- Recover the disinfected components and rinse them thoroughly with lukewarm drinking water.
- After disinfecting the accessories, shake them vigorously and place them on a paper towel, or alternatively dry them with a hot air jet (e.g. hair dryer).
- Dispose of the solution according to the disinfectant manufacturer's instructions.

Sterilisation

The sterilisation procedure described in this section is only effective on treated components if it is adhered to in all points and only if the components to be treated are sanitised beforehand, and is validated in accordance with ISO 17665-1.

Equipment: Steam steriliser with fractionated vacuum and overpressure compliant with EN 13060.

Procedure: Package each individual component to be treated in sterile barrier system or packaging in accordance with EN 11607. Insert the packed components into the steam steriliser, making sure to keep the vessel (9) upright.

Perform the sterilisation cycle in accordance with the instructions for use of the equipment by selecting a temperature of:

(method A): 134°C and a time of 10 minutes first.

(method B): 121°C and a time of 20 minutes first.

Storage:

After sanitising, disinfection or sterilisation, reassemble the vessel and connecting pipes according to the instructions given in the 'CONNECTION DIAGRAM'.

At the end of each use, store the medical device complete with accessories in a dry and dust-free place.

- The manual aspirated flow control (15) is a sterile disposable product and must be replaced each time it is used.
- The ANTIBACTERIAL HYDROPHOBIC FILTER (7) is a single-use product and must be replaced each time it is used.

Table of planned methods / patient accessories								
	6a	6b	6c	8	9	10	10a	11
HYGIENIC PREPARATION IN THE HOME								
<u>Sanitisation</u>								
method A	✓	✓	✓	✓	✓	✓	✓	✓
method B	✓	✓	✓	✓	✓	✓	✓	✓
method C	✓	✓	✓	✓	✓	✓	✓	✓
<u>Disinfection</u>								
	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	\	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES
HYGIENIC PREPARATION IN A CLINIC OR HOSPITAL								
<u>Sterilisation</u>								
method A	✓ MAX 50 TIMES	✓ MAX 50 TIMES	✓ MAX 50 TIMES	\	\	\	\	\
method B	\	\	\	\	✓ MAX 50 TIMES	✓ MAX 50 TIMES	✓ MAX 50 TIMES	✓ MAX 50 TIMES
✓: planned \: not planned								

VESSEL HOLDER

For ease of use, the vessel holder can be disassembled and reassembled on the side (RIGHT or LEFT) considered most suitable (Figure 7).

Using a tool, unhook the vessel holder from the bottom of the appliance and rotate it to disengage it. Reattach it, proceeding in reverse order, on the chosen side. Make sure that the vessel holder is properly hooked to the bottom of the appliance. This system allows up to two vessel holders to be connected (Figure 8).

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FLAEM NUOVA S.p.A. non risponde di ulteriori estensioni del periodo di garanzia assicurate da terzi.