
# Delphinus



IT > MANUALE ISTRUZIONI  
EN > USER MANUAL  
FR > MODE D'EMPLOI  
NL > GEBRUIKSAANWIJZING  
EL > ΕΓΧΕΙΡΙΔΙΟ ΟΔΗΓΙΩΝ ΧΡΗΣΗΣ  
PL > INSTRUKCJA OBSŁUGI

A = # DELPHINUS F1000

B = # RF7-2

Schema di collegamento - Connection diagram -- Schéma de branchement -  
Verbindungsschema - Διαγράμμα σύνδεσης - Schemat połączeń

**INDICE**

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
## INDEX

Intended use	15
Indications for use	15
Contraindications	15
Intended users	15
Target group of patients	15
Operating environment	15
Warnings concerning possible malfunctions	15
Warnings	15
Warnings on interference risks during use in diagnostic investigations	17
Case history of faults and their resolution	17
Disposal	18
Notification of serious events	18
Symbols on device or packaging	18
Information on restrictions or incompatibilities with certain substances	19
Compressor unit specifications	19
Nebulizer Technical Specifications	19
Technical specifications of device	20
Environmental conditions	20
Duration	20
Equipment and material information	21
Operating Instructions	21
Hygienic preparation	22
Air filtration	25

# Aerosol therapy device

These user instructions are provided for DELPHINUS F1000 and RF7-2 model devices. The FLAEM Aerosol Therapy device consists of a compressor unit (A), a nebuliser and some accessories (B).

## INTENDED USE

Medical device for the administration of medication through inhalation, with inhalation therapy and medication prescribed by a Doctor.

## INDICATIONS FOR USE

Treatment of respiratory diseases. Medications must be prescribed by a doctor who has assessed the patient's general condition.

### **CONTRAINDICATIONS**

- The medical device should NOT be used for patients who are unable to breathe on their own or who are unconscious.
- Do not use the device in anaesthetic or assisted ventilation circuits.

## INTENDED USERS

The devices are intended for use by legally authorised medical personnel/health workers (doctors, nurses, therapists, etc.). The device can be used directly by the patient.

### **TARGET GROUP OF PATIENTS**

Adults, children of all ages, infants. Before using the device, the user manual must be read carefully and an adult responsible for safety must be present if the device is used by infants, children of any age or persons with limited abilities (e.g. physical, mental or sensory). It is up to the medical personnel to assess the patient's condition and capabilities in order to determine, when prescribing therapy, whether the patient is able to operate the aerosol safely on their own or whether the therapy should be administered by a responsible person.

Please refer to medical personnel for the evaluation of the use of the device on particular types of patients such as pregnant women, lactating women, infants, incapacitated persons or persons with limited physical capabilities.

## OPERATING ENVIRONMENT

This device can be used in healthcare facilities, such as hospitals, outpatient clinics, etc., or even at home.

### **WARNINGS CONCERNING POSSIBLE MALFUNCTIONS**

- Should your appliance fail to perform, please contact the authorised service centre for clarification.
- The manufacturer should be contacted to report problems and/or unexpected events relating to operation and if necessary for clarification of use and/or maintenance/hygienic preparation.
- Please also refer to the case history of faults and their resolution.

## WARNINGS

- Use the device only as a therapeutic inhaler. This medical device is not intended as a life-saving device. Any other use is considered improper and can be dangerous. The manufacturer is not liable for any improper use.
- Always consult your general practitioner for identification of treatment.
- Follow the instructions of your doctor or respiratory rehabilitation therapist regarding the type of medicine, dosage and treatment indications.
- If you should experience allergic reactions or other problems while using the device, stop using it immediately and consult your doctor.
- Keep this manual safe for further reference.
- If the packaging is damaged or opened, contact the distributor or service centre.
- Do not expose the device to particularly extreme temperatures.
- Do not place the device near heat sources, in sunlight or in an excessively hot environment.

- The time required to switch from storage to operating conditions is about 2 hours.
- It is forbidden to access the compressor unit opening in any way. Repairs may only be carried out by personnel authorised by the manufacturer. Unauthorised repairs invalidate the warranty and may pose a danger to the user.
- Suffocation risk:
  - Some components of the appliance are small enough to be swallowed by children, so keep the appliance out of the reach of children.
- Strangulation risk :
  - Do not use the supplied connecting hose and cables outside of their intended use, they could cause a strangulation hazard, take special care for children and people with special difficulties, often these people are not able to correctly assess the dangers.
- Fire risk:
  - THIS device is not suitable for use in the presence of an anaesthetic mixture that is flammable with air, or with oxygen or nitrous oxide.
- Risk of electrocution:
  - Before first use, and periodically during the life of the product, check the integrity of the structure of the appliance and the power cord to ensure that there is no damage; if it is damaged, do not plug it in and immediately take the product to an authorised service centre or your dealer.
  - Keep the power cable away from animals (e.g. rodents), otherwise such animals could damage the insulation of the power cable.
  - Always keep the power cable away from hot surfaces.
  - Never obstruct the ventilation slots located on both sides of the compressor unit.
  - Do not handle the compressor unit with wet hands. Do not use the compressor unit in humid environments (e.g. while bathing or showering). Do not immerse the compressor unit in water; if this happens, pull the plug immediately. Do not pull out or touch the compressor unit while immersed in water, unplug it first. Take it immediately to an authorised FLAEM service centre or your dealer.
- Risk of ineffectiveness of therapy:
  - Performance may vary with particular types of drugs (e.g. those with high viscosity or in suspension). For further information, please refer to the package leaflet provided by the drug manufacturer.
  - Use the nebuliser in the correct position, as upright as possible; do not tilt the nebuliser beyond an angle of 30 degrees, in any direction, to prevent the drug from spilling into the mouth or being overdispersed, reducing the effectiveness of the treatment.
  - Pay attention to the indications provided with the drug and avoid using the devices with substances and dilutions other than those recommended.
  - Only use the device in a dust-free environment, otherwise therapy may be impaired.
  - Do not obstruct or insert objects into the filter and its housing in the unit.
  - In the case of substances that are too dense, dilution with suitable saline solution may be necessary, as prescribed by a doctor.
  - Only use original Flaem accessories or spare parts, no liability is accepted if non-original parts or accessories are used.
- Risk of infection :
  - We recommend personal use of the accessories to avoid any risk of infection.
  - Follow hygienic preparation before each use. Ensure that tubing and accessories are not stored in the vicinity of other accessories or devices for different therapies (e.g. infusions).
  - At the end of therapy, do not leave the medicine inside the nebuliser and proceed with the hygienic preparation operations.
  - If the nebuliser is used for several types of medicine, the residues must be completely removed. Therefore, perform hygienic preparation after each inhalation, also to achieve the highest degree of hygiene and to optimise the service life and operation of the device.
- Risk of injury:
  - Do not place the device on a soft supporting surface such as a sofa, a bed or a tablecloth.
  - Always operate it on a hard surface clear of any objects.

## WARNINGS ON INTERFERENCE RISKS DURING USE IN DIAGNOSTIC INVESTIGATIONS

This device is designed to meet the current requirements for electromagnetic compatibility. As far as EMC requirements are concerned, electromedical devices require special care when being installed and used, and are therefore required to be installed and/or used in accordance with the manufacturer's specifications. Risk of potential electromagnetic interference with other devices. Mobile or portable RF radio and telecommunication devices (mobile phones or wireless connections) may interfere with the operation of electromedical devices. The device may be susceptible to electromagnetic interference in the presence of other devices used for specific diagnosis or treatment. For more information visit [www.flaemnuova.it](http://www.flaemnuova.it).

## CASE HISTORY OF FAULTS AND THEIR RESOLUTION

Before performing any operation, switch off the appliance and disconnect it from the mains.

Problem	Cause	Remedy
The device is not working	Power cable not correctly plugged into the appliance socket or mains power socket	Plug the power cable correctly into the sockets
The device does not nebulise or nebulises poorly	The medicine was not inserted into the nebuliser	Pour the right amount of medicine into the nebuliser
	The nebuliser has not been assembled correctly.	Disassemble and reassemble the nebuliser correctly according to the connection diagram on the cover.
	The nebuliser is obstructed	Hygienic preparation of the nebuliser. Deposits of medicine due to lack of hygienic preparation of the nebuliser impair its efficiency and function. Strictly follow the instructions in the HYGIENIC PREPARATION chapter.
	The accessories are not properly connected to the appliance	Check the correct connection between the air inlet of the appliance and the accessories (see connection diagram on cover).
	The pipe is bent or damaged or twisted	Unwind the pipe and check it for crushing or punctures. Replace it if necessary.
	The air filter is dirty	Replace the filter
The device is louder than usual	The filter is inserted	Insert the filter properly into the housing

If, after checking the conditions described above, the device still is not working properly, we recommend that you contact your trusted dealer or an authorised FLAEM service centre nearest to you. You can find a list of all Service Centres at <http://www.flaemnuova.it/it/info/assistenza>

## DISPOSAL

### Compressor unit

 In accordance with Directive 2012/19/EC, the symbol on the equipment indicates that the equipment to be disposed of (excluding accessories) is considered as waste and must therefore be subject to 'separate collection'. Therefore, the user must deliver this waste (or have it delivered) to the separate collection centres set up by the local authorities, or hand it over to the retailer when purchasing a new appliance of an equivalent type. Separate waste collection and subsequent treatment, recovery and disposal operations promote the production of equipment from recycled materials and limit the negative environmental and health effects caused by improper waste management. Unauthorised disposal of the product by the user entails the application of the administrative sanctions provided for in the laws transposing Directive 2012/19/EC of the member state or country where the product is disposed of.

### Nebuliser and accessories

To be disposed of as general waste after a sanitisation cycle.

### Packaging

Product box

Product packaging bag and tube

Heat-shrinkable nebuliser film and accessories

## NOTIFICATION OF SERIOUS EVENTS

Serious events occurring in connection with this product must be reported immediately to the manufacturer or the competent authority.

An event is considered serious if it causes or may cause, directly or indirectly, death or an unforeseen serious deterioration in a person's state of health.

COUNTRY	AUTHORITY
Ireland	<b>Health Products Regulatory Authority</b> Kevin O'Malley House, Earlsfort Centre, Earlsfort Terrace, IE - Dublin 2 E-mail: <a href="mailto:devicesafety@hpra.ie">devicesafety@hpra.ie</a>
Malta	<b>Medical Devices Unit Medicines Authority</b> Sir Temi Żammit Buildings, Malta Life Sciences Park San Ġwann SĠN 3000 E-mail : <a href="mailto:devices.medicinesauthority@gov.mt">devices.medicinesauthority@gov.mt</a>

## SYMBOLS ON DEVICE OR PACKAGING

 Medical CE marking ref. regulation 2017/745 EU and subsequent updates

 Class II device

 Before use: Caution check instructions for use

 Function switch off

 Function switch on

 Phthalate- and bisphenol-free

 Model number

 Temperature limits

 Atmospheric pressure limits

 Batch Code

Quality mark

 Serial number of the device

 Manufacturer

 Type BF applied part

 Alternating current

 Attention

 See instructions for use

 Medical device

 Moisture limits

 Production date

 Unique device identifier

Protection rating of the envelope: IP21.  
(Protected against solid bodies larger than 12 mm. Protected against access with a finger; Protected against vertically falling drops of water).

**IP21**

## INFORMATION ON RESTRICTIONS OR INCOMPATIBILITIES WITH CERTAIN SUBSTANCES

- Interactions: The materials used in the device are biocompatible materials and comply with statutory regulations, however possible allergic reactions cannot be completely excluded.
- Use the medicine as soon as possible once it has been opened and avoid leaving it in the nebuliser; once therapy has ended, do not leave the medicine inside the nebuliser and proceed with hygienic preparation.

## COMPRESSOR UNIT SPECIFICATIONS

**Model:** DELPHINUS F1000

Power supply:	230V~ 50Hz 140VA
Max. pressure:	2.6 ± 0.4 bar
Air flow to compressor:	10 l/min approx
Noise level (at 1 m):	56 dB (A) approx
Operation:	Continued
Dimension:	18(L) x 30(W) x 10(H) cm

## NEBULIZER TECHNICAL SPECIFICATIONS

**Model:** RF7-2

### RF7 Dual Speed Plus nebuliser

Minimum drug capacity: 2 ml

Maximum drug capacity: 8 ml

## TECHNICAL SPECIFICATIONS OF DEVICE

**Model:** DELPHINUS F1000 combined RF7-2

Operating pressure (with neb.): 0.95 bar.

	Speed selector B2.4	
	in pos Max	in pos Min with valve
<sup>(1)</sup> Release:	0.53 ml/min approx.	0.23 ml/min approx.
<sup>(2)</sup> MMAD:	2.91 $\mu\text{m}$	2.76 $\mu\text{m}$
<sup>(2)</sup> Respirable fraction < 5 $\mu\text{m}$ (FPF):	76%	78.5%

(1) Data collected according to Flaem's internal procedure.  
(2) In vitro characterisation performed by TÜV Rheinland Italia S.r.l. in cooperation with the University of Parma. More details are available on request.

### APPLIED PARTS

Type BF applied parts are: patient accessories (B3, B3.2, B4, B5)

Weight: 2.100 kg

## ENVIRONMENTAL CONDITIONS

### Operating conditions:

Ambient temperature: Between +10°C and +40°C

Relative air humidity: Between 10% and 95%

Atmospheric pressure: Between 69 KPa and 106 KPa

### Storage and transport conditions:

Ambient temperature: Between -25°C and +70°C

Relative air humidity: Between 10% and 95%

Atmospheric pressure: Between 69 KPa and 106 KPa

## DURATION

**Model: DELPHINUS F1000** Service life 1000 hours.  
(Compressor unit)

**Model: RF7-2**  
(Nebuliser and accessories) The expected average life is 1 year, however it is advisable to replace the nebuliser every 6 months during intensive use (or sooner if the nebuliser is clogged) to ensure maximum therapeutic effectiveness.

EQUIPMENT AND MATERIAL INFORMATION		
The equipment includes:	Information on materials	
<b>A -</b>	<b>Compressor unit - Model: DELPHINUS F1000</b> A1 - Switch A2 - Air intake A3 - Air filter A4 - Nebuliser holder A5 - Carrying handle A6 - Power cable A7 - Cable compartment	
<b>B -</b>	<b>Nebuliser and accessories - Model: RF7-2</b> B1 - Connecting tube (compressor/nebuliser unit) B2 - RF7 Dual Speed Plus Nebuliser B2.1 - Lower part B2.2 - Nozzle B2.3 - Upper part B2.4 - Speed selector with valve	
	B3 - Mouthpiece with valve B3.1 - Expiratory valve B3.2 - Non-invasive nasal piece	Polypropylene
	B4 - Adult SoftTouch mask	Polypropylene +
	B5 - Paediatric SoftTouch mask	Thermoplastic Elastomers
	B9 - Manual nebulisation control	Polypropylene
IMPORTANT NOTE: There is an identification label on the packaging, remove it and affix it in the spaces provided on page 2. When replacing the Nebuliser and accessories, perform the same procedure.		

## OPERATING INSTRUCTIONS


**Before each use, wash your hands thoroughly and clean your appliance as described in the section "HYGIENE PREPARATION". During application, it is recommended to adequately protect oneself from drips. This appliance is suitable for the administration of medicinal substances (solutions and suspensions), for which aerosol administration is envisaged; such substances must in any case be prescribed by a doctor. In the case of substances that are too dense, dilution with suitable saline solution may be necessary, as prescribed by a doctor.**

1. Insert the plug of the power cable (A7) into a mains socket corresponding to the voltage of the appliance. It must be positioned in such a way that disconnection from the power grid is not difficult.
2. Insert the nozzle (B2.2) into the upper part (B2.3) by pressing as indicated by the 2 arrows in the 'Connection diagram' in point B2. Insert the Speed selector with valve (B2.4) into the upper part (B2.3) as shown in the 'Connection diagram' in point B2. Pour the medicine prescribed by the doctor into the lower part (B2.1). Close the nebuliser by turning the upper part (B2.3) clockwise.
3. Connect the accessories as shown in the 'Connection Diagram' on the cover.
4. Sit comfortably while holding the nebuliser in your hand, place the mouthpiece to your mouth or use a nasal piece (if provided) or mask. If you use the mask accessory, place it on your face as shown in the figure (with or without the use of the elastic band).
5. Switch on the device by pressing the switch (A1) and inhale and exhale deeply; after inhaling, it is recommended to hold your breath for a moment so that the inhaled aerosol droplets can settle. Then exhale slowly.
6. When the application is finished, switch off the device and disconnect it from the mains.

**ATTENTION:** If, after the therapy session, an obvious deposit of moisture forms inside the tube (B1), detach the tube from the nebuliser and dry it with the ventilation of the compressor itself; this action prevents possible mould growth inside the tube.

To facilitate the connection of the connection tube (B1) to the compressor unit, act on the end of it by simultaneously rotating and inserting it, and for removal by rotating and extracting it.


### WAYS OF USING THE 'RF7 DUAL SPEED PLUS' NEBULISER WITH SPEED SELECTOR AND VALVE SYSTEM

It is professional, fast and suitable for the administration of all types of drugs, including the most expensive ones, even in patients with chronic diseases. Thanks to the geometries of the internal ducts of the RF7 Dual Speed Plus nebuliser, a suitable and active particle size was obtained for treatment down to the lower respiratory tract.


To speed up inhalation therapy, position the speed selector with valve (B2.4) by pressing the MAX marking with your finger.

To make the inhalation therapy more effective, position the speed selector with valve (B2.4) by pressing with a finger on the side opposite the Max marking, in this case you will obtain an optimal administration of the drug while minimising its dispersion in the surrounding environment, thanks to the valve system with which the ampoule, mouthpiece and mask are equipped.


### USE OF MANUAL NEBULISATION CONTROL

In order to achieve continuous nebulisation, it is advisable not to apply the manual nebulisation control (B9), especially in the case of children or care-dependent persons. The manual nebulisation control helps limit the dispersion of the drug into the surrounding environment.


To activate nebulisation, plug the hole of the manual nebulisation control (B9) with a finger and inhale gently and deeply; it is recommended that after inhaling you hold your breath for a moment so that the inhaled aerosol droplets can settle,


In the meantime, to stop nebulisation, remove your finger from the hole of the manual nebulisation control, thus avoiding the waste of medication and optimising its administration. Then exhale slowly.


### SoftTouch masks

**SoftTouch** masks have an outer edge made of **soft, biocompatible material** that ensures an optimal fit on the face, and are also equipped with the **innovative Dispersion Limiter**. These characteristic features allow for greater sedimentation of the drug in the patient, and again limit its dispersion.

Upon inhalation, the tab, which acts as a **Dispersion Limiter**, bends inwards in the mask.

Upon exhalation, the tab, which acts as a **Dispersion Limiter**, bends outwards from the mask.

## HYGIENIC PREPARATION

Switch off the appliance before each hygienic preparation operation and disconnect it from the mains.

### **Compressor unit (A) and pipe exterior (B1)**

Use only a cloth moistened with antibacterial detergent (non-abrasive and free of solvents of any kind).

### **Nebuliser and accessories**

Open the nebuliser by turning the upper part (B2.3) anticlockwise, detach the nozzle (B2.2) and the speed selector with valve (B2.4) from the upper part (B2.3) as shown in the "Connection diagram" in B2.

Then proceed according to the instructions below.

### **Sanitisation**

Before and after each use, sanitise the nebuliser and accessories by choosing one of the methods provided in the table and described below.

**method A:** Sanitise the accessories under warm drinking water (approx. 40°C) with mild dishwashing detergent (non-abrasive).

**method B:** Sanitise the accessories in the dishwasher with the hot cycle (70°C).

**method C:** Sanitise the accessories by soaking in a solution of 50% water and 50% white vinegar, then rinse thoroughly with warm drinking water (approx. 40°C).

After sanitising the accessories, shake them vigorously and place them on a paper towel, or alternatively dry them with a hot air jet (e.g. hair dryer).

### **Disinfection**

After sanitising the nebuliser and accessories, disinfect them using one of the methods provided in the table and described below. Each method can be performed for a limited number of times (see figure in table).

**method A:** Obtain a disinfectant of the electrolytic chloride type (active ingredient: sodium hypochlorite), specifically for disinfection, available in all pharmacies.

#### **Execution:**

- Fill a container of a suitable size to hold all the individual components to be disinfected with a solution of drinking water and disinfectant, observing the proportions indicated on the packaging of the disinfectant.
- Completely immerse each individual component in the solution, taking care to avoid the formation of air bubbles in contact with the components. Leave the components immersed for the period of time indicated on the disinfectant packaging, and associated with the concentration chosen for the preparation of the solution.
- Recover the disinfected components and rinse them thoroughly with lukewarm drinking water.
- Dispose of the solution according to the disinfectant manufacturer's instructions.

**method B:** Disinfect accessories by boiling them in water for 10 minutes; use demineralised or distilled water to avoid hard water deposits.

**method C:** Disinfect the accessories in a hot bottle steamer (not microwave). Carry out the process strictly according to the hot bottle steamer's instructions. For disinfection to be effective, choose a steriliser with an operating cycle of at least 6 minutes.

If you also want to carry out sterilisation, skip to the **Sterilisation** section.

After disinfecting the accessories, shake them vigorously and place them on a paper towel, or alternatively dry them with a hot air jet (e.g. hair dryer).

**At the end of each use, store the device complete with its accessories in a dry and dust-free place.**

### **Sterilisation**

Device: Steam steriliser with fractionated vacuum and overpressure compliant with EN 13060.

**Execution:** Package each individual component to be treated in sterile barrier system or packaging in accordance with EN 11607. Place the packaged components in the steam steriliser. Perform the sterilisation cycle in accordance with the instructions for use of the equipment by selecting a temperature of 134°C and a time of 4 minutes first.

**Storage:** Store sterilised components according to the instructions for use of the system or sterile barrier packaging chosen.

The sterilisation procedure was validated in accordance with ISO 17665-1.

**At the end of each use, store the device complete with its accessories in a dry and dust-free place.**

Table of planned methods / patient accessories

Accessory patient									
Method	B2.1	B2.2	B2.3	B2.4	B3	B3.2	B4	B5	B9
<b>HYGIENIC PREPARATION IN THE HOME</b>									
<b>Sanitisation</b>									
method A	✓	✓	✓	✓	✓	✓	✓	✓	✓
method B	✓	✓	✓	✓	✓	✓	✓	✓	✓
method C	✓	✓	✓	✓	✓	✓	✓	✓	✓
<b>Disinfection</b>									
method A	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES
method B	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES
method C	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES
<b>HYGIENIC PREPARATION IN A CLINIC OR HOSPITAL</b>									
<b>Disinfection</b>									
method A	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES
method B	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES
<b>Sterilisation</b>									
	✓ MAX 300 TIMES Assembled nebuliser				✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES	✓ MAX 300 TIMES
✓: planned    ✗: not planned									


## AIR FILTRATION

The appliance is equipped with an extraction filter (A3) to be replaced with a replacement filter when it is dirty or changes colour. Do not wash or re-use the same filter. Regular replacement of the filter is necessary to help ensure proper compressor performance. The filter must be checked on a regular basis. Contact your dealer or authorised service centre for replacement filters.

**To replace the filter, pull it out as shown in the figure.**

The filter is designed so that it is always fitted in its housing. Do not replace the filter during use. **Only use original Flaem accessories or spare parts, no liability is accepted if non-**

**original parts or accessories are used.**

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**FLAEM NUOVA S.p.A. non risponde di ulteriori estensioni del periodo di garanzia assicurate da terzi.**