
SINCE 1966
MADE IN ITALY

EasyNeb 3

by **FLAEM**®

- I** MANUALE ISTRUZIONI D'USO - Apparecchio per aerosolterapia ad ultrasuoni
- GB** INSTRUCTIONS FOR USE MANUAL - Ultrasound aerosol therapy unit
- F** MODE D'EMPLOI - Appareil pour aérosol thérapie à ultrasons
- NL** GEBRUIKSAANWIJZING - Inhalator voor aerosoltherapie met ultrasonen
- D** GEBRAUCHSANWEISUNG - Ultraschallgerät für Inhalationstherapie
- E** MANUAL DE INSTRUCCIONES - Nebulizador de ultrasonido para aerosolterapia
- P** INSTRUÇÕES DO MANUAL - Aparelho electrónico de tratamento com aerosol

ITALIANO
Pag. 1

ENGLISH
pag. 11

FRANÇAIS
Pag. 21

NEDERLANDS
Pag. 31

DEUTSCH
Pag. 41

ESPAÑOL
PÁG. 51

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Schema di collegamento - Assembly diagram - Schéma de montage - Verbindungsschema
Anschlusschema - Croquis de conexiones - Esquema de ligação

EasyNeb 3

Mod. EASYNEB 3

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Mod. EASYNEB 3

Ultrasound aerosol therapy unit

We are pleased you have purchased our product and we thank you for your trust in us. We aim at fully satisfying our customers by offering them state-of-the-art products for the treatment of respiratory diseases.

Read these instructions carefully and retain them for future reference. Only use the accessory as described in this manual. This is a home medical device to nebulise and administer medication prescribed or recommended by your doctor upon assessing the patient's general conditions. Please note that the full range of Flaem products is visible on the website www.flaem.it.

THE UNIT STANDARD FITTINGS INCLUDE:

- | | |
|--------------------------|---|
| A Nebuliser | B Mouthpiece |
| A1) operation button | C Adult nasal prong |
| A2) indicator light | D Child nasal prong |
| A3) flow adjustment knob | E Cap for cleaning |
| A4) respiratory mask | F Power supply unit |
| A5) air filter | G Transport Case |
| A6) transparent cover | H Spare parts kit
(O-rings for mask, filters and disks) |

IMPORTANT SAFEGUARDS

- This device is also intended for direct use by the patient.
- Before using the unit for the first time, and periodically throughout its life, check the integrity of its structure and the power cable, making sure that there is no damage; if it is damaged, do not insert the plug and take the product immediately to an authorized FLAEM technical assistance centre or to your local dealer. Take the appliance immediately to an authorised service centre or your dealer.
- Should your device fail to provide the expected performance, contact the authorised service centre for clarifications.
- This unit is especially suitable for administering drugs in solution form and can also work with drugs in suspension form. However, in the case of some glycol containing drugs the instructions for use advise against using ultrasound devices, since these drugs cannot be nebulised with them.
- Interactions: the materials used that come into contact with drugs have been tested with a wide range of medications. However, due to the variety and continuous evolution of medication, the possibility of a chemical interaction cannot be excluded. Once opened, the medication should be used as soon as possible, and prolonged contact between drug and nebuliser chamber should be avoided. The materials used in the device are biocompatible in accordance with the provisions of Directive 93/42 EC and subsequent amendments. However, the possibility of occurrence of allergic reactions cannot be entirely excluded.
- Use the unit with a dose of at least 4 ml to be nebulised; if necessary, dilute the drug with physiological solution.
- It is normal for there to be a non-nebulised drug residual of around 0.5 – 0.8 ml at the end of treatment; therefore, when nebulisation becomes discontinuous, stop the treatment as it is pointless to insist on nebulising the remainder, thus causing the unit to overheat.
- Never press the nebuliser button in the absence of medication in order not to overheat the transducer.
- Keep the power cable away from animals (for example, rodents) which could damage the

insulation.


- Use the device only in dust-free conditions, otherwise treatment could be compromised.
- Always keep the unit in an upright position to prevent the medication from spilling out.
- Remove the unit from the carrying bag (G) before use!
- When the adjustment knob is in the "STOP" position, no aerosol is emitted.
- Never operate the device without solution (drug).
- Do not use the tubes and the cords of the device for purposes other than those specified as they may cause strangulation hazard; utmost attention must be paid to children and persons with special needs; usually, such persons are not capable of assessing the hazard correctly.
- The average service life of the device is about 600 hours.
- Children and disabled persons should always use the unit under strict supervision of an adult who has read this manual.
- Several parts of the device are small enough to be swallowed by children; therefore keep the device out of reach of children.
- The unit should not be used in the presence of anaesthetic mixture inflammable with air, oxygen or nitrous oxide.
- The casing of the device and the power supply unit are not protected against liquid penetration. Do not wash the device under running water or submerge in water.
- Do not expose the device to high temperatures.
- Do not put the device and power adaptor near heat sources, direct sunlight, or in hot environments.
- Do not place any objects inside the air slots.
- The power unit provided has been specifically designed for use with the EasyNeb 3 unit. Never use the power unit with other appliances or for uses other than those indicated in this manual and never use EasyNeb 3 with other power units.
- Never use adapters for mains voltages different from those shown on the power supply unit (plate details).
- Always keep the power supply and power cord away from heated surfaces.
- Do not handle the unit with wet hands. Do not use the device in humid environments (e.g. while taking a bath or shower). Do not immerse the device in water; should this accidentally occur, unplug the device immediately. Do not remove or touch the device immersed in water before pulling out the plug. Take the device immediately to an authorized service center or to your local dealer.
- Repairs must be performed by authorized only, by complying with the information provided by the manufacturer. Any unauthorized repairs will void the warranty and may pose a safety hazard for the user.

WARNING: Do not modify this equipment without authorization of the manufacturer.


- The Manufacturer, the Retailer and the Importer are responsible for the safety, reliability and performance of the device only if: a) the device is used in accordance with the operating instructions b) the electrical system of the premises where the device is used is in compliance with current laws.
- The manufacturer must be contacted about any problems and/or unexpected events concerning operation and for any clarifications on use, maintenance/cleaning.
- The amount of time required between preservation conditions and use is approximately 2 hours.

INSTRUCTIONS FOR USE


Before each use, clean hands thoroughly and clean the device as described in the section on "CLEANING SANITISATION DISINFECTION". During use, it is advisable to protect yourself from any dripping. It is recommended that each person use their own nebulizer cup and accessories to prevent risk of infection due to contamination. This device is suitable for the administration of medical substances and not, for which the administration via aerosol is foreseen; these substances are to be in any case prescribed by the Doctor. In case of too thick substances, the dilution with a suitable physiological solution could be needed, according to the medical prescription.


1. Choose a place to comfortably operate the unit. Remove the respiratory mask from the unit. (FIG. 1)


2. Pour the medication into its container; if the medication volume is 2 ml, dilution of the medication with saline solution is recommended to obtain a volume of at least 4 ml. (FIG. 2)


3. Put the respiratory mask back in place (FIG. 3). Make certain that the adjustment knob is in the desired position; (-) for treatment of the lower airways, (+) for treatment of the upper airways.

4. Connect the nebuliser to the power mains using the special power unit as indicated in figure 4. This must be positioned so that it is not difficult to make the disconnection from the mains. After using the device unplug it; the socket must be positioned so that it is not difficult to disconnect from the mains.

5. Sit comfortably holding nebuliser in a vertical position, bring mask near your face and press button (A1) (FIG. 5); the nebuliser led (A2) will show a steady green light and the nebulised medication will start coming out. Breathe smoothly and deeply. To increase the effect of the therapy, after inspiration, hold your breath for a second so the drops of inhaled aerosol are able to deposit along the airways. In order to obtain full medication consumption, release button (A1) each time after breathing in. Then slowly breathe out.

6. Accessories: the unit can be used with a respiratory mask, mouthpiece or a nasal prong by placing these in their specific housings. When using a nasal prong inhale deeply in order to avoid condensation which could lead to an incorrect nebulisation.

ATTENTION: When the medication is about to run out, nebulisation becomes intermittent. Switch off the application by pressing the button (A1). Despite the unit having a sensor indicating the presence of medication, it is always best not to use the appliance dry (without medication).

AUTOMATIC SWITCH-OFF: The appliance switches off automatically after about 10 minutes; nonetheless, if you wish to continue with the therapy, simply press the button again and the appliance will start operating.


CLEANING, SANITIZATION AND DISINFECTION

When treatment is completed remove the power cord from the handset keeping it away from wet or damp areas. Turn the adjustment knob to the "STOP" position.


CLEANING THE OUTSIDE OF THE DEVICE

Only use a cloth dampened with antibacterial detergent (non-abrasive and free of solvents of any kind).


CLEANING THE TRANSDUCER


1. Remove the mask and, if necessary, the nasal prong or the mouthpiece.

2. Pour a spoon (4 ml) of demineralised or drinking water in the transducer area.


3. Put the respiratory mask back in place.


4. Use the supplied cap to close the nebulisation hole located inside the mask (as shown in the picture). The cap should be used only for cleaning purposes.


5. Insert again the power cord and operate the unit for 1-2 minutes.

6. Remove the power cord and the respiratory mask

7. Pour the water off by inclining the nebuliser as shown in the picture.

8. Dry the transducer area with a sterile gauze without pressing the transducer. If the transparent cover (A6) looks wet or dirty, take it out and clean it along with its disk.

SANITIZATION

Before and after each use, sanitize the accessories choosing one of the methods described below

(Method A): The accessories that can be sanitized with the method A are (A4-A6-B-C-D-E)

Sanitize the accessories in warm (approximately 40° C) drinking water with a gentle, non abrasive dish detergent or in the dishwasher with hot cycle.

(Method B): The accessories that can be sanitized with the method B are (A4-A6-B-C-D-E)

Sanitize the accessories by immersion in a solution consisting of 50% water and 50% white vinegar. Finally rinse thoroughly with warm (approximately 40° C) drinking water.

(Method C): The accessories that can be sanitized with the method C are (A4-A6-B-C-D)

Sanitize the accessories by boiling in water for 10 minutes; it is preferable the use of demineralized or distilled water to avoid scale deposits.

After sanitization shake the accessories vigorously and place them on a paper towel, or alternatively, dry them with a jet of hot air (e.g. using a hair dryer).

DISINFECTION

Transducer disinfection

Repeat all the steps previously described, that is from point 1 to point 8 of Cleaning transducer section, replacing water with a chemical disinfectant electrolytic oxidizing chlorine type, (active ingredient: sodium hypochlorite) specifically designed for disinfection and available in all pharmacies, carefully following the relevant instructions. Then rinse with clean water and finally dry the transducer area with a sterile gauze, without exerting unnecessary pressure on the transducer itself.

Accessory disinfection

The accessories that can be disinfected are (A4-A6-B-C-D-E)

The disinfection procedure described in this paragraph is to be performed before using the accessories and is effective only if all the prescriptions herein are met and all the components involved are previously sanitized.

The disinfectant to be used must be of electrolytic oxidizing chlorine type, (active ingredient: sodium hypochlorite), specific for disinfection and is available in all pharmacies.

Procedure:

- Fill a container, suitable to contain all the individual components to be disinfected, with a solution of drinking water and disinfectant, while respecting the proportions indicated on the packaging of the disinfectant itself.


- Completely immerse each individual component in the solution, taking care to avoid the formation of air bubbles in contact with the components. Leave the components immersed for the period of time indicated on the packaging of the disinfectant, and associated with the concentration chosen to prepare the solution.
- Remove the components now disinfected and rinse thoroughly with lukewarm drinking water.

At the end of each use store the device complete with accessories in a dry place away from dust.

AIR FILTERING


On average every 4 months or 70 applications replace the air purification filter (A5) located inside the unit. (Available as spare part, code ACO47). Do not replace the filter during use.

Only use original accessories and spare parts by Flaem, we disclaim any liability in the event of using non original spare parts or accessories.

OPTIONAL ACCESSORIES (TO BE PURCHASED SEPARATELY)


The following accessories can be purchased from your dealer or an authorised after-sales centre.

Cigarette-lighter cord, 12V – Code 04232

Instructions for use of the cigarette-lighter cord in cars, caravans, boats

Plug the unit into the cigarette-lighter socket (12VDC) through the special power supply cord (as shown in the picture).

Paediatric kit – Code ACO85

consisting of: connection tube, child mask and operating instructions

SPARE PARTS

Description

- Mouthpiece
- Respiratory mask
- Adult nasal prong
- Child nasal prong
- Set of 2 air filters

Code

ACO55
ACO63
ACO54
ACO56
ACO47

SYMBOLS

NEBULISER

CE Medical Marking ref. Dir. 93/42 EEC and subsequent amendments. In compliance with EN 60601-1-11

Class II unit

Important: check the operating instructions

PHTHALATES & BPA FREE

In compliance with European standard EN 10993-1 Biological evaluation of medical devices and 93/42/EEC Medical Devices Directive"

Minimum and maximum room temperature

Minimum and maximum atmospheric pressure

Direct current

Serial number of the device

Manufacturer

BF-type fixed part

Protect the device against moisture

Minimum and maximum air moisture

Hazard: electrocution.
Consequence: Death.
Do not use the unit while bathing or showering.

POWER SUPPLY UNIT

Compliant with Directives

Power supply unit for home use

Plug for very low safety voltage

NOTICE (ref. A Directive 2007/47/EC)

Flaem Nuova SpA has been manufacturing top-quality electro-medical appliances for the past 50 years or more, using carefully selected and tested materials in compliance with the requirements of reference directives and standards, with patient/user safety as number one priority. Because of this, we are able to state that the use of the Flaem device and original accessories purchased by you does not involve any risk whatsoever for the user's health.

The specific tests conducted at accredited laboratories (as per UNI EN ISO 10993-17: 2004 rule on Biological Assessment Part 17 Definition of allowable limits for leachable substances) actually demonstrate that the BPA used is released in amounts much lower than the maximum allowed for the safety threshold.

TROUBLESHOOTING

PROBLEM	REMEDY
1. While pressing the handset button the LED is not lit.	<ul style="list-style-type: none"> - When operating through mains supply make sure that the handset and the power supply unit are properly plugged in. - When operating through a cigarette-lighter socket (optional accessory), check that the cord is properly plugged in and the fuse is not burnt out.
2. In the case of poor nebulisation	<ul style="list-style-type: none"> - Make sure that the drug is not too cold – drug temperature should be close to room temperature.
3. If the unit does not nebulise	<ul style="list-style-type: none"> - Make sure the drug is not glycol-based and always remember that the use of an ultrasonic nebulizer is not recommended. - Check that the nebulising chamber contains medication. - Make sure that quantity of medication is between 4 and 8 ml max. - Check that the medication is diluted sufficiently with saline solution 0.9%. - Make sure that the flow adjustment knob is not in the "STOP" position. - Check that the transparent cover (A6) and the disk positioned below are clean.
4. If the appliance interrupts nebulisation and the green LED flashes, this means there is no more medication	<ul style="list-style-type: none"> - Put some more medication into the tray or proceed to clean and disinfect.
5. Failure to clean the nebuliser from drug deposits may endanger proper operation of the unit.	<ul style="list-style-type: none"> - Carefully follow the instructions for CLEANING, SANITIZATION AND DISINFECTION.

Should the unit again fail to operate after following the above instructions, please contact your local dealer or authorised after-sales centre.

TECHNICAL SPECIFICATIONS

NEBULISER

Mod. EASYNEB 3

Power	15V $\overline{=}$
Absorption	850 mA max
Frequency	2,4 MHz
Automatic switch-off:	after 10 minutes
OPERATION	30 MINUTES ON 30 MINUTES OFF
Max chamber capacity	8 ml.
Max. operating speed ⁽¹⁾	0,7ml/min approx
MMAD ⁽²⁾	5,7 μ m
Breathable fractions < 5 μ m ⁽²⁾	40.9%
Noise (at 10 cm) ⁽¹⁾	<50 dB(A)
Dimensions	7x4x21h cm
Weight	280 g
Carrying bag dimensions	20x13x13h cm
Weight of bag (complete with appliance and accessories)	1 Kg

⁽¹⁾Tests performed with phys. sol. 0.9% NaCl according to Flaem internal procedure..

⁽²⁾In vitro testing certified by TÜV Rheinland LGA Products GmbH - Germany in compliance with the European Standard EN 13544-1 for nebulising systems. Further details are available upon request.

Operating conditions:	Temperature . min 10°C; max 40°C Air humidity: min 10%; max 95% RH
Storage conditions:	Temperature . -25°C; max 70°C Air humidity: min 10%; max 95% RH
Operating/storage atmospheric pressure:	min. 690 hPa; max. 1060 hPa P

POWER SUPPLY UNIT

Mod. UME315-1510

Use only power supply unit code 13476

Primary	100-240V ~ 50/60Hz
Secondary	15V $\overline{=}$ 1A
Dimensions:	7x5x7 cm
Weight:	115 g

APPLIED PARTS

Type BF applied parts are: (A4,B,C,D)

ELECTROMAGNETIC COMPATIBILITY

This device was designed to satisfy the currently required requisites for electromagnetic compatibility (EN 60601-1-2:2015). Electro-medical devices require particular care during installation and use relative to EMC requirements. Users are therefore requested to install and/or use these devices following the manufacturer's specifications. There is a risk of potential electromagnetic interference with other devices. RF mobile or portable radio and telecommunications devices (mobile telephones or wireless connections) can interfere with the functioning of electro-medical devices. For further information visit our website www.flaemnuova.it. The Medical Device may be subject to electromagnetic interference if other devices are used for specific diagnosis or treatments. Flaem reserves the right to make technical and functional modifications to the product with no prior warning.

DISPOSAL OF DEVICE:

In conformity with Directive 2012/19/EC, the symbol shown on the device to be disposed of indicates that it is considered as waste and is therefore subject to "sorted waste collection". The user must therefore take (or have taken) the above waste to a pre-sorted waste collection centre set up by the local authorities, or else give it back to the dealer when purchasing a new appliance of the same type. Pre-sorted waste collection and the subsequent treatment, recovery and disposal operations favor the production of appliances made of recycled materials and limit the negative effects of any incorrect waste management on the environment and public health. The unlawful disposal of the product by the user could result in administrative fines as provided by the laws transposing Directive 2012/19/EC of the European member state or of the country in which the product is disposed of.

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© Copyright by FLAEM NUOVA S.p.A., 2018