Medical Devices

# Universal Plus

- IT > MANUALE ISTRUZIONI D'USO
  - EN > INSTRUCTIONS FOR USE MANUAL
  - FR > MODE D'EMPLOI
  - NL > GEBRUIKSAANWIJZING
  - DE > BEDIENUNGSANLEITUNG
  - ES > MANUAL DE INSTRUCCIONES
  - PL > INSTRUKCJA OBSŁUGI
  - EL > ΕΓΧΕΙΡΙΔΙΟ ΟΔΗΓΙΩΝ ΧΡΗΣΗΣ
- < AR كتيب الاستخدام

— 50 —

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MOD. UNIVERSAL PLUS

Ultrasound aerosoltherapy unit

Congratulations on your purchase and thank you for choosing our unit.

Our goal is to fully satisfy customers by offering them cutting-edge systems for the treatment of respiratory tract ailments. Carefully read these instructions and keep them for future reference. Use the unit only as described in this instruction manual. **This is a home medical device to nebulise and administer medication prescribed or recommended by your doctor upon assessing the patient's general conditions.** Please note that the full range of Flaem products is visible on the website [www.flaem.it](http://www.flaem.it).

## THE UNIT STANDARD FITTINGS INCLUDE:

### A - Nebulizer

- A1 - ON/OFF button
- A2 - ON/OFF operation LED
- A3 - Upper part
- A4 - Flow adjustment knob
- A5 - Table support
- A6 - DC plug socket
- A7 - Spare parts kit (nebulizer upper part  
O ring and drug anti dispersion valve)

### C - Accessories

- C1 - Mouthpiece
- C2 - Non-invasive nasal prong
- C3 - Paediatric SoftTouch mask
- C4 - Adult SoftTouch mask
- C5 - Extension tube

### B - Power supply unit

- B1 - Jack

### D - Carry bag

## IMPORTANT WARNINGS


- This device is also intended for direct use by the patient.
- Before using the unit for the first time, and periodically throughout its life, check the integrity of its structure and the power cable, making sure that there is no damage; if it is damaged, do not insert the plug and take the product immediately to an authorized FLAEM technical assistance centre or to your local dealer. Take the appliance immediately to an authorised service centre or your dealer.
- Should your device fail to provide the expected performance, contact the authorised service centre for clarifications.
- This unit is especially suitable for administering drugs in solution form and can also work with drugs in suspension form. However, in the case of some glycol containing drugs the instructions for use advise against using ultrasound devices, since these drugs cannot be nebulised with them.
- Interactions: the materials used that come into contact with drugs have been tested with a wide range of medications. However, due to the variety and continuous evolution of medication, the possibility of a chemical interaction cannot be excluded. Once opened, the medication should be used as soon as possible, and prolonged contact between drug and nebuliser chamber should be avoided.
- The appliance is not suitable for essential oils or balsamic substances.
- Do not run the appliance if it does not contain any solution.
- If the upper part is not correctly fitted in, the appliance does not work; If the upper part is removed during the operation, the appliance switches off automatically.  
To start it again place the upper part in its housing and press button (A1).
- Use the unit with a dose of at least 4 ml to be nebulised; if necessary, dilute the drug with physiological solution.
- It is normal for there to be a non-nebulised drug residual of around 0.5 – 0.8 ml at the end of treatment; therefore, when nebulisation becomes discontinuous, stop the treatment as it is pointless to insist on nebulising the remainder, thus causing the unit to overheat.
- Always keep the unit in an upright position to prevent the medication from spilling out.
- Do not use the tubes and the cords of the device for purposes other than those specified as they may cause strangulation hazard; utmost attention must be paid to children and persons with special needs; usually, such persons are not capable of assessing the hazard correctly.
- The average service life of the device is about 600 hours.
- The Medical Device is maintenance-free, aside from the replacement of accessories.
- Children and disabled persons should always use the unit under strict supervision of an adult who has read this manual.

- Several parts of the device are small enough to be swallowed by children; therefore keep the device out of reach of children.
- The unit should not be used in the presence of anaesthetic mixture inflammable with air, oxygen or nitrous oxide.
- Keep the power unit away from any liquid substance.
- Do not expose the device to high temperatures.
- Do not put the device and power adaptor near heat sources, direct sunlight, or in hot environments.
- The power supply unit included can be connected to the power mains in all countries by means of suitable adapters. Always keep the power supply unit and power cord away from hot surfaces.
- Keep the power cable away from animals (for example, rodents) which could damage the insulation.
- The power unit provided has been specifically designed for use with the unit. Never use the power unit with other appliances or for uses other than those indicated in this manual and never use with other power units.
- Never use adapters for mains voltages different from those shown on the power supply unit (plate details).
- Always keep the power supply and power cord away from heated surfaces.
- Do not handle the unit with wet hands. Do not use the device in humid environments (e.g. while taking a bath or shower). Do not immerse the device in water; should this accidentally occur, unplug the device immediately. Do not remove or touch the device immersed in water before pulling out the plug. Take the device immediately to an authorized service center or to your local dealer.
- Repairs must be performed by authorized FLAEM personnel only by complying with the information provided by the manufacturer. Unauthorized repairs shall void the warranty and may pose a safety hazard for the user.
- ATTENTION: Do not make any changes to this device without manufacturer authorization
- The Manufacturer, the Retailer and the Importer are responsible for the safety, reliability and performance of the device only if: a) the device is used in accordance with the operating instructions b) the electrical system of the premises where the device is used is in compliance with current laws.
- The manufacturer should be contacted for reporting problems and/or unexpected events related to device operation.
- Interactions: the materials used that come into contact with medications have been tested with a wide range of medications. However, due to the variety and continuous evolution of medication, the possibility of interactions cannot be excluded. We recommend that once the medication has been opened you should use it as soon as possible and avoid prolonged contact of the medication with the nebulizer cup. The materials used in the medical device are biocompatible in accordance with the provisions of Directive 93/42 EC and subsequent amendments. However, the possibility of occurrence of allergic reactions cannot be entirely excluded.
- The manufacturer must be contacted about any problems and/or unexpected events concerning operation and for any clarifications on use, maintenance/cleaning.
- The amount of time required between preservation conditions and use is approximately 2 hours.


## INSTRUCTIONS FOR USE

**Before each use, clean hands thoroughly and clean the device as described in the section on "SANITISATION DISINFECTION". During use, it is advisable to protect yourself from any dripping. It is recommended that each person use their own nebulizer cup and accessories to prevent risk of infection due to contamination. This device is suitable for the administration of medical substances and not, for which the administration via aerosol is foreseen; these substances are to be in any case prescribed by the Doctor. In case of too thick substances, the dilution with a suitable physiological solution could be needed, according to the medical prescription.**


1. Remove the upper part (A3) from the unit body (fig. 1)

2. Insert the medication dose (max 6 ml) in the relevant tray (if a syringe is used, remove the needle) (fig. 2); The medication is best diluted to at least 4 ml using saline solution to optimise nebulisation.


3. Reinsert the unit upper part and insert the jack (B1) of the power supply cord into socket (A6) and then connect it to a mains outlet adequate to the device's voltage. The position of the socket must be such that the device can be easily unplugged from the mains network.


4. Connect the accessories as indicated in the "Assembly diagram".

5. Sit comfortably, holding the nebulizer (A) in your hands, place the mouthpiece (C1) over your mouth or use the nasal prong (C2) or mask (C3/C4).


While using the nasal prong (C2) we recommend breathing in strongly to prevent the formation of condensation which would prevent correct nebulisation. Should you opt for the mask (C3/C4) accessory, place it on your face as shown in the picture (fig. 4).

6. To give the dose of medicine to the patient while he/she is lying down (fig. 5), use the extension tube (C5), one end of which must be connected to the mask (C3/C4) while the other must be connected to the appliance suitably housed in the special table stand (A5).

7. Press the ON/OFF (A1) button for a moment (fig. 6); the nebulizer LED (A2) will show a steady green light and the nebulised medication will start to come out.
8. Regulate the supply flow by means of the special knob (A4); turn it clockwise (position +) to treat the upper airways, or anticlockwise (position -) for the lower airways (fig. 6).  
To increase the effect of the therapy, after inspiration, hold your breath for a second so the drops of inhaled aerosol are able to deposit along the airways. Then expire slowly.
9. To stop nebulisation press button (A1).
10. When the medication is about to run out, nebulisation becomes intermittent. Switch off the application by pressing the button (A1). Despite the unit having a sensor indicating the presence of medication, it is always best not to use the appliance dry (without medication).

- AUTOMATIC SWITCH-OFF:** The appliance switches off automatically after about 10 minutes; nonetheless, if you wish to continue with the therapy, simply press the button again and the appliance will start operating.
11. Upon completing application, switch off the device and unplug it.

## SANITIZATION AND DISINFECTION

- Switch off the device before all SANITIZATION AND DISINFECTION procedures and unplug the power cable from the socket.
- Separate the unit (A) from the power, disconnecting the power cable (fig.3).
- Remove the top part (A3) from the body (Fig. 1) and the connected accessory. .

### SANITIZATION

Before and after each use, sanitise the components **A-C1-C2-C3-C4-C5** by choosing one of the methods described below.

**method A:** Sanitise accessories **A-C1-C2-C3-C4-C5** under potable hot water (approximately 40°C) with a gentle, non abrasive dish detergent.

**method B:** Sanitise accessories **C1-C2-C3-C4-C5 (excluding appliance A)** in the dishwasher with a hot cycle.

**method C:** Sanitise the unit and its accessories **A-C1-C2-C3-C4-C5** inside a plastic, porcelain or glass container (not metal); pour a 50% water and 50% white vinegar solution directly on top of these parts. Finally, rinse thoroughly with drinking water.

If you want to also perform the cleaning for DISINFECTION, jump to the DISINFECTION paragraph.

### DISINFECTION

After sanitising the nebuliser cup and the accessories, disinfect them choosing one of the methods described below.

**method A:** Accessories **C1-C2-C3-C4-C5** can be disinfected.

The disinfectant must be an electrolytic chloroxidizer (active principle: sodium hypochlorite) specific for disinfecting, which is available in any pharmacy.

#### Procedure:

- Fill a container, suitable to contain all the individual components to be disinfected, with a solution of drinking water and disinfectant, while respecting the proportions indicated on the packaging of the disinfectant itself.
- Completely immerse each individual component in the solution, taking care to avoid the formation of air bubbles in contact with the components. Leave the components immersed for the period of time indicated on the packaging of the disinfectant, and associated with the concentration chosen to prepare the solution.
- Remove the components now disinfected and rinse thoroughly with lukewarm drinking water.
- Dispose of the solution as recommended by the disinfectant manufacturer.

**method B:** Sanitise the accessories **C1-C2-C3-C4-C5** by boiling them in water for 10 minutes; use demineralised or distilled water to prevent calcium deposits.

**method C:** Sanitise the accessories **C1-C2-C3-C4-C5** with a hot steam steriliser for baby-bottle (not the microwave type). Perform the process faithfully following the instructions of the steriliser. To ensure that the disinfection is effective, choose a steriliser with an operating cycle of at least 6 minutes.

After the CLEANING AND / OR DISINFECTION operations dry the area of the transducer, without

exerting unnecessary pressure on the transducer itself, dry the accessories and the outside of the appliance with a paper towel or even better with a jet of hot air (e.g. hair dryer). Finally, check that the "power plug" area of the device is completely dry.

**At the end of each use store the device complete with accessories in a dry place away from dust.**

## SPARE PARTS

Description	Code
- Extension tube	ACO104T
- Mouthpiece and Nasal prong	ACO044B
- Adult SoftTouch mask	ACO611P
- Paediatric SoftTouch mask	ACO612P

## SYMBOLS

	EC Marking medical ref. Dir. 93/42 EEC and subsequent updates		Device serial number
	Class II device		Manufacturer
	Attention: check the instructions for use		Type BF applied part
	In compliance with: European Standard EN 10993-1 "Biological Evaluation of medical devices" and the European Directive 93/42/EEC "Medical Devices". Phthalate-free. In compliance with: Reg. (EC) no. 1907/2006		Enclosure protection rating: IP21. (Protected against solid bodies over 12 mm. Protected against access with a finger; Protected against vertically falling water drops.)
	Minimum and maximum room temperature		Minimum and maximum air moisture
	Minimum and maximum atmospheric pressure		Power supply for home use
	Direct current		Approvals
	European conformity		Plug for very low safety voltage

### NOTICE (ref. A Directive 2007/47/EC)

Flaem Nuova SpA has been manufacturing top-quality electro-medical appliances for the past 50 years or more, using carefully selected and tested materials in compliance with the requirements of reference directives and standards, with patient/user safety as number one priority. Because of this, we are able to state that the use of the Flaem device and original accessories purchased by you does not involve any risk whatsoever for the user's health.

The specific tests conducted at accredited laboratories (as per UNI EN ISO 10993-17 rule on Biological Assessment Part 17 Definition of allowable limits for leachable substances) actually demonstrate that the BPA used is released in amounts much lower than the maximum allowed for the safety threshold.

## TROUBLESHOOTING

PROBLEM	REMEDY
1. While pressing the handset button the green led is not lit.	- Make sure that the handset and the power supply unit are properly plugged in.
2. In the case of poor nebulisation	- Make sure that the drug is not too cold – drug temperature should be close to room temperature.
3. If the unit does not nebulise	<ul style="list-style-type: none"> <li>- Make sure that the upper part is correctly fitted in.</li> <li>- Check that the nebulizer contains medication.</li> <li>- Make sure that quantity of medication is between 4 and 6 ml max.</li> <li>- Check that the medication is diluted sufficiently with saline solution 0.9%.</li> <li>- Make sure the adjustment knob and the disc underneath are perfectly clean.</li> <li>- Make sure the medication is not glycol based.</li> </ul>
4. If the appliance interrupts nebulisation and the green LED flashes, this means there is no more medication	- Put some more medication into the tray or proceed to clean and disinfect.
5. Failure to clean the nebulizer from drug deposits may endanger proper operation of the unit.	- Carefully follow the instructions for cleaning and disinfection.

Should the unit again fail to operate after following the above instructions, please contact your local dealer or authorised after-sales centre.

## DISPOSAL OF DEVICE:

In conformity with Directive 2012/19/EC, the symbol shown on the device to be disposed of indicates that it is considered as waste and is therefore subject to “sorted waste collection”. The user must therefore take (or have taken) the above waste to a pre-sorted waste collection centre set up by the local authorities, or else give it back to the dealer when purchasing a new appliance of the same type. Pre-sorted waste collection and the subsequent treatment, recovery and disposal operations favor the production of appliances made of recycled materials and limit the negative effects of any incorrect waste management on the environment and public health. The unlawful disposal of the product by the user could result in administrative fines as provided by the laws transposing Directive 2012/19/EC of the European member state or of the country in which the product is disposed of.

## PACKAGING DISPOSAL

Product box

Corrugated tube packing bag and table stand packing bag (where provided)

Power supply packaging box

Heat shrink film accessories

Product packaging bag

## TECHNICAL SPECIFICATIONS

### Mod. UNIVERSAL PLUS NEBULIZER

Power:	15V 
Absorption:	850 mA
Frequency:	2.4 MHz
Operation:	30 min ON/ 30 min OFF
Automatic switch-off :	after 10 minutes
Max tray capacity :	6 ml.
Maximum Delivery rate:	$\geq 0.7 \text{ ml} / \text{min}^{(1)}$
MMAD:	4.15 $\mu\text{m}$ approx <sup>(2)</sup>
Breathable fraction < 5 $\mu\text{m}$ (FPF)	69.61% approx <sup>(2)</sup>

(1) Data shown is in accordance with Flaem Nuova internal procedure I29-P07.5. Values shown for the aerosol delivery may vary depending on the patient's respiratory capacity.

(2) Data calculated based on the values established by Malvern Spraytec and compared with the values of the tests performed on other devices in compliance with the European Standard for aerosol therapy devices Standard EN 13544-1.

Noise (at 10 cm):	<50 dB(A)
Dimensions:	9x5x12h cm
Weight:	200 g
In conformity with:	 0051

### APPLIED PARTS

The BF type applied parts are:	patient accessories (C1, C2, C3, C4)
Operating conditions:	Temperature: min 10°C; max 40°C Air humidity: min 10%; max 95% RH
Storage conditions:	Temperature: -25°C; max 70°C Air humidity: min 10%; max 95% RH
Operating/storage atmospheric pressure:	min. 690 hPa; max. 1060 hPa P

### POWER SUPPLY UNIT

#### Use only power supply unit code 17532

Primary:	100-240V ~ 50/60Hz
Secondary:	15V  1A
Dimensions:	10x7x9 cm
Weight:	510 g

### ELECTROMAGNETIC COMPATIBILITY

This device was designed to satisfy the currently required requisites for electromagnetic compatibility (EN 60601-1-2). Electro-medical devices require particular care during installation and use relative to EMC requirements. Users are therefore requested to install and/or use these devices following the manufacturer's specifications. There is a risk of potential electromagnetic interference with other devices. RF mobile or portable radio and telecommunications devices (mobile telephones or wireless connections) can interfere with the functioning of electro-medical devices. For further information visit our website [www.flaemnuova.it](http://www.flaemnuova.it). The Medical Device may be subject to electromagnetic interference if other devices are used for specific diagnosis or treatments. Flaem reserves the right to make technical and functional modifications to the product with no prior warning.

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© 2020 FLAEM NUOVA®
**FLAEM NUOVA S.p.A. non risponde di ulteriori estensioni del periodo di garanzia assicurate da terzi.**
FLAEM NUOVA S.p.A.