FLAEM®

Koala

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Aerosol therapy device

These user instructions are provided for P0615EM and RF7-2 model devices. The FLAEM Aerosol Therapy device consists of a compressor unit (A), a nebuliser and some accessories (B) and spare parts (D).

INTENDED USE. Medical Device for the administration of medication through inhalation, with inhalation therapy and medication prescribed by a doctor.

OPERATING INSTRUCTIONS. Treatment of respiratory diseases. Medications must be prescribed by a doctor who has assessed the patient's general condition.

⚠ CONTRAINDICATIONS. • The medical device should NOT be used for patients who are unable to breathe on their own or who are unconscious. • Do not use the device in anaesthesia or assisted ventilation circuits.

INTENDED USERS. The devices are intended for use by legally authorised medical personnel/health workers (doctors, nurses, therapists, etc.). The device can be used directly by the patient.

⚠ TARGET GROUP OF PATIENTS. Adults, children of all ages, infants. Before using the device, the user manual must be read carefully and an adult responsible for safety must be present if the device is used by infants, children of any age or persons with limited abilities (e.g. physical, mental or sensory). It is up to the medical personnel to assess the patient's condition and capabilities in order to determine, when prescribing therapy, whether the patient is able to operate the aerosol safely on their own or whether the therapy should be administered by a responsible person.

Medical personnel should be responsible for assessing the use of the device on specific types of patients such as pregnant women, lactating women, infants, incapacitated persons or persons with limited physical capabilities.

OPERATING ENVIRONMENT. This device can be used in healthcare facilities, such as hospitals, outpatient clinics, etc., or even at home.

⚠ WARNINGS REGARDING POSSIBLE MALFUNCTIONS. • Should your appliance fail to perform, contact the authorised service centre for clarification. • The manufacturer should be contacted to report problems and/or unexpected occurrences relating to operation and if necessary for clarification of use and/or maintenance/hygienic preparation. • Please also refer to the case history of faults and their resolution.

WARNINGS. • Use the device only as a therapeutic inhaler. This medical device is not intended as a life-saving device. Any other use is considered misuse and may be dangerous. The manufacturer is not liable for any misuse. • Always consult your doctor for identification of treatment. • Follow the instructions of your doctor or respiratory rehabilitation therapist regarding the type of medication, the dosage and treatment indications. • If you should experience allergic reactions or other problems while using the device, stop using it immediately and consult your doctor. • Keep this manual in a safe place for further reference. • If the packaging is damaged or opened, contact the distributor or service centre. • Do not expose the device to extreme temperatures. • Do not place the device near heat sources, in sunlight or in an excessively hot environment. • The time required to switch from storage to operating conditions is approx. 2 hours. • It is forbidden to access the opening of the compressor unit in any way. Repairs may only be carried out by personnel authorised by the manufacturer. Unauthorised repairs invalidate the warranty and may pose a danger to the user.

WARNINGS FOR USE WITH POWER SUPPLY UNIT. The device must only be used with one of the following power supply units - Model: UME318-1215, Model: MP12M-120150- AG

Risk of suffocation: Some components of the appliance are small enough to be swallowed by children, so keep the appliance out of the reach of children.

Risk of strangulation: Do not use the supplied connecting hose and cables outside of their intended use, they could cause a strangulation hazard, take special care for children and people with special difficulties, often these people are not able to correctly assess the dangers.

Risk of fire: This device is not suitable for use in the presence of an anaesthetic mixture that is flammable with air, or with oxygen or nitrous oxide.

Risk of electric shock: - Before first use, and periodically during the life of the product, check the integrity of the structure of the appliance and the power supply unit cord to ensure that there is no damage; if it is damaged, do not plug it in and immediately take the product to an authorised service centre or your dealer. - Keep the power cable of the power supply unit away from animals (e.g. rodents), otherwise such animals could damage the insulation of the power cable. - Always keep the power supply cable away from hot surfaces. - Do not handle the device with wet hands. Do not use the device in humid environments (e.g. while bathing or showering). Do not immerse the device in water.

Risk of ineffectiveness:- Performance may vary with particular types of drugs (e.g. those with high viscosity or in suspension). For more information, refer to the package leaflet provided by the drug manufacturer. - Use the nebuliser in the correct upright position if possible, do not tilt the nebuliser more than 30 degrees, in any direction, to prevent the medication from spilling into the mouth or being overly dispersed, reducing the effectiveness of the treatment. - Pay attention to the directions supplied with the drug and avoid using the devices with substances and dilutions other than those recommended. - Only use the device in a dust-free environment, otherwise the therapy may be negatively affected. - Do not obstruct or insert objects into the filter and its housing in the device. - In the case of too dense substances, dilution with suitable saline solution may be necessary, as prescribed by a doctor. - Only use original Flaem accessories or spare parts; no liability is accepted if non-original parts or accessories are used.

Risk of infection:- We recommend personal use of the accessories to avoid any risk of infection due to contact.

- Follow hygienic preparation before each use. Ensure that the connecting tube and the accessories are not stored near other accessories or devices for other therapies (e.g. infusions). - At the end of the treatment, do not leave the medicine inside the nebuliser and proceed with hygienic preparation.

- If the nebuliser is used for several types of medicine, the residues must be completely removed. Therefore, perform hygienic preparation after each inhalation session, also to achieve the highest degree of hygiene and to optimise the service life and operation of the device.

Risk of injury:- Do not place the appliance on a soft supporting surface, such as a sofa, bed or tablecloth.

- Always operate it on a hard, unobstructed surface.

WARNINGS ON INTERFERENCE RISKS DURING USE IN DIAGNOSTIC INVESTIGATIONS

This device is designed to meet the current requirements for electromagnetic compatibility. As far as EMC requirements are concerned, electromedical devices require special care when being installed and used, and are therefore required to be installed and/or used in accordance with the manufacturer's specifications. Risk of potential electromagnetic interference with other devices. Mobile or portable RF radio and telecommunication devices (mobile phones or wireless connections) may interfere with the operation of electromedical devices. The device may be susceptible to electromagnetic interference in the presence of other devices used for specific diagnosis or treatment. For more information visit www.flaemnuova.it.

TROUBLESHOOTING. Before performing any operation, switch off the appliance and disconnect it from the mains.

Problem	Cause	Remedy
The device is not working.	The power supply connections were not made correctly.	Make the power supply connections correctly.
The device does not nebulise or nebulises poorly.	The medicine was not inserted into the nebuliser.	Pour the medicine into the nebuliser.
	The nebuliser has not been assembled correctly.	Disassemble and reassemble the nebuliser correctly according to the connection diagram on the cover.
	The nebuliser is obstructed.	Perform the hygienic preparation of the nebuliser. Deposits of medicine due to lack of hygienic preparation of the nebuliser impair its efficiency and function. Strictly follow the instructions in the HYGIENIC PREPARATION chapter.
		Check the correct connection between the air inlet of the appliance and the accessories (see connection diagram on cover).
	The accessories are not properly connected to the appliance.	Unwind the pipe and check it for crushing or punctures. Replace it if necessary.
	The pipe is bent or damaged or twisted.	Replace the filter.
The device is louder than usual.	The air filter is dirty.	Insert the filter properly into the housing.
	The filter is not inserted.	

If, after checking the conditions described above, the device still is not working properly, we recommend that you contact your trusted dealer or an authorised FLAEM service centre nearest to you. You can find a list of all Service Centres at <http://www.flaemnuova.it/it/info/assistenza>

DISPOSAL

Compressor unit

 In accordance with Directive 2012/19/EC, the symbol on the equipment indicates that the equipment to be disposed of (excluding accessories) is considered as waste and must therefore be subject to 'separate collection'. Therefore, the user must deliver this waste (or have it delivered) to the separate collection centres set up by the local authorities, or hand it over to the retailer when purchasing a new appliance of an equivalent type. Separate waste collection and subsequent treatment, recovery and disposal operations promote the production of equipment from recycled materials and limit the negative environmental and health effects caused by improper waste management. Unauthorised disposal of the product by the user entails the application of the administrative sanctions provided for in the laws transposing Directive 2012/19/EC of the member state or country where the product is disposed of.

Nebuliser and accessories. To be disposed of as general waste after a sanitisation cycle.

Packaging		Product box		Product and tube packaging bag
		Heat-shrinkable nebuliser film and accessories		Sack packaging bag

REPORTING SERIOUS EVENTS. Please report any serious incidents occurring in connection with this device to the manufacturer and to the competent authority of the Member State where you reside. An event is considered serious if it causes or may cause, directly or indirectly, death or an unforeseen serious deterioration in a person's state of health.

INFORMATION ON RESTRICTIONS OR INCOMPATIBILITIES WITH CERTAIN SUBSTANCES

• Interactions: The materials used in the device are biocompatible materials and comply with mandatory regulations; however, possible allergic reactions cannot be completely excluded. • Use the medicine as soon as possible once it has been opened and avoid leaving it in the nebuliser; at the end of the treatment, do not leave the medicine inside the nebuliser and proceed with hygienic preparation.

SYMBOLS ON DEVICE OR PACKAGING

 Medical CE marking ref. regulation 0051 2017/745 EU and subsequent updates		IP21 Protection rating of the envelope: IP21. (Protected against solid bodies larger than 12 mm. Protected against access with a finger; Protected against vertically falling drops of water).	
 Before use: Caution check instructions for use		 Extra-low voltage safety plug	
 Class II device	 Manufacturer	 Serial number of the device	
 Type BF applied part	 Phthalate- and bisphenol-free	 See instructions for use	
 Switched "ON" When switching off the unit, the switch stops compressor operation on only one of the two power phases.	 Alternating current	 Model number	
 Switched "OFF"	 Medical device	 Attention	
 Temperature limits	 Moisture limits	 Atmospheric pressure limits	
 Production date	 Batch Code	 Unique device identifier	
 Direct Current	 Plug for safety extra-low voltage	 Power supply unit for domestic use	

NEBULISER TECHNICAL SPECIFICATIONS. Model: RF7-2

RF7 Dual Speed Child Nebuliser

Minimum drug capacity: 2 ml

Maximum drug capacity: 8 ml

COMPRESSOR UNIT TECHNICAL SPECIFICATIONS. Model: P0615EM

Power supply: 12V DC 1.25A

Noise level (at 1 m): 49,5 dB (A) approx

Max. pressure: 1.7 bar approx.

Operation: 30 minutes ON - 30 minutes OFF

Air flow to compressor: 14 l/min approx

Dimension: 18(L) x 16.5(D) x 13(H) cm

TECHNICAL SPECIFICATIONS OF DEVICE. Model: P0615EM combined with RF7-2

	Speed selector B2.4	
	in pos Max	in pos Min
Operating pressure (with neb.): [bar]	0.75	
Delivery rate: [ml/min] approx	0.42 ⁽¹⁾	0.20 ⁽¹⁾
(Aerosol Output) AO: [ml]	0.482 ⁽³⁾	0.311 ⁽²⁾
Percentage of fill volume emitted per minute: [%]	5.62 ⁽³⁾	1.79 ⁽²⁾
(Aerosol output rate) AOR: [ml/min]	0.112 ⁽³⁾	0.036 ⁽²⁾
(Mass Median Aerodynamic diameter) MMAD: [µm]	3.72 ⁽³⁾	1.64 ⁽²⁾
(Geometrical Standard Deviation) GSD:	2.36 ⁽³⁾	1.88 ⁽²⁾
Respirable Fraction (< 5 µm): [%]	63.5 ⁽³⁾	95.8 ⁽²⁾
Fraction of aerosol particles < 2 µm: [%]	23.6 ⁽³⁾	63.3 ⁽²⁾
Fraction of aerosol particles > 2 µm ; < 5 µm: [%]	39.8 ⁽³⁾	32.5 ⁽²⁾
Fraction of aerosol particles > 5 µm: [%]	36.5 ⁽³⁾	4.2 ⁽²⁾
Residual volume: [ml]	1.06 ⁽³⁾	1.16 ⁽²⁾

(1) The delivery rate was measured with a 0.9% Na Cl saline solution at 23°C by weight loss according to Flaem internal procedure. (2)(3) In vitro characterization according to EN ISO 27427 with 4ml, 0.1% Salbutamol in 0.9% NaCl solution with adult ventilatory pattern. Performance are likely to be different for paediatric or infant populations.

(2) Characterization performed by TÜV Rheinland Italia S.r.l. in collaboration with the University of Parma.

(3) Characterization performed by University of Parma.

Further details are available on request.

APPLIED PARTS

Type BF applied parts are:

patient accessories (B3, B4, B5, B6, M1)

Weight:

1.1 Kg

ELECTRICAL CHARACTERISTICS OF THE POWER SUPPLY UNIT. Model: UME318-1215 / Model: MP12M-120150- AG

Input (primary): 100-240V 50-60Hz - Output (secondary): DC 12V 1.5A

ENVIRONMENTAL CONDITIONS

Operating conditions:

Ambient temperature: Between +10°C and +40°C

Relative air humidity: Between 10% and 95%

Atmospheric pressure: Between 69 KPa and 106 KPa

Storage and transport conditions:

Ambient temperature: Between -25°C and +70°C

Relative air humidity: Between 10% and 95%

Atmospheric pressure: Between 69 KPa and 106 KPa

DURATION

Model: P0615EM

(Compressor unit)

Service life 500 hours

Model: RF7-2

(Nebuliser and accessories)

The expected average life is 1 year, however it is advisable to replace the nebuliser every 6 months during intensive use (or sooner if the nebuliser is clogged) to ensure maximum therapeutic effectiveness.

APPLIANCE AND MATERIAL INFORMATION		
The equipment includes:		Information on materials
A -	Compressor unit - Model: P0615EM	
	A1 - Switch, A2 - Air intake, A3 - Air filter A4 - Nebuliser holder, A8 - Power supply unit	
B -	Nebuliser and accessories - Model: RF7-2	
	B1 - Connecting tube (compressor/nebuliser unit)	
	B2 - RF7 Dual Speed Child Nebuliser	
	B2.1 - Lower part, B2.2 - Nozzle	
	B2.3 - Upper part, B2.4 - Speed selector	
	B3 - Mouthpiece	Polypropylene
	B4 - Adult SoftTouch mask	Polypropylene + Thermoplastic Elastomers
C -	B5 - Paediatric SoftTouch mask	
	B6 - Paediatric nasal piece	
	R - Adjustable connection	Polypropylene
	M1 - B015 size 1 (1 to 11 months)	Silicone
	Bag for accessories	
	IMPORTANT NOTE: There is an identification label on the packaging, remove it and affix it in the spaces provided on page 2. When replacing the Nebuliser and accessories, perform the same procedure.	
D -	Replacement filters	

OPERATING INSTRUCTIONS

Before each use, wash your hands thoroughly and clean your appliance as described in the section "HYGIENE PREPARATION". During application, it is recommended to adequately protect oneself from drips. This appliance is suitable for the administration of medicinal substances (solutions and suspensions), for which aerosol administration is envisaged; such substances must in any case be prescribed by a doctor. In the case of substances that are too dense, dilution with suitable saline solution may be necessary, as prescribed by a doctor.

1. Plug the power supply unit (A8) into the mains socket. It must be positioned in such a way that disconnection from the power grid is not difficult. Insert the "L" plug of the power supply cable into the DC plug of the device (FIG. 1).

2. Insert the nozzle (B2.2) in the upper part (B2.3) pressing as shown by the 2 arrows in the "Connection diagram" in point B2. Insert the Speed selector (B2.4) in the upper part (B2.3) as shown in the "Connection diagram" in point B2. Pour the medication prescribed by the doctor into the lower part (B2.1). Close the nebuliser by turning the upper part (B2.3) clockwise.
3. Connect the accessories as shown in the 'Connection Diagram' on the cover.
4. Sit comfortably while holding the nebuliser in your hand, place the mouthpiece to your mouth or use a nasal piece (if provided) or mask. If you use the mask accessory, place it on your face as shown in the (Fig. 2) (with or without the use of the elastic band).
5. Switch on the device by pressing the switch (A1) and inhale and exhale deeply; after inhaling, it is recommended to hold your breath for a moment so that the inhaled aerosol droplets can settle. Then exhale slowly.
6. When the application is finished, switch off the device and disconnect it from the mains.

ATTENTION: If, after the therapy session, an obvious deposit of moisture forms inside the tube (B1), detach the tube from the nebuliser and dry it with the ventilation of the compressor itself; this action prevents possible mould growth inside the tube.

To facilitate the connection of the connection tube (B1) to the compressor unit, act on the end of it by simultaneously rotating and inserting it, and for removal by rotating and extracting it.


HOW TO USE THE "RF7 DUAL SPEED CHILD" NEBULIZER WITH THE SPEED SELECTOR

Professional and fast, this device is suitable for administering all types of medications, including more expensive ones, even in patients with chronic diseases. The geometry of the internal lines of the RF7 Dual Speed Child nebulizer ensures the ideal granulometry for active treatment all the way down to the lower respiratory tract.



To speed up the inhalation therapy, move the speed selector button (B2.4) by pressing on MAX with your finger (Fig. 3).

To increase the effectiveness of the inhalation therapy, press with your finger on the wording MIN of the speed selector (B2.4). In this position, the speed selector acts as a valve and allows nebulizing the optimal amount of drug to the lower respiratory tract, reducing its dispersion in the environment (Fig. 4).

SoftTouch masks

SoftTouch masks have an outer edge made of **soft, biocompatible material** that ensures an optimal fit on the face, and are also equipped with the **innovative Dispersion Limiter**. These characteristic features allow for greater sedimentation of the drug in the patient, and again limit its **dispersion**.

HYGIENIC PREPARATION

Switch off the appliance before each hygienic preparation operation and disconnect it from the mains.

Compressor unit (A) and pipe exterior (B1)

Use only a cloth moistened with antibacterial detergent (non-abrasive and free of solvents of any kind).

Nebuliser and accessories

Open the nebuliser by turning the upper part (B2.3) anticlockwise, detach the nozzle (B2.2) and the speed selector (B2.4) from the upper part (B2.3) as shown in the "Connection diagram" in B2.

Then proceed according to the instructions below.

Sanitisation

Before and after each use, sanitise the nebuliser and accessories by choosing one of the methods provided in the table and described below.

method A: Sanitise the accessories under warm drinking water (approx. 40°C) with mild dishwashing detergent (non-abrasive).

method B: Sanitise the accessories in the dishwasher with the hot cycle (70°C).

method C: Sanitise the accessories by soaking in a solution of 50% water and 50% white vinegar, then rinse thoroughly with warm drinking water (approx. 40°C).

After sanitising the accessories, shake them vigorously and place them on a paper towel, or alternatively dry them with a hot air jet (e.g. hair dryer).

Disinfection

After sanitising the nebuliser and accessories, disinfect them using one of the methods provided in the table and described below. Each method can be performed for a limited number of times (see figure in table).

method A: Obtain a disinfectant of the electrolytic chloride type (active ingredient: sodium hypochlorite), specifically for disinfection, available in all pharmacies.

Execution:

- Fill a container of a suitable size to hold all the individual components to be disinfected with a solution of drinking water and disinfectant, observing the proportions indicated on the packaging of the disinfectant.
- Completely immerse each individual component in the solution, taking care to avoid the formation of air bubbles in contact with the components. Leave the components immersed for the period of time indicated on the disinfectant packaging, and associated with the concentration chosen for the preparation of the solution.
- Recover the disinfected components and rinse them thoroughly with lukewarm drinking water.
- Dispose of the solution according to the disinfectant manufacturer's instructions.

method B: Disinfect accessories by boiling them in water for 10 minutes; use demineralised or distilled water to avoid hard water deposits.

method C: Disinfect the accessories in a hot bottle steamer (not microwave). Carry out the process strictly according to the hot bottle steamer's instructions. For disinfection to be effective, choose a steriliser with an operating cycle of at least 6 minutes.

After disinfecting the accessories, shake them vigorously and place them on a paper towel, or alternatively dry them with a hot air jet (e.g. hair dryer).

At the end of each use, store the device complete with its accessories in a dry and dust-free place.

Table of planned methods / patient accessories						
Accessorio paziente	Sanitisation			Disinfection		
	method A	method B	method C	method A	method B	method C
B2.1	✓	✓	✓	300	300	300
B2.2	✓	✓	✓	300	300	300
B2.3	✓	✓	✓	300	300	300
B2.4	✓	✓	✓	300	300	300
B3	✓	✓	✓	300	300	300
B4	✓	✓	✓	300	300	300
B5	✓	✓	✓	300	300	300
B6	✓	✓	✓	300	300	300
R	✓	✓	✓	300	300	300
M1	✓	✓	✓	300	300	300
✓: planned \: not planned 300: ✓ MAX 300 TIMES						

AIR FILTRATION. The appliance is equipped with an extraction filter (A3) to be replaced with a replacement filter when it is dirty or changes colour. Do not wash or re-use the same filter. Regular replacement of the filter is necessary to help ensure proper compressor performance. The filter must be checked on a regular basis. Contact your dealer or authorised service centre for replacement filters. **To replace the filter, pull it out as shown in the (Fig. 5).** The filter is designed so that it is always fitted in its housing. Do not replace the filter during use. **Only use original Flaem accessories or spare parts, no liability is accepted if non-original parts or accessories are used.**

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FLAEM NUOVA S.p.A. non risponde di ulteriori estensioni del periodo di garanzia assicurate da terzi.