
FLAEM[®] ProfiNeb

Medical Devices

Schema di collegamento - Connection diagram - Schéma de raccordement - Anschlussplan - Aansluitschema - Esquema de conexión - Схема соединения - مخطط التوصيل - Σχεδιάγραμμα σύνδεσης - Schemat połączeniowy - Bağlantı şeması - 연결도식 - 连接图

MANUALE ISTRUZIONI D'USO

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USER INSTRUCTIONS MANUAL

ENGLISH page 10

MODE D'EMPLOI

FRAANÇAIS page 19

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NEDERLANDS pagina 28

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GERMAN Seite 37

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Mod. P0915EM

DEVICE FOR AEROSOL THERAPY

We are pleased you have purchased our product and we thank you for your trust in us. Our goal is our customers' complete satisfaction, offering them cutting edge products to treat respiratory disorders.

Read these instructions carefully and retain them for future reference. Only use the accessory as described in this manual. This is a medical device for home use to nebulise and administer medication prescribed or recommended by your doctor upon assessing the patient's general conditions. Please note that the full range of products is available on the website www.flaemnuova.it.

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- E1 - "L" jack

C - Accessories

- C1 - RF7 Dual Speed Plus Nebuliser
 - C1.1 - Lower part
 - C1.2 - Nozzle
 - C1.3 - Upper part
 - C1.4 - Speed selector with valve
- C2 - Mouthpiece
 - C2.1 - Exhalation valve
- C3 - Paediatric SoftTouch mask
- C4 - Adult SoftTouch mask
- C5 - Clip

IMPORTANT NOTICES

- This device is also intended for direct use by the patient
- Before using the product for the first time and periodically during its service life, check the power supply cable to make sure there is no damage. Should there be damage, do not plug in the cable and immediately bring the product to a FLAEM authorised service centre or to your trusted dealer.
- Should your device fail to provide the expected performance, contact the authorised service centre for clarifications.
- The expected medical life of the accessories is 1 year. It is, however, advisable to replace the nebuliser cup every 6 months in the event of intense use (or earlier if the cup is obstructed) to always guarantee maximum therapeutic efficacy.
- In the presence of children and non self-sufficient individuals, the device must be used under the close supervision of an adult who has read this manual.
- Some parts of the device are small enough to be swallowed by children; keep the device out of the reach

of children.

- Do not use the supplied tubing and cables for anything other than their intended use. These parts could cause a strangling hazard: pay close attention to children and persons with particular difficulties as they are often unable to accurately evaluate danger.
- The apparatus is unsuitable for use in presence of flammable anaesthetic mixture with air, oxygen or nitrous oxide.
- The power supply provided with your device can be connected to the mains power supply in all countries around the world by using suitable adapters. Always keep the power supply and the power cable away from hot surfaces.
- Keep the power supply cable away from animals (for example, rodents) which could damage the insulation.
- Do not handle the power supply with wet hands. Do not use the device in damp environments (for example, while taking a bath or shower). Do not immerse the device in water; in the event of immersion immediately disconnect the plug. Do not remove or touch the immersed device; unplug the power cable first. Immediately bring the device to an authorised FLAEM service centre or to your trusted dealer.
- Use the device only in dust-free conditions, otherwise treatment could be compromised.
- Do not wash the device under running water or by immersion and keep it safe from being sprayed by water or other liquids.
- Do not expose the device to particularly extreme temperatures.
- Keep the device away from sources of heat, direct sunlight or excessively warm locations.
- Always use it on a rigid surface that is clear of obstacles.
- Do not obstruct or introduce objects in the filter and its relevant housing in the device.
- Repairs must be done only by authorised FLAEM personnel only, by complying with the information provided by the manufacturer. Unauthorised repairs void the warranty and may be hazardous for the user.
- The average expected duration for the compressor is 400 hours.
- ATTENTION: Do not modify this device without authorisation from the manufacturer.
- The Manufacturer, the Vendor and the Importer shall be held responsible for safety, reliability and performance only if: a) the device is used in compliance with the instructions for use b) the wiring where the device is being used is in compliance with safety regulations and current laws.
- Interactions: the materials used in contact with medication have been tested with a vast range of medications. However, in view of the variety and continuous evolution of pharmaceuticals, interactions cannot be ruled out. We recommend using the medication as soon as possible once it has been opened and preventing prolonged exposure in the nebuliser cup. The materials used in the device are biocompatible in accordance with the provisions of Directive 93/42 EC and subsequent amendments. However, the possibility of occurrence of allergic reactions cannot be entirely excluded.
- The manufacturer must be contacted about any problems and/or unexpected events concerning operation and for any clarifications on use, maintenance/cleaning.
- The amount of time required between preservation conditions and use is approximately 2 hours.

INSTRUCTIONS USE

Before each use, wash your hands accurately and clean the device as described in the paragraph "CLEANING, SANITISATION AND DISINFECTION". The nebuliser cup and accessories are for personal use only to avoid potential risks of transmitting infection. This device is suitable for the administration of medical substances and not, for which the administration via aerosol is foreseen; these substances are to be in any case prescribed by the Doctor. In case of too thick substances, the dilution with a suitable physiological solution could be needed, according to the medical prescription.

During application, it is advisable to suitably protect yourself from any dripping.


1. Insert the "L" jack (E1) of the power supply cord into socket (A5) and then connect it to a mains outlet adequate to the device's voltage. After using the device unplug it, the position of the socket must be such that the device can be easily unplugged from the mains network.
2. Insert the nozzle (C1.2) in the upper part (C1.3) pressing as shown by the 2 arrows in the "Connection diagram" in point C1. Insert the Speed selector with valve (C1.4) in the upper part (C1.3) as shown in the "Connection diagram" in point C1.
3. Pour the medication prescribed by the doctor into the lower part (C1.1). Close the nebuliser by turning the upper part (C1.3) clockwise.
4. Connect the accessories as indicated in the "Assembly diagram."
5. Sit comfortably holding the nebuliser in your hand, place the mouthpiece onto your mouth or alternatively use the mask. Should you opt for the mask accessory, place it on your face as shown in the picture (with or without using the elasticated strap).
6. Start the device by means of switch (A1) and breathe deeply in and out. After inhaling, we recommend holding your breath for an instant to allow the inhaled drops of aerosol to deposit. Then exhale slowly.
7. Upon completing application, switch off the device and unplug it.

WARNING: If after the therapy session an evident deposit of moisture forms within the pipe (B), detach the pipe from the nebuliser and dry it with the very ventilation from the compressor; this operation prevents possible blooms of mould inside the pipe.


USE METHODS OF THE "RF7 DUAL SPEED PLUS" NEBULISER WITH SPEED SELECTION AND VALVULAR SYSTEM

It is professional, quick, recommended to administer all types of medications, including the costlier ones, even in patients with chronic diseases. Thanks to the geometry of the RF7 Dual Speed Plus nebuliser cup internal ducts, we have obtained a recommended, active particle size to treat even the lower respiratory tract.


For faster inhalation therapy, position the valve speed selector (C1.4) pressing on the word MAX with your finger.


For more effective inhalation therapy, position the valve speed selector (C1.4) pressing on the opposite side of the word Max with your finger. In this case, you have optimal assumption of the medication, minimising leaks into the surrounding environment, thanks to the valvular system with which the nebuliser cup, mouthpiece and mask are equipped.

SOFTTOUCH MASKS

The **SoftTouch** masks have an external edge made in **soft, biocompatible material** that guarantees optimal adherence to the face and are also equipped with the **innovative Dispersion Limiting device**. These distinctive features allow increased sedimentation of the medication within the patient's body, also **limiting its dispersion**.

During inhalation, the tab that works as a **Dispersion Limiting device** bends toward the inside of the mask.

During exhalation, the tab that works as a **Dispersion Limiting device** bends toward the outside of the mask.

CLEANING, SANITISING, DISINFECTING

Switch off the device before any cleaning procedure and unplug the power cable from the socket.

DEVICE (A) AND TUBING EXTERIOR (B)

Use only a damp cloth with antibacterial soap (non-abrasive and with no solvents of any sort).

ACCESSORIES

Open the nebuliser by turning the upper part (C1.3) anticlockwise, remove the nozzle (C1.2) from the upper part (C1.3) as shown in the "Assembly diagram" in point C1.

Then proceed according to the following instructions.

SANITISATION

Before and after each use, sanitise the nebuliser cup and the accessories, choosing one of the methods described below.

method A: Sanitise accessories **C1.1-C1.2-C1.3-C1.4-C2-C3-C4** under potable hot water (approximately 40°C) with a gentle, non abrasive dish detergent.

method B: Sanitise accessories **C1.1-C1.2-C1.3-C1.4-C2-C3-C4** in the dishwasher with a hot cycle.

method C: Sanitise accessories **C1.1-C1.2-C1.3-C1.4-C2-C3-C4** by immersing them in a solution of 50% water and 50% white vinegar, then rinse thoroughly under potable hot water (approximately 40°C).

If you also wish to perform cleaning by DISINFECTING, jump to the DISINFECTING paragraph.

After having sanitised the accessories, shake them vigorously and place them on a paper towel. Alternatively, dry them with a jet of hot air (for example, a hair dryer).

DISINFECTION

After sanitising the nebuliser cup and the accessories, disinfect them choosing one of the methods described below.

method A: Accessories **C1.1-C1.2-C1.3-C1.4-C2-C3-C4** can be disinfected.

The disinfectant must be an electrolytic chlorine oxidiser (active principle: sodium hypochlorite) specific for disinfecting, which is available in any pharmacy.

Directions:

- Fill a container of suitable size to hold all of the parts to be disinfected with a solution of potable water and disinfectant, according to the proportions indicated on the packaging of the disinfectant.
- Completely immerse each part in the solution, taking care to avoid the formation of air bubbles on the parts. Leave the parts immersed for the amount of time indicated on the packaging of the disinfectant associated with the concentration chosen for the solution.
- Remove the disinfected parts and rinse with plenty of warm potable water.
- Dispose of the solution following the instructions provided by the disinfectant manufacturer.

method B: Sanitise the accessories **C1.1-C1.2-C1.3-C1.4-C2-C3-C4** by boiling them in water for 10 minutes; use demineralised or distilled water to prevent calcium deposits.

method C: Sanitise the accessories **C1.1-C1.2-C1.3-C1.4-C2-C3-C4** with a hot steam steriliser for baby-bottle (not the microwave type). Perform the process faithfully following the instructions of the steriliser. To ensure that the disinfection is effective, choose a steriliser with an operating cycle of at least 6 minutes.

After having disinfected the accessories, shake them vigorously and place them on a paper towel. Alternatively, dry them with a jet of hot air (for example, from a hair dryer).

At the end of each use store the device complete with accessories in a dry place away from dust.

AIR FILTERING


The device is quipped with an extraction filter (A3) that must be replaced when it is dirty or changes colour. Do not wash or reuse the same filter. Regular filter replacement is necessary to help ensure proper compressor performance. The filter must be checked regularly.

Contact your retailer or authorised service centre for replacement filters.

Remove the filter as shown in the picture to replace it.

The filter was designed to always remain secured in its seat. Do not replace the filter during use.

Only use original accessories and spare parts by Flaem, we disclaim any liability in the event of using non original spare parts or accessories.

SYMBOLS

EC Marking medical ref. Dir. 93/42 EEC and subsequent updates

Class II device

Attention: check the instructions for use

"ON"

When the device is switched off, the on/off switch stops compressor function in only one phase of the two-phase power supply.

"OFF"

In compliance with: European Standard EN 10993-1 "Biological Evaluation of medical devices" and the European Directive 93/42/EEC "Medical Devices". Phthalate-free. In compliance with: Reg. (EC) no. 1907/2006

Minimum and maximum room temperature

Minimum and maximum atmospheric pressure

Socket for low safety voltage

Direct current

Device serial number

Manufacturer

Type BF applied part

Alternating current

IP21

Degree of protection of the casing: IP21. (Protected against solid foreign objects larger than 12mm. Protected against access with a finger; Protected against vertically dripping water).

Minimum and maximum air moisture

Power supply for home use

Plug for low safety voltage

TECHNICAL SPECIFICATIONS

Mod. P0915EM

Power source:	12V DC 1.25A
Max pressure:	2 ± 0.3 bar
Compressor air delivery:	approx 9 l/min
Noise (at 1 m):	approx 51 dB (A)
Operation:	30 minutes ON - 30 minutes OFF

Operating Conditions:

Temperature:	min 10°C; max 40°C
Air humidity:	min 10%; max 95%
Atmospheric pressure:	min 69KPa; max 106KPa

Storage conditions:

Temperature:	min -25°C; max 70°C
Air humidity:	min 10%; max 95%
Atmospheric pressure:	min 69KPa; max 106KPa

Dimensions:	12(L) x 12(P) x 6(H) cm
Weight:	0.58 Kg

Power supply

Only use power supply mod. UME318-1215 (cod. 17111)	
Primary:	100/240V~ 50/60Hz
Secondary:	12V DC 1.5 A

APPLIED PARTS

Type BF applied parts are:	patient accessories (C2, C3, C4)
----------------------------	----------------------------------

RF7 Dual Speed Plus ampoule

Minimum medication capacity:	2 ml
Maximum medication capacity:	8 ml
Operating pressure (with neb.):	0.75 bar

	C1.4 speed selector	
	Max	Min With valve
⁽¹⁾ Delivery:	approx. 0.42 ml/min	approx. 0.18 ml/min
⁽²⁾ MMAD:	2.44 µm	2.95 µm
⁽²⁾ Breathable fraction < 5 µm (FPF):	80.7 %	74.7 %

(1) data recorded according to internal Flaem procedure I29-P07.5.
 (2) In vitro testing performed by TÜV Rheinland Italia S.r.l. in collaboration with the University of Parma and in compliance with the EN 13544-1: 2007 + A1 European Standard for aerosol therapy devices Standard. More details are available on request.

TROUBLE-SHOTING

Switch off the device before any procedure and unplug the power cable from the socket.

PROBLEM	CAUSE	SOLUTION
The device does not work	The power adaptor has not been correctly inserted into the socket of the device or the power socket	Correctly insert the power adaptor in sockets.
The device does not nebulize or nebulizes insufficiently.	The medication has not been inserted in the nebulizer	Pour the right amount of medication in the nebulizer
	The nebulizer has not been properly fitted.	Disassemble the nebulizer and reassemble it correctly as shown in the connection diagram on the cover
	The nebulizer nozzle is clogged	Disassemble the nebulizer, remove the nozzle and perform cleaning operations Failure to remove medication deposits from the nebulizer affects its efficiency and operation. Strictly comply with the instructions contained in the CLEANING, SANITISATION AND DISINFECTION chapter
	The air pipe is not correctly connected to the device	Ensure the device's accessories are properly connected to the air intake (see connection diagram on the cover).
	The air pipe is bent, damaged or twisted	Unwind the pipe and ensure it is not crushed or punctured. If necessary replace it.
	The air filter is dirty	Replace the filter
The apparatus is noisier than usual	The Filter is not correctly inserted in its housing	Fully and properly insert the Filter in its housing

If after verifying the above mentioned conditions the device should not operate properly, contact your trusted retailer or the nearest authorized FLAEM service center. You can find the list of all the Service Centres at <http://www.flaemnuova.it/it/info/assistenza>

DEVICE DISPOSAL

In compliance with the Directive 2012/19/EC, the symbol printed on the device shows that the device to be disposed of is considered waste and must therefore be an item of "differentiated collection". Consequently, the user must take it (or have it taken) to the designated collection sites provided by the local authorities, or turn it in to the dealer when purchasing an equivalent new device. Differentiated waste collection and the subsequent treatment, recycling and disposal procedures promote the production of devices made with recycled materials and limit the negative effects on the environment and on health caused by potential improper waste management. If the user disposes of the product in an unauthorised fashion, administrative penalties will be applied according to the transposition laws of directive 2012/19/EC of the member state or country where the product will be disposed of.

ELECTROMAGNETIC COMPATABILITY

This device was designed to satisfy the currently required requisites for electromagnetic compatibility (EN 60601-1-2). Electro-medical devices require particular care during installation and use relative to EMC requirements. Users are therefore requested to install and/or use these devices following the manufacturer's specifications. There is a risk of potential electromagnetic interference with other devices. RF mobile or portable radio and telecommunications devices (mobile telephones or wireless connections) can interfere with the functioning of electro-medical devices. For further information visit our website www.flaemnuova.it. The Device may be subject to electromagnetic interference if other devices are used for specific diagnosis or treatments. Flaem reserves the right to make technical and functional modifications to the product with no prior warning.

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FLAEM NUOVA S.p.A. non risponde di ulteriori estensioni del periodo di garanzia assicurate da terzi.