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I

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Neo Spire

Device for respiratory physiotherapy

We are pleased you have purchased our product and we thank you for your trust in us. It is very important that the health care worker and/or the patient read and understand the device use and maintenance information.

INTENDED USE

Neo spire is a device for respiratory physiotherapy specifically designed for the removal of tracheobronchial secretions in a physiological manner without the need to stimulate cough, without irritating the patient's airway and without requiring the need for cooperation from the patient.

It can be used in health care structures like hospitals and in home care. The product should be used by legally authorized medical personnel (doctors, nurses and therapists), by the home therapist, or by the patients themselves under the guidance of qualified medical personnel. It is highly recommended, to ensure proper functioning and to extend the life of the device, that the health professional and/or the patient read and understand the device use and maintenance information.

The device does not require lubrication, is manageable, simple to use, reliable, resistant and quiet.

Note: Use only original FLAEM accessories

OPERATING INSTRUCTIONS

Patients with reduced expectorant capacity

Cystic Fibrosis

Pre and post surgery

- Pre-surgical preparation of COPD patients
- Post-surgical rehabilitation: thoracoabdominal
- Post-surgical rehabilitation: heart surgery
- Post-surgical rehabilitation: ENT

Other groups of patients with secretions:

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CONTRAINDICATIONS

Patients with marked arterial hypotension

DEVICE NOMENCLATURE

Ⓐ - Physiotherapy device (main unit)

- A1 - On/Off switch
- A2 - Air intake
- A3 - Filter Lid
- A4 - Dispenser support
- A5 - Handle for transport
- A6 - Power supply cord socket
- A7 - Power cable

Ⓑ - Connection tube (main unit/dispenser)

Accessories (single-patient)

- C1 - Dispenser
 - C1.1 - Dispenser body
 - C1.2 - Fitting
- C2 - Mouthpiece

⚠ IMPORTANT WARNINGS

The manufacturer makes every effort to ensure that every product is of the highest quality and safety; however, as for any electrical device, basic safety regulations must always be followed in order to avoid harming persons and things.

- Before using the product for the first time, and periodically during its lifetime, check the integrity of the device structure and of the power cable to make sure there is no damage. In the event of damage, do not plug in the cable and immediately take the product to an authorised FLAEM service centre or to your local dealer.
- The average estimated life of accessories is one year. However, in order to assure continued maximum therapeutic effectiveness of **accessories kit, single-patient** use is advisable and it should be replaced every six months during intensive utilisation periods.
- Children and people who are not self-sufficient must use the device under the close supervision of an adult who has read this manual.

- Some parts of the device are small enough to be swallowed by children; therefore, keep the device out of the reach of children.
- Do not use the provided tubing and cables for any other purpose than those specified, as they can cause risk of strangulation. Be particularly attentive with children and individuals with disabilities because they are often unable to correctly assess risk.
- The apparatus is unsuitable for use in presence of flammable anaesthetic mixture with air, oxygen or nitrous oxide.
- Always keep the power cable away from hot surfaces.
- Always keep the power cable away from animals (for example, rodents) which could damage the insulation.
- Do not handle the plug with wet hands. Do not use the device in damp conditions (for example when taking a bath or shower). Do not immerse the device in water; in the event of immersion immediately disconnect the plug. Do not remove or touch the immersed device; unplug the power cable first. Immediately bring the device to an authorised FLAEM service centre or to your trusted dealer.
- The device casing is not waterproof. Do not wash the device under running water or by immersion and keep it safe from being sprayed by water or other liquids.
- Do not expose the device to particularly extreme temperatures.
- Do not place the device near sources of heat, in direct sunlight or in excessively hot rooms.
- Do not obstruct or put objects into the filter or its related housing in the device.
- Never obstruct the air inlets positioned on both sides of the device.
- Always use it on a rigid surface that is clear of obstacles.
- Make sure there is no material obstructing the air vents before each use.
- Do not insert any objects into the ventilation slits.
- Always unplug the power cable immediately after use.
- The device is equipped with an easy to inspect safety fuse in case of malfunction. Before doing so, unplug the power cable. Repairs must be done only by authorised FLAEM personnel. Unauthorised repairs void the warranty and may be hazardous for the user.
- The average operating life of the compressor is 1500 hours approximately.

ATTENTION: Do not modify this device without authorisation from the manufacturer.


- The Manufacturer, the Vendor and the Importer shall be held responsible for safety, reliability and performance only if: a) the device is used in compliance with the instructions for use; b) the wiring where the device is being used is in compliance with safety regulations and current laws.
- The Manufacturer must be contacted to report problems and/or unexpected events related to operation.

INSTRUCTIONS FOR USE


Before each utilisation, wash your hands accurately and clean the device as described in the paragraph “CLEANING, SANITISATION, DISINFECTION AND STERILISATION”. During application it is advisable to suitably protect yourself from any dripping. The dispenser and accessories are for personal use only to avoid potential risks of transmitting infection.

It is recommended to replace them every six months.

1. Insert power supply cord (A7) into socket (A6) and then connect it to a mains outlet adequate to the device voltage. The position of the socket must be such that the device can be easily unplugged from the mains network.

2. Connect one end of the connection pipe (B) to the device. Connect the other end of the pipe (B) to the dispenser (C1).

3. Connect accessories as indicated in the "Connection diagram" on the cover. Sit comfortably holding the dispenser in your hand and put the mouthpiece on your mouth
4. Start the device by means of switch (A1) and breathe deeply in and out.
5. Upon completing application, switch off the device and unplug it. NEO SPIRE can be repeatedly used throughout the day in sessions of about 15 minutes, as prescribed by your doctor.
6. You can insert the dispenser in its support (A4) if a break is needed during the therapy session; this prevents contaminating the mouthpiece by placing it on dirty surfaces.
7. If after the therapy session an evident deposit of moisture forms within the pipe (B), detach the pipe from the nebuliser and dry it with the very ventilation from the compressor; this operation prevents possible blooms of mould inside the pipe.

CLEANING, SANITISATION, DISINFECTION AND STERILISATION.

Switch off the device before any cleaning procedures and always unplug the power cable from the socket.

CLEANING

Device and tubing exterior

Use only a damp cloth with antibacterial soap (non-abrasive and with no solvents of any sort).

Dispenser and accessories

Use of the dispenser and accessories must be personal to avoid any risk of contagious infection. Disassemble the dispenser as shown in the "Connection diagram", sec. C1.

SANITISATION

Before and after each use, sanitise the dispenser and accessories **C1.1-C1.2** and **C2** with potable hot water (approximately 40°C) with a gentle dish detergent or in the dishwasher with a hot cycle.

After having sanitised the accessories, shake them vigorously and place them on a paper towel. Alternatively, dry them with a jet of hot air (for example, a hair dryer).

DISINFECTION

Accessories **C1.1-C1.2** and **C2** can be disinfected

(method A) The disinfection procedure described in this paragraph is to be followed before using the accessories. It is effective on the parts that undergo this treatment only if each step is followed and only if the parts to be treated have previously been sanitised.

The disinfectant must be an electrolytic chloroxidizer (active principle: sodium hypochlorite) specific for disinfecting, which is available in any pharmacy.

Procedure:

- Fill a container big enough to hold all of the parts to disinfect with a solution of potable water and disinfectant, according to the proportions indicated on the packaging of the disinfectant.
- Completely immerse each part in the solution, taking care to avoid the formation of air bubbles on the parts. Leave the parts immersed for the amount of time indicated on the packaging of the disinfectant associated with the concentration chosen for the solution.
- Remove the disinfected parts and rinse abundantly with warm potable water.
- After having disinfected the accessories, shake them vigorously and place them on a paper towel. Alternatively, dry them with a jet of hot air (for example, a hair dryer).
- Dispose of the solution following the instructions provided by the disinfectant manufacturer.

(method B) Disinfect the accessories by boiling in water for 20 minutes. It is required to use de-mineralised or distilled water to prevent calcium deposits.

- After having disinfected the accessories, shake them vigorously and place them on a paper towel. Alternatively, dry them with a jet of hot air (for example, a hair dryer).

STERILISATION (for hospital use)

Accessories **C1.1-C1.2 and C2** can be sterilised

The sterilisation procedure described in this paragraph is effective on the parts that undergo this treatment only if every step is followed and only if the parts to be treated have previously been sanitised. The procedure is validated in its conformity to ISO 17665-1.

Device: Fractionated vacuum overpressure steam steriliser in accordance with EN 13060.

Procedure: Wrap every single part to be treated with a sterile barrier system or packaging in accordance with Norm EN 11607. Place the packed components in the steam steriliser.


Run the sterilisation cycle according to the operating instructions of the device by selecting a temperature of 134°C and a time of 10 minutes first.

Storage: Store the sterilised parts as per the instructions for use of either the sterile barrier system or packaging.

AIR FILTERING

The device is fitted with a suction filter, located under the Filter Lid (A3), which must be replaced when it is dirty or discoloured. To replace it, lift the Filter Lid (A3) and extract the Filter.

Do not wash or reuse the same Filter. Regular filter replacement is necessary to assure proper compressor performance. The Filter must be regularly checked. Contact your retailer or authorised service centre for replacement filters.

TROUBLE-SHOOTING

Switch off the device before any procedure and unplug the power cable from the socket.

PROBLEM	CAUSE	SOLUTION
The device does not switch on.	The power cable has not been correctly inserted into the socket of the device or the power socket	Correctly insert the power cable in sockets.
The device does not work	The dispenser was not assembled correctly.	Disassemble and reassemble the dispenser correctly as shown in the connection diagram on the cover.
	The air pipe is not correctly connected to the device	Check correct connection between the device air inlet and the accessories (see connection diagram on the cover).
	The air pipe is bent, damaged or twisted	Unwind the pipe and ensure it is not crushed or punctured. If necessary replace it.
	The air filter is dirty	Replace the filter
The apparatus is noisier than usual	Filter Lid not correctly inserted in its housing	Fully and properly insert the Filter Lid in its housing
If after verifying the above mentioned conditions the device should not operate properly, contact your trusted retailer or the nearest authorised FLAEM service centre.		

DEVICE SPECIFICATIONS

Mod. P0713EM

Power source:	230 V ~ 50 Hz / 210 VA
Fuse:	1 x T2A 250V
Compressor air delivery:	approx 17 l/min
Noise (at 1 m):	56 dB (A) (approx.)
Operation:	Continuous
Operating conditions:	Temperature min. 10°C; max. 40°C Relative humidity min. 10%; max 95%
Storage conditions:	Temperature min. -25°C; max 70°C Relative humidity min. 10%; max. 95%

Operating/storage atmospheric pressure:	min. 690 hPa; max. 1060 hPa
Device dimensions	26 (L) x 12 (W) x 23.5 (H) cm
Weight	2.2 Kg

APPLIED PARTS

Type BF applied parts are: patient accessories (C2)

SYMBOLS

Class II device

Type BF applied part

Attention: check the instructions for use

Alternating current

IP21 Degree of protection of the casing: IP21.
(Protected against solid foreign objects larger than 12mm. Protected against access with a finger; Protected against vertically dripping water.)

CE Marking medical ref. Dir EEC 93/42 and subsequent updates

"ON" When the device is turned off, the switch interrupts compressor operation only on one of the two power supply phases
"OFF"

Risk: electrocution.
Consequence: death.
Do not use the device while taking a bath or a shower

In compliance with: European Standard EN 10993-1 "Biological Evaluation of medical devices" and the European Directive 93/42/EEC "Medical Devices". Phthalate-free. In compliance with: Reg. (CE) no. 1907/2006

Manufacturer

Device serial number

DEVICE DISPOSAL

In compliance with the Directive 2012/19/EC, the symbol printed on the device shows that the device to be disposed of is considered waste and must therefore be a "separate collection" item. Consequently, the user must take it (or have it taken) to the differentiated collection sites provided by the local authorities, or turn it in to the dealer when purchasing an equivalent new device. Differentiated waste collection and the subsequent treatment, recycling and disposal procedures promote the production of devices made with recycled materials and limit the negative effects on the environment and on health caused by potential improper waste management. The unlawful disposal of the product by the user could result in administrative fines as provided by the laws transposing Directive 2012/19/EC of the European member state or of the country in which the product is disposed of.

ELECTROMAGNETIC COMPATIBILITY

This device was designed to satisfy the currently required requisites for electromagnetic compatibility (EN 60601-1-2:2007). Electro-medical devices require particular care during installation and use relative to EMC requirements. Users are therefore requested to install and/or use these devices following the manufacturer's specifications. There is a risk of potential electromagnetic interference with other devices, in particular with other analysis and treatment devices. RF mobile or portable radio and telecommunications devices (mobile telephones or wireless connections) can interfere with the functioning of electro-medical devices. For further information visit our website www.flaemnuova.it or www.flaem.it

Flaem reserves the right to make technical and functional modifications to the product with no prior warning.

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© Copyright by FLAEM NUOVA S.p.A., 2016