OxyLast[®] Med Easy

Apparecchio per terapia mediante ossigeno e pressione negativa
Therapy device by oxygen and negative pressure

Mod. P0613EM

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Mod. P0613EM

IMPORTANT NOTICES

- The medical device must only be used by Medical personnel. Before using the product for the first time, and periodically during its service life, check the integrity of the device structure and of the power cable to make sure there is no damage. In the event of damage, do not plug in the cable and immediately take the product to an authorised service centre or to your trusted dealer.
- Some parts of the device are small enough to be swallowed by children; keep the device out of the reach of children.
- Do not use the supplied tubing and cables for anything other than their intended use. These parts could cause a strangling hazard: pay close attention to children and persons with particular difficulties as they are often unable to accurately evaluate danger.
- The power supply cable must be positioned so that the device can be easily unplugged from the power mains.
- Always keep the power cable away from hot surfaces.
- Keep the power cable away from animals (for example, rodents) which could damage the insulation.
- To avoid the risk of shock, the device must be connected to a system with adequate ground connection.
- Do not handle the device with wet hands. Do not use the device in damp environments. Do not immerse the device in water; in the event of immersion immediately disconnect the plug. Do not remove or touch the immersed device; unplug the power cable first. Immediately bring the device to an authorised FLAEM service centre or to your trusted dealer.
- The device casing is not waterproof. Do not wash the device under running water or by immersion and keep it safe from being sprayed by water or other liquids.
- **ATTENTION: OxyLast® Med Easy must be used by connecting to a cylinder or to a centralised system or to an O2 concentrator providing hyperbaric O2 (O2 Input: flow MAX 5L/min' and pressure MAX 50Kpa) with a 95 ± 3 % concentration. ALWAYS CLOSE OXYGEN DELIVERY at the end of each treatment.**
- **When the ALARM indicator light turns on (red LED flashing), decrease the suction value using the Suction Adjustment knob. If the problem persists, close the O2 , sitch-off the device and take it to an authorized service center or to your trusted representative.**
- Do not expose the device to particularly extreme temperatures.
- Keep the device away from sources of heat, direct sunlight or excessively warm locations.
- Never obstruct the air vents located on the sides of the device.
- Always use it on a rigid surface that is clear of obstacles.
- Repairs must be done only by authorised FLAEM personnel. Unauthorised repairs void the warranty and may be hazardous for the user.
- The average expected duration for the product is 1000 hours.
- **ATTENTION: Do not modify this device without authorisation from the manufacturer.**
- The Manufacturer, the Vendor and the Importer shall be held responsible for safety, reliability and performance only if: a) the device is used in compliance with the instructions for use b) the wiring where the device is being used is in compliance with safety regulations and current laws.
- The manufacturer must be contacted to report problems and/or unexpected events related to operation.

INTENDED USE:

Device for the treatment, with the combined effect of TPN and hyper-oxygenation, of wounds / injuries and adjuvant to postoperative healing (see operating instructions).

OPERATING ENVIRONMENT:


The medical device must be used in healthcare facilities, such as hospital outpatient clinics, cosmetic medicine centres, NOT at home. It can only be used by Medical personnel.

It must be used with O₂ with features indicated in the user manual of this medical device and generated or loaded from certified cylinders or certified centralised systems which supply oxygen or a O₂ concentrator.

OPERATION PRINCIPLE:


The device has the purpose of simulating a miniature hyperbaric chamber which may be positioned on a limited skin surface (buttocks, abdomen, thighs, etc.). It acts via a cap placed on the injured part to be treated, on which a given vacuum time is generated and immediately afterwards oxygen is injected, via the same tubes.

OxyLast® Med Easy amplifies the benefits of oxygen therapy combined with negative pressure therapy. OxyLast® Med Easy acts on the area to be treated in two STAGES as shown below:



STEP 1

OxyLast® Med Easy creates a negative pressure zone in the area to be treated, as well as acting as negative pressure therapy, promoting optimal oxygen absorption.

STEP 2

OxyLast® Med Easy enriches the internal atmosphere of the cap with oxygen.

In the presence of an abundant exudate, it is therefore necessary to have a gauze, which should not protrude from the edges of the cap.

INDICATIONS FOR USE:

- Acute injuries: trauma, burns;
- Chronic injuries: pressure, decubitus, leg, diabetic ulcers;
- Dehiscid wounds, infected wounds, infected wounds after sternotomy;
- Surgical wounds: skin grafting, flaps, wound bed preparation.

CONTRAINDICATIONS:

The medical device must NOT be used on patients with:

- Bleeding;
- Difficult haemostasis;
- In anticoagulant therapy;
- Bandage near a blood vessel;
- Open abdomen after laparotomy without closure of the abdominal wall.

The medical device must NOT be used in the following conditions:

- Untreated osteomyelitis;
- Pre-abscess lesions;
- Dry eschar, if debridement of the dry necrotic tissue is not performed;
- Coagulation disorders/haematologic diseases;
- Malignant lesions (TPN may lead to cellular proliferation);
- Exposure of vital organs (except cases of open abdomen after laparotomy without abdominal closure);
- Non-stabilised bone fragments;
- Allergy to certain components required for the procedure.

BASIC INFORMATION:


The device is configured with a structure designed to contain the components that generate and manage the aspiration flow and the oxygen flow needed for the therapeutic treatment.

The device is connected to a cylinder or to a centralised system or to an Oxygen concentrator providing hyperbaric O₂ (Oxygen Input: flow MAX 5L/min' and pressure MAX 50Kpa) with a 95 ± 3 % concentration (For the operating instructions with the cylinder or with the centralised system or with an Oxygen concentrator, refer to the user manual or the procedures in place in the facility using the medical device).

The device is connected to the patient via tubes connected into caps or silicone jigs that are placed onto the injured parts thus forming the miniature hyperbaric chambers.

A keypad makes it possible to select cycles, start and switch off the device.

CONTROL KEYPAD (1) AND FUNCTIONS


- 1.1 - **POWER:** green LED on indicates that the appliance is powered on
- 1.2 - **OXYGEN:** blue LED on indicates that oxygen is delivered
- 1.3 - **VACUUM:** blue LED on indicates that the appliance is sucking air (for 30 sec)
- 1.4 - **START:** Switch on key
- 1.5 - **TIME:** Key to set the oxygen delivery time
 - 1.5.1 - blue LED indicates oxygen delivery time
 - MIN - minimum time (10 sec.)

MED - medium time (20 sec.)
MAX - maximum time (30 sec.)

- 1.6 - **STOP**: Switch off key
- 1.7 - **ALARM**: red LED flashing, indicates a suction problem
- 1.8 - **OUT 1**: Air suction / oxygen delivery 1st pipe connection
- 1.9 - **OUT 2**: Air suction / oxygen delivery 2nd pipe connection
- 1.10 - Suction vacuum meter
- 1.11 - **VACUUM**: Suction adjustment knob (Vacuum from 0 to - 125 mbar MAX approx)

STANDARD EQUIPMENT INCLUDES:

- 1 - Control keypad
- 2 - Power cable socket
- 3 - On/off switch
- 4 - Connection of O₂ at 95 ±3%: MAX 5 L/min' - 50Kpa
- 5 - Connection tube :
 - (1x) OxyLast® Med Easy → O₂ cylinder / O₂ centralised system)
 - (1x) Filter → cap
- 6- (a-b-c-d-e-f) (1x) caps in various sizes

- 7 - Connection Tube (1x) : OxyLast® Med Easy **OUT 1**(1.8) /**OUT 2** (1.9) → filter
- 8 - Anti-bacterial filter (1x)
- 9 - Cap tube connection reducer (1x)
- 10 - Suction control connector (1x)
- 11 - Strap L. 420 mm (x 1)
- 12- Strap L. 660 mm (x 1)
- 13 - Strap L. 1200 mm (x1)
- 14 - Power cable

OPERATING INSTRUCTIONS


Prior to any use, accurately wash your hands and wear gloves if required. Clean the device as described in the paragraph "CLEANING, SANITISATION, DISINFECTION AND STERILISATION".

The accessories (5-6a/6b/6c/6d/6e/6f-7-9-10) must be replaced for each new patient.

The anti-bacterial filter (8) is single use.

1 - Insert the power cable (14) in the socket (2) and connect it to a power socket matching the appliance's voltage (see rated data) in a safety-compliant system, press the ON/OFF switch to "I" and the Green POWER light will switch on.

2 - Connect the accessories as indicated in the "Assembly diagram"

3 - Insert tube (6) into socket 1.8 (OUT1) or 1.9 (OUT2).

You may choose to work with one or two caps simultaneously.

When one outlet only is used, e.g. OUT1, the tube of the other one, OUT2, must be disconnected.

4 - Insert the tube (5) to the oxygen connection (4), then connect it to your oxygen dispenser (see the oxygen manual). Enable oxygen delivery. ATTENTION: the oxygen inlet flow must be set at MAX 5 L/min' and the pressure must be MAX 50 Kpa. The flow and oxygen concentration at $95 \pm 3\%$ must be regulated by qualified personnel.

5 - Press the START (1.4) key to switch on the appliance. The time is set to MED by default.

6 - Press the TIME(1.5) key to choose a different oxygen delivery time.

The LEDs (1.5.1) show the selected time.

MIN - minimum time (10 sec.)

MED - medium time (20 sec.)

MAX - maximum time (30 sec.)

7 - Turn the VACUUM regulator (1.11) to change the vacuum setting to be created within the cap. It increas

creases by turning it clockwise, it decreases anti-clockwise.

8 - Take the cap, close the hole of the suction Control Connector by means of the appropriate plug (see DWG 1) and place it on the area to be treated, then start treatment. The straps supplied may be used to hold the cap in place on the area to be treated (see DWG 4).

During (step 1) the cap goes into a vacuum and remains attached to the user, whereas in (step 2) a light pressure must be applied to prevent the cap from detaching.



If the cap needs to be removed in the suction stage, open the plug as shown in DWG 2.

There is the option of performing treatment using the suction control connector in manual mode by closing or opening the hole with one finger (see DWG 3).

WARNING. At any time, you may

- stop the appliance by pressing the STOP key (1.6)
- adjust the oxygen time (see point 6)
- modify the vacuum setting within the cap (see pint 7).

9 - Upon completing therapy switch off the appliance by pressing the STOP key (1.6), unplug it. SHUT OFF OXYGEN DELIVERY, then disconnect the tube (5) from the oxygen fitting (4)

DWG 4

CLEANING, SANITATION, DISINFECTION AND STERILISATION

Switch off the device before any cleaning procedure and unplug the power cable from the socket. Disconnect tube (5) from fitting (4).

The accessories 5-6a/6b/6c/6d/6e/6f-7-9-10 are for single use only.

The anti-bacterial filter (8) is single use.

CLEANING THE DEVICE AND TUBES 5-7

Use only a damp cloth with antibacterial soap (non-abrasive and with no solvents of any sort).

SANITISATION, DISINFECTION AND STERILISATION OF THE ACCESSORIES

SANITISATION

before and after each use, sanitise accessories (6a/6b/6c/6d/6e/6f-9-10) under potable hot water (approximately 40°C) with a gentle, non abrasive dish detergent or in the dishwasher with a hot cycle. After having sanitised the accessories, shake them vigorously and place them on a paper towel, or dry them with a jet of hot air (for example, a hair dryer).

DISINFECTION

Accessories (6a/6b/6c/6d/6e/6f-9-10) can be disinfected

The disinfection procedure described in this paragraph is to be followed before using the accessories. It is effective on the parts that undergo this treatment only if each step is followed and only if the parts to be treated have previously been sanitised. The disinfectant must be an electrolytic chloroxidiser (active principle: sodium hypochlorite) specific for disinfecting, which is available in any pharmacy.

Directions:

- Fill a container big enough to hold all of the parts to disinfect with a solution of potable water and disin-

fectant, according to the proportions indicated on the packaging of the disinfectant.

- Completely immerse each part in the solution, taking care to avoid the formation of air bubbles on the parts. Leave the parts immersed for the amount of time indicated on the packaging of the disinfectant associated with the concentration chosen for the solution.
- Remove the disinfected parts and rinse abundantly with warm potable water.
- After having disinfected the accessories, shake them vigorously and place them on a paper towel, or dry them with a jet of hot air (for example, from a hair dryer).
- Dispose of the solution following the instructions provided by the disinfectant manufacturer.

STERILISATION

Accessories (**6a/6b/6c/6d/6e/6f**) can be sterilised.

The tubes, filter, connectors and all fittings CANNOT be sterilised as the polymers they are made of cannot withstand high temperatures.

The sterilisation procedure described in this paragraph is effective on the parts that undergo this treatment only if every step is followed and only if the parts to be treated have previously been sanitised. The procedure is validated in its conformity to ISO 17665-1.

Device: Steam steriliser with fractionated vacuum and overpressure in accordance with EN 13060.

Directions: Wrap every single part to be treated with a sterile barrier system or packaging in accordance with standard EN 11607. Place the packed components in the steam steriliser. Run the sterilisation cycle according to the operating instructions of the device by selecting a temperature of 134°C and a time of 10 minutes first.

Storage: Store the sterilised parts as per the instructions for use of either the sterile barrier system or packaging.

TROUBLE-SHOOTING

Problem	Cause	Solution
LED 1.1 POWER does not switch on	The power supply cable is not correctly attached.	Insert the power cable.
	Fuses cut-off	Take it to an authorised centre.
The LED 1.1 POWER is on but the device does not work	Circuit breaker tripped	Allow the medical device to cool for at least 40 minutes. Try and switch it back on and if it still does not work take it to an authorised centre.
The LED ALARM 1.7 is on	Over-vacuum protection tripped	When the ALARM indicator light turns on (red LED flashing), decrease the suction value using the Suction Adjustment knob. If the problem persists, close the O2, sitch-off the device and take it to an authorized service center or to your trusted representative.

SYMBOLS

EC Marking medical ref. Dir. 93/42 EEC and subsequent updates

Risk: electrocution.
Consequence: Death.
Using the device while taking a bath or shower is prohibited

Attention: check the instructions for use

Type BF applied part

Do not smoke

Device serial number

Do not use naked flames

Manufacturer

Oxygen

Protective earth

"ON"

"OFF"

TECHNICAL SPECIFICATIONS

Mod. P0613EM

Power source:	220V-240V ~ 50Hz 100W
Fuses:	2 x T2A 250V
Max suction negative pressure:	-200 mbar/-150 mmHg approx
Noise (at 1 m):	approx 56 dB (A)
Operation:	Continuous
Hyperbaric oxygen inlet concentrated at 95 ± 3 %	MAX 5 L/min' MAX 50 Kpa
Dimensions:	(L) 40.5x (D) 22 x (H) 17 cm
Weight:	5.5 Kg approx

Operating Conditions:	
Temperature:	min. 10°C; max 40°C
Air humidity:	min. 10%; max 95%
Atmospheric pressure:	min. 69KPa; max 106KPa
Transport and storage conditions:	
Temperature:	min. -25°C; max 70°C
Air humidity:	min. 10%; max 95%
Atmospheric pressure:	min. 69KPa; max 106KPa

Safety approvals:

Type BF applied parts are: patient accessories (6a,6b,6c,6d,6e,6f)

Disposal: In compliance with the Directive 2012/19/EC, the symbol printed on the device shows that the device to be disposed of is considered waste and must therefore be an item of "differentiated collection". Consequently, the user must take it (or have it taken) to the designated collection sites provided by the local authorities, or turn it in to the dealer when purchasing an equivalent new device. Differentiated waste collection and the subsequent treatment, recycling and disposal procedures promote the production of devices made with recycled materials and limit the negative effects on the environment and on health caused by potential improper waste management. If the user disposes of the product in an unauthorised manner, administrative penalties shall be applied according to the transposition laws of directive 2012/19/EC of the member state or country where the product is disposed of.

Packaging disposal

Product box

Polystyrene inserts

Accessories packaging bag
and cable packaging bag

Electromagnetic compatibility:

This device was designed to satisfy the currently required requisites for electromagnetic compatibility (EN 60601-1-2). Electro-medical devices require particular care during installation and use relative to EMC requirements. Users are therefore requested to install and/or use these devices following the manufacturer's specifications. There is a risk of potential electromagnetic interference with other devices. RF mobile or portable radio and telecommunications devices (mobile telephones or wireless connections) can interfere with the functioning of electro-medical devices. For further information visit our website www.flaemnuova.it. The Medical Device may be subject to electromagnetic interference if other devices are used for specific diagnosis or treatments. Flaem reserves the right to make technical and functional modifications to the product with no prior warning.

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FLAEM NUOVA S.p.A. non risponde di ulteriori estensioni del periodo di garanzia assicurate da terzi.
FLAEM NUOVA S.p.A.