**SfigmoLife Evo**

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## Product Overview

**OPERATING PRINCIPLES.** The blood pressure monitor is a convenient and practical medical device developed by using modern electronic technology and the principle of indirect blood pressure measurement. The blood pressure is measured indirectly through the oscillometric method. During the process of deflation or inflation, the monitor senses pulse wave information and obtains systolic and diastolic blood pressure data through a series of complex conversions and calculations. The main function of the blood pressure monitor is to measure and display the blood pressure and pulse rate of the human body.

**INTENDED USE.** The blood pressure monitor is intended for use by individuals older than twelve (12) years to measure the systolic and diastolic blood pressure and pulse rate at home and in medical institutions.

**INDICATIONS FOR USE.** Measure blood pressure (systolic and diastolic) and pulse rate.

**CONTRAINDICATION.** - Do not use this product on pregnant women, children, or individuals unable to express their own consciousness. - Do not use it on patients with a pacemaker, defibrillator, or implanted metallic or electronic device.

- Do not use it on patients with severe arrhythmia.

## Product Profile

**PRODUCT CHARACTERISTICS.** - Large font display - One-button operation - Indication of blood pressure range - Body movement indicator - Cuff air-tightness indicator - Voice broadcast function (optional) - Data Recording for multiple users - USB Type-C port

**DEVICE DESCRIPTION.** The blood pressure monitor uses an inflatable cuff wrapped around the patient's upper arm. After pressing the "START/STOP" button, the cuff is automatically inflated by an internal pump. Systolic and diastolic pressures are measured using the oscillometric method. The deflation rate is controlled by an internal exhaust valve. A quick-release mechanism allows the cuff pressure to be fully released in case of an emergency. The device has a maximum pressure safety limit of 300 mmHg and will not inflate beyond this value. If an irregular heartbeat is detected during measurement, an indicator will be activated. After the measurement, systolic and diastolic pressures along with pulse rate are displayed on the screen and stored in memory. The cuff then deflates automatically to 0 mmHg.

**STRUCTURE COMPOSITION.** The Blood Pressure Monitor is composed of a main unit, cuff, battery and USB Type-C cable.

## PRODUCT APPEARANCE


1. Power Adapter Port (Img. 1)
2. Air Hole (Img. 2)

**CAUTION.** Please use the cuff specified by the manufacturer; otherwise, it may

cause measurement errors. For the power adapter socket, please use the USB cable specified by the manufacturer.

| PACKING LIST                                                             |          |
|--------------------------------------------------------------------------|----------|
| Name                                                                     | Quantity |
| Main device                                                              | 1 set    |
| Large-size cuff (optional)<br>Suitable for arm circumference:<br>22-42cm | 1 piece  |
| 1.5V AA-size alkaline battery (optional)                                 | 4 pieces |
| USB Type-C cable                                                         | 1 piece  |
| User manual                                                              | 1 piece  |

## INTRODUCTION TO DISPLAY SCREEN

According to the WHO Hypertension Classification, the following values can be used as indicators for high blood pressure in home measurements.

When the systolic value is greater than or equal to 140 mmHg (18.7 kPa) or the diastolic value is greater than or equal to 90 mmHg (12 kPa), the “▲” will light up. This indicator applies to home blood pressure measurements.

- Do not self-diagnose or self-treat based solely on the measurement results.
- Consult a doctor if necessary.

### Safety Instructions

**ELECTRICAL SAFETY.** - The design of the blood pressure monitor complies with the relevant national safety requirements for medical electrical equipment, such as IEC 60601-1 and other standards. - According to the General Requirements for Safety of Medical Electrical Equipment (IEC 60601-1), the classification is as follows:

- a. Classified by the type of protection against electric shock: internal power supply, class II equipment.
- b. Classified by the degree of protection against electric shock: BF type application part.
- c. Classified by the degree of protection against harmful liquid: not applicable.
- d. Classified by the degree of safety when using flammable anesthetic gas mixed with air or flammable anesthetic gas mixed with oxygen or nitrous oxide: non-AP/ APG type.
- e. Classified by the electromagnetic compatibility: Group I, Class B.

**⚠Warning.** The information you should know about how to avoid harm and injury to users and medical personnel.

**⚠Caution.** The emphasized important information.

**⚠WARNING.** - Do not perform high temperature disinfection or use electron beam and gamma-rays for sterilization. - Make sure that the environment where the device is used is not interfered by strong electromagnetic interference sources, such as wireless transmitters, cell phones, microwave ovens, etc. - This device cannot be used in places where there are high-frequency electric surgical knives, microwave ovens and mobile phones. - This device should not be used close to or stacked with other equipment. If it must be used close to or stacked with other equipment, it should be observed and verified that it can operate normally under its intended configuration. - This product is only used for blood pressure measurement. Accidents may occur when used for purposes other than measuring blood pressure. - This product is intended only for daily monitoring of blood pressure and is not suitable for the diagnosis of hypertension. DO NOT adjust medication based on measurement results from this blood pressure monitor. Take medication as prescribed by your physician. Only a physician is qualified to diagnose and treat High Blood Pressure. It is very dangerous for the user to interpret the results and begin treatment on their own. Please give the test results to a doctor who knows

your health condition for diagnosis. If the user thinks that there is a problem with the measurement result, he/she should consult professionals or after-sales personnel. - For those with peripheral circulatory disorders caused by certain diseases, the blood pressure values measured at the wrist and upper arm may be obvious differences. - Do not use the device on the injured arm or the arm under medical treatment. - Do not apply the arm cuff to an arm receiving an intravenous drip or blood transfusion. - Do not use the device in the case of flammable anesthetic gases mixed with air or flammable anesthetic gases mixed with oxygen or nitrous oxide. - Consult your physician before using the device on the arm with an arteriovenous (A-V) shunt. - The air-hose or the USB Type-C cable may cause accidental strangulation in children. - Contains small parts that may cause a choking hazard if swallowed by children. - Please do not disassemble, repair, or modify the device without permission. - The airbag may be at risk in a prolonged over-inflated state. - Consult your physician before using the device for any of the following conditions: • If you have had a mastectomy. • Common arrhythmias such as atrial or ventricular premature beats or atrial fibrillation, arterial sclerosis, poor perfusion, diabetes, age, pregnancy, pre-eclampsia, renal diseases. • People with severe blood flow problems or blood disorders as cuff inflation can cause bruising. Acute internal bleeding may be caused by compression of the arm during measurement • Do not take measurements more than necessary. It may cause bruising due to blood flow restriction. - The information in this manual may be modified without prior notice. The service life of the device is 5 years. When the device exceeds its effective service life, please dispose of it according to relevant local regulations. - Do not sterilize the device with low-temperature steam or other similar methods. - This manual is intended for general users without specialized training. - Users are expected to be familiar with the basic operation and functions of the device. - The device contains precision electronic components. Please do not open or tamper with it. - This product may be transported by air, sea, or land after packaging, but it should not be stored outdoors during transportation. - Please note that patient motion, trembling, or shivering may affect the measurement reading. - Remove the arm cuff if it fails to deflate during measurement. - Use only the special arm cuff provided by the manufacturer for this device. Use of other arm cuffs may result in incorrect measurement results. - Do not use the device in a moving vehicle (car, airplane).

- This device contains precision components. Avoid exposure to extreme temperatures, humidity, direct sunlight, shock, or dust. - Do not fold the cuff tightly or store the tube in a twisted state for extended periods, as this may shorten the lifespan of the components. - Used equipment, components, and batteries should not be treated as regular household waste and must be disposed of in accordance with applicable regulations. - If the appearance of the cuff is found to be damaged, do not take blood pressure measurements. It may damage your skin. - Do not press

the air hose while taking a measurement. - Do not expose it to devices or subject it to strong shocks or vibrations. - Do not inflate the arm cuff when it is not wrapped around your arm. - Do not drop the device or expose it to strong shocks or vibrations. - The accuracy of the product has been validated. It is generally recommended to have the device calibrated once a year to ensure proper and accurate performance. Do not attempt to calibrate the device yourself. If calibration is needed, please contact your dealer or our after-sales service team. - In case of pressure discomfort or if the pressure does not release, immediately turn off the unit or disconnect the tube plug to quickly deflate the cuff. - Avoid eating, drinking tea/cafeine/alcohol, exercising, bathing, and smoking for 30 minutes before taking a measurement. Rest for at least 5 minutes before taking the measurement. - Before performing another blood pressure measurement, please rest quietly for 2-3 minutes, allowing the artery to return to its normal state between measurements. - This device is intended for upper arm measurements only and is not suitable for use on the wrist or other parts of the body.

**USE OF ADAPTER WITH USB TYPE-C CABLE.** - Do not use the USB Type-C cable or Adapter if the device or the cable is damaged. Turn off the power and unplug the power cord immediately. - Plug the adapter with a USB Type-C cable into the appropriate voltage outlet.

- Never plug in or unplug the device from the electrical outlet with wet hands. - Securely insert the power plug into the outlet. - When disconnecting the power plug from the outlet, be sure to grasp and pull the power plug itself, not the cord.

- Unplug the adapter when not in use. - Disconnect the power plug before cleaning.

8. The adapter is an optional configuration supplied by the manufacturer. If the package does not include one, users should purchase a DC 5V/1A adapter that meets the safety requirements of IEC 60601-1 and IEC 60601-1-11. Use of an adapter that does not meet these requirements may damage the device and may be hazardous to users.

### **COMMON FAULTS AND TROUBLESHOOTING**

1. The table below lists common faults that may occur during use and provides corresponding troubleshooting methods. If you encounter a similar issue, please refer to the table for guidance.

| Common faults                                                                                                                      | Troubleshooting                                                                                                                                                                                                                                                                                                  |
|------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Although equipped with batteries, but it does not work after pressing the [Start/Stop] button.                                     | Check whether the positive and negative polarity of the battery is reversed.                                                                                                                                                                                                                                     |
|                                                                                                                                    | Please replace with new batteries.                                                                                                                                                                                                                                                                               |
| Repeated measurement failures, or the measurement value is seriously low (or high).                                                | Check whether the cuff is connected and wrapped correctly.                                                                                                                                                                                                                                                       |
|                                                                                                                                    | Check whether the cuff is wrapped too tightly or too loosely. Check whether the clothes rolled up are too tight to press the measurement position on the arm, if so, please take off the clothes and measure again.                                                                                              |
|                                                                                                                                    | Make sure to measure blood pressure in a quiet and relaxed state. It is best to take a few deep breaths before the measurement to relax yourself.                                                                                                                                                                |
| Although the machine is in good working condition, but the results of each measurement are different.                              | Please read the various precautions listed in "Factors affecting blood pressure measurement" carefully. Note: Blood pressure is dynamic, so there must be some differences between each measurement.                                                                                                             |
| The measured blood pressure value is different from the value that the doctor measured for yourself, or even seriously undersized. | Write down the blood pressure value measured by you every day and consult your doctor. Note: Frequent consultation with the doctor will cause psychological stress, which will lead to the blood pressure value measured in front of the doctor to be higher than the value measured in a relaxed state at home. |

|                                                                  |                                                                        |
|------------------------------------------------------------------|------------------------------------------------------------------------|
| Although the air pump works, but the air pressure does not rise. | Check whether the cuff connection is good and there is no air leakage. |
|                                                                  | If it is damaged due to long usage, please buy a new product.          |

1. The following table shows error indicators and their possible causes that may occur during the measurement process. Please use the correct method and try the measurement again.

| Error display | Cause of failure                                     |
|---------------|------------------------------------------------------|
| Err1          | Self-check error.                                    |
| Err2          | Cuff air-tightness issue.                            |
| Err3          | Air leakage error.                                   |
| Err4          | Air pressure error; the valve may not open properly. |
| Err5          | Weak signal detected.                                |
| Err6          | Pressure exceeds the measurement range.              |
| Err7          | Improper pressurization due to arm or body movement. |
| Err8          | Pressure exceeds hardware protection limit.          |
| Err9          | Signal saturation error.                             |
| Err10         | System error.                                        |
| Err11         | Measurement timeout error.                           |
| Err12         | Incorrect cuff type.                                 |
| Err13         | Device not calibrated.                               |
| Err14         | Unknown error.                                       |
| Err15         | Measurement range exceeded.                          |

|                              |                                                                                                                                                                                                                                                                                               |
|------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Err16                        | Insufficient valid pulse was detected during measurement. This error code appears to ensure measurement accuracy. It may result from tachycardia, bradycardia, geriatric, excessive pressure variation, or from low pulse pressure (the difference between systolic and diastolic pressures). |
| The low battery icon flashes | Battery power is low or depleted.                                                                                                                                                                                                                                                             |

**⚠WARNING.** If the issues above persist or if any other malfunction occurs, please contact your local distributor.

#### SYMBOLS PRESENT ON THE DEVICE OR PACKAGING

|                                                                                     |                                                       |                                                                                     |                                                                                                        |
|-------------------------------------------------------------------------------------|-------------------------------------------------------|-------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------|
|    | Serial number                                         | <b>IP20</b>                                                                         | Protected against solid foreign objects of 12,5 mmφ and greater; Non-protected for falling water drops |
|    | Manufacturer                                          |    | Caution                                                                                                |
|    | Type BF applied part Type BF applied part             |    | Direct current                                                                                         |
|    | Class II equipment                                    |    | Date of Manufacture                                                                                    |
|    | Unique device identifier                              |    | Medical device                                                                                         |
|    | CE Marking                                            |    | Model number                                                                                           |
|   | Refer to instruction manual                           |   | Importer and distributor                                                                               |
|  | MR unsafe                                             |  | Country and date of manufacture                                                                        |
|  | Temperature limits for storage and transport          |  | Humidity limit for storage and transport                                                               |
|  | Atmospheric pressure limits for storage and transport |  | Authorized representative                                                                              |

## DISPOSAL

In accordance with Directive 2012/19/EU, the symbol shown on the equipment indicates that the device to be disposed of (excluding accessories) is considered waste and must therefore be subject to “separate collection.” Accordingly, the user must take (or have taken) the waste to designated separate collection centers set up by local authorities, or hand it over to the retailer when purchasing a new device of equivalent type.

Separate collection of the waste and subsequent treatment, recovery, and disposal operations promote the production of equipment using recycled materials and help to limit the negative effects on the environment and human health caused by improper waste management.

Improper disposal of the product by the user may result in the application of administrative sanctions as provided for by the national legislation transposing Directive 2012/19/EU in the member state or country where the product is disposed of.

### Accessories

They must be disposed of as general waste after undergoing a sanitization process

### Packaging

Product box

Packaging bag

## INSTRUCTIONS FOR USE

**PREPARATION FOR USE.** Open the battery cover on the back of the device. Insert four AA-size batteries into the battery compartment, ensuring that the polarity (+ and –) matches the markings inside the compartment. After inserting the batteries, close the battery cover until a “click” sound is heard, indicating it is securely in place.

1. Open the battery cover and insert the battery (Img. 3).

2. Close the battery cover (Img. 4).

During use, if the “ ” appears, it indicates that the battery level is low. During use, if the “ ” starts flashing, it means the battery has run out of power. Please replace all 4 batteries at the same time. Do not mix the old and new batteries. It is recommended to use alkaline batteries for longer service life. If the device will not be used for an extended period (more than three months), please remove the batteries.

**⚠Caution.** The machine can be connected to an external power adapter. The power supply should meet the following conditions Class II equipment; Output voltage: DC 5V; Output current: ≥1A; Compliance with IEC 60601-1 or IEC 62368-1 standards; Must provide at least two MOOP insulation layers between the AC input and DC output

## HOW TO USE THE CUFF

1. Insert the air-hose plug on the cuff into the air-hose socket on the left side of the machine, as shown in the figure below (Img. 5).
2. Remove any tight-fitting clothing from your upper arm. Take off jackets, sweaters, or other bulky garments to ensure that the cuff is in direct contact with your skin. Clothing may interfere with pulse detection and result in inaccurate measurements.
3. Wrap the cuff snugly around your upper arm. Ensure that the lower edge of the cuff is positioned approximately 2–3 cm above the elbow. The cuff should fit comfortably and be able to rotate slightly around the arm (Img. 6).
4. Rest your lower arm flat on a table with your palm facing upward. Sit upright and ensure that the center of the cuff is at the same level as your heart.


⚠5. Make sure that the rubber air hose is not twisted, folded, or compressed by your arm.

**Caution.** • Constriction of the upper arm caused by rolling up your sleeve may result in inaccurate readings. Do not roll up tight sleeves during measurement; instead, ensure the cuff has direct contact with your bare arm. • Do not rest your arm on the air hose.

• Blood pressure can be measured on either the left or right arm. However, readings may differ between arms, so it is recommended to consistently use the same arm for measurements. • If there is a significant difference in readings between the two arms, please consult a physician for advice.

## CORRECT POSTURE FOR MEASUREMENT

1. Unfold the cuff and place it around the upper arm, with the lower edge of the cuff 2–3 cm away from the elbow joint. Adjust the cuff to ensure it is snug, but not overly tight, and that it can be rotated to fit the arm in any position (Img. 7).
2. Place the lower arm flat on the table, palms up. Tuck in the shoulder and elbow, so

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that the arm is level with the body. Position and align cuff at the same location (just above the inside of the elbow joint) as the heart, which should be held at chest level (Img. 8).

3. Make sure the cuff is not twisted. Be sure the airway is free, and the hose and the pressure bulb are correctly connected before inflating the cuff. Be sure not to pinch the patient's arm when positioning the cuff (Img. 9).

## MEASUREMENT OF BLOOD PRESSURE

In order to obtain accurate measurement results, please do not eat, smoke, drink alcohol, bathe, or engage in any strenuous physical activity within 30 minutes before taking a measurement. Measurements should be taken at the same time each day, as much as possible. Please sit quietly for 5 minutes before measuring.

1. Press the "START/STOP" button. The deflation icon on the display screen will flash, indicating that the system is performing a zero-reset test. After a few seconds, the inflation icon will appear, showing that the system has reset and the air pump is inflating the cuff.

2. Once the inflation is complete, the device will automatically begin measuring. The heartbeat icon will flash, the deflation valve will slowly release air, and the pressure will gradually decrease. Remain still and silent throughout the measurement. Do not talk or move your body or hands.

3. After the measurement is finished, the device will display your systolic pressure, diastolic pressure, and pulse rate. The cuff will automatically and quickly release any remaining air in the cuff.

4. Take off the cuff.

5. Press the "START/STOP" button to turn off the device. Even if you forget to turn it off, the device will automatically turn off within 1 minute after the measurement is complete. The user can use the "START/STOP" button to quickly deflate the cuff.

**Statement 1.** The blood pressure value measured by the device is equivalent to the value measured by the auscultation method, and its error meets the requirements of ISO 81060-2. For related information, please call the customer service hotline.

**Statement 2.** If the device is stored or used outside of the temperature or humidity range specified by the manufacturer, the system may not meet the generated performance specifications.

**Caution.** • Too frequent measurements may cause injury to the patient due to blood flow interference. • Pressurization of the cuff can temporarily cause loss of function of simultaneously used monitoring ME equipment on the same limb.

## AUXILIARY FUNCTIONS FOR CORRECT MEASUREMENT

1. Self-checking function of cuff: During the measurement process, the system automatically checks the tightness of the cuff and provides an indication: When the icon  is displayed, the tightness of the cuff is normal; when the icon  is flashing, Please press the "START/STOP" button to stop the measurement, and then adjust it

to the correct positioning and tightness before remeasurement. If you need further support, please contact the manufacturer or the local distributor.

2. Body movement: If the machine detects body movements during the measurement process, the icon  will flash, indicating that the movements may cause incorrect measurement results. Please measure again.

**USER SWITCHING.** To facilitate blood pressure management, this product supports 3 groups of users: User A and User B. and Other. Press the "USER" button to display the current user. Continue to press the "USER" button to switch between users. Press the "SET" key to confirm the ideal user for further measurements.

**USE OF MEMORY FUNCTION.** This product can automatically record 45 groups of measurement data under each user ID. (including time, date, blood pressure value, and pulse rate). When you need to check history records, you can access the recorded data by pressing the "MEM" button.

1. Pressing the "MEM" key once enters the memory review mode and displays the average value of the most recent three measurements. Pressing the "MEM" key again displays the first set of stored data.

2. Each time the "MEM" button is pressed, the measurement results will be displayed in the order from 1st to 45th.

3. Press the "MEM" button to browse the stored measurement records. Each time you press the "MEM" button, the device will first display the sequence number of the  measurement, followed by the corresponding measurement result.

4. Delete memory records. (Please refer to 4.13 "Deletion of Memory Records").

**Caution.** If the storage capacity is full, the measurement data from the 1st to the 45th reading will be sequentially replaced and stored by the most recent measurement data.

**TIME DISPLAY.** Press the "SET" button, and the current time will be displayed on the screen.

**SETTINGS OF TIME AND DATE.** The blood pressure value varies at different times of the day. This machine can automatically record each blood pressure reading and pulse rate to facilitate better management of your blood pressure. After installing the battery or powering on the machine, the time and date need to be reset. Follow the steps below to set the time and date (for example, if the date is 2019-11-15 and the time is 08:28 a.m.):

1. Press the "SET" button once to turn on the machine then hold down the "SET" button for about 3 seconds to enter the time setting mode. The two digits indicating the year will flash, then short press the "MEM" button to increase the year and the "USER" button to decrease the year, press the "SET" button to confirm the year setting.

2. After setting the year, continue to press the "SET" button to switch to the month setting, the digits indicating the month start flashing, then short press the "MEM"

button to increase, and the “USER” button to decrease the month.

3. In the same way, set the month, day, hour, and minute in sequence.

**UNIT SETTING OF MMHG OR KPA.** After completing the voice function settings, continue to press the “SET” button to enter the mmHg/kPa unit setting. Then, press the “USER” button to decrease or the “MEM” button to increase to switch between mmHg and kPa. Unit mmHg is set by default for the result displayed.

**DELETION OF MEMORY VALUE.** After completing the mmHg/kPa unit setting, continue to press the “SET” button to enter the memory deletion mode for the current user. The memory icon and user icon will be displayed. Then, press the “USER” button to decrease or the “MEM” button to increase to toggle between YES (delete) and NO (do not delete). Press the “SET” button again to confirm the selected option. By default, “Do not delete” is displayed.

| TECHNICAL PARAMETERS             |                                                                                                                                                                                                                                   |
|----------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Name                             | Blood pressure monitor                                                                                                                                                                                                            |
| Display mode                     | LCD digital display                                                                                                                                                                                                               |
| Method of measurement            | Oscillometric method                                                                                                                                                                                                              |
| Measurement range                | Pressure: 0mmHg ~ 290mmHg<br>0kPa ~ 38.6kPa<br>DIA: 20 mmHg (2.7 kPa) ~ 200 mmHg (26.6 kPa)<br>SYS: 35 mmHg (4.7 kPa) ~ 255 mmHg (33.9kPa)<br>Pulse rate: 30bpm ~ 200bpm                                                          |
| Accuracy                         | Static Pressure: $\pm 0.4$ kPa ( $\pm 3$ mmHg)<br>Pulse Rate: $\pm 5$ bpm or $\pm 5\%$ of the reading, whichever is greater.                                                                                                      |
| Operation mode classification    | Continuous operation                                                                                                                                                                                                              |
| Electrical safety classification | Internal power supply equipment<br>(When using only the batteries),<br><input checked="" type="checkbox"/> Class II equipment (When using the adapter and cable),<br><input checked="" type="checkbox"/> BF type application part |

| TECHNICAL PARAMETERS                            |                                                                                                                                                                                                                    |
|-------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Types of equipment                              | Non-AP/APG type equipment (the equipment that cannot be used in the presence of flammable anesthetic gas mixed with air or flammable anesthetic gas mixed with oxygen or nitrous oxide)                            |
| The degree of protection against harmful liquid | IP20                                                                                                                                                                                                               |
| Pressure sensor                                 | Semiconductor pressure sensor                                                                                                                                                                                      |
| Pressurization mode                             | Automatic pressurization by using pressure pump                                                                                                                                                                    |
| Deflation mode                                  | Automatic quick deflation                                                                                                                                                                                          |
| Power supply                                    | Four AA-size dry batteries (DC 6V) or Power adapter: Input AC100V~240V, 50/60Hz, 0.5A; Output DC 5V 1A                                                                                                             |
| Battery life                                    | Approximately 330 measurements can be performed by using four AA-size high-performance alkaline dry batteries [Under the condition of room temperature of 23 °C and pressurization to 192mmHg (25.3kPa) each time] |
| Suitable for arm circumference                  | 22-42cm (large)                                                                                                                                                                                                    |
| Overall dimension                               | 108*108*55mm                                                                                                                                                                                                       |
| Weight                                          | Approx. 380g (device only excluding batteries)                                                                                                                                                                     |
| Memory storage                                  | Up to 45 records for each of 3 users.                                                                                                                                                                              |
| Automatic shutdown                              | Automatic shutdown after no operation within 60s                                                                                                                                                                   |

|                                     |                                          |
|-------------------------------------|------------------------------------------|
| Operating environment               | Temperature: 5°C ~ 40°C                  |
|                                     | Humidity: ≤85%RH, no condensation        |
|                                     | Atmospheric pressure: 80.0kPa ~ 105.0kPa |
| Conditions of storage and transport | Temperature: -20°C ~ 55°C                |
|                                     | Humidity: ≤93%RH, no condensation        |
|                                     | Atmospheric pressure: 50.0kPa ~ 106.0kPa |

### Electromagnetic Compatibility

**⚠ CAUTION.** • This device complies with the electromagnetic compatibility requirements of IEC 60601-1-2 and IEC 80601-2-30. • Users should install and operate the device in accordance with the electromagnetic compatibility information provided.

• Please refer to the annex for detailed guidelines and the manufacturer's declarations.

**⚠ WARNING.** • Do not use the device near active HF surgical equipment or inside the RF shielded room of a magnetic resonance imaging (MRI) system, where the intensity of electromagnetic disturbances is high. • Avoid using this device adjacent to or stacked with other equipment, as it may cause improper operation. If such use is necessary, ensure both devices are monitored to confirm they function correctly.

• The use of accessories and cables other than those specified or provided by the manufacturer may increase electromagnetic emissions or reduce electromagnetic immunity, potentially resulting in malfunction. • Portable RF communication equipment should not be used within 30 cm (12 inches) of any part of the Upper Arm Blood Pressure Monitor, including specified cables. Otherwise, device performance may be degraded. • Operating below the minimum specified amplitude or value in the user manual may lead to inaccurate results. The minimum amplitude of the patient's physiological signal is 30 bpm. • Using the device in environments with strong electromagnetic sources (e.g., MRI or CT imaging rooms) may negatively affect its performance. Avoid use in such settings. **⚠** • Essential Performance is defined as performance, the absence or degradation of which would present an unacceptable risk. The essential performance functions of this device comply with Annex I of IEC 60601-1-2, EN 60601-1-2, IEC 80601-2-30 (Subclauses 201.17 & 202), and EN IEC 80601-2-30 (Subclauses 201.17 & 202). • Even if other equipment complies with relevant emission standards, this device may still experience interference from it.

**CABLE**

| No | Cable Name       | Cable Length (m) | Shielded or not | Remarks |
|----|------------------|------------------|-----------------|---------|
| 1  | USB Type-C cable | 0.7              | no              | \       |

**Guidelines and Manufacturer's Statements – Electromagnetic Emission**

The Blood Pressure Monitor is intended to be used in the following specified electromagnetic environment, and the buyers or users of the device shall ensure that it is used in such electromagnetic environment:

| Emission Test                                             | Conformity | Electromagnetic Environment - Guideline                                                                                                                                               |
|-----------------------------------------------------------|------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| RF Emissions<br>CISPR 11                                  | Group 1    | It uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.               |
| RF Emissions<br>CISPR 11                                  | Category B | It is suitable for used in domestic establishments and in establishment directly connected to a low voltage power supply network which supplies buildings used for domestic purposes. |
| Harmonic Emissions<br>IEC 61000-3-2                       | Class A    |                                                                                                                                                                                       |
| Voltage Fluctuations / Flicker Emissions<br>IEC 61000-3-2 | Complies   |                                                                                                                                                                                       |

### Guidelines and Manufacturer's Statements – Electromagnetic Immunity

The Blood Pressure Monitor is intended to be used in the following specified electromagnetic environment, and the buyers or users of the device shall ensure that it is used in such electromagnetic environment:

| Immunity Test                                     | IEC 60601 Test Level                                                                    | Coincidence Level                                                               | Electromagnetic Environment - Guideline                                                                                                               |
|---------------------------------------------------|-----------------------------------------------------------------------------------------|---------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------|
| Electrostatic Discharge (ESD)<br>IEC 61000-4-2    | $\pm 8$ kV contact<br>$\pm 2$ kV, $\pm 4$ kV, $\pm 8$ kV,<br>$\pm 15$ kV air            | $\pm 8$ kV contact<br>$\pm 2$ kV, $\pm 4$ kV,<br>$\pm 8$ kV, $\pm 15$ kV<br>air | The floor should be wood, concrete or ceramic tile.<br>If the floor is covered with synthetic material, the relative humidity should be at least 30%. |
| Electrical Fast Transient/ Burst<br>IEC 61000-4-4 | $\pm 2$ kV for power supply lines<br>$\pm 1$ kV for Input/output lines                  | $\pm 2$ kV for power supply lines                                               | Mains power quality should be that of a typical commercial or hospital environment.                                                                   |
| Surge<br>IEC 61000-4-5                            | $\pm 0.5$ , $\pm 1$ kV line to line<br>$\pm 0.5$ , $\pm 1$ kV, $\pm 2$ kV line to earth | $\pm 0.5$ , $\pm 1$ kV line to line                                             | Mains power quality should be that of a typical commercial or hospital environment.                                                                   |

| Immunity Test                                                                                                   | IEC 60601 Test Level                                                                                                                                                                                        | Coincidence Level                                                                                                                                                                                                               | Electromagnetic Environment - Guideline                                                                                                                                                                                                                                                                  |
|-----------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <p>Voltage dips, short interruptions, and voltage variations on power supply input lines<br/>IEC 61000-4-11</p> | <p>&lt;5 % UT(&gt;95% dip in UT.)<br/>for 0.5 cycle<br/>&lt;5 % UT(&gt;95% dip in UT.)<br/>for 1 cycle<br/>70% UT(30% dip in UT)<br/>for 25/30 cycles<br/>&lt;5% UT(&gt;95 % dip in UT)<br/>for 5/6 sec</p> | <p>&lt;5 % UT<br/>(&gt;95% dip in UT.)<br/>for 0.5 cycle<br/>&lt;5 % UT<br/>(&gt;95% dip in UT.)<br/>for 1 cycle<br/>70% UT<br/>(30% dip in UT)<br/>for 25/30 cycles<br/>&lt;5% UT<br/>(&gt;95 % dip in UT)<br/>for 5/6 sec</p> | <p>Mains power quality should be that of a typical commercial or hospital environment. If the user of the electronic sphygmomanometer requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.</p> |
| <p>Power Frequency Magnetic Field (50/60Hz)<br/>IEC 61000-4-8</p>                                               | <p>30A/m</p>                                                                                                                                                                                                | <p>30A/m</p>                                                                                                                                                                                                                    | <p>The power frequency magnetic field should have the horizontal characteristics of a typical site in a typical commercial or hospital environment.</p>                                                                                                                                                  |
| <p>Note: UT refers to AC network voltage before the test voltage is applied.</p>                                |                                                                                                                                                                                                             |                                                                                                                                                                                                                                 |                                                                                                                                                                                                                                                                                                          |

### Guidelines and Manufacturer's Statements – Electromagnetic Immunity

The Blood Pressure Monitor is intended to be used in the following specified electromagnetic environment, and the buyers or users of the device shall ensure that it is used in such electromagnetic environment:

| Immunity test                                    | IEC 60601 Test Level                                   | Coincidence level                                          | Electromagnetic Environment - Guideline |
|--------------------------------------------------|--------------------------------------------------------|------------------------------------------------------------|-----------------------------------------|
| Conducted<br>RF transmission<br>IEC 61000-4-6    | 3 Vrms<br>150 kHz ~<br>80 MHz<br>6 Vrms in<br>ISM bans | 3 Vrms<br>150 kHz ~<br>80 MHz<br>6 Vrms nelle<br>bande ISM | See the box below                       |
| Trasmissione<br>RF irradiata<br>IEC<br>61000-4-3 | 10 V/m<br>80 MHz ~<br>2.7 GHz                          | 10 V/m<br>80 MHz ~<br>2.7 GHz                              |                                         |

The portable and mobile RF communication equipment should not be used closer than the recommended isolation distance to any part of the Blood Pressure Monitor, including the cables. The distance shall be calculated by a formula corresponding to the transmitter frequency. The recommended isolation distance:

$$d = 1.2 \sqrt{P}$$

$$d = 1.2 \sqrt{P} \quad 80 \text{ MHz} \sim 800 \text{ MHz}$$

$$d = 2.3 \sqrt{P} \quad 800 \text{ MHz} \sim 2.7 \text{ GHz}$$

In the formula, P is the maximum rated output power of the transmitter (in watts W) provided by the transmitter manufacturer; D is the recommended isolation distance in meters (m). The field intensity of a stationary RF transmitter is determined by surveying the electromagnetic field a, and at each frequency range B should be lower than the coincidence level. The interference may occur near devices marked with the following symbol. (V)

Note 1: The formula of higher frequency band should be adopted at 80MHz and 800MHz.

Note 2: This guideline may not be applicable to all cases, as the electromagnetic propagation is influenced by the absorption and reflection of buildings, objects, and human bodies.

| Immunity test                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  | IEC 60601 Test Level | Coincidence level | Electromagnetic Environment - Guideline |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------|-------------------|-----------------------------------------|
| <p>A. The field strength of fixed transmitters, such as base stations for wireless (cellular/cordless telephones and mobile ground radios, amateur radios, AM and FM radios, and television broadcasts etc., cannot be theoretically predicted with any accuracy. In order to evaluate the electromagnetic environment of fixed RF transmitter, the electromagnetic field survey should be conducted. If the measured field intensity of the site where the LT-P Blood Pressure Monitor is used is higher than the RF conformance level of the above application, the device should be observed to verify its normal operation. If abnormal performance of the device is observed, the additional measures may be necessary, such as reorienting or re-positioning the LT-P Blood Pressure Monitor.</p> <p>B. The field intensity should be lower than 3 V/m in the whole frequency range of 150KHz~80MHz.</p> |                      |                   |                                         |

#### Guidance and manufacturer's declaration - Electromagnetic Immunity

Radiated RF IEC 61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment)

| Test Frequency (MHz) | Band (MHz) | Service           | Modulation                          | Modulation (W) | Distance (m) | IMMUNITY TEST LEVEL (V/m) |
|----------------------|------------|-------------------|-------------------------------------|----------------|--------------|---------------------------|
| 385                  | 380        | TETRA 400         | Pulse modulation 18 Hz              | 1.8            | 0.3          | 27                        |
|                      | -          |                   |                                     |                |              |                           |
|                      | 390        |                   |                                     |                |              |                           |
| 450                  | 430        | GMRS 460, FRS 460 | FM $\pm$ 5 kHz deviation 1 kHz sine | 2              | 0.3          | 28                        |
|                      | -          |                   |                                     |                |              |                           |
|                      | 470        |                   |                                     |                |              |                           |
| 710                  | 704        | LTE Band 13, 17   | Pulse modulation 217 Hz             | 0.2            | 0.3          | 9                         |
| 745                  | -          |                   |                                     |                |              |                           |
| 780                  | 787        |                   |                                     |                |              |                           |

|      |      |                                                              |                         |     |     |    |
|------|------|--------------------------------------------------------------|-------------------------|-----|-----|----|
| 810  | 800  | GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5       | Pulse modulation 18 Hz  | 2   | 0.3 | 28 |
| 870  | –    |                                                              |                         |     |     |    |
| 930  | 960  |                                                              |                         |     |     |    |
| 1720 | 1700 | GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1,3,4,25; UMTS | Pulse modulation 217 Hz | 2   | 0.3 | 28 |
| 1845 | –    |                                                              |                         |     |     |    |
| 1970 | 1990 |                                                              |                         |     |     |    |
| 2450 | 2400 | Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7         | Pulse modulation 217 Hz | 2   | 0.3 | 28 |
|      | –    |                                                              |                         |     |     |    |
|      | 2570 |                                                              |                         |     |     |    |
| 5240 | 5100 | WLAN 802.11 a/n                                              | Pulse modulation 217 Hz | 0.2 | 0.3 | 9  |
| 5500 | –    |                                                              |                         |     |     |    |
| 5785 | 5800 |                                                              |                         |     |     |    |

## Hygienic Preparation

To protect the blood pressure monitor from damage and to ensure accurate measurements, please follow these guidelines:

1. After use or during storage, properly organize and store the device and its accessories to prevent exposure to strong impacts or vibrations.
2. Do not expose the device and its accessories to high temperatures, high humidity, dust, or direct sunlight.
3. Since the cuff contains an air-tight inner bladder, please handle it carefully and avoid folding, pulling, or twisting it.
4. Do not disassemble or attempt to repair this device without authorization.
5. Do not replace internal components without authorization.
6. If the machine will not be used for an extended period (more than 3 months), remove the batteries.
7. If the machine becomes dirty, clean it with a soft, dry cloth. If the machine is particularly soiled, use a soft cloth dampened with water, wring it out, and then wipe it dry. If necessary, gently disinfect with absorbent cotton moistened with alcohol, avoiding the control panel and nameplate. Do not use detergent or cleaning powder.
8. Do not spill water or other liquids into the machine.
9. It is recommended to perform maintenance every 3 months. Users should not attempt any maintenance other than battery replacement.

**Note.** It is recommended that the device be calibrated at least once a year by the national metrology and testing department, the manufacturer, or an authorized repair center. If a static pressure manometer test mode is required, please contact the manufacturer or distributor for further information.

**Note.** The manufacturer, importer, and distributor shall not be held responsible for any device malfunction resulting from failure to comply with the above storage and usage precautions.

## Common Questions and Answers about Blood Pressure Measurement

1. What factors affect blood pressure measurement?
  - Blood pressure should be measured under consistent conditions for accurate comparison. Ideally, measurements should be taken in a quiet environment at the same time each day.
  - Factors such as poor arm support, mental stress, or tension can increase blood pressure. Ensure you are in a comfortable and relaxed state before measuring, and avoid moving your arm during the measurement process.
  - The cuff should be positioned at the same height as your heart. If it is too high or too low, the results may be inaccurate.
  - Repeated measurements can cause blood to accumulate in the arm, and this may lead to numbness and incorrect readings. Therefore, it is recommended to wait 2-3 minutes (or longer) between consecutive measurements.
  - Other factors that can affect blood pressure measurement results include consu-

ming alcohol, tea, coffee, smoking, bathing, taking medication, or eating within 30 minutes of measurement, as well as the body movements used during the measurement process.

2. Why is the blood pressure measured at home lower than at the hospital?

- Blood pressure measured at home may be 20 mmHg- 30 mmHg (2.7kPa-4.0kPa) lower than at the hospital because hospital visits can cause anxiety, leading to higher readings. At home, you are usually more relaxed. Therefore, it is important to know your normal blood pressure at home.
- If the cuff is positioned higher than the heart during measurement, the reading may be lower. Please ensure the correct posture is used during measurement.

3. Why is the blood pressure measured at home higher than at the hospital?

- In people taking antihypertensive medication, readings may be higher when the medication wears off.
- Incorrect cuff positioning can prevent accurate detection of the pulse signal, resulting in elevated readings. Ensure the arm is properly positioned.
- Poor posture, such as bending forward, crossing legs, or sitting on a low surface, can also lead to high readings if the arm is below heart level. Always maintain proper posture.

4. Why do I feel pain or numbness from the cuff during measurement?

- The cuff inflates to temporarily stop blood flow in the arteries, which may cause mild pain or numbness. This is normal and not harmful. Continue the measurement without worry.

5. When is the best time to measure blood pressure?

- In the morning, measure your blood pressure after getting up, using the bathroom, and before breakfast.
- Measure at night before going to bed. If you take medication for high blood pressure, measure before taking your medication.
- At other times, ensure your body and mind are calm. Try to measure at the same time each day.

6. Why does my blood pressure reading vary each time?

- Blood pressure changes with every heartbeat. For example, a person with a pulse rate of 70 beats per minute experiences approximately 100, 800 blood pressure changes per day.
- Since blood pressure fluctuates constantly, a single measurement may not affect the full picture. It is advisable to take 2-3 readings in a row, each measurement should gap of at least 2-3 minutes.
- When measuring continuously, note the following:
  - Arm pressure may restrict blood flow to the fingers, causing congestion. If this happens, loosen the cuff, raise the arm, and perform some hand-clenching and stretching before remeasuring.
  - Improper cuff placement, especially if too close to the elbow, may result in inaccurate readings. Always ensure proper posture and cuff positioning.

7. Why is managing blood pressure at home important?

- Tracking your daily blood pressure helps identify trends and is useful for managing

your health. These records are valuable for medical diagnosis. For best accuracy, record the measurement conditions (e.g., time, medications, lifestyle factors) along with the readings.

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**SETTINGS OF TIME AND DATE.** The blood pressure value varies at different times of the day. This machine can automatically record each blood pressure reading and pulse rate to facilitate better management of your blood pressure.

After installing the battery or powering on the machine, the time and date need to be reset. Follow the steps below to set the time and date (for example, if the date is 2019-11-15 and the time is 08:28 a.m.):

1. Press the "SET" button once to turn on the machine then hold down the "SET" button for about 3 seconds to enter the time setting mode. The two digits indicating the year will flash, then short press the "MEM" button to increase the year and the "USER" button to decrease the year, press the "SET" button to confirm the year setting.
2. After setting the year, continue to press the "SET" button to switch to the month setting, the digits indicating the month start flashing, then short press the "MEM" button to increase, and the "USER" button to decrease the month.
3. In the same way, set the month, day, hour, and minute in sequence.

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© 2025 FLAEM NUOVA®
**FLAEM NUOVA S.p.A.**