# FLAEM<sup>®</sup>

Medical Devices

# Zeffiro

Code: SPC04

**REF** DL-08B

IT > MANUALE ISTRUZIONI  
EN > USER MANUAL  
DE > BETRIEBSANLEITUNG  
FR > MANUEL D'UTILISATION  
NL > GEBRUIKERSHANDLEIDING  
EL > ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ  
PL > INSTRUKCJA OBSŁUGI

ITALIANO

ENGLISH

DEUTSCH


FRANÇAIS

NEDERLANDS

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POLSKI

Schema di collegamento - Connection diagram - Anschlussplan - Schéma de connexion - Aansluitschema - Διάγραμμα σύνδεσης - Schemat połączeń

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### 1.1 Generale

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## 1 IMPORTANT INFORMATION

### 1.1 General

Please read these instructions for use carefully and completely. Do not discard, so you can consult them at a later date. If you fail to comply with the instructions for use, injury or damage to the device cannot be ruled out.

If symptoms persist or if your condition worsens, discontinue the application and contact your doctor immediately.

### 1.2 Structure of safety instructions

Safety-critical warnings are categorized according to hazard levels in these instructions for use:

- The signal word **WARNING** is used to indicate hazards which, without precautionary measures, can result in serious injury or even death.
- The signal word **CAUTION** is used to indicate hazards which, without precautionary measures, can result in minor to moderate injury or impair treatment.
- The signal word **NOTICE** is used to indicate general precautionary measures which are to be observed to avoid damaging the product during use.

### 1.3 Treatment of anyone who requires assistance

Anyone who requires assistance must be supervised constantly by an adult during inhalation therapy. This is the only way to ensure safe and effective treatment.

Individuals in this group often underestimate the hazards involved, thus resulting in a risk of injury. The product contains small parts. Small parts can block the airways and lead to a choking hazard. Therefore, make sure that you always keep all components of the product out of the reach of babies and infants. This product is only designed for patients who are able to breathe by themselves and are conscious. This is the only way to ensure effective treatment and avoid a choking hazard.

### 1.4 Hygiene

Only use the SPACER and accessories that have been cleaned, disinfected and dried as described in chapter 4 for inhalation therapy. Contamination and residual moisture encourage the growth of bacteria, increasing the risk of infection. Therefore, please observe the following hygiene instructions:

- **Each spacer must only be used at home by a single patient for hygiene reasons.**
- Make absolutely sure you also carry out cleaning and drying before using the device for the first time.
- For cleaning and disinfection, always use fresh tap water of drinking water quality.
- Make sure all components are dried properly whenever you have performed cleaning or disinfection.
- Do not keep the SPACER and accessories in a damp environment or together with damp objects.

## 2 PRODUCT DESCRIPTION

### 2.1 Components

Check that all components of your Zeffiro product are contained in your package. If anything is missing, please notify the dealer from whom you purchased the Zeffiro product immediately.

No.	Component	No.	Component
1(L).	Adult Mask	4.	Storage body
1(C).	Child Mask	5.	Back piece
2.	Dust cap	6.	Valves
3.	Mouthpiece	7.	Whistle

Check the "Assembly scheme" at page 2 to identify each component and its assembly way.

### 2.2 Intended purpose

The SPACER is a holding chamber, and is used in conjunction with medication sprays or "metered dose inhalers" in the treatment of diseases of the airways.

### 2.3 Intended use

Spacer for aerosol is intended for use by treating chronic obstructive pulmonary disease and asthma with medicine-assisted inhalation device.

The SPACER is a holding chamber, and is used in conjunction with medication sprays or "metered dose inhalers" in the treatment of diseases of the airways.

### 2.4 Combination of medical devices

The Spacer for aerosol should be used in conjunction with "metered dose inhalers" in the treatment of diseases of the airways. The metered dose inhalers are approximate rectangular,  $23\pm 2\text{mm}$  long  $\times$   $16.5\pm 2\text{mm}$  wide.

### 2.5 Contraindications

There are no contraindications known.

### 2.6 Product combinations

The Spacer for Aerosol has three specifications: 'C' is child mask, 'L' is adult mask, and 'N' is without mask. Patients can choose the appropriate specifications based on the size of the face.

### 2.7 Description of function

The SPACER helps to: minimize coordination mistakes when using metered dose inhalers avoid deposition of large quantities of medication in the oropharyngeal area, and the undesirable side effects associated therewith. With a flexible connecting adaptor it can be used with all conventional metered dose inhalers.

### 2.8 Material information

PETG/PP/ABS	Inhalation chamber, Mouthpiece
Silicone	Valve in the chamber, mask
Thermoplastic rubber	Back piece, dust cap

All materials are latex-free.

### 2.9 Shelf life

The shelf life of the Spacer for Aerosol is 3years. And the product is recommended be replaced after 12 months or be replaced after using 1460 times.

## 3 INHALATION

### 3.1 Function check

Clear pathway through the valve

After receiving the SPACER, check that the path through the valve in the mouthpiece is unobstructed.

If the valve is obstructed, the SPACER must not be used.

#### Position of the valve

**Before every use**, check that the valve is in the correct position inside the SPACER mouthpiece:

- Check that the valve is in the position shown inside the mouthpiece.

If necessary, use a stream of water to correct its position.

- Dry the SPACER completely before using it again.

### 3.2 Preparation for inhalation

#### WARNING

Since the inhalation chamber of the SPACER is not completely closed, small particles can get in, and might be breathed in during an inhalation session (choking hazard). Therefore, it is imperative to make sure there are no foreign bodies inside the SPACER before every use.

#### CAUTION

Check all product components and the accessories before each use. Replace any broken, misshapen or seriously discoloured parts. Please also follow the instructions for assembly given below. Damaged components and an incorrectly assembled holding chamber may impair the function of the holding chamber and thus the treatment as well.

#### CAUTION

If the a mask is being used, make sure that the expiratory valve shim is pressed out. If the valve shim is deflected inward, the dosage delivered may be insufficient when inhaling.

### 3.3 Performing the inhalation

- 1) Carefully examine the product for damage, missing parts or foreign objects (Fig.1). Remove any foreign objects prior to use.
- 2) Remove cap from the chamber and place the mask connector over chamber's mouthpiece until it fits securely (Fig. 2).

#### 3) CAUTION

Before the patient starts inhaling, ensure that all parts are firmly connected to each other. Otherwise, insufficient dosage may be delivered.

- 4) Remove cap from the MDI, shake MDI immediately before use (Fig. 3) and insert MDI into the backpiece of the chamber (Fig. 4).
- 5) Try to exhale, and hold mask against the face, press MDI and inhale slowly. Hold mask in place and breath in and out through the chamber for 5-6 breaths (Fig. 5). If there is not a mask, place lips around mouthpiece, press MDI and inhale slowly (Fig. 6). Slow down inhalation if a tone is heard. Hold breath for 5-10 seconds.

#### CAUTION

If you use a mask, make sure that the mask completely covers both corners of the mouth and the nose. Otherwise, the treatment can be less effective due to escaping aerosol, because of underdosage, for example. For possible side effects of escaping aerosol, please read the instructions for use of the medication in question

 **CAUTION**

Children younger than 18 months should inhale four to six times; children older than 18 months should inhale deeply two to four times. Otherwise, insufficient dosage may be delivered. For adults, one breath is enough to inhale the medication from the spacer.

6) If the doctor asks you to take medicine for more than one spray, please take a rest for one minute, and then repeat steps 4- 5.

#### **4 HYGIENIC PREPARATION**

Regularly clean all product components and the accessories used, if they are visibly dirty immediately after use. The chamber is usually cleaned once a week.

##### **4.1 Preparation**

- Take the metered dose inhaler out of the connecting ring on the inhalation chamber and close it with the protective cap provided for this purpose.
- If applicable, detach the mask from the mouthpiece.
- Disconnect the mouthpiece from the inhalation chamber (if applicable).

##### **4.2 Cleaning**

- Remove the back piece and disassemble the other parts of the spacer (Fig. 7).
- Soak both parts for 15 minutes in lukewarm water with liquid detergent. Agitate gently (Fig. 8).
- Rinse in clean water (Fig. 9).
- Shake out excess water. Do not rub dry (Fig. 10).
- Let air dry in vertical position (Fig. 11).
- Replace back piece when unit is completely dry and ready for use (Fig. 12).

##### **4.3 Visual inspection**

Inspect all product components after each cleaning and disinfection. Replace any broken, misshapen or seriously discoloured parts.

##### **4.4 Drying and storage**

- Dry all individual parts completely.
- Because of the antistatic properties of the inhalation chamber, you can dry them with a clean tea towel.
- Place the protective cap over the mouthpiece (if applicable). Info: The protective cap protects the mouthpiece from damage during transport and storage.
- Store the SPACER and its accessories in a dry, dust-free place (safely out of continuous

direct sunlight).


## 5 MISCELLANEOUS

### 5.1 Disposal

All product components can be disposed of with domestic waste unless this is prohibited by the disposal regulations prevailing in the respective member countries.

### 5.2 Explanation of symbols

The following symbols are displayed on the packaging of your Zeffiro spacer product and packaging:

	Manufacturer		Medical device
	Please follow the instructions for use.		Caution
	Medical CE marking ref. regulation 2017/745 EU and subsequent updates		Use-by date
	Lot no.		Date of manufacture
	Distributor		Unique Device Identifier
	Importer		Model / Reference
	Multiple use – single patient		Temperature limits
	Atmospheric pressure limits		Moisture limits
 Authorized representative in the European Community / European Union			

### 5.3 ENVIRONMENTAL CONDITIONS

#### Storage and Operating conditions:

Ambient temperature between +5°C and +40°C  
Relative air humidity between 15% and 95%  
Atmospheric pressure between 60KPa and 106KPa

#### Transport conditions:

Ambient temperature between -10°C and +40°C  
Relative air humidity  $\leq$  80%  
Atmospheric pressure between 60KPa and 106KPa

### 5.3 Contact

For all product information and in the event of defects or questions about usage, please contact your trusted dealer or an authorised FLAEM service centre nearest to you. You can find a list of all Service Centres at <http://www.flaemnuova.it/it/info/assistenza>

### 5.4 Notification of serious events

Serious events occurring in connection with this product must be reported immediately to the manufacturer or the competent authority.

An event is considered serious if it causes or may cause, directly or indirectly, death or an unforeseen serious deterioration in a person's state of health.

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**Taian Dalu Medical Instrument Co., Ltd**
**Lotus NL B.V.**
**Flaem Nuova S.p.A.**