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## Clinical Thermometer

ECOLOGICAL MERCURY FREE

**A NEW NONTOXIC THERMOMETER LIQUID.** An alloy of gallium, indium and tin is used as the thermometer liquid. It is completely non-toxic. The thermometer must be stored at a temperature no lower than  $-10^{\circ}\text{C}$  and no higher than  $+42^{\circ}\text{C}$ . Mercury-free gallium-in-glass clinical thermometers are analogue maximum thermometers for measuring body temperature. They are manufactured and tested in accordance with the requirements Standard EN 12470 Clinical thermometers. Part 1: Metallic liquid-in-glass thermometers with maximum device.

**APPLICATI ON.** Before taking a temperature, make sure that the indicator, liquid is below the first, numbered mark  $35^{\circ}\text{C}$ . Position the thermometer according to the type of measurement site (oral, axillary). After about 4 min the thermometer will have measured the body temperature and can be read. At the end of the measurement the thermometer liquid must be brought back down below the first numbered mark. Given the light weight of the liquid the new case helps the operations and the reading thanks to the lens integrated into the transparent case (see figure). The thermometer has a hygienic outer glass sheath and can be disinfected an unlimited number of times. The temperature of the cleaning and disinfectant solution must not exceed  $42^{\circ}\text{C}$ .

### TAKING TEMPERATURES

**Oral measurement (mouth).** Oral measurement is the method preferred by doctors and provides fast, accurate read-

ings. The thermometer bulb must be placed under the tongue on the left or right side, near one of the two back molars. DO NOT USE THIS METHOD FOR INFANTS OR YOUNG CHILDREN!

**Axillary measurement (underarm).** Underarm measurements are used fairly often. The thermometer bulb must be positioned in the armpit and held in place by exerting pressure with the arm.

#### SPECIAL RULES FOR TAKING TEMPERATURES

1. To avoid inaccurate readings, allow at least 4 minutes for the measurement.
2. Delay taking temperature for half an hour after showering or bathing.
3. If there has been a sudden drop or rise in temperature, retake the temperature more frequently.
4. Always shake the metallic liquid down to 35 °C or below before each use.
5. Clean the thermometer after every use by washing gently in lukewarm water and detergent. Rinse well.
6. Keep a written record of readings and times taken to give to your doctor.
7. The thermometer is not designed for use in premature babies or as ovulation thermometer!

**INFANTS AND CHILDREN.** Take the temperature under the arm or in the groin as explained above, and not in the mouth! Ensure that the temperature is held in position gently but firmly for at least 4 minutes. Infants and children should sit or lie down to ensure correct measuring. The manufacturer has verified the compliance of all thermometers with the prescribed margin of error of (1 /10) °C.

#### DISPOSAL

Do not dispose of in municipal waste but give platforms collection.

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Lotus Global Co., Ltd.
FLAEM NUOVA S.p.A.