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PulmoWaves 2

Device for respiratory rehabilitation physiotherapy

We are pleased you have purchased our product and we thank you for your trust in us. It is very important that the health care worker and/or the patient read and understand the device use and maintenance information.

INTENDED USE

Pulmowaves 2 is a medical device is for respiratory rehabilitation physiotherapy and in particular for the removal of tracheobronchial secretions, including a nebulising device on demand.

The Medical Device should be prepared and used in the presence of or by qualified personnel for respiratory rehabilitation physiotherapy who have read and understood this instruction manual provided, or adults who have received specific information from the aforementioned qualified personnel and who have read and understood this user manual provided.

The medical device can also be used for home-care treatment.

The therapeutic session may be set to (5-10-15-20 minutes) according to the doctor's assessment and the physical condition of the patient.

The device does not require lubrication, is manageable, simple to use, reliable, resistant and quiet.

Note: Use only original FLAEM accessories

OPERATING INSTRUCTIONS

- All drugs can be nebulised with the device.
- COPD chronic obstructive pulmonary disease (Stage 0 to Stage IV in accordance with the GOLD guidelines).
- Bronchiectasis.
- Cystic Fibrosis.
- Neuromuscular, orthopedic or other disorders involving the respiratory system (bronchial hypersecretion, ineffective cough) or who needs a quick effective removal of secretions.

CONTRAINDICATIONS FOR USE

- Children younger than 3 years old.
- Patients with cardiovascular instability.

DEVICE NOMENCLATURE

A Respiratory physiotherapy device

- A1 Switch
- A2 Digital pressure gauge
- A3 Flow regulator with filter
- A4 Pulsed air output
- A5 12Vdc output power supply socket for Nebulise
- A6 Pressure switch input
- A7 Filter tube
- A8 Filter
- A9 Pressure switch tube
- A10 Pulsed air tube
- A11 Power supply socket
- A12 Power supply cable
- A13 Tube connector clip
- A14 Protective gasket

B Nebuliser Device

- B1 Device switch
- B2 Nebuliser
 - B2.1 Lower part
 - B2.2 Upper part
 - B2.3 – Support with valve
- B3 Nebuliser power supply socket
- B4 Power supply cable for nebuliser device
- B5 Cable plug B4

C ACCESSORIES

- C1 PEP adjustment ring
- C2 Mouthpiece
- C3 PEP Fitting
- C4 Nose clip
- C5 Valve for use without nebuliser device

COMPOSITION

The device is a system composed of a pneumatic compressor, an ultrasonic nebulising device used to nebulise saline solution only, and a dispenser equipped with a PEP (Positive Expiratory Pressure) device.

FUNCTIONS OF THE MEDICAL DEVICE

The medical device has the following functions:

- PEP system (Positive Expiratory Pressure), which removes phlegm.
- Vibrations of the inhaled air that removes and loosens secretions.
- Aerosol therapy to humidify and treat the respiratory tract.

OPERATION

The operation is based on the generation of vibrations by the pneumatic compressor that is automatically activated by the inhalation of the patient. These vibrations, whose amplitude is adjustable via the flow regulator on the device console, enables the removal of mucus that clogs the airways. To avoid drying the airway during the treatment, the flow inhaled by the patient draws the aerosol of the saline solution, produced by the nebuliser.

During the treatment, the presence of a 5-position adjustable PEP device allows for the management of the exhalation phase of the patient, the flow intensity of which is approximately visible by the pressure gauge located in the central area of the console.

IMPORTANT WARNINGS

The manufacturer makes every effort to ensure that every product is of the highest quality and safety; however, as for any electrical device, basic safety regulations must always be followed in order to avoid harming persons and things.

- Warning: The device should be carried in the bag provided and enclosed in its protective bag
- The patient is also considered as the operator.
- Before using the product for the first time, and periodically during its lifetime, check the integrity of the device structure and of the power cable to make sure there is no damage. In the event of damage, do not plug in the cable and immediately take the product to an authorised FLAEM service centre or to your local dealer.
- The average life span of the accessories is approximately 6 months. The filter (A8) is for single patient use and must be replaced within the period of use indicated on the packaging.
- Children and people who are not self-sufficient must use the device under the close supervision of an adult who has read this manual.
- Some parts of the device are small enough to be swallowed by children; therefore, keep the device out of the reach of children.
- Do not use the provided tubing and cables for any other purpose than those specified, as they can cause risk of strangulation. Be particularly attentive with children and individuals with disabilities because they are often unable to correctly assess risk.
- The apparatus is unsuitable for use in presence of flammable anaesthetic mixture with air, oxygen or nitrous oxide.
- Always keep the power supply cable away from hot surfaces.
- Always keep the power cable away from animals (for example, rodents) which could damage the insulation.
- Do not handle the plug with wet hands. Do not use the device in damp conditions (for example when taking a bath or shower). Do not immerse the device in water; in the event of immersion immediately disconnect the plug. Do not remove or touch the immersed device; unplug the power cable first. Immediately bring the device to an authorised FLAEM service centre or to your trusted dealer.
- Use the device only in dust-free conditions, otherwise treatment could be compromised.
- The device casing is not waterproof. Do not wash the device under running water or by immersion and keep it safe from being sprayed by water or other liquids.
- Do not expose the device to particularly extreme temperatures.
- Do not place the device near sources of heat, in direct sunlight or in excessively hot rooms.
- Do not obstruct or put objects into the filter or its related housing in the device.
- Never obstruct the air inlets positioned on both sides of the device.
- Always use it on a rigid surface that is clear of obstacles.
- Make sure there is no material obstructing the air vents before each use.
- Do not insert any objects into the ventilation slits.
- Always unplug the power cable immediately after use.
- The device is equipped with an easy to inspect safety fuse in case of malfunction. Before doing so, unplug the power cable. Repairs must be done only by authorised FLAEM personnel. Unauthorised repairs void the warranty and may be hazardous for the user.
- The average duration of use intended for the device is 1000 h for the compressor and 400 h for the P0112EM nebuliser.

ATTENTION: Do not modify this device without authorisation from the manufacturer.


- The Manufacturer, the Vendor and the Importer shall be held responsible for safety, reliability and performance only if: a) the device is used in compliance with the instructions for use; b) the wiring where the device is being used is in compliance with safety regulations and current laws.

- The Manufacturer must be contacted to report problems and/or unexpected events related to operation.
- No maintenance and/or repairs must be carried out while the device is in use.
- The Manufacturer, the Vendor and the Importer shall be held responsible for safety, reliability and performance only if: a) the device is used in compliance with the instructions for use; b) the wiring where the device is being used is in compliance with safety regulations and current laws.
- The Manufacturer must be contacted to report problems and/or unexpected events related to operation.
- The manufacturer must be contacted about any problems and/or unexpected events concerning operation and for any clarifications on use, maintenance or cleaning.
- The amount of time required to switch from preservation to operating conditions is approximately 2 hours
- The components (A7, A9, A10, C1, C2, C3, C4, C5) are for personal use and should be cleaned before each use, as described in the "Cleaning, sanitation, disinfecting and sterilising" section.
- Interactions: the materials used in contact with medications have been tested with a vast range of medications. However, in view of the variety and continuous evolution of medications, interactions cannot be ruled out. We recommend the medication is used as soon as possible after opening and its prolonged exposure to the nebuliser is avoided. The materials used in the medical device are biocompatible in accordance with the provisions of Directive 93/42 EC as amended. However, the possibility of occurrence of allergic reactions cannot be entirely ruled out.
- Should your device fail to perform as expected, contact the authorised service centre for clarifications.
- The supplied cable is fitted with a protective gasket (A14) to prevent the penetration of liquids into the device. Do not separate this gasket from the cable.

INSTRUCTIONS FOR USE

Before each utilisation, wash your hands accurately and clean the device as described in the paragraph "CLEANING, SANITISATION, DISINFECTATION AND STERILISATION (for hospital use)". The parts (A7-A9-A10-C2-C3) are for personal use only to avoid potential risks of transmitting infection. It is recommended to replace them every six months.

1. Insert power supply cord (A12) into socket (A11) and then connect it to a mains outlet adequate to the device voltage. The position of the socket must be such that the device can be easily unplugged from the mains network.


Tube assembly

- Connect tube (A10) to the pulsed air outlet (A4) and tube (A7) to the pressure switch inlet (A6).
- Insert the plug (B5) of the cable (B4) into the socket (A5).
- Screw the ring nut by turning it clockwise.
- Connect the filter (A8) to the tube (A7).

WARNING


The protection cap of the 12Vdc socket must be removed to insert the plug and must be replaced after use to protect the socket.

Mouthpiece


- Connect mouthpiece (C2) to PEP connector (C3).
- Connect the two tubes (A9-A10) coming out of (A4-A8) to PEP connector C3.

Device connection


- Insert the coaxial cable plug (B4) coming out of (A5) into socket (B3) of the nebuliser device.
- Insert the protective gasket (A14) into the power supply cable plug (A12).
- Place the plug on power supply cable (A12) into socket (A11) on the device A.

Nebuliser device


- Open the upper part (B2.2) of the nebuliser (B2) counterclockwise and pour the saline solution into (B2.1).
- Apply the valve support (B2.3) to the lower part (B2.2).
- Close the nebuliser (B2) of the Nebuliser Device B and connect the connector (C3) and mouthpiece (C2) as shown in the connection diagram.

Use without nebuliser device


It is possible to use the Pulmowaves device without the nebuliser, by closing the PEP fitting (C3) and with the valve (C5) as shown in the picture.

Device activation


- Turn on the device by pressing power switch A1, then press button B1 of the nebuliser.
- The compressor will start only with the patient's inhalation through mouthpiece C2.

Flow regulator


Adjust the ring nut **C1** of the connector **C3** and the flow regulator **A3** of the device for the desired respiratory physiotherapy according to the values reported in the following table.

Connector adjustment ring PEP	Flow regulator on MAX		Flow regulator set on MIN	
	Diam. Hole	Val. Flow	Val. Press	Val. Press
Ø 2	7,5 l/min	7 mbar	1,2 l/min	2 mbar
Ø 3	6,4 l/min	6 mbar	1,1 l/min	1 mbar
Ø 4	3,9 l/min	3 mbar	1,0 l/min	0 mbar
Ø 5	3,3 l/min	2 mbar	0,8 l/min	0 mbar
Ø 6	3,2 l/min	1 mbar	0,6 l/min	0 mbar

Digital pressure gauge


- During the patient's exhalation phase, the pressure gauge (A2) displays the positive pressure generated according to the hole selected on the ring nut (C1) of the PEP connector (C3).
- The pressure gauge is only designed to indicate the correct patient exhalation. It is not an indicator of the precise value of exhalation flow.

Use the "+" and "-" keys to set the therapy duration. The therapy may last 5, 10, 15 or 20 minutes. At the end of the set therapy, the 4 LEDs flash and an acoustic warning beeps 5 times.


Tube connector clip

Attach tubes (A9-A10) and cable B4 to the proper hooks of the Clip A13 as shown in the detail.

Nose clip

- For the correct use of the PEP connector, the patient must put on the nose clip (C4), as shown in the figure.

CLEANING, SANITISATION, DISINFECTION AND STERILISATION

Witch off the device before any cleaning procedures and always unplug the power cable from the socket.
DEVICE (A) and (B)

Use only a damp cloth with antibacterial soap (non-abrasive and with no solvents of any sort).

ACCESSORIES

Accessories (A8-C2-C3-C4-C5-B2) and connection tubes (A7-A9-A10) are for personal, single-patient use to avoid any risks of transmitting infection.

Open the nebuliser (B2) by turning the upper part anticlockwise, remove the upper part as shown in the figure.

NOTE: For sanitisation, disinfection and sterilisation operations the PEP C3 connector must not be disassembled.

ANTI-BACTERIAL FILTER

WARNING

- The filter (A8) is for single patient use and must be replaced as indicated on the package label.
- It cannot be washed, disinfected or sterilised and is not recyclable.

CLEANING AT HOME - SANITISATION AND DISINFECTION

SANITISATION

Before and after each use, sanitise the accessories and the nebuliser by choosing one of the methods described below.

method A: Sanitise accessories (C1-C2-C3-C4-C5-A7-A9-A10-B2) under potable hot water (approximately 40°C) with a gentle, non abrasive dish detergent.

method B: Sanitise accessories (C1-C2-C3-C4-C5) in the dishwasher with a hot cycle.

method C: Sanitise accessories (C1-C2-C3-C4-C5-A7-A9-A10-B2) by immersing them in a solution of 50% water and 50% white vinegar, then rinse thoroughly under potable hot water (approximately 40 °C).

If you want to also perform the cleaning for DISINFECTION, jump to the DISINFECTION paragraph.

DISINFECTION

After sanitising the accessories, disinfect them choosing one of the methods described below

method A: Accessories (B2-C1-C2-C3-C4-C5-A7-A9-A10) can be disinfected.

The disinfectant must be an electrolytic chloroxidizer (active principle: sodium hypochlorite) specific for disinfecting, which is available in any pharmacy.

Implementation:

- Fill a container big enough to hold all of the parts to disinfect with a solution of potable water and disinfectant, according to the proportions indicated on the packaging of the disinfectant.
- Completely immerse each part in the solution, taking care to avoid the formation of air bubbles on the parts. Leave the parts immersed for the amount of time indicated on the packaging of the disinfectant associated with the concentration chosen for the solution.
- Remove the disinfected parts and rinse abundantly with warm potable water.
- Dispose of the solution following the instructions provided by the disinfectant manufacturer.

method B: Disinfect the accessories (C1-C2-C3-C4-C5-A7-A9-A10-B2) by boiling them in water for 10 minutes; use demineralised or distilled water to prevent calcium deposits.

method C: Disinfect the accessories (C1-C2-C3-C4-C5-A7-A9-A10-B2) with a hot steam steriliser for baby-bottle (not the microwave type). Perform the process faithfully following the instructions of the steriliser. To ensure that the disinfection is effective, choose a steriliser with an operating cycle of at least 6 minutes.

After the CLEANING AND / OR DISINFECTION operations dry the area of the transducer, without exerting unnecessary pressure on the transducer itself, dry the accessories and the outside of the appliance with a paper towel or even better with a jet of hot air (e.g. hair dryer). Finally, check that the “power plug” area of the device is completely dry.

CLEANING IN A CLINICAL OR HOSPITAL SETTING - DISINFECTION AND STERILISATION

Before disinfection or sterilisation, sanitise the nebuliser cup and the accessories, choosing one of the methods described below.

method A: Sanitise accessories (C1-C2-C3-C4-C5-A7-A9-A10-B2) under potable hot water (approximately 40°C) with a gentle, non abrasive dish detergent

method B: Sanitise accessories (C1-C2-C3-C4-C5) in the dishwasher with a hot cycle.

DISINFECTION

Accessories (C1-C2-C3-C4-C5-A7-A9-A10-B2) can be disinfected.

The disinfectant must be an electrolytic chloroxidizer (active principle: sodium hypochlorite) specific for disinfecting, which is available in any pharmacy.

Implementation:

- Fill a container big enough to hold all of the parts to disinfect with a solution of potable water and disinfectant, according to the proportions indicated on the packaging of the disinfectant.
- Completely immerse each part in the solution, taking care to avoid the formation of air bubbles on the parts. Leave the parts immersed for the amount of time indicated on the packaging of the disinfectant associated with the concentration chosen for the solution.
- Remove the disinfected parts and rinse abundantly with warm potable water.
- Dispose of the solution following the instructions provided by the disinfectant manufacturer.

If you want to also perform the STERILISATION, jump to the STERILISATION paragraph.

After the CLEANING AND / OR DISINFECTION operations dry the area of the transducer, without exerting unnecessary pressure on the transducer itself, dry the accessories and the outside of the appliance with a paper towel or even better with a jet of hot air (e.g. hair dryer). Finally, check that the “power plug” area of the device is completely dry.

STERILISATION

Accessories (C1-C2-C3-C5-A7-A9-A10-B2 Assembled) can be sterilised.

Device: Fractionated vacuum overpressure steam steriliser in accordance with EN 13060.


Implementation: Wrap every single part to be treated with a sterile barrier system or packaging in accordance with Norm EN 11607. Place the packed components in the steam steriliser. Run the sterilisation cycle according to the operating instructions of the device by selecting a temperature of 134 °C and a time of 10 minutes first.

For the SILICONE TUBING (A7-A9-A10): temperature of 134°C and a duration of 10 minutes for a maximum of 30 times.

Storage: Store the sterilised parts as per the instructions for use of either the sterile barrier system or packaging.

The sterilisation procedure is validated in its conformity to ISO 17665-1.

AIR FILTERING

The device is fitted with a suction filter, located inside the flow regulator (A3), which must be replaced when it is dirty or discoloured. To replace it, lift the flow regulator (A3) and extract the Filter.

Do not wash or reuse the same Filter. Regular replacement is necessary to assure proper compressor performance. The Filter must be regularly checked. Contact your retailer or authorised service centre for replacement filters.

SPARE PARTS

Description of parts	KIT code				
	ACO662P Six-month pack consumables	ACO663P Patient pack individ. cons.	ACO664P Outpatient pack cons.	ACO665P Outpatient pack daily	ACO077B Patient individ.
Set of tubes shq L.800, L.690, L.90	1	1	-	3	-
Anti-bacterial filter	24	1	1	-	-
Mouthpiece	1	1	1	-	-
PEP dispenser	1	1	1	-	-
Nose clip	1	1	1	-	-
Clip for external tubes	5	5	-	15	-
Manual	1	1	1	1	-
Porex filter for compressor	1	-	-	-	-
RF8 Plus	-	-	-	-	1

Use only with original Flaem accessories.

DEVICE SPECIFICATIONS

Compressor Model P0112EM F1000		Ultrasonic Nebuliser Device Model P0112EM	
Power source:	230 V ~ 50 Hz / 140 VA	Power supply	12 V
Fuse	1x T2AL 250V	Nebulisation ml/min ⁽¹⁾	0.32 ml/min'
Noise (at 1 m):	55 dB (A) (approx.)	MMAD ⁽²⁾	4.26 µm
		Respirable fraction <5µm ⁽²⁾	51.7 % (approx.)
		Medication cup capacity	8 ml
		⁽¹⁾ Data recorded according to internal Flaem procedure I29-P07.5.	
		⁽²⁾ Data calculated on the basis of the values obtained with Malvern Spraytec, and compared with the values of the tests performed at the University of Parma.	
Appliance dimensions	22 (L) x 11 (P) x 23 (H) cm	Weight	2.200 Kg
Bag dimensions	27 (L) x 17 (P) x 25 (H) cm	guarantee	2 year

Operating conditions:

Min temperature 5°C; max. 40 °C

Air humidity min. 10%; max. 95%

Storage and transport conditions:

Min temperature -25°C; max. 70°C

Air humidity min. 10%; max. 95%

Operating / Storage Atmospheric Pressure:

min. 69 KPa; max. 106 KPa

APPLIED PARTS

Type BF applied parts are:
patient accessories (C2)

SYMBOLS

CE Marking medical ref. Dir EEC 93/42 and subsequent updates

Class II device

Attention: check the instructions for use

"ON"

When the device is turned off, the switch interrupts compressor operation only on one of the two power supply phases

"OFF"

In compliance with: European Standard EN 10993-1 "Biological Evaluation of medical devices" and the European Directive 93/42/EEC "Medical Devices". Phthalate-free. In compliance with: Reg. (CE) no. 1907/2006

Output power supply for Nebulizer Device. (Mod. P0112EM).

Temperature limit

Device serial number

Manufacturer

Type BF applied part

IP21

Degree of protection of the casing: IP21. (Protected against solid foreign objects larger than 12mm. Protected against access with a finger; Protected against vertically dripping water).

Keep dry

Alternating current

Humidity limit

Atmospheric pressure limit

PACKAGING DISPOSAL

Heat shrink film accessories

Product packaging bag and tube packaging bag

DEVICE DISPOSAL

In compliance with the Directive 2012/19/EC, the symbol printed on the device shows that the device to be disposed of is considered waste and must therefore be a "separate collection" item. Consequently, the user must take it (or have it taken) to the differentiated collection sites provided by the local authorities, or turn it in to the dealer when purchasing an equivalent new device. Differentiated waste collection and the subsequent treatment, recycling and disposal procedures promote the production of devices made with recycled materials and limit the negative effects on the environment and on health caused by potential improper waste management. If the user disposes of the product in an unauthorised fashion, administrative penalties will be applied according to the transposition laws of directive 2012/19/EC of the member state or country where the product will be disposed of.

ELECTROMAGNETIC COMPATIBILITY

This device was designed to satisfy the currently required requisites for electromagnetic compatibility (EN 60601-1-2:2007). Electro-medical devices require particular care during installation and use relative to EMC requirements. Users are therefore requested to install and/or use these devices following the manufacturer's specifications. There is a risk of potential electromagnetic interference with other devices, in particular with other analysis and treatment devices. RF mobile or portable radio and telecommunications devices (mobile telephones or wireless connections) can interfere with the functioning of electro-medical devices. For further information visit our website www.flaemnuova.it or www.flaem.it

Flaem reserves the right to make technical and functional modifications to the product with no prior warning.

FAULT IDENTIFICATION

Switch off the device before any procedure and unplug the power cable from the socket.

PROBLEM	CAUSE	SOLUTION
The device does not work	The power cable (A12) has not been inserted properly into the device power socket (A11) or into the mains socket.	Insert the power cable correctly into the socket (A11) and into the relative mains socket.
The nebuliser (B2) does not nebulise and/or dispenses the drug correctly.	The switch (B1) on the device (B) is set to "OFF". The cable connector (B4) has not been inserted into the 12VDC socket (A5). The device power cable (B4) has not been inserted.	Set the switch (A1) of the device (A) to "ON". Insert the cable connector (B4) into the socket (A5). Insert the device cable (B4).
	The medication or saline solution has not been inserted into the nebuliser (B2) (bowl).	Pour the correct amount of medication or saline solution into the nebuliser (B2).
	The nebuliser (B2) has not been properly fitted	Disassemble the nebuliser (B2) and reassemble it correctly as shown in the diagram on the cover.
	The nebuliser nozzle (B2) is clogged	Disassemble the nebuliser (B2), remove the nozzle and perform cleaning operations. Failure to remove medication deposits from the nebuliser affects its efficiency and proper operation. Strictly comply with the instructions contained in the cleaning, sanitisation and disinfection chapter.
The digital pressure gauge (C) does not light up the LEDs and/or does not work properly.	The power cable (A12) has not been inserted properly into the device power socket (B4) or into the mains socket.	Insert the power cable correctly into the socket (A11) and into the relative mains socket.
	The air tubes (A9) (A10) are not properly connected to the unit (A) and/or the PEP connector (C3).	Check that the tubes are correctly connected to the relative sockets of the device (A). See the connection diagram on the cover.
	Air tubes (A9) (A10) are bent and/or damaged and/or kinked. The ring nut (C1) has not been positioned correctly.	Unroll the tubes (A9) (A10) and check that they are not crushed, punctured or bent, if so, replace them. Adjust the PEP ring nut (C1) correctly by matching the holes.
	The flow regulator filter (A3) is dirty.	Replace the filter (A3).
	Anti-bacterial filter (A8) clogged or dirty, bent tube (A9) and/or faulty circuit.	Restore the tube (A9), replace the filter (A8) and/or contact an authorised service centre to have the electronic circuit replaced.
	The electronic circuit does not work	Contact a FLAEM service centre to replace the electronic circuit.

When I inhale, the compressor does not activate and no pulsed air comes out of the mouthpiece.	The power cable (A12) has not been inserted properly into the device power socket (A11) or into the mains socket.	Insert the power cable correctly into the socket (A11) and into the relative mains socket.
	The switch (A1) is set to "OFF". The pressure switch is faulty.	Set switch (A1) to "ON" or contact a FLAEM service centre to replace the pressure switch.
	The air tubes (A9) (A10) are not properly connected to the unit (A) and/or the PEP connector (C3)	Check that the tubes are correctly connected to the relative sockets of the device (A). See the connection diagram on the cover.
	The flow regulator filter (A3) is dirty.	Replace the filter (A3).
The device (A) is very noisy.	The pressure regulator (A3) is not correctly inserted in its housing	Properly insert the pressure regulator (A3) in its housing.

If, after checking the aforementioned conditions, the device still does not work properly, we recommend that you contact your local dealer or an authorised FLAEM service centre.

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