
# FLAEM<sup>®</sup> Wi.Neb

Medical Devices

IT > MANUALE ISTRUZIONI D'USO - EN > INSTRUCTIONS FOR USE

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# Wi.Neb

Mod. P0209EM

## Portable and Ultra-Compact Aerosol Therapy Device

Congratulations on your purchase of Up-Wi.Neb Go and thank you for confidence in us. Our goal is to fully satisfy our customers by offering them cutting-edge systems for treating ailments of the respiratory system. Please read these instructions carefully and keep them handy for future reference. Use the device only as described in this manual. This is a home medical device designed to administer nebulizer medications prescribed or recommended by your physician upon assessing the patient's general conditions. Please visit [www.flaem.it](http://www.flaem.it) to see the entire range of Flaem products.


### IMPORTANT SAFEGUARDS


- This device is also intended for direct use by the patient
- Before using the device for the first time, and periodically throughout its life, check the power supply adaptor cable for any damage; if you detect any damage, do not plug it in. Take the device immediately to an authorized FLAEM service center or to your local dealer.
- Should your device fail to provide the expected performance, contact the authorised service centre for clarifications.
- The average duration provided for compressor is: 400 Hours.
- Children and persons with reduced physical, sensory, or mental capabilities may only use this device under the strict supervision of an adult who has read this manual.
- Several parts of the device are small enough to be swallowed by children; therefore keep the device out of reach of children.
- Do not use the tubes and the cords of the device for purposes other than those specified as they may cause strangulation hazard; utmost attention must be paid to children and persons with special needs; usually, such persons are not capable of assessing the hazard correctly.
- This device is not suited for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrogen protoxide.
- The power adaptor of your Wi-Neb can be plugged into different types of wall outlets worldwide by using the appropriate adaptors. Always keep the power adapter and the power supply cord away from hot surfaces.
- Keep the power cable away from animals (for example, rodents) which could damage the insulation.
- Do not handle the power adaptor with wet hands. Do not use the device in humid environments (e.g. while taking a bath or shower). Do not immerse the device in water; should this accidentally occur, unplug the device immediately. Do not remove or touch the device immersed in water before pulling out the plug. Take the device immediately to an authorized FLAEM service center or to your local dealer.
- Use the device only in dust-free conditions, otherwise treatment could be compromised.
- Do not wash the device under running water or submerge in water. Keep away from splashing water or other liquids.
- Do not expose the device to high temperatures.
- Do not put the device and power adaptor near heat sources, direct sunlight, or in hot environments.
- Do not obstruct or insert any objects in the filter or in its housing in the device.
- Repairs must be performed by authorized FLAEM personnel only, by complying with the information provided by the manufacturer. Any unauthorized repairs will void the warranty and may pose a safety hazard for the user.
- The Manufacturer, the Retailer and the Importer are responsible for the safety, reliability and performance of the device only if:
  - a) the device is used in accordance with the operating instructions
  - b) the electrical system of the premises where the device is used is in compliance with current laws.
- Interactions: The materials used that come into contact with medications have been tested with a wide range of medications. However, due to the variety and continuous evolution of medication, the possibility of a chemical interaction cannot be excluded. We recommend that once the medication has been opened you should use it as soon as possible and avoid prolonged contact of the medication with the nebulizer. The materials used in the device are biocompatible in accordance with the provisions of Directive 93/42 EC and subsequent amendments. However, the possibility of occurrence of allergic reactions cannot be entirely excluded.
- The manufacturer must be contacted about any problems and/or unexpected events concerning operation and for any clarifications on use, maintenance/cleaning.
- The amount of time required between preservation conditions and use is approximately 2 hours.

## OPERATION WITH EXTERNAL POWER SUPPLY

### OPERATION WITH POWER ADAPTOR


Plug power adaptor in the electrical receptacle. The socket must be placed so as to allow for convenient disconnection from the power mains.


Insert the "L" jack of the power adaptor cable in the DC plug of the device.

### OPERATION WITH CIGARETTE LIGHTER RECEPTACLE


Insert the cigarette lighter plug in the 12V DC cigarette lighter receptacle of the vehicle/ camper/ boat. Check the User Manual of the vehicle to make sure that the voltage supplied is 12V DC.

Insert the "L" jack of the cigarette lighter cord in the DC plug of the device.

### HOW TO CONNECT THE ACCESSORIES

Prior to each use, carefully wash your hands and clean your Wi Neb according to the instructions provided in the “CLEANING, SANITIZATION AND DISINFECTION” section. During the application you should use adequate protection against any drips. The average lifespan of accessories is 1 year; however, we recommend you replace the nebulizer every 6 months in intensive use (or earlier, if the nebulizer is blocked) to always ensure maximum therapeutic effect. The nebulizer and accessories are for personal use only to prevent possible risk of infection

due to contamination.


**Use original FLAEM accessories only.**

This device is suitable for the administration of medical substances and not, for which the administration via aerosol is foreseen; these substances are to be in any case prescribed by the Doctor. In case of too thick substances, the dilution with a suitable physiological solution could be needed, according to the medical prescription.

Connect one end of the connection hose (E) to the device.

Connect the other end of the hose (E) to the L-fitting (F).

Connect the L-fitting (F) to the nebulizer (G).

Insert the nebulizer (G) in its housing on the device (A).


Open the nebulizer by turning the lid anticlockwise.

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Close the nebulizer by simply turning the lid clockwise.


Connect the accessory desired.

If you use the face mask, fit it to the face as shown in the picture.


RF7 Dual Speed Plus nebulizer provides 2 different ways to administer medication; consult your doctor for information on the most appropriate method to suit your needs.

To speed up the inhalation therapy, press on the word MAX of the speed selector button.

For a more effective inhalation therapy, press on the opposite side of the word Max located on the speed selector button; this allows you to inhale an optimal amount of the medication and reduces dispersion of the medication into the surrounding environment to a minimum thanks to the valve system provided in the nebulizer, mouthpiece, and mask.


Turn on the device by pressing the on/off switch to the ( I ) position; the blue LED will turn on.

## CLEANING SANITISATION DISINFECTION

Switch off the device before any cleaning procedure and unplug the power cable from the socket.

### DEVICE AND TUBING EXTERIOR

Use only a damp cloth with antibacterial soap (non-abrasive and with no solvents of any sort).

### ACCESSORIES

Open the nebuliser by turning the upper part (G1) anticlockwise, remove the nozzle (G2) from the upper part (G1) by pressing as shown by the 2 arrows in the "Connection diagram" in point G.

Then proceed according to the following instructions.

### SANITISATION

Before and after each use, sanitise the nebuliser cup and the accessories, choosing one of the methods described below.

**method A:** Sanitise accessories **F-G.1-G.2-G.3-G.4-H-I-L-M** under potable hot water (approximately 40°C) with a gentle, non abrasive dish detergent.

**method B:** Sanitise accessories **F-G.1-G.2-G.3-G.4-H-I-L-M** in the dishwasher with a hot cycle.

**method C:** Sanitise accessories **F-G.1-G.2-G.3-G.4-H-I-L-M** by immersing them in a solution of 50% water and 50% white vinegar, then rinse thoroughly under potable hot water (approximately 40 °C).

If you want to also perform the cleaning for DISINFECTION, jump to the DISINFECTION paragraph.

### DISINFECTION

After sanitising the nebuliser cup and the accessories, disinfect them choosing one of the methods described below.

**method A:** Accessories **G.1-G.2-G.3-G.4-H-I-L-M** can be disinfected.

The disinfectant must be an electrolytic chloroxidizer (active principle: sodium hypochlorite) specific for disinfecting, which is available in any pharmacy.


#### Implementation:

- Fill a container big enough to hold all of the parts to disinfect with a solution of potable water and disinfectant, according to the proportions indicated on the packaging of the disinfectant.
- Completely immerse each part in the solution, taking care to avoid the formation of air bubbles on the parts. Leave the parts immersed for the amount of time indicated on the packaging of the disinfectant associated with the concentration chosen for the solution.
- Remove the disinfected parts and rinse abundantly with warm potable water.
- Dispose of the solution following the instructions provided by the disinfectant manufacturer.

**method B:** Sanitise the accessories **G.1-G.2-G.3-G.4-H-I-L-M** by boiling them in water for 10 minutes; use demineralised or distilled water to prevent calcium deposits.

**method C:** Sanitise the accessories **G.1-G.2-G.3-G.4-H-I-L-M** with a hot steam steriliser for baby-bottle (not the microwave type). Perform the process faithfully following the instructions of the steriliser. To ensure that the disinfection is effective, choose a steriliser with an operating cycle of at least 6 minutes.

At the end of each use store the device complete with accessories in a dry place away from dust.


Drying : After sanitization/disinfection shake the accessories vigorously and place them on a paper towel, or alternatively, dry them with a jet of hot air (e.g. using a hair dryer).


### FILTER REPLACEMENT

- Under normal conditions of use, the filter should be replaced after every 100 applications, or at least once a year.
- You must check the filter on a regular basis (every 10-15 applications).
- Replace the filter if it is dirty, a brownish or grayish color, or clogged.
- Replace the filter if it is damp or wet.
- Do not wash or dry the filter.
- Do not replace the filter during use.

**Only use original accessories and spare parts by Flaem, we disclaim any liability in the event of using non original spare parts or accessories**

To remove the "AIR FILTER" cover, insert a small coin in the rear slot and gently pry open.

Remove the dirty filter using household tweezers and insert a new filter making sure to position it all the way back in its seat.

Replace the "AIR FILTER" cover by 1 positioning the front part in the seat and 2 gently pressing on the rear part until it snaps securely into place.

## TROUBLESHOOTING

<p><b><u>If the device does not work</u></b> when operated using power adaptor or cigarette lighter cable)</p>	<ul style="list-style-type: none"> <li>• Check that the power adaptor plug is properly plugged in. Check that the electrical receptacle functions by plugging in a different device to see if it works.</li> <li>• Check if the power cord and/or plug are faulty or damaged. In this case do not use the device.</li> <li>• Check the plug for connecting the cable to the 12 V DC cigarette lighter receptacle is properly inserted and that the fuse, found in the connection plug is not blown.</li> </ul>
<p><b><u>If the device provides low mist</u></b></p>	<ul style="list-style-type: none"> <li>• Check that the medicine is not too cold; it should be close to room temperature.</li> <li>• Check that there is at least 2ml of medicine in the nebulizer.</li> </ul>
<p><b><u>If the device provides no mist</u></b></p>	<ul style="list-style-type: none"> <li>• Check for medicine in the nebulizer.</li> <li>• Check that the nebulizer nozzle is not clogged.</li> <li>• Check that the air inlet of the device and the accessories are properly connected (see connection diagram in fig. 1).</li> <li>• Make sure that the device is connected to the power adaptor and that the power adaptor is connected to the receptacle.</li> </ul>
<p><b><u>The nebulizer seems inefficient</u></b></p>	<ul style="list-style-type: none"> <li>• Failure to clean the nebulizer from medicine deposits compromises the efficiency and functioning of the device.</li> <li>• Carefully read and follow the instructions included in the <b>CLEANING, SANITIZATION AND DISINFECTION</b> section.</li> </ul>

If after checking the conditions listed above the nebulizing function of the device still does not work contact your local dealer or authorized service center for assistance.

## SPARE PARTS AVAILABLE

Description	Code
RF7 PLUS nebulizer kit composed of: nebulizer, mouthpiece with valve and nasal prong	ACO436P
1m connection tube	ACO35
Adult mask	ACO437P
Pediatric mask	ACO438P
Air filter replacement kit for device (2 pcs)	ACO439P
Multi-volt universal power adaptor	ACO458P
	ACO460P with USA plug
	ACO214P
Power cord for 12V DC cigarette lighter receptacle	ACO 211P

## TECHNICAL SPECIFICATIONS

Device mod. P0209EM	
Power supply:	12V DC
Absorption:	1.25 A
Max Pressure:	2 ± 0,3 bar
Compressor air output:	9 l/min approx
Noise level (at 1 m):	approx. 52 dB (A)
Operation:	30 minutes ON – 30 minutes OFF
Operating conditions:	Temperature: min 10°C; max 40°C Air humidity: min 10%; max 95% Atmospheric pressure: min 69KPa; max 106KPa
Storage conditions:	Temperature: min -25°C; max 70°C Air humidity: min 10%; max 95% Atmospheric pressure: min 69KPa; max 106KPa
In conformity with Directive 93/42 EEC	
Safety approvals:	
Dimensions (W)x(D)x(H)	13x9.5x5.5 cm
Weight	0.55 Kg

## APPLIED PARTS


Type BF applied parts are:	patient accessories (H, I, L, M)
----------------------------	----------------------------------

<b>RF7 PLUS Nebulizer</b>	
Medication minimum capacity:	2 ml
Medication maximum capacity:	8 ml
Operating pressure (with nebulizer):	0.75 bar
(1) Delivery:	Min With valve: approx. 0.18 ml/min
	Max: approx. 0.42 ml/min
(2) MMAD:	Min With valve: 2.95 $\mu\text{m}$
	Max: 2.44 $\mu\text{m}$
(2) Breathable fractions < 5 $\mu\text{m}$ (FPF):	Min With valve: 74.7%
	Max: 80.7%
(1) Data detected according to Flaem I29-P07.5 internal procedure. (2) In vitro testing performed by TÜV Rheinland Italia S.r.l. in collaboration with the University of Parma and in compliance with the EN 13544-1: 2007 + A1 European Standard for aerosol therapy devices Standard. More details are available on request.	
Use only power adaptor mod. UME318-1215 (cod. 13014)	
Primary:	100/240V~ 50/60Hz
Secondary:	12V DC 1.5 A
<b>Universal power adaptor (USA alternative)</b>	
Use only power adaptor mod. UML318-1215 (cod. 13373)	
Primary:	100/240V~ 50/60Hz
Secondary:	12V DC 1.5 A
<b>Power cord for 12V DC cigarette lighter receptacle</b>	
If the fuse inside the connection plug blows, replace with 3A.	
<b>Disposal of device</b>	
 In conformity with Directive 2012/19/CE, the symbol shown on the device to be disposed of indicates that it is considered as waste and is therefore subject to "sorted waste collection". The user must therefore take (or have taken) the above waste to a pre-sorted waste collection centre set up by the local authorities, or else give it back to the dealer when purchasing a new appliance of the same type. Pre-sorted waste collection and the subsequent treatment, recovery and disposal operations favor the production of appliances made of recycled materials and limit the negative effects of any incorrect waste management on the environment and public health. The unlawful disposal of the product by the user could result in administrative fines as provided by art. 50 and subsequent amendments of Legislative Decree no. 22/1997.	

SYMBOLS					
	Class II device		Manufacturer		Temperature limit
	Type BF applied part		Alternate current		Atmospheric pressure limit
	Important: check the operating instructions		Power adaptor for household use		Humidity limit
	Direct current		Plug for safety extra low voltage		<b>DISPOSAL PACKAGING</b>
	Switch "ON"		CE mark for medical devices ref. Dir. 93/42 EEC and subsequent amendments		PBox product and power supply
	Switch "OFF"				TÜV approval
	Outlet for safety extra low voltage		Serial number		Heat shrink film accessories
	In compliance with European standard EN 10993-1 Biological evaluation of medical devices and 93/42/EEC Medical Devices Directive	<b>IP21</b>	<b>Enclosure protection rating: IP21.</b> (Protected against solid bodies over 12 mm. Protected against access with a finger; Protected against vertically falling water drops).		Bag packing bag

### ELECTROMAGNET COMPATIBILITY

This device was designed to satisfy the currently required requisites for electromagnetic compatibility (EN 60601-1-2). Electro-medical devices require particular care during installation and use relative to EMC requirements. Users are therefore requested to install and/or use these devices following the manufacturer's specifications. Risk of potential electromagnetic interference with other devices. Mobile or portable RF radio and telecommunication devices (mobile phones or wireless connections) may interfere with the operation of medical devices. For more information, visit the website [www.flaemnuova.it](http://www.flaemnuova.it). The appliance may be susceptible to electromagnetic interference in the presence of other devices used for specific diagnoses or treatments. Flaem reserves the right to make technical and functional changes to the product without prior notice.

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© Copyright by FLAEM NUOVA S.p.A., 2022
**FLAEM NUOVA S.p.A. non risponde di ulteriori estensioni del periodo di garanzia assicurate da terzi.**