

AirFeel[®] IMT

The new inspiratory muscles trainer.



1

IMT TECHNOLOGY. WHAT IS AND HOW IT IS MADE

The AirFeel® IMT is an inspiratory muscle training device. It is recommended for expand pulmonary capacity, increase thickness of the diaphragm and for improve pelvis base muscle's strength.



AirFeel® IMT | Cod. ACO798P


Dimensions: Ø42 x 110 mm
Weight: 55 g
Materials: PP, SUS, Silicone
Master carton: 24 pz



Operative range
6-60 cm H₂O



Accessible cap for
optimal cleaning



Easy reading
of the set value



Easy pressure
setting

Thanks to IMT technology you can reduce the work of the respiratory muscles, improve respiratory efficiency, obtain a decrease in the perception of effort and less sign of fatigue, improve the attenuated metaboreflex of the respiratory muscles.

AirFeel® IMT is indicated for clinical asthma, heart failure, occupational applications, sports applications, and other activities that require good breathing capacity.



SPORT APPLICATION
(Cycling, running and swimming)



CLINICAL APPLICATION
(Cystic fibrosis, COPD, asthma, heart failure stroke)



OCCUPATIONAL APPLICATION
(Hypoxia, weightlifting)



AirFeel® IMT includes a mouthpiece, a nose clip and the instruction manual.

The device is made of medical grade materials: polipropilene, silicone and stainless steel are suitable for cleaning, chemical and thermal disinfection, as well as being autoclavable. A practical nose clip is supplied to complete the equipment, useful for maximising breathing exercises.

2 HOW AIRFEEL® IMT TECHNOLOGY WORKS

You can use AirFeel® IMT in two positions:

When in resting position 1 the spring pushes the valve cap avoiding the air passage.

When in position 2 patient breaths the pressure applied to the valve (blue arrow) makes the cap to slide down and let the air (green arrow) to pass through. The pressure can be set by rotating the central shaft and compressing/releasing the spring.



AirFeel® IMT



Flaem Nuova Spa
Via Colli Storici 221/223/225
25015 Desenzano del Garda (BS) Italy
Tel. (+39) 030 9910168
info@gruppoflaem.it
www.gruppoflaem.it



www.fluem.eu



COD. 20815AO REV. DATE 10/2024 © Copyright by FLAEM NUOVA S.p.A.
FLAEM NUOVA reserves the right to make technical and functional alterations to this product without prior notice.