
# Rhino Clear Mobile

by **FLAEM**<sup>®</sup>

over 40 years  
MADE IN ITALY

- Doccia Nasale Portatile Ricaricabile  
Manuale Istruzioni D'uso
- Rechargeable Portable Nasal Irrigator  
Instructions for use manual
- Douche nasale Portative Rechargeable  
Mode D'emploi
- Draagbare en oplaadbare Neusdouche  
Gebruiksaanwijzing
- Mobile Nasendusche mit Akku  
Bedienungsanleitung
- Nebulizador Nasal Portátil Recargable  
Manual De Instrucciones
- Переносной аккумуляторный  
ирригатор для носа  
инструкция по эксплуатации
- 可充电便携式鼻腔冲洗器  
使用说明书
- دش أنفي محمول قابل لإعادة الشحن  
دليل تعليمات الاستخدام

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ENGLISH  
page 9

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NEDERLANDS  
pagina 25

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# Rhino Clear Mobile

## Rechargeable Portable Nasal Irrigator

**IMPORTANT: Before using the portable nasal irrigator, allow it to charge for at least 24 hours. Subsequently, we recommend keeping the unit inserted on the recharger base connected to the mains.**

Congratulations on your purchase and thank you for choosing our unit. Our goal is to fully satisfy consumers by offering them cutting-edge systems for preventing and treating ailments of the upper airways. Visit our Internet website [www.flaem.it](http://www.flaem.it) to view the full range of Flaem products.

**Read these instructions carefully before using the nasal douche in order to use it correctly. Furthermore, we advise you to store this manual with care should you need it for further referencing. Rhino Clear Mobile is a medical device for home use.**

### What is it?

**Rhino Clear Mobile** is a nasal irrigator specifically designed for flushing and treating the nasal cavities of both adults and children with saline solutions or medications for aerosol therapy. Once charged, the nasal irrigator can be transported and used without being connected to the power supply and charged using the recharger base provided. Its use must be prescribed by a Doctor who has assessed the patient's general condition.

### How does it work?

**Rhino Clear Mobile** is designed to nebulise an aqueous solutions for flushing the nasal cavities and to collect the solution, as it falls, in a dedicated collection chamber inside the nasal irrigator. It can be used with saline solutions (isotonic and hypertonic solutions, thermal waters) or medicated solutions for aerosol therapy, as prescribed by a doctor. The ergonomic irrigator compressor body constitutes a comfortable grip fitted with a button for convenient activation of nebulisation and is also equipped with an inlet valve that allows the inflow of air and restricts noise levels. The nebuliser body has two adjacent chambers, one for holding the solution to be nebulised and the other for collecting the solution. The nasal adaptors, which come in 3 different sizes to suit the user's nostrils, are fitted onto the nebuliser body. The front part of the recharger base houses a special nasal adaptor storage compartment with protective lid. A second recess on the bottom of the base houses the power supply cord.

### When should it be used?

**Rhino Clear Mobile** can be used for:

1. Physiological flushing of nasal cavities (saline solution or thermal water). This helps maintain optimum nasal hygiene in children. In adults, it aids physiological flushing of the nose and can be used to prepare the respiratory tract for subsequent administration of medication, favouring even distribution.
2. Therapeutic flushing: Use **Rhino Clear Mobile** to administer medication only when prescribed by your doctor or pharmacist, respecting the type of medication, doses and duration of treatment indicated.


### WARNING:

The pull-out nebuliser chamber available as an accessory (code ACO409P), allows full personalisation of the appliance for each member of the family. Personalised use of the nebulisation chamber allows optimum hygiene and prevents contagion. The average life span of the accessories is approximately 1 year.

**Use original FLAEM accessories only.**


## This nasal irrigator set includes:

- (1a) "A" nasal adaptor with soft rubber tip, for children under 3
- (1b) "B" nasal adaptor with soft rubber tip, for children aged 3 to 10
- (1c) "C" nasal adaptor with soft rubber tip, suitable from 10 years
- (2) Separator
- (3) Atomiser
- (4) Pull-out, washable nebuliser body made up of 2 separate chambers:  
(4a) 15 ml nebulisation chamber  
(4b) Collection chamber
- (5) Irrigator compressor body with ergonomic grip  
(5a) Dispensing button  
(5b) Recharge led  
(5c) Inlet valve with lid
- (6) Recharger base with nasal adaptor compartment  
(6a) Adaptor compartment with lid  
(6b) Power supply cord compartment
- (7) Fun child's colour sticker
- (8) Handy carry bag

## 1. Connection and charging

- 1.1 - Always make sure that the mains voltage matches the information on the rating plate of the recharger base.
- 1.2 - Place the nasal irrigator complete with recharger base on a stable surface.
- 1.3 - Connect the power supply cord to the power socket (FIG. A).
- 1.4 - When the irrigator compressor body is placed in the recharger base, the green indicator light (5b) on the button will indicate the recharging status. If, after using the nasal irrigator, it is replaced on the recharger base, the indicator light will come on again. You cannot use the nasal wash while it is charging.
- 1.5 - Charging for first use takes about 24 hours and assures the utmost autonomy (6/7 applications of 10 ml). At the end of charging, disconnect the base from the power mains by unplugging the device from the socket, the position should ensure disconnection from the power mains can be carried out conveniently.
- 1.6 - After long periods of inactivity, recharge the unit by standing it on the recharger base for 24 hours, at least once every three months, to avoid damage to the battery inside. The internal batteries have an average life of 1 year and can be recharged about 200 times.

### IMPORTANT INFORMATION TO KEEP THE BATTERY EFFICIENT

- We recommend keeping the unit on the charger base when it is not in use
- After long periods of inactivity, recharge the unit by standing it on the charger base for 24 hours, at least once every three months, to avoid damage to the battery

## 2. Operating instructions

- Before each use, thoroughly wash your hands and clean the nasal irrigator and its accessories as directed in section 3. Cleaning, Sanitization, Disinfection and Sterilization.

- Take care during use as the solution may drip from the nose.

2.1 - Before using for the first time, charge the batteries as indicated in the previous section.

2.2 - To use the nasal irrigator correctly, the man symbol (♂) (4c), on the nebulisation chamber (4), must always face the patient.

The nebulisation chamber can either be fitted so that this symbol is on the same side as the power button (5a), so the nasal irrigator can be activated using the thumb, or with the symbol on the opposite side to the button, to be activated using the index finger.

To obtain one of these conditions, when starting to screw the nebulisation chamber (4) clockwise onto the compressor body (5) (figure B), the man symbol must be on the opposite side to the position desired (figure C).

2.3 - Pour the solution (max 15 ml) into the nebulisation chamber (4a) as shown in figure D. Check that the atomiser (3) is in place.

2.4 - Apply the separator (2), using the ball (2a) next to the man symbol (♂) (4c) on the outside of the nebulisation chamber (4) to guide it, as indicated in Figure E. Then insert the nasal adaptor that best fits your nostrils (1a/1b/1c).

2.5 - Place the hole on the nasal adaptor inside a nostril (FIG. F).

2.6 - Start the nasal irrigator by holding the dispensing button (5a) pressed down while breathing in. Release the button before breathing out; during this stage, to help clear your nose of mucous, gently blow through the nostril being treated with the nasal irrigator, keeping the other nostril closed.






2.7 - Also during nebulisation in one nostril, press your finger on the other nostril to keep it closed.

2.8 - Repeat steps from 2.5 to 2.7 for your other nostril, trying to divide the solution evenly between the two.

2.9 - When treating children or persons with disabilities, who are not able to coordinate the various manoeuvres of step 2.6, we suggest switching nostrils approximately every 15 seconds, pausing in between to allow the discharge of mucous.

2.10 - Before administering aerosol treatments, it is a good idea to wash your nose with saline solution. This will allow the medication to be distributed more effectively over the mucous membrane.

2.11 - After completing the treatment, switch off the unit. Clean and/or sanitise the unit as directed in section 3. Cleaning, Sanitization, Disinfection and Sterilization. The average life of the device includes about 1,000 applications.

### **3.Cleaning, Sanitization, Disinfection and Sterilization. (For clinical or hospital use only)**

Careful cleaning of the unit is extremely important to ensure long life and proper operation.

#### **3.1 - Cleaning - Irrigator compressor body**

To clean use only a damp cloth and an antibacterial detergent (non-abrasive and solvent free).

**WARNING:** The irrigator compressor body is not protected against the penetration of liquids; never wash it under running water or by soaking.

#### **3.2 - NASAL ADAPTORS, SEPARATOR, ATOMISER and NEBULISER BODY, WITH THE EXCEPTION OF THE COMPRESSOR BODY**

After use, remove the nasal adaptor (1a-1b-1c), the separator(2),and the atomiser (3)and unscrew the nebuliser body(4). Then clean the above parts performing the following operations:

3.3 - **Sanitization:** The above parts can be sanitised choosing one of these methods: a) in warm (approximately 40° C) drinking water with mild dishwashing liquid (non abrasive) or in the dishwasher with hot cycle. b) by soaking in a 50% water to 50% white vinegar solution. Finish by rinsing thoroughly with drinking water.

3.4 - **Drying:** After sanitization shake the accessories vigorously and place them on a paper towel, or alternatively, dry them with a jet of hot air (e.g. using a hair dryer).

Reassemble the nasal irrigator as shown in the drawing on page 10.

3.5 - **Disinfection:** The accessories that can be disinfected are (1a-1b-1c-2-3-4).

**(Method A):** The disinfection procedure described in this paragraph is to be performed before using the accessories and is effective only if all the prescriptions herein are met and all the components involved are previously sanitized.

The disinfectant to be used must be of electrolytic oxidizing chlorine type, (active ingredient: sodium hypochlorite), specific for disinfection and is available in all pharmacies.

##### **Procedure:**

- Fill a container, suitable to contain all the individual components to be disinfected, with a solution of drinking water and disinfectant, while respecting the proportions indicated on the packaging of the disinfectant itself.
- Completely immerse each individual component in the solution, taking care to avoid the formation of air bubbles in contact with the components. Leave the components immersed for the period of time indicated on the packaging of the disinfectant, and associated with the concentration chosen to prepare the solution.
- Remove the components now disinfected and rinse thoroughly with lukewarm drinking water.
- Dispose of the solution as recommended by the disinfectant manufacturer.

**(Method B):** By boiling in water for 10 minutes; preferably use demineralised or distilled water to avoid scale deposition.

After disinfection shake the accessories vigorously and place them on a paper towel, or alternatively, dry them with a jet of hot air (e.g. using a hair dryer).

#### **3.6 - Sterilization:**

Sterilization intended for clinical or hospital use is the responsibility of the same structure that it performs and must be performed according to the instructions in this manual.

The accessories that can be sterilized are (2-3-4).


The sterilization procedure described in this paragraph is effective only if all the prescriptions herein are met, all the components involved are previously sanitized and the procedure itself is validated according to ISO 17665-1 rule. **Unit:** Steam sterilizer with fractionated vacuum and overpressure according to EN 13060 rule.

**Procedure:** Package each component to be treated in sterile barrier system or packaging complying with EN 11607 rule. Place the packed components in the steam sterilizer.

Run the sterilization cycle according to the instructions manual of the device, selecting a temperature of 134° C and a duration of 10 minutes.

**Storage:** Store sterilized components as per instructions of the sterile barrier package or system chosen.

- 3.7 - **INLET VALVE.** The nasal irrigator features a inlet valve (5c) which should be cleaned when dirty; it is housed under the lid located on the bottom of the irrigator compressor body (FIG. G). To clean the valve first remove the lid (by levering with a biro and holding down the other side of the lid) and then the valve; finally wash the valve in water and detergent. Before reassembling, make sure that the valve is completely dry.

#### **ADDITIONAL CLEANING INFORMATION of the atomization chamber (with the exception of the compressor body)**

In order to avoid deposits that may clog the atomizer and the atomization chamber and consequently affect the operation of the nasal irrigator, it is recommended that cleaning and sanitizing is performed after each use following the instructions for use contained in this instruction manual included with the nasal irrigator.

After rinsing accurately the parts under warm tap water (40°C approx.), shake the wet parts vigorously in order to eliminate any remaining droplets of water. Then, place the parts on a paper towel or dry them using a jet of warm air.

**WARNING:** Should your nasal irrigator not atomize properly after cleaning and sanitizing is performed, do not use sharp tools (i.e. screwdrivers, scissors ends, needles, hairpins, etc.) to try to clean or unclog the atomizer or atomization chamber holes as this could damage these parts and impair the correct operation of the unit. Instead, perform additional cleaning and sanitizing cycles as described in the instruction manual and above recommended, until the nasal irrigator starts to operate properly again.

## **4. Important cautionary notes**

Like all electrical appliances, especially where children are concerned, Rhino Clear Mobile requires a number of basic precautions.

- 4.0 - This M.D. is also intended for direct use by the patient.
- 4.1 - Before using the device for the first time and periodically during the product's life, check the integrity of its structure and of the power cable, making sure that there is no damage; if it is damaged, do not insert the plug and take the product immediately to an authorised technical assistance centre or to your local dealer.
- 4.2 - Do not allow the nasal irrigator to be used as a toy. Careful attention is necessary when used by or near children.
- 4.3 - Children and disabled persons should always use the nasal irrigator under strict supervision of an adult who has read this manual.
- 4.4 - Some components of the device are small enough to be swallowed by children and can cause suffocation; store the device out of the reach of children.
- 4.5 - Never use adaptors for power supplies different from the one indicated on the recharger base of the nasal irrigator (plate details). Keep the cord away from heated surfaces
- 4.6 - Keep the power supply cord away from animals (for example rodents) as they may damage the insulation of the power supply cord.
- 4.7 - The unit should not be used in the presence of inflammable anaesthetic mixtures containing oxygen or nitrous oxide.
- 4.8 - Do not handle the recharger base plug with wet hands. Never submerge the recharger base in water; should this occur accidentally, unplug it immediately. Do not touch the unit under water before pulling out the plug. Take the unit immediately to an authorised service centre or your dealer.
- 4.9 - Never wash the irrigator compressor body under running water or by soaking. Do not use the nasal irrigator when taking a bath or shower.
- 4.10 - Repairs may only be performed by authorised personnel, following the information in the restricted area of

the website [www.flaem.it](http://www.flaem.it). Unauthorised repairs void the warranty and may pose a danger for the user.

- 4.11 - Warning: if batteries need replacing take the unit to an authorised service centre.
- 4.12 - The Manufacturer, the Retailer and the Importer are responsible for safety, reliability and performance only if: a) the unit is used in accordance with the operating instructions; b) the electrical system of the premises where the unit is used is in conformity with current laws.
- 4.13 - The manufacturer must be contacted about any problems and/or unexpected events concerning operation and for any clarifications on the use, maintenance or cleaning.
- 4.14 - Microbial contamination: nasal irrigator and the accessories are only for personal use to avoid any risk of infection.
- 4.15 - Interactions: the materials used have undergone extensive biocompatibility tests (ISO 10993-5 and ISO 10993-10) in compliance with the essential requirements set forth in the European Medical Device Directive 93/42/EEC. Phthalates free. In conformity with Reg. (EC) no. 1907/2006.
- 4.16 - Carry out the cleaning procedures after every application for maximum hygiene, to obtain optimal performance and extend the life of your unit.
- 4.17 - Do not obstruct the inlet valve on the bottom of the irrigator compressor body; never insert any object into its lid.
- 4.18 - WARNING: Do not modify this device without authorisation from the manufacturer.**
- 4.19 - Do not use the power cord outside the intended use, as it may cause dangers of strangulation, pay particular attention to children and people with particular difficulties, often these people are not able to properly assess the dangers.
- 4.20 - The time it takes to move from the storage condition to the operating condition is about 4 hours.

 **Disposal:** In conformity with Directive 2012/19/EC, the symbol shown on the appliance to be disposed of indicates that this is considered as waste and must therefore undergo "sorted waste collection". The user must therefore take (or have taken) the above waste to a pre-sorted waste collection centre set up by the local authorities, or else give it back to the dealer when purchasing a new appliance of the same type. Pre-sorted waste collection and the subsequent treatment, recovery and disposal operations favour the production of appliances made of recycled materials and limit the negative effects of any incorrect waste management on the environment and public health. The unlawful disposal of the product by the user could result in administrative fines as provided by the laws transposing Directive 2012/19/EC of the European member state in which the product is disposed of.

## 5. Troubleshooting

PROBLEM	REMEDY
<p><b>5.1</b> - If the green indicator light (5b) does not come on when the irrigator compressor body is placed on the re-charger base</p>	<p>Check that the power cord is properly plugged into the power socket. Check the electric socket by plugging another appliance into it. Check for faults in the power cord or the plug. In the event of a fault, always contact your retailer or an authorised service centre.</p>
<p><b>5.2</b> - Nebulisation is poor or too short</p>	<p>The batteries could be almost flat, recharge the nasal irrigator for at least 24 hours. Check that the atomiser (3) is set properly in place. Check that the nasal irrigator is clean and, if not, carry out the cleaning operations; carefully follow the <b>Cleaning, Sanitization, Disinfection and Sterilization</b> instructions contained in the manual.</p>
<p><b>5.3</b> - The nasal irrigator does not nebulise</p>	<p>The batteries could be almost flat, recharge the nasal irrigator for at least 24 hours. Check that the solution has been properly poured into the nebulisation chamber (4a). Make sure that the quantity of liquid does not exceed 15 ml. When using medication per aerosol therapy, make sure it is diluted correctly with saline solution 0.9% (as prescribed by your doctor). Check that the nasal irrigator is clean and, if not, carry out the cleaning operations; carefully follow the <b>Cleaning, Sanitization, Disinfection and Sterilization</b> instructions contained in the manual.</p>
<p><b>5.4</b> - The nasal irrigator stops nebulisation (turns off) during operation.</p>	<p>This means that the thermal protection device inside the nasal irrigator has tripped, probably due to an overload. Allow the irrigator to cool down for a few minutes then push the button to turn the irrigator back on. Allow to cool for at least 10 minutes. Check that the nasal irrigator is clean and, if not, carry out the cleaning operations; carefully follow the <b>Cleaning, Sanitization, Disinfection and Sterilization</b> instructions contained in the manual.</p>
<p><b>NOTE: Should the unit again fail to operate after following the above instructions, please contact your local dealer or an authorised service centre.</b></p>	

## 6. Technical specifications

Model:	<b>P0907EM</b>
Recharger base voltage:	220-230V~ 50-60Hz / 18VA
Nebulisation chamber max. capacity:	15 ml
Nebulisation time for 10 ml <sup>(2)</sup> :	approx. 2 minutes (1 ml in approx. 12 seconds)
Cumulative mass percentage oversize > 10 micron:	between 95 and 96 % <sup>(1)</sup>

(Fraction not inhaled useful for the treatment of the upper airways). <sup>(1)</sup> Values measured with a computerised laser diffractometer Malvern Mastersizer, using the procedure devised by the Pharmaceutical Department of Parma University.  
<sup>(2)</sup> Average values measured with 0.9% NaCl saline solution in compliance with Flaem internal procedure.

Nasal irrigator dimensions:	7x6x15.5h cm
Nasal irrigator weight:	210 g
In conformity with EEC directive 93/42:	 0051

### Approvals:

Recharger base dimensions:	15x9x8h cm
Total weight (recharger base with unit and accessories):	480 g
Nasal irrigator internal batteries:	Ni-MH, rechargeable, 1.2V 620/650mAh
Operating conditions:	Temperature: min +5°C; max 40°C Air humidity: min 15%; max 90% <b>RH</b>
Storage conditions:	Temperature: min -10°C; max 45°C Air humidity: min 10%; max 75% <b>RH</b>
Operating/storage atmospheric pressure:	min. 690 hPa; max. 1060 hPa <b>P</b>

The type BF applied parts are the 3 nasal adaptors (1a, 1b, 1c)

## 7. Symbols

Class II unit

Medical CE marking in compliance with Directive 93/42 and subsequent amendments. In accordance with EN 60601-1-11

Alternating current

Follow operating instructions

Type BF applied part

Operating instructions

**IP45**

IP45 Protected against solid foreign bodies with  $\varnothing$  1 mm or larger and against jets of water

TÜV approval

### Electromagnetic compatibility:

This device was designed to meet the current requirements for electromagnetic compatibility (EN 60 601-1-2:2015). For additional information please go to [www.flaem.it](http://www.flaem.it). The Medical Device may be subject to electromagnetic interference if other devices are used for specific diagnosis or treatments.

Flaem Nuova reserves the right to make technical and functional modifications to the product with no prior notice.

# Rhino Clear Mobile

## Douche nasale Portative Rechargeable

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FLAEM NUOVA S.p.A.