Rhino Clear Sprint

by FLAEM®

SINCE 1966
MADE IN ITALY

**PHTHALATES
& BPA FREE**

- Istruzioni d'Uso
- Instructions for Use
- Mode d'Emploi
- Gebruiksaanwijzingen
- Bedienungsanleitung
- Instrucciones de Uso
- Руководство по использованию
- Instrukcja Obsługi
- Instruções de Utilização
- Kullanım talimatları
- 使用指南
- تعليمات الاستخدام

**PATENTED AND
PATENTED PENDING**

ITA
2

ENG
7

FRA
12

NL
17

DE
22

ES
27

PYC
32


POL
37









POR
42

TÜR
47

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Rhino Clear Sprint

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Rhino Clear Sprint

Nasal douche for the whole family

Thank you for your purchase. Our aim is the complete satisfaction of our customers by offering cutting-edge products for the prevention and therapy of the respiratory tract. Visit the website **www.flaem.it**

ENG

Read these instructions carefully before using the nasal douche in order to use it correctly. In addition, we advise you to keep this manual with care should you need it for further consultation. Rhino Clear Sprint is a medical device for home use.

What is it?

Rhino Clear Sprint is a nasal douche specifically designed to irrigate and treat adult and child nasal cavities, either with saline solution or with medicated solutions for aerosol use. The nasal douche is portable and does not need to be plugged in to be used. Its use must be prescribed by a Doctor who has assessed the patient's general conditions.

How does it work?

Rhino Clear Sprint is intended for the washing of nasal cavities, collecting the outflowing solution in a special retrieval chamber present in the nasal douche. It can be used with saline solutions, (isotonic, hypertonic and thermal waters) aerosol medicines, as prescribed by a doctor. The ergonomic compressor body of the douche makes a comfortable handle provided with a button which can easily activate nebulisation and it is also provided with an aspiration valve which allows air to enter and decreases noise. The flask has been designed to accommodate two adjacent chambers, one for the nebulising solution and the other for the outflowing solution. Nasal adaptors, available in 3 different sizes according to the size of the nasal opening of the subject, are to be applied to the flask.

When should it be used?

Rhino Clear Sprint can be used for:

1. Physiological washing of the nasal cavities (physiological solution or thermal water). Useful for children to help them with nasal hygiene. In adults it favours physiological washing of the nose and it can be used before an aerosol therapy of the respiratory tract in order to help the distribution of the medicine.
2. Therapeutic washing: the administration of medicines through **Rhino Clear Sprint** has to be prescribed by a doctor or by a pharmacist who will decide the type of medicine, the administration dose and the treatment duration.

ATTENTION:

The flask, available as an accessory (cod. ACO409P), allows for a complete personalization of the device for each member of the family. For correct hygiene and in order to avoid the spread of infection, a personal use of the device is advisable. The average life span of the accessories is approximately 1 year. Use original FLAEM accessories only.

The standard kit of the device includes (fig A):

- (1) Nasal adaptor 'A' with a soft rubber tip, suitable for children up to 3 years of age
Nasal adaptor 'B' with a soft rubber tip, suitable for children from 3 to 10 years of age
Nasal adaptor 'C', with soft rubber tip, suitable for 10 years of age and over
- (2) Separator
- (3) Atomizer
- (4) Removable and washable flask consisting of 2 separate chambers:
 - (4a) 15 ml nebulisation chamber
 - (4b) Retrieval chamber
- (5) Compressor body with ergonomic handle
 - (5a) Activation button
 - (5b) Battery compartment
 - (5c) Aspiration valve with lid

1. Inserting batteries

1. Open the battery compartment (5b) by pressing and lifting the lid in the direction of the arrows (fig B).

2. Place 2 batteries in the compartment, paying attention to the polarities written inside the compartment. Use two alkaline batteries (AA).
3. Place the lid back on the compartment by lowering the upper part of the compartment until it is completely fastened.

Remove the batteries when the device is not going to be used for a prolonged period of time. Leaving the batteries in the compartment for too long can cause leaking which can damage the device.

ATTENTION:

- Does not contain parts which can be repaired by the user. The warranty does not cover batteries or damage caused by old batteries.
- Only use good quality batteries. Always use new batteries and do not use old and new batteries together. Use batteries of the same brand and of the same type.

Battery disposal: Exhausted batteries are to be disposed of in the appropriate recipient.

2. Instructions for use

- **Before each usage, wash hands thoroughly and clean the nasal douche and its accessories as indicated in paragraph 3. Cleaning, Sanitization, Disinfection and Sterilization.**
 - **During application it is advisable to be careful as the solution can drip from the nose.**
- 2.1 - Before usage place batteries as indicated in paragraph 1.
 - 2.2 - For a correct use of the nasal douche, the symbol of the little man (♂) (4c), placed on the flask (4), always has to be directed towards the patient. The flask can be fixed on the position of the little man on the side of the activation button (5a) allowing activation of the nasal douche with the thumb, or with the little man on the opposite side of the button in order to allow activation with the index finger.
In order to obtain one of the aforementioned conditions, when screwing the flask clockwise (4) on the douche compressor body (5) (fig C.), the symbol of the little man needs to be located on the side opposite to the position you would like to obtain (fig. C).
 - 2.3 - Pour the solution (max 15 ml) in the nebulisation slot (4a) of the flask, as shown in (fig. E). Verify the presence of the atomizer (3).
 - 2.4 - Apply the separator (2) directing it with the dot (2a) in correspondence with the little man (♂) (4c) located on the external part of the flask (4), as shown in (fig F). Then insert the most suitable nasal adaptor for the size of your nasal openings (1A-1B-1C)
 - 2.5 - Insert the hole of the nasal adaptor in one nasal opening (fig. G).
 - 2.6 - Activate the nasal douche by keeping the activation button pressed (5a) during the inhalation phase and release it before exhaling; in this phase, in order to help the discharge of the mucus, gently blow through the nasal opening that has been treated with the nasal douche, keeping the other nasal opening closed.
 - 2.7 - During nebulisation it is advisable to close the other nostril where the douche is not being applied, with a finger.
 - 2.8 - Repeat procedures **2.5** to **2.7** for the other nasal opening as well, trying to evenly distribute the solution in the two nasal openings.
 - 2.9 - In the case of children or people who are not self-sufficient, or not able to perform the various phases of point **2.6**, it is advisable to carry out the application alternating the two nasal openings after about 15 seconds, with breaks necessary to allow for the discharge of mucus.
 - 2.10 - In case of medicine administration through aerosol, it is advisable to carry out a wash with a saline solution first, in order to obtain a more efficient distribution of the medicine on the nasal mucous membrane.
 - 2.11 - At the end of the application, turn off the device and carry out cleaning and/or sanitizing following the instructions in paragraph 3. Cleaning, Sanitization, Disinfection and Sterilization.

3. Cleaning, Sanitization, Disinfection and Sterilization. (For clinical or hospital use only)

A thorough cleaning is extremely important for the lifespan and well-functioning of the device. After usage, unassemble the nasal douche by following the drawing (see fig A), remove the batteries and proceed with cleaning and/or sanitizing of the components performing the following operations:

3.1 - **Cleaning: Irrigator compressor body**

For component (5) only use a damp cloth with antibacterial detergent (non-abrasive and solvents free).

The nasal douche has an aspiration valve (5c) which is to be cleaned when dirty. To access the valve lift the lid of the battery compartment and remove the valve by lifting it with a ballpoint pen (fig H); then wash the valve with water and detergent. Before reassembling, make sure that the valve is completely dry.

ATTENTION:

The compressor body of the douche is not protected against the penetration of liquids; therefore do not wash it with running water or immersion. The compressor body of the douche is not protected against the penetration of liquids; therefore do not wash it with running water or immersion.

3.2 - **Sanification: Clean the components (1A-1B-1C-2,3,4) choosing one of the following methods.**

Method A: Wash the components with hot drinking water (about 40°C) with a delicate dish detergent (non-abrasive).

Method B: through immersion in a solution composed of 50% water and 50 % white vinegar. Then rinse well with hot drinking water (about 40°C);

Drying: After washing or sanitizing the components, shake them well and place them on a paper napkin or dry them with hot air (for example with a blow dryer). Reassemble the nasal douche as shown in the picture (see fig A).

3.3 - **Disinfection: The accessories that can be disinfected are (1A-1B-1C-2,3,4)**

Method A: The disinfection procedure described in this paragraph is to be performed before using the accessories and is effective only if all the prescriptions herein are met and all the components involved are previously sanitized.

The disinfectant to be used must be of electrolytic oxidizing chlorine type, (active ingredient: sodium hypochlorite), specific for disinfection and is available in all pharmacies.

Procedure:

- Fill a container, suitable to contain all the individual components to be disinfected, with a solution of drinking water and disinfectant, while respecting the proportions indicated on the packaging of the disinfectant itself.
- Completely immerse each individual component in the solution, taking care to avoid the formation of air bubbles in contact with the components. Leave the components immersed for the period of time indicated on the packaging of the disinfectant, and associated with the concentration chosen to prepare the solution.
- Remove the components now disinfected and rinse thoroughly with lukewarm drinking water.
- Dispose of the solution as recommended by the disinfectant manufacturer.

Method B: boil in water for 10 minutes; it is advisable to use demineralized water or distilled water in order to avoid chalky deposits.

After disinfection shake the accessories vigorously and place them on a paper towel, or alternatively, dry them with a jet of hot air (e.g. using a hair dryer).

3.4 - **Sterilization:**

Sterilization intended for clinical or hospital use is the responsibility of the same structure that it performs and must be performed according to the instructions in this manual.

The accessories that can be sterilized are (2-3-4)

The sterilization procedure described in this paragraph is effective only if all the prescriptions herein are met, all the components involved are previously sanitized and the procedure itself is validated according to ISO 17665-1 rule.

Unit: Steam sterilizer with fractionated vacuum and overpressure according to EN 13060 rule.

Procedure: Package each component to be treated in sterile barrier system or packaging complying with EN 11607 rule. Place the packed components in the steam sterilizer. Run the sterilization cycle according to the instructions manual of the device, selecting a temperature of 134° C and a duration of 10 minutes.

Storage: Store sterilized components as per instructions of the sterile barrier package or system chosen.

ADDITIONAL CLEANING INFORMATION of the atomization chamber (with the exception of the compressor body).

In order to avoid deposits that may clog the atomizer and the atomization chamber and consequently affect the operation of the nasal irrigator, it is recommended that cleaning and sanitizing is performed after each use following the instructions for use contained in this instruction manual included with the nasal irrigator. After rinsing accurately the parts under warm tap water (40°C approx.), shake the wet

parts vigorously in order to eliminate any remaining droplets of water. Then, place the parts on a paper towel or dry them using a jet of warm air. **WARNING:** Should your nasal irrigator not atomize properly after cleaning and sanitizing is performed, do not use sharp tools (i.e. screwdrivers, scissors ends, needles, hairpins, etc.) to try to clean or unplug the atomizer or atomization chamber holes as this could damage these parts and impair the correct operation of the unit. Instead, perform additional cleaning and sanitizing cycles as described in the instruction manual and above recommended, until the nasal irrigator starts to operate properly again.

4. Important warning

Just as any electrical device, especially in the presence of children, Rhino Clear Sprint needs to be used carefully.

- The nasal douche is not to be used as a toy. The utmost care and attention is required when the douche is used by children or when children are present.
- In the presence of children or of people who are not self-sufficient, use the nasal douche under the careful supervision of an adult who has read the enclosed manual.
- Do not use batteries with a different alimentation other than those listed on the plate.
- Remove the batteries before storing Rhino Clear Sprint.
- It is a device which cannot be used in the presence of flammable anaesthetics with oxygen or with nitrogen protoxide.
- The casing of the device is not protected against the penetration of liquids. Do not wash the device under running water or by immersing it and keep it away from splashes of water or of other liquids.
- Do not use the nasal douche when showering or bathing.
- Do not expose the device or the batteries to extreme temperatures. Do not place the device and the battery near sources of heat, sunlight or in extremely hot places.
- Repairs may only be performed by authorised personnel following the information in the restricted area of the website www.flaem.it. Unauthorised repairs void the warranty and may pose a danger for the user.
- **ATTENTION:** Do not modify this device without authorisation from the manufacturer.
- The Manufacturer, Importer and Seller shall be liable for the safety, reliability and performance only if the device is used in accordance with the operating instructions.
- The manufacturer must be contacted about any problems and/or unexpected events concerning operation and for any clarifications on use, maintenance or cleaning.
- Microbial contamination: you must foresee a personal use of the nasal shower and accessories to avoid any risk of infection.
- Interactions: the materials used have been tested according to the bio-compatibility standards (ISO 10993) in compliance with the essential requirements of the Medical Device Directive 93/42 EEC.
- Clean the device after each use, to assure optimal hygiene and to improve the service life and operation of the device.
- Do not obstruct the suction valve on the bottom of the douche compressor body; never insert anything into the cover.
- This Medical Device is intended for direct use by the patient.
- The time it takes to switch from operating condition is about 4 hours.

SYMBOLS

 Batteries

 0051

Medical CE marking in compliance with Directive 93/42 and subsequent amendments. In accordance with EN 60601-1-11

BF type applied part

Follow usage instructions

IP45

Protected against solid foreign bodies with \varnothing 1 mm or larger and against jets of water

TÜV approval

5. Troubleshooting

| PROBLEMS | SOLUTIONS |
|---|--|
| 5.1 – If the nasal douche nebulises poorly or for a short time. | Batteries might be partly exhausted. Remove and replace batteries. Check that the atomizer (3) has been positioned correctly. Check that the nasal douche is clean and, if not, clean it again by carefully following the instructions reported in the manual regarding Cleaning, Sanitization, Disinfection and Sterilization. |
| 5.2 – If the nasal douche does not nebulise | Batteries might be partly exhausted or positioned wrongly. Remove and replace the batteries. Check that the solution has been poured in correctly in the nebulisation slot (4a) of the flask (4). Check that the quantity of liquid is not above 15 ml. If using an aerosol medicine, check that it has been sufficiently diluted with physiological solution 0,9% (strictly follow the doctor's prescription). Check that the nasal douche is clean and, if not, clean it again by carefully following the instructions reported in the manual regarding Cleaning, Sanitization, Disinfection and Sterilization. |
| NOTE: If after verifying the abovementioned conditions the device still does not nebulise, we advise you to contact your seller or the authorized assistance centre. | |

6. Technical characteristics

| | |
|--|---|
| Model: | P0709EM |
| Input voltage | 3,0 Vdc – 4,5 W |
| Max capacity of the nebulisation slot | 15 ml |
| Nebulisation time for 10 ml of physiological solution NaCl 0,9% about 2 minutes (1 ml about 12 seconds) | |
| Cumulative Percentage of Mass Oversize > 10 microns between 95 and 96 % ⁽¹⁾ (Non-respirable fraction useful for the treatment of the upper respiratory tract) - ⁽¹⁾ Values measured with a computerized laser diffractometer Malvern Mastersizer according to the procedure of the Pharmaceutical Department of the University of Padova. | |
| Nasal douche dimensions | 7x6x15,5h cm |
| Nasal douche weight (no batteries) | 185 g |
| Complying with directive 93/42 ECC |  |
| Approvals: | |
| Nasal douche batteries | 2 1,5 V alkaline batteries (AA type) |
| Conditions of use: | Temperature: min +5°C; max 40°C Air humidity: min 15%; max 90% RH |
| Storing conditions: | Temperature: min -10°C; max +45°C Air humidity: min 10%; max 75% RH |
| Atmospheric pressure in operating/storage conditions: | min. 690 hPa; max. 1060 hPa P |
| IP45 | |
| Applied BF-type parts are: | (1A-1B-1C) |

Disposal: For an adequate disposal, remember that the device contains batteries. In compliance with Directive 2012/19/EC, the symbol on the device indicates that the device to dispose of is considered as waste, and should therefore be disposed of accordingly. The user shall take the abovementioned waste (or have it taken) to a proper separate collection of waste centre created by local administrations, or hand it in to the seller when buying a new device of the same kind. The separate collection of waste and the following operations of treating, recycling and disposal, favour the production of devices with recycled materials and limit the negative effects on the environment and on health, caused by inappropriate waste handling. The unlawful disposal of the product by the user could result in administrative fines as provided by the laws transposing Directive 2012/19/EC of the European member state or of the country in which the product is disposed of.

Electromagnetic compatibility

This device was designed to meet the current requirements for electromagnetic compatibility (EN 60 601-1-2:2015). For additional information please go to www.flaem.it. The Medical device may be subject to electromagnetic interference if other devices are used for specific diagnosis or treatments.

Flaem Nuova reserves the right to make technical and functional modifications to the product with no prior notice.

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© Copyright by FLAEM NUOVA S.p.A.,
FLAEM NUOVA S.p.A.