SINCE 1966
MADE IN ITALY

FLAEM[®]

a healthy lifestyle

- I** MANUALE ISTRUZIONI D'USO
- GB** INSTRUCTIONS FOR USE MANUAL
- F** MODE D'EMPLOI
- NL** GEBRUIKSAANWIJZING
- D** BEDIENUNGSANLEITUNG
- E** MANUAL DE INSTRUCCIONES
- RUS** РУКОВОДСТВО ПО ЭКСПЛУАТАЦИИ

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Lella the Ladybug

COMPLETE SYSTEM FOR AEROSOL THERAPY

Mod. P0406EM F700

We are pleased you have purchased our product and we thank you for your trust in us. We aim at fully satisfying our customers by offering them state-of-the-art products for the treatment of respiratory diseases.

Read these instructions carefully and retain them for future reference. Only use the accessory as described in this manual. This is a home medical device to nebulise and administer medication prescribed or recommended by your doctor upon assessing the patient's general conditions. Please note that the full range of Flaem products is visible on the website www.flaem.it.

THIS DEVICE COMES EQUIPPED WITH:

A - Aerosol Therapy Device

- A1 - Switch
- A2 - Air outlet port
- A3 - Air filter
- A4 - Nebulizer port
- A5 - Socket for power cord
- A6 - Power cord

B - Connection tube (main unit / nebulizer)

C - Accessories

- C1 - RF6 Basic² nebuliser
 - C1.1 - Lower part
 - C1.2 - Complete nozzle
 - C1.3 - Upper part
 - C1.4 - Cap
- C2 - Mouthpiece
- C3 - Adult mask (PVC)
- C4 - Child mask (PVC)

C5 - RHINO CLEAR® Nasal shower

- C5.1 - Nasal transparent adapter "A", suitable for children up to 3 years of age
 - Nasal white adapter "B", suitable for children 3 to 10 years old
 - Nasal blue adapter "C", suitable for children 10 years old and up
- C5.2 - Separator
- C5.3 - Atomizer
- C5.4 - Shower body with ergonomic handle
 - C5.4.1 - 10 ml nebulization chamber
 - C5.4.2 - Collection chamber
- C5.5 - Dispensing button
 - C5.5.1 - Button
 - C5.5.2 - Rubber cap
- C6 - Nasal shower-tube connect fitting

E - Carrying bag

F - "Lella the Ladybug" fairy tale

IMPORTANT SAFEGUARDS

- This device is also intended for direct use by the patient.
- Before using the product for the first time, and periodically during its lifetime, check the integrity of the device structure and of the power cable to make sure there is no damage. In the event of damage, do not plug in the cable and immediately take the product to an authorised FLAEM service centre or to your trusted dealer.
- Should your device fail to provide the expected performance, contact the authorised service centre for clarifications.
- The expected medical life of the accessories is 1 year. It is, however, advisable to replace the nebuliser cup every 6 months in the event of intense use (or earlier if the cup is obstructed) to always guarantee maximum therapeutic efficacy.
- Children and people who are not self-sufficient must use the device under the close supervision of an adult who has read this manual.
- Some parts of the device are small enough to be swallowed by children; therefore, keep the device out of the reach of children.
- Do not use the supplied tubing and cables for anything other than their intended use. These parts could cause a strangling hazard: pay close attention to children and persons with particular difficulties as they are often unable to accurately evaluate danger.
- The apparatus is unsuitable for use in presence of flammable anaesthetic mixture with air, oxygen or nitrous oxide.
- Always keep the power supply cable away from hot surfaces.
- Keep the power cable away from animals (for example, rodents) which could damage the insulation.
- Do not handle the device with wet hands. Do not use the device in damp environments (for example, while taking a bath or shower). Do not immerse the device in water; in the event of immersion immediately disconnect the plug. Do not remove or touch the immersed device; unplug the power cable first. Immediately bring the device to an authorised FLAEM service centre or to your trusted dealer.

- Use the device only in dust-free conditions, otherwise treatment could be compromised.
- The casing of this unit is not protected against the penetration of liquids. Do not wash the device under running water or submerge in water. Keep away from splashing water or other liquids.
- Do not expose the device to particularly extreme temperatures.
- Do not place the device near sources of heat, in direct sunlight or in excessively hot rooms.
- Do not obstruct or put objects into the filter or its related housing in the device.
- Never obstruct the air vents located on both sides of the device.
- Always use it on a rigid surface that is clear of obstacles.
- Make sure there is no material obstructing the air vents before each use.
- Do not put any objects in the air vents.
- Repairs, including the replacement of the supply cord, are to be carried out by FLAEM authorised personnel only, by complying with the information provided by the manufacturer. Unauthorized repairs will void the warranty and may represent a danger to the user.
- The average expected duration for the compressor series are: F400: 400 hours, F700: 700 hours, F1000: 1000 hours, F2000: 2000 hours.
- **WARNING:** Do not modify this device without authorisation from the manufacturer.
- The Manufacturer, the Vendor and the Importer shall be held responsible for safety, reliability and performance only if: a) the device is used in compliance with the instructions for use b) the wiring where the device is being used is in compliance with safety regulations and current laws.
- Interactions: the materials used in contact with medication have been tested with a vast range of medications. However, in view of the variety and continuous evolution of pharmaceuticals, interactions cannot be ruled out. We recommend using the medication as soon as possible once it has been opened and preventing prolonged exposure in the nebuliser cup.
- The manufacturer must be contacted about any problems and/or unexpected events concerning operation and for any clarifications on use, maintenance/cleaning.
- Interactions: The materials used in the medical device are biocompatible in accordance with the provisions of Directive 93/42 EC and subsequent amendments. However, the possibility of occurrence of allergic reactions cannot be entirely excluded.
- The amount of time required between preservation conditions and use is approximately 2 hours.

INSTRUCTIONS FOR USE

RF6 Basic² NEBULISER


Before each use, clean hands thoroughly and clean the device as described in the section on “CLEANING SANITISATION DISINFECTION”. During use, it is advisable to protect yourself from any dripping. It is recommended that each person use their own nebulizer cup and accessories to prevent risk of infection due to contamination.

This device is suitable for the administration of medical substances and not, for which the administration via aerosol is foreseen; these substances are to be in any case prescribed by the Doctor. In case of too thick substances, the dilution with a suitable physiological solution could be needed, according to the medical prescription.


1. Plug the power cord (A6) into the socket (A5) and plug the power cord into a power socket corresponding to the voltage of the unit. This must be positioned so that it is not difficult to make the disconnection from the mains.
2. Assemble the nozzle as indicated in the “Connection diagram” in point C1.2 Insert the complete nozzle in the bottom part (C1.1). Insert the Upper ventilation cover (C1.4) in the upper part (C1.3) as shown in the “Connection diagram” in point C1.
3. Pour the medication prescribed by the physician into the lower part (C1.1). Close the nebuliser by turning the upper part (C1.3) clockwise.
4. Connect accessories as indicated in the “Connection diagram” on the cover.
5. Sit comfortably holding the nebuliser in your hand, place the mouthpiece onto your mouth or alternatively use the nose piece or mask. Should you opt for the mask accessory, place it on your face as shown in the picture (with or without using the elasticated strap).

The cap (C1.4) provided with the nebuliser RF6 Basic² makes it possible to deliver the medication in 2 modes:

“standard mode”
for optimum delivery and minimum waste of medication

“high speed mode” for a quicker treatment

6. Turn on the device using the switch (A1) and then breathe in and exhale deeply.

7. When the treatment is finished, shut off the device and unplug the power cord.

WARNING: If after the therapy session moisture accumulates in the tube (B), remove the tube from the nebulizer and dry it using the device's compressor fan; this operation will prevent mold from forming inside the tube.

NASAL WASH RHINO CLEAR®

Before each use, wash your hands thoroughly and sanitize RHINO CLEAR® as described in section "CLEANING SANITIZATION DISINFECTION". The NASAL WASH SYSTEM and accessories are for personal use only to avoid potential risks of transmitting infection. During the application it is recommended to be adequately protected against dripping.

- 1) Plug the power cord (A6) into the socket (A5) and plug the power cord into a power socket corresponding to the voltage of the unit. This must be positioned so that it is not difficult to make the disconnection from the mains.
- 2) Insert the atomizer (C5.3) into the nebulization chamber (C5.4.1). Place the solution (max 10 ml) in the nebulizing chamber (C5.4.1).
- 3) Apply the separator (C5.2) verifying that the arrow on the shower body is aligned with the separator (Pic. B). Then, insert the correct adapter for your age, or the most appropriate to fit the size of your nostrils (C5.1A/C5.1B/C5.1C).
- 4) Connect the RHINO CLEAR® device through the tube (B) and the fitting (C6).
- 5) Turn on the device using the switch (A1).
- 6) Place the device close to the nose and gently insert the adapter nosepiece into one nostril.
- 7) During treatment it is preferable to hold the other nostril closed, where the shower is not applied, by applying light finger pressure.
- 8) Run the shower holding down the manual dispensing button (C5.5) during the inhalation phase and release it before exhaling; at this stage, to help expel the mucus, blow gently through the nostril treated directly into the RHINO CLEAR®.
- 9) Repeat from Step 6 to Step 8 for the other nostril, trying to evenly distribute the solution between the two nostrils.
- 10) When used on children or other dependent people, or in any case when there are problems in coordinating the various phases of paragraph 8), we recommend to run the application alternating the nostrils about 20 seconds each, with pauses required to expel the mucus.
- 11) The operations described in the preceding paragraphs is useful for nasal cleaning using physiological solutions. The treatment with drugs is advisable after a first cleaning with physiological solution in order to achieve a more effective distribution of the drug on the nasal mucosa.
- 12) When the treatment is finished, shut off the device and unplug the power cord. Clean it following the instructions provided in the "CLEANING SANITISATION DISINFECTION" paragraph.

CLEANING SANITISATION DISINFECTION

Switch the device off and disconnect it from the power supply before proceeding with any cleaning operation.

CLEANING THE DEVICE (A) AND THE TUBE'S OUTER SURFACE (B)

Only use a cloth with antibacterial detergent (non-abrasive and with no solvents of any kind)

ACCESSORIES

We recommend personal use of accessories to avoid any risk of spreading the infection.

Disassemble the accessories as described below:

C1 - Open the nebulizer by turning the upper part (C1.3) anticlockwise, remove the nozzle from the upper part (C1.3) and then disassemble it as shown in the "Assembly diagram" in section C1.2.

C5 - Disassemble it as shown in the "Assembly diagram" in section C5.

SANITISATION AND DISINFECTION

SANITISATION

Before and after each use, sanitise the nebuliser cup and the accessories, choosing one of the methods described below.

method A: Sanitise accessories **C1.1-C1.2-C1.3-C1.4-C2-C3-C4-C5.1A-C5.1B-C5.1C-C5.2-C5.3-C5.4-C5.5-C6** under potable hot water (approximately 40°C) with a gentle, non abrasive dish detergent.

method B: Sanitise accessories **C1.1-C1.2-C1.3-C1.4-C2-C3-C4-C5.1A-C5.1B-C5.1C-C5.2-C5.3-C5.4-C5.5-C6** in the dishwasher with a hot cycle.

method C: Sanitise accessories **C1.1-C1.2-C1.3-C1.4-C2-C3-C4-C5.1A-C5.1B-C5.1C-C5.2-C5.3-C5.4-C5.5-C6** by immersing them in a solution of 50% water and 50% white vinegar, then rinse thoroughly under potable hot water (approximately 40 °C).

If you want to also perform the cleaning for DISINFECTION, jump to the DISINFECTION paragraph.

After having sanitised the accessories, shake them vigorously and place them on a paper towel. Alternatively, dry them with a jet of hot air (for example, a hair dryer).

DISINFECTION

After sanitising the nebuliser cup and the accessories, disinfect them choosing one of the methods described below.

method A: Accessories **C1.1-C1.2-C1.3-C1.4-C2-C3-C4-C5.1A-C5.1B-C5.1C-C5.2-C5.3-C5.4-C5.5-C6** can be disinfected.

The disinfectant must be an electrolytic chloroxidizer (active principle: sodium hypochlorite) specific for disinfecting, which is available in any pharmacy.

Implementation:

- Fill a container big enough to hold all of the parts to disinfect with a solution of potable water and disinfectant, according to the proportions indicated on the packaging of the disinfectant.
- Completely immerse each part in the solution, taking care to avoid the formation of air bubbles on the parts. Leave the parts immersed for the amount of time indicated on the packaging of the disinfectant associated with the concentration chosen for the solution.
- Remove the disinfected parts and rinse abundantly with warm potable water.
- Dispose of the solution following the instructions provided by the disinfectant manufacturer.

method B: Disinfect the accessories **C1.1-C1.2-C1.3-C1.4-C2-C5.1A-C5.1B-C5.1C-C5.2-C5.3-C5.4-C5.5-C6** by boiling them in water for 10 minutes; use demineralised or distilled water to prevent calcium deposits.

method C: Disinfect the accessories **1C1.1-C1.2-C1.3-C1.4-C2-C5.1A-C5.1B-C5.1C-C5.2-C5.3-C5.4-C5.5-C6** with a hot steam steriliser for baby-bottle (not the microwave type). Perform the process faithfully following the instructions of the steriliser. To ensure that the disinfection is effective, choose a steriliser with an operating cycle of at least 6 minutes.

After having disinfected the accessories, shake them vigorously and place them on a paper towel. Alternatively, dry them with a jet of hot air (for example, from a hair dryer).

At the end of each use store the device complete with accessories in a dry place away from dust.


AIR FILTERING

The unit is equipped with an air filter (A3) which should be replaced when dirty or when its colour changes. Do not wash or re-use the used filter. The filter must be regularly replaced to help ensure efficient compressor performance. The filter should be regularly checked. For spare filters contact your dealer or an authorised service centre.

To replace the filter:

Lift the lid and remove the filter using a screwdriver (see fig.). The filter is designed to remain fixed in its housing.

Only use original accessories and spare parts by Flaem, we disclaim any liability in the event of using non original spare parts or accessories.

TROUBLESHOOTING

Switch off the device before any procedure and unplug the power cable from the socket.

PROBLEM	CAUSE	SOLUTION
The device does not work	The power cable has not been correctly inserted into the socket of the device or the power socket	Correctly insert the power cable in sockets.
The device does not nebulize or nebulizes insufficiently.	The medication has not been inserted in the nebulizer	Pour the right amount of medication in the nebulizer
	The nebulizer has not been properly fitted.	Disassemble the nebulizer and reassemble it correctly as shown in the connection diagram on the cover.
	The nebulizer nozzle is clogged	Disassemble the nebulizer, remove the nozzle and perform cleaning operations. Failure to remove medication deposits from the nebulizer affects its efficiency and operation. Strictly comply with the instructions contained in the "CLEANING, SANITISATION, DISINFECTION" chapter
	The air pipe is not correctly connected to the device	Ensure the device's accessories are properly connected to the air intake (see connection diagram on the cover).
	The air pipe is bent, damaged or twisted	Unwind the pipe and ensure it is not crushed or punctured. If necessary replace it.
	The air filter is dirty	Replace the filter
The apparatus is noisier than usual	The Filter is not correctly inserted in its housing	Fully and properly insert the Filter in its housing

If after verifying the above mentioned conditions the device should not operate properly, contact your trusted retailer or the nearest authorized FLAEM service center.

SYMBOLS

TÜV approval

EC Marking medical ref. Dir 93/42 EEC and subsequent updates

Class II device

Attention: check the instructions for use

"OFF" for part of Equipment

"ON" for part of equipment

In compliance with: European Standard EN 10993-1 "Biological Evaluation of medical devices" and European Directive 93/42/EEC "Medical Devices." Phthalate-free. In compliance with: Reg. (EC) no. 1907/2006

Minimum and maximum room temperature

Minimum and maximum atmospheric pressure

Device serial number

Manufacturer

Type BF applied part

Alternating current

Keep dry

Minimum and maximum air moisture

TECHNICAL DATA

Mod. P0406EM F700

Voltage:	230V~ 50Hz 140VA	115V ~ 50Hz	220V ~ 50Hz
Safety certifications:			

Max Pressure:	2.6 ± 0.4 bar
Compressor air output:	10 l/min approx
Dimensions (L)x(P)x(H)	21x16x13 cm
Weight	1.600 Kg
Noise level (at 1 m):	55 dB (A) approx
Operation:	Continuous use

APPLIED PARTS

Type BF applied parts include: patient accessories (C2,C3,C4,C5.1A-C5.1B-C5.1C)

RF6 Basic² Nebulizer

Medication minimum capacity	2 ml
Medication maximum capacity:	8 ml
Operating pressure (with neb.)	0.8 bar

	Mod. High speed	Mod. standard
⁽¹⁾ Delivery:	0.42 ml/min approx.	0.32 ml/min approx.
⁽²⁾ MMAD:	3.6 µm	3.6 µm
⁽²⁾ Breathable fractions < 5 µm (FPF):	67%	66%

(1) Data detected according to Flaem I29-P07.5 internal procedure. (2) In vitro characterization carried out at the Inamed Research GmbH & Co. KG on behalf of TÜV Rheinland Product Safety GmbH – Germany in compliance with the European Standard for aerosol therapy appliances, Standard EN 13544-1. More details are available on request.

NASAL SHOWER SPECIFICATIONS

Mod.: Rhino Clear[®]

Capacity:	drug 10 ml maximum
10 ml dispensing time ⁽¹⁾	2'
MMD ⁽²⁾	48 µm
% <10 µm ⁽²⁾	3.2%

⁽¹⁾ Tests carried out with 0,9% NaCl physiological solution according to Flaem internal procedure

⁽²⁾ Data gathered with Laser Malvern Mastersizer system (in a certified laboratory)

Operating conditions:

Temperature: min 10°C; max 40°C
Air humidity: min 10%; max 95%
Atmospheric pressure: min 69KPa; max 106KPa

Storage conditions:

Temperature: min -25°C; max 70°C
Air humidity: min 10%; max 95%
Atmospheric pressure: min 69KPa; max 106KPa

ELECTROMAGNETIC COMPATIBILITY

This device was designed to satisfy the currently required requisites for electromagnetic compatibility (EN 60601-1-2). Electro-medical devices require particular care during installation and use relative to EMC requirements. Users are therefore requested to install and/or use these devices following the manufacturer's specifications. There is a risk of potential electromagnetic interference with other devices. RF mobile or portable radio and telecommunications devices (mobile telephones or wireless connections) can interfere with the functioning of electro-medical devices. For further information visit our website www.flaemnuova.it. The Medical Device may be subject to electromagnetic interference if other devices are used for specific diagnosis or treatments. Flaem reserves the right to make technical and functional modifications to the product with no prior warning.

DEVICE DISPOSAL

In compliance with the Directive 2012/19/EC, the symbol printed on the device shows that the device to be disposed of is considered waste and must therefore be an item of "differentiated collection". Consequently, the user must take it (or have it taken) to the designated collection sites provided by the local authorities, or turn it in to the dealer when purchasing an equivalent new device. Differentiated waste collection and the subsequent treatment, recycling and disposal procedures promote the production of devices made with recycled materials and limit the negative effects on the environment and on health caused by potential improper waste management. The unlawful disposal of the product by the user could result in administrative fines as provided by the laws transposing Directive 2012/19/EC of the European member state or of the country in which the product is disposed of.

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© Copyright by FLAEM NUOVA S.p.A., 2020
FLAEM NUOVA S.p.A. non risponde di ulteriori estensioni del periodo di garanzia assicurate da terzi.