FLAEM[®]

Medical Devices

SpaceNeb

IT > MANUALE ISTRUZIONI D'USO
EN > OPERATING INSTRUCTIONS MANUAL
DE > BEDIENUNGSANLEITUNG
FR > MODE D'EMPLOI
NL > GEBRUIKSAANWIJZING
EL > ΕΓΧΕΙΡΙΔΙΟ ΟΔΗΓΙΩΝ ΧΡΗΣΗΣ
PL > INSTRUKCJA OBSŁUGI I KONSERWACJI
ES > MANUAL DE INSTRUCCIONES DE USO
SV > BRUKSANVISNING
HR > PRIRUČNIK S UPUTAMA ZA UPORABU
DA > BRUGSANVISNING
CS > NÁVOD K POUŽITÍ
SK > NÁVOD NA POUŽITIE
RO > MANUAL DE INSTRUȚIUNI DE UTILIZARE
FI > KÄYTTÖOHJE
PT > MANUAL DE INSTRUÇÕES DE UTILIZAÇÃO
BG > РЪКОВОДСТВО ЗА УПОТРЕБА
AR < دليل تعليمات الاستخدام

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Rhino Clear
Space
(ACO820P)

Aerosol therapy device

These operating instructions are provided for the P0905EM REF 1 and RF8-2E model devices. The FLAEM Aerosol Therapy device consists of a compressor unit (A), a nebuliser and some accessories (B).

INTENDED PURPOSE. Medical Device for the administration of medication through inhalation, with inhalation therapy and medication prescribed by a doctor.

OPERATING INSTRUCTIONS. Treatment of respiratory diseases. Medications must be prescribed by a doctor who has assessed the patient's general condition.

⚠️ CONTRAINDICATIONS. • The medical device should NOT be used for patients who are unable to breathe on their own or who are unconscious • Do not use the device in anaesthetic or assisted ventilation circuits.

INTENDED USERS. The devices are intended for use by legally authorised medical personnel/health workers (doctors, nurses, therapists, etc.). The device can be used directly by the patient.

⚠️ TARGET PATIENT GROUP. Adults, children of all ages, infants. Before using the device, the user manual must be read carefully and an adult responsible for safety must be present if the device is used by infants, children of any age or persons with limited abilities (e.g. physical, mental, or sensory). It is the medical personnel's responsibility to assess the patient's condition and capabilities in order to determine, when prescribing use of the device, whether the patient is able to operate the aerosol device safely on their own or whether the therapy should be administered by a responsible person.

Medical personnel should be responsible for assessing the use of the device on specific types of patients such as pregnant women, lactating women, infants, incapacitated persons or persons with limited physical capabilities.

OPERATING ENVIRONMENT. This device can be used in healthcare facilities, such as hospitals, outpatient clinics, etc., or even at home.

⚠️ WARNINGS ABOUT POSSIBLE MALFUNCTIONS. • Should your appliance fail to perform, please contact the authorised service centre for clarification. • The manufacturer should be contacted to report problems and/or unexpected events relating to operation and if necessary for clarification of use and/or maintenance/hygienic preparation. • Please also refer to the case history of faults and their resolution.

WARNINGS. • Use the device only as a therapeutic inhaler. This medical device is not intended as a life-saving device. Any other use is considered misuse and there may be hazards. The manufacturer is not liable for any improper use.

• Always consult your general practitioner for identification of treatment. • Follow the instructions of your doctor or respiratory rehabilitation therapist regarding the type of medicine, the dosage and the treatment indications. • If you should experience allergic reactions or other problems while using the device, stop using it immediately and consult your doctor. • Keep this manual safe for further reference. • Do not expose the device to extreme temperatures. • If the packaging is damaged or opened, contact the distributor or service centre. • The time required to switch from storage to operating conditions is about 2 hours. • The appliance must not be used like a toy. The utmost attention should be paid when used by children. • In the event of damage to the device after being dropped (e.g. detached parts), do not use the medical device but contact an authorised service centre. • Do not use the appliance while driving, or in any other situation where you should not be distracted in order not to endanger yourself, people or pets nearby and the objects around you. • Before using it for the first time and periodically during the product's life, check the integrity of the device's structure to make sure that it is not damaged; should there be damage, immediately take the product to an authorised service centre (see Service sheet) or to your usual dealer. • It is forbidden to access the device opening in any way. Repairs may only be carried out by personnel authorised by the manufacturer. Unauthorised repairs invalidate the warranty and may pose a hazard to the user.

Risk of suffocation: Some parts of the appliance are small enough to be swallowed by children; therefore, keep the appliance out of the reach of children.

Risk of strangulation: Do not use the supplied connecting hose and cables outside of their intended use, they could cause a strangulation hazard. Take special care for children and people with special difficulties, often these people are not able to correctly assess the dangers.

Fire hazard: Do not place the appliance near sources of heat, light or in excessively hot and humid premises. Store as described in the chapter on ENVIRONMENTAL CONDITIONS. - THIS APPLIANCE IS NOT suitable for use in the presence of a flammable anaesthetic mixture with oxygen or nitrous oxide. - The device must be recharged via the supplied cable. The device must be recharged from external power sources as specified in the section "Technical Specifications External Power Source".

Risk of electrocution: - Do not use the appliance in the shower or bath and do not use it with wet hands. - The appliance does not contain parts that can be repaired by the user. - Do not wash the appliance under running water or by immersion. - If you use a power supply unit for battery charging, it must be certified according to EN 60601-1. - Keep the power cable away from animals (e.g. rodents), otherwise such animals could damage the insulation of the power cable.

Risk of ineffectiveness of the treatment: - Performance may vary with particular types of drugs (e.g. those with high viscosity or in suspension). For further information, please refer to the package leaflet provided by the drug manufacturer. - Use the nebuliser in the correct position, as upright as possible; do not tilt it beyond an angle of 30 degrees, in any direction, to prevent the drug from spilling into the mouth or being excessively dispersed, reducing the effectiveness of the treatment. - Pay attention to the instructions provided with the drug and avoid using the devices with substances and dilutions other than those recommended. - Only use original Flaem accessories or spare parts, no liability is accepted if non-original parts or accessories are used

Risk of infection: - We recommend personal use of the accessories to avoid any risk of infection through transmission. - Follow hygienic preparation before each use. - At the end of therapy, do not leave the medicine inside the medicine Tray and carry out the hygienic preparation operations. - If the nebuliser is used for several types of medicine, the residues must be completely removed. Therefore, perform hygienic preparation after each inhalation, also to achieve the highest degree of hygiene and to optimise the service life and operation of the device.

Risk of injury: Improper use of the pressure-generating device at the fitting can cause serious injury. Do not place the appliance on a soft supporting surface such as a sofa, a bed or a tablecloth. - Always operate it on a hard surface clear of any objects.

WARNINGS ON RISKS OF INTERFERENCE DURING USE IN DIAGNOSTIC INVESTIGATIONS. This appliance device is designed to meet the current requirements for electromagnetic compatibility. As far as EMC requirements are concerned, electromedical devices require special care when being installed and used, and are therefore required to be installed and/or used in accordance with the manufacturer's specifications. Risk of potential electromagnetic interference with other devices. Mobile or portable RF radio and telecommunication devices (mobile phones or wireless connections) may interfere with the operation of electromedical devices. The device may be susceptible to electromagnetic interference in the presence of other devices used for specific diagnosis or treatment. For more information visit www.flaemuova.it.

TROUBLESHOOTING

Problem	Cause	Remedy
The device switches on but does not nebulise or does it poorly.	The appliance is not set to maximum speed (the only speed compatible with the nebuliser chamber).	Set the appliance to maximum speed.
	The medication has not been inserted in the nebuliser.	Pour the right amount of medication in the nebuliser.
	The nebuliser has not been assembled correctly.	Disassemble and reassemble the nebuliser correctly according to the connection diagram on the cover.
	The nebuliser is obstructed.	Perform the hygienic preparation of the nebuliser. Deposits of medicine due to lack of hygienic preparation of the nebuliser impair its efficiency and function. Strictly follow the instructions in the HYGIENIC PREPARATION chapter.
The device does not switch on.	The air filter is dirty.	Replace the filter.
	Fully discharged battery.	Charge battery.

If, after checking the conditions described above, the device still is not working properly, we recommend that you contact your trusted dealer or an authorised FLAEM service centre nearest to you. You can find the list of all the Service Centres at <http://www.flaemnuova.it/it/info/assistenza>

DISPOSAL

Appliance

 In accordance with Directive 2012/19/EC, the symbol on the equipment indicates that the equipment to be disposed of (excluding accessories) is considered as waste and must therefore be subject to 'separate collection'. Therefore, the user must deliver this waste (or have it delivered) to the separate collection centres set up by the local authorities, or hand it over to the retailer when purchasing a new appliance of an equivalent type. Separate waste collection and subsequent treatment, recovery and disposal operations promote the production of equipment from recycled materials and limit the negative environmental and health effects caused by improper waste management. Unauthorised disposal of the product by the user entails the application of the administrative sanctions provided for in the laws transposing Directive 2012/19/EC of the member state or country where the product is disposed of.

Accessories. Accessories such as masks and mouthpiece are to be disposed of as general waste after a disinfection cycle.

Packaging

 Product box	 Battery heat shrink film	 Product packaging bag	 Accessory heat shrink film
---	--	---	---

REPORTING SERIOUS EVENTS. Please report any serious incidents occurring in connection with this device to the manufacturer and to the competent authority of the Member State where you reside. An event is considered serious if it causes or may cause, directly or indirectly, death or an unforeseen serious deterioration in a person's state of health.

INFORMATION ON RESTRICTIONS OR INCOMPATIBILITIES WITH CERTAIN SUBSTANCES. • Interactions: The materials used in the device are biocompatible materials and comply with statutory regulations, however possible allergic reactions cannot be completely excluded.

SYMBOLS ON DEVICE OR PACKAGING

 Unique device identifier	 Serial number of the device	 Class II device	 Production date
 Manufacturer	 Type BF applied part	 Direct current	 Moisture limits
 Model number	 Phthalate- and bisphenol-free	 Attention	 Identifies the model variant
 Medical device	 Temperature limits	 Atmospheric pressure limits	 Before use: Caution check instructions for use
 Speed level selector button	 ON/OFF button	 See instructions for use	 Batch code
IP21 Protection rating of the envelope: IP21. (Protected against solid bodies larger than 12 mm. Protected against access with a finger)	 Medical CE Marking ref. Regulation 2017/745 EU as amended	 Quality mark	

COMPRESSOR UNIT TECHNICAL SPECIFICATIONS. Model: P0905EM REF 1

Power supply: 5V DC 2A Max. pressure: 1.5 ± 0.3 | Air flow to compressor: 6 l/min approx
bar | Noise level (at 1 m): 49 dB (A) approx

Operation: Use 30 minutes ON/ 40 minutes OFF
Dimension: 7,5(L) x 7,5(D) x 15(H) cm

NEBULISER TECHNICAL SPECIFICATIONS. Model: RF8-2E

RF8 Evo nebuliser

Minimum drug capacity: 2 ml

Maximum drug capacity: 8 ml

TECHNICAL CHARACTERISTICS POWER SOURCE

External: 5V DC 2A

OPERATION AND BATTERY AUTONOMY. The appliance runs on internal lithium batteries. The battery maintains maximum efficiency for 300 charging cycles. Thereafter, the runtime will progressively decrease and battery replacement is recommended; this procedure must only be carried out by an authorised service centre. Autonomy with new battery at full charge is approx. 50 minutes. When not in use for long periods of time, however, recharge the device's battery at least every 3 months in order to avoid possible damage to the battery.

INTERNAL BATTERY CHARGING TIME

Approximately 2 hours The battery must be recharged at an ambient temperature between 5 and 40°C.

TECHNICAL SPECIFICATIONS OF DEVICE. Model: P0905EM REF 1 combined RF8-2E

Operating pressure at maximum speed (with neb.) [bar]: **0.60**

	Without particle selector B2.2	With particle selector B2.2
Dispensing speed: [ml/min]	0.32 approx ⁽¹⁾	0.23 approx ⁽¹⁾
(Aerosol output) AO: [ml]	0.222 ⁽³⁾	0.285 ⁽²⁾
Percentage of filling volume emitted per minute: [%].	3.12 ⁽³⁾	2.18 ⁽²⁾
(Aerosol output rate) AOR: [ml/min]	0.062 ⁽³⁾	0.044 ⁽²⁾
(Mass Median Aerodynamic Diameter) MMAD: [µm]	4.31 ⁽³⁾	2.75 ⁽²⁾
(Geometric standard deviation) GSD:	2.39 ⁽³⁾	2.87 ⁽²⁾
Breathable fraction (< 5 µm): [%]	56.7 ⁽³⁾	71.6 ⁽²⁾
Aerosol particle fraction < 2 µm: [%].	19.0 ⁽³⁾	37.9 ⁽²⁾
Aerosol particle fraction > 2 µm ; < 5 µm: [%]	37.7 ⁽³⁾	33.7 ⁽²⁾
Aerosol particle fraction > 5 µm: [%]	43.3 ⁽³⁾	28.4 ⁽²⁾
Residual volume: [ml]	1.094 ⁽³⁾	0.951 ⁽²⁾

(1) The delivery rate was measured using 0.9% saline solution at 23°C by weight loss, in accordance with Flaem internal procedures.

Note: if the data are marked with both (2) and (3), proceed; if the data are marked only with (2), skip to the section below the line.

(2)(3) In vitro characterisation in accordance with EN ISO 27427 using 4 ml of 0.1% salbutamol in 0.9% NaCl solution and adult respiratory parameters. The performance may vary for the paediatric or infant population.

(2) Characterisation performed by TÜV Rheinland Italia S.r.l. in collaboration with the University of Parma.

(3) Characterisation performed by the University of Parma.

More details are available on request

APPLIED PARTS

Type BF applied parts are: patient accessories (B3, B4, B5)

Weight: 0.400 kg

ENVIRONMENTAL CONDITIONS

Operating conditions:

Ambient temperature: Between +5°C and +40°C

Relative air humidity: Between 15% and 90%

Atmospheric pressure: Between 69KPa and 106 KPa

Storage and transport conditions:

Ambient temperature: Between 15°C and +35°C

Relative air humidity: Between 10% and 75%

Atmospheric pressure: Between 69KPa and 106 KPa

DURATION

Model: P0905EM REF 1

(Compressor unit)

Service life 300 hours.

Model: RF8-2E

(Nebuliser and accessories)

The expected average life is 1 year, however it is advisable to replace the nebuliser every 6 months during intensive use (or sooner if the nebuliser is clogged) to ensure maximum therapeutic effectiveness.

APPLIANCE AND MATERIAL INFORMATION

The equipment includes:		Information on materials
A -	Compressor unit - Model:P0905EM REF 1	
	A1 - Switch, A2 - Air Inlet	
	A3 - Air filter	
	A3.1 - Filter, A3.2 - Filter housing cap	
	A4 - Speed selector button	
	A4.1 - Speed level indicator LED	
B -	A5 - USB-C charging socket	
	A5.1 - Socket protection cap, A5.2 - Battery status indication LED	
	A6 - Power cable (USB-C)	
	Nebuliser and accessories - Model: RF8-2E	
	B1 - Connecting hose (compressor/nebuliser unit) 2 m	
	B1.1 - Connection	
B -	B2 - RF8 Evo nebuliser	
	B2.1 - Lower part, B2.2 - Particle selector, B2.3 - Upper part, B2.4 - Drugsaver	
	B3 - Mouthpiece	Polypropylene
	B4 - Adult SoftTouch mask	Polypropylene +
C -	B5 - Paediatric SoftTouch mask	Thermoplastic Elastomers
	IMPORTANT NOTE: Remove the identification label on the packaging and affix it in the appropriate spaces on page 4. When replacing the Nebuliser and accessories, perform the same procedure.	
Carry bag		

OPERATING INSTRUCTIONS

Before each use, wash your hands thoroughly and clean your appliance as described in the section "HYGIENE PREPARATION". During application, it is recommended to adequately protect oneself from any dripping. This appliance is suitable for the administration of medicinal substances (solutions and suspensions), for which aerosol administration is envisaged; such substances must in any case be prescribed by a doctor. In the case of substances that are too dense, dilution with suitable saline solution may be necessary, as prescribed by a doctor.

MODE OF USE

Mode 1 - Direct use (Appliance held by the user) (p. 2.). The device is held directly by the user during drug administration.

Mode 2 - Use with connecting hose (p. 3). The compressor is placed on a flat, stable surface (e.g. a table). The ampoule is connected to the main unit via the specific flexible hose, allowing the nebuliser to be brought closer to the patient's face.

This configuration is especially ideal for: paediatric patients, elderly patients or those requiring assistance who may have difficulty holding the device or cases where greater freedom of movement is required.

1. Insert the lower part (B2.1) into the compressor (A) (Fig. 1). Rotate the lower part by 90° so that the drug quantity window is visible.
2. Insert the prescribed drug into the lower part (B2.1) and return the nebuliser to its original position (Fig.1). Assemble the nebuliser according to the instructions in the section "INSTRUCTIONS FOR USE OF THE RF8 Evo NEBULIZER".
3. Connect the accessories as shown in the 'Connection Diagram' on the cover.
4. If you want to use the hose, connect the various devices as shown in the "Connection diagram" on the cover (page 3).
5. Sit comfortably while holding the nebuliser in your hand, place the mouthpiece to your mouth or use a nasal piece (if provided) or mask. If you use the mask accessory, place it on your face as shown in (Fig. 2).
6. Switch on the device by operating the switch (A1) and inhale and exhale deeply. After breathing in, it is recommended to hold your breath for a moment so that the inhaled aerosol droplets can settle. Then exhale slowly.
7. When the application is finished, switch off the device by pressing the button (A1).

ATTENTION: Once the appliance is switched on, LED A5.2 indicates the battery status see section "BATTERY STATUS".

*When switched on, the device is at maximum speed see (LED A4.1), the only speed compatible with aerosol therapy ampoule use. If you accidentally change the speed, reset the maximum speed by pressing the speed selector key (A4).

Speeds other than maximum speed are only to be used when the medical device is used in combination with the Flaem nasal douche.

INSTRUCTIONS FOR USING THE RF8 Evo NEBULIZER

The RF8 Evo nebuliser can be used in two different configurations that can be selected according to the required performance. The specifications of each configuration can be found in the section "APPLIANCE TECHNICAL SPECIFICATIONS. Model: P0905EM REF 1 combined RF8-2E"

Configuration without particle selector B2.2 (Fig. 6)	Configuration with particle selector B2.2 (Fig. 3-4)
Nebulisation suitable for the mid-to-upper respiratory tract (tracheobronchial).	Nebulisation suitable for the lower respiratory tract (alveoli).
- Insert the upper part of B2.3 into the lower part of B2.1. Turn clockwise (Fig. 6). - Apply the Drugsaver (B2.4) (Fig. 7).	- Set the particle selector switch B2.2 on the upper part of B2.3 (Fig. 3). - Insert the upper part (B2.3 + B2.2) into the lower part of B2.1. Turn clockwise (Fig. 4). - The "DrugSaver" cap limits drug leakage during the treatment and reduces the risk of the drug entering the eyes during exhalation (Fig. 5).

BATTERY STATUS

LED A5.2 status	Meaning
Steady green light	When connected to USB and with battery charge completed (device switched off)
Steady blue light	Device in operation (standalone or while charging)
Flashing yellow light	Battery charging, when connected to USB-C power supply (device switched off)
Steady blue light and flashing red light 3 times every 2 minutes	Low battery warning, still allows use (standalone or while charging)
Steady red light for 10 seconds, immediate switch-off.	Low battery ALARM, device shuts down (standalone and while charging)

CHARGING THE BATTERY

Locate the socket: Find the charging socket (A5) on the body of the device.


Remove the protection: Gently lift the protective cap (A5.1).

Connect the cable: Plug the USB-C cable into the socket (A5) and connect the other end to an external power source (5V DC 2A).

Waiting time: A full charge takes about 2 hours.

The operation of the device is guaranteed even while charging the battery.

SoftTouch masks

The **SoftTouch** masks have an external edge made in **soft, biocompatible material** that guarantees optimal adherence to the face and are also equipped with the **Leaking Limiting device**. These characteristic features allow for greater sedimentation of the drug in the patient, and again **limit its dispersion**.

Upon inhalation, the tab, which acts as a **Dispersion Limiting Device**, **bends inwards in the mask**.

Upon exhalation, the tab, which acts as a **Dispersion Limiting Device** **bends outwards from the mask**.

HYGIENIC PREPARATION

Switch off the appliance before each hygienic preparation operation and disconnect it from the mains.

Compressor unit (A), hose exterior (B1) and connector (B1.1).

Use only a cloth moistened with antibacterial detergent (non-abrasive and free of solvents of any kind).

Nebuliser and accessories

Disassemble the nebuliser as shown in the "Connection diagram" under B2.

Then proceed according to the instructions below.

Sanitisation

Before and after each use, sanitise the nebuliser and accessories by choosing one of the methods provided in the table and described below.

method A: Sanitise the accessories under warm drinking water (approx. 40°C) with mild dishwashing liquid (non-abrasive).

method B: Sanitise the accessories in the dishwasher with the hot cycle (70°C).

method C: Sanitise the accessories by soaking in a solution of 50% water and 50% white vinegar, then rinse thoroughly with warm drinking water (approx. 40°C).

After sanitising the accessories, shake them vigorously and place them on a paper towel, or alternatively dry them with a hot air jet (e.g. hair dryer).

Disinfection

After sanitising the nebuliser and accessories, disinfect them using one of the methods provided in the table and described below. Each method can be performed for a limited number of times (see figure in table).

method A: Obtain a disinfectant of the electrolytic chloride type (active ingredient: sodium hypochlorite), specifically for disinfection, available in all pharmacies.

Execution: - Fill a container big enough to hold all of the parts being disinfected with a solution of potable water and disinfectant, according to the proportions indicated on the packaging of the disinfectant. - Completely immerse each part in the solution, taking care to avoid the formation of air bubbles on the parts. Leave the components immersed for the period of time indicated on the disinfectant packaging, and associated with the concentration chosen for the preparation of the solution. - Remove the disinfected parts and rinse abundantly with warm potable water. - Dispose of the solution according to the instructions provided by the disinfectant manufacturer.

method B: Disinfect accessories by boiling in water for 10 minutes; use demineralised or distilled water to avoid lime deposits.

method C: Disinfect the accessories in a hot bottle steamer (not microwave). Carry out the process strictly according to the hot bottle steamer's instructions. For disinfection to be effective, choose a steriliser with an operating cycle of at least 6 minutes.

If you also want to carry out sterilisation, skip to the **Sterilisation** section.

After disinfecting the accessories, shake them vigorously and place them on a paper towel, or alternatively dry them with a hot air jet (e.g. hair dryer).

At the end of each use, store the device complete with its accessories in a dry and dust-free place.

Sterilisation: Equipment: Steam steriliser with fractionated vacuum and overpressure compliant with EN 13060.

Execution: Package each individual component to be treated in sterile barrier system or packaging in accordance with EN 11607. Place the packaged components in the steam steriliser. Perform the sterilisation cycle in accordance with the instructions for use of the equipment by selecting a temperature of 134°C and a time of 4 minutes first.

Storage: Store sterilised components according to the instructions for use of the system or sterile barrier packaging chosen.

The sterilisation procedure was validated in accordance with ISO 17665-1.

At the end of each use, store the device complete with its accessories in a dry and dust-free place.

Table of planned methods / patient accessories								
Method	B2.1	B2.2	B2.3	B2.4	B3	B3.2	B4	B5
HYGIENIC PREPARATION AT HOME								
Sanitisation								
method A	✓	✓	✓	✓	✓	✓	✓	✓
method B	✓	✓	✓	✓	✓	✓	✓	✓
method C	✓	✓	✓	✓	✓	✓	✓	✓
Disinfection								
method A	300	300	300	300	300	300	300	300
method B	300	300	300	300	300	300	300	300
method C	300	300	300	300	300	300	300	300
HYGIENIC PREPARATION IN THE CLINIC OR HOSPITAL								
Disinfection								
method A	300	300	300	300	300	300	300	300
method B	300	300	300	300	300	300	300	300
Sterilisation								
	300 Assembled nebuliser				300	300	300	300
✓: foreseen \: not foreseen 300: ✓ MAX 300 TIMES								

AIR FILTRATION

The appliance is equipped with an extraction filter (A3.1) to be replaced when it is dirty or changes colour. Do not wash or re-use the same filter. Regular replacement of the filter is necessary to help ensure proper compressor performance. The filter must be checked on a regular basis. Contact your dealer or authorised service centre for replacement filters.

To replace the filter, pull it out as shown in (Fig.8).

The filter was designed to always remain secured in its housing. Do not replace the filter during use. **Only use original Flaem accessories or spare parts, no liability is accepted if non-original parts or accessories are used.**

NASAL DOUCHE

For complete care of the entire respiratory system, we recommend using it in combination with the Rhino Clear Space nasal douche (Model RC3-1). While aerosols act on the lower airways (lungs and bronchi), Rhino Clear Space is specifically designed to treat the upper respiratory tract, providing deep, targeted therapy for the nose and throat.

Rhino Clear Space is available for sale separately (Code ACO820P) or inside the package in the full version.

Please refer to the nasal douche user manual for proper use.

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FLAEM NUOVA S.p.A. non risponde di ulteriori estensioni del periodo di garanzia assicurate da terzi.