



close to your baby

 SINCE 1966 
MADE IN ITALY

Baby Nasal Aspirator

Hygienic


Practical and easy to use

Comfortable soft
nasal adaptor



**PHTHALATES
& BPA FREE**

www.flaem.it




Nasal Aspirator

Baby Flaem Nasal Aspirator does not require a person to create suction. The patented system makes it easy to use by simply pushing in the plunger same as you would with a syringe and slowly releasing to suck out the mucus from the child's nostril.

In addition the nozzle can stay safely in the baby's nostril because no excessive pressure is created when pushing the plunger due to the special configuration of the device.

The nasal aspirator can simply be disassembled and assembled for easy cleaning and sanitizing.



6 - unit DISPLAY BOX		
SINGLE UNIT BOX including 1 Aspirator and 1 spare soft nozzle	ACO423P	 8 018294 010904

Patented and Patent Pending Product 14733A0 Copyright by FLAEM NUOVA S.p.A., 2011 - FLAEM NUOVA reserves the right to make technical and functional alterations to this product without prior notice.



Manufacturer:
 FLAEM NUOVA S.p.A.
 Via Colli Storici, 221
 25015 S. Martino della Battaglia (BS) Italy
 Tel. +39 030 9910168 - Fax +39 030 9910287
 info@flaemnuova.it
www.fluem.it



Export Department of FLAEM NUOVA

medexpo international s.r.l.

3, Via del Passero - 20147 Milano - Italy
 Phone (39) 02 4120769
 Fax (39) 02 48376295
 E-mail: info@medexpointernational.com