



*close to your baby*

 SINCE 1966   
MADE IN ITALY

## ***Baby Nasal Aspirator***

Hygienic


Practical and easy to use

Comfortable soft nasal adaptor



**PHTHALATES  
& BPA FREE**

[www.flaem.it](http://www.flaem.it)



## Nasal Aspirator

Baby Flaem Nasal Aspirator does not require a person to create suction. The patented system makes it easy to use by simply pushing in the plunger same as you would with a syringe and slowly releasing to suck out the mucus from the child's nostril.

In addition the nozzle can stay safely in the baby's nostril because no excessive pressure is created when pushing the plunger due to the special configuration of the device.

The nasal aspirator can simply be disassembled and assembled for easy cleaning and sanitizing.



	6 - unit DISPLAY BOX		
	SINGLE UNIT BOX including 1 Aspirator and 1 spare soft nozzle	<b>ACO423P</b>	<p>8 018294 010904</p>

Patented and Patent Pending Product 14733A0 Copyright by FLAEM NUOVA S.p.A., 2011 - FLAEM NUOVA reserves the right to make technical and functional alterations to this product without prior notice.



Manufacturer:  
FLAEM NUOVA S.p.A.  
Via Colli Storicci, 221  
25015 S. Martino della Battaglia (BS) Italy  
Tel. +39 030 9910168 - Fax +39 030 9910287  
info@flaemnuova.it  
**www.fluem.it**



Export Department of FLAEM NUOVA

**medexpo international s.r.l.**

3, Via del Passero - 20147 Milano - Italy  
Phone (39) 02 4120769  
Fax (39) 02 48376295  
E-mail: info@medexpointernational.com