



RespirAir

Reliable and compact
nebulizer system
for frequent use



2-speed Nebulizer
with valves



Comfortable
adult and paediatric
Soft Touch masks



High efficiency
air filter



The respiratory
therapy

MADE
IN ITALY

www.flaem.it

RespirAir

Medical Device for the administration of drugs by inhalation, whose inhalation therapy and drugs must be prescribed by the doctor.



Technical Features

Use	Continuous	
Power supply/Power	230V ~ 50Hz / 140VA	
Air filtering	The device is equipped with a professional inlet air filter, which blocks impurities and preserves the integrity of the compressor, helping to guarantee the useful life of the appliance.	
Compressor air output	10 l/min (approx)	
Max pressure	2.6 ± 0.4 bar	
Sound level (at 1 m)	54 dB (A) (approx)	
Dimensions (length x width x height)	16 x 16 x 10 cm	
Weight	1,45 Kg	
Compliance	In compliance with EU regulation 2017/745 Medical Device Directive and following edits	CE 0051
	Phthalates free. In conformity with Reg. (EC) no. 1907/2006	PHthalates & BPA FREE
	Quality mark	TUV Rheinland CERTIFIED

RF7 Dual Speed Plus nebulizer technical features

2-speed Nebulizer with valves: effectiveness and speed of treatment are the main features of this nebulizer with two dispensing speeds and a system of valves, which, thanks to its new internal geometry, allows the nebulisation of the medication with optimal particle size suitable for the treatment of the entire respiratory tract, reducing the medication's dispersion and final residual volume



	POS-1 (MAX)	POS-2
Delivery ml/min	0.53 approx. (1)	0.32 approx. (1)
MMAD	2.91 µm (2)	2.76 µm (2)
Breathable fraction < 5 µm (FPF)	76 % (2)	78.5% (2)
Maximum fill volume	8 ml	

(1) Data shown is in accordance with Flaem Nuova internal procedure. Values shown for the aerosol delivery may vary depending on the patients respiratory capacity.

(2) In vitro testing performed by TÜV Rheinland Italia S.r.l. in collaboration with the University of Parma. More details are available on request.

Accessories

RF7 Dual Speed Plus Nebulizer
Mouthpiece with valve
Soft Touch Adult mask
Soft Touch Child mask
1 m connection tube

RespirAir

Code	FU02P00
Model	P0611EM F1000 + RF7-2
Master carton	8 pieces



8 018294 012038

To know more

