

Wi.Neb go
 Portable and silent nebulizer,
 with rechargeable Li-Ion
 battery pack



**Speed, Silence,
 Comfort**



**2-speed Nebulizer
 with valve system**



**Comfortable adult and
 paediatric Soft Touch
 masks**



**The respiratory
 therapy**

— 50 —
 years of italian excellence

www.flaem.it

Wi.Neb go

Portable and silent nebulizer,
with rechargeable Li-Ion
battery pack

Made in Italy



Technical Features

12 Vdc powered rotary piston compressor (does not require lubrication and maintenance)

Air filtering	For better quality, the compressed air is sucked in from the external environment and filtered with a high efficiency filter system (density= 230 mg/m ³)		
Compressor air output	9 l/min (approx.)		
Max Pressure	2 ± 0.3 bar		
Noise level (at 1 m)	52 dB (A) (approx.)		
Dimensions	13 (W) x 9.5 (D) x 5 (H) cm		
Weight of device	570 g with battery - 430 g without battery		
Power supply/Power	Multivoltage PSU 100/240 V ~ 50/60 Hz/12V DC		
Approvals	0051 In compliance with 93/42/EEC Medical Device Directive (for the 230 V - 50 Hz unit only)	PHTHALATES & BPA FREE	Complies with: European standard EN 10993-1 "Biological Evaluation of Medical Devices" and European Directive 93/42/EEC "Medical Devices" Phthalates free. In conformity with Reg. (EC) no. 1907/2006
		TÜV Rheinland CERTIFIED Type Approved Safety Respirator Production Surveillance www.tuv.com ID: 021000020	In compliance with std. EN 60601-1 and EN 60601-1-1

RF7 Dual Speed Plus nebulizer technical features

Delivery ml/min	POS-1 0.18 approx. ⁽¹⁾ 2.95 µm ⁽²⁾ 74.7 % ⁽²⁾	POS-2 (MAX) 0.42 approx. ⁽¹⁾ 2.44 µm ⁽²⁾ 80.7% ⁽²⁾
MMAD		
Breathable fraction < 5 µm (FPF)		
Maximum fill volume	8 ml ⁽²⁾	

⁽¹⁾ Data shown is in accordance with Flaem Nuova internal procedure I29-P07.5. Values shown for the aerosol delivery may vary depending on the patients respiratory capacity – ⁽²⁾ In vitro testing performed by TÜV Rheinland Italia S.r.l. in collaboration with the University of Parma and in compliance with the EN 13544-1: 2007 + A1 European Standard for aerosol therapy devices Standard. More details are available on request.



Included accessories

- RF7 Dual Speed Plus Nebulizer Kit composed of:
Nebulizer, mouthpiece, nasal prong and Drug saving device
- 1 m connection tube
- Soft Touch Adult Mask
- Soft Touch Pediatric Mask
- Multivoltage power supply unit
- 12VDC cigarette lighter power cord
- Carrying bag
- Li-Ion battery pack

Technical features of battery code ACO377P

Type	Li-ion 7.4V
Battery autonomy	About 45 minutes (at 25° C)
Charging time	About 2 hours
Operating conditions	Temp.: min 10°C; max 40°C RH: min 10%; max 75% Atm. Press.: P min 86KPa; max 106KPa
Storage conditions	Temp.: min -5°C; max 35°C RH: min 10%; max 75% Atm. Press.: P min 86KPa; max 106KPa



Wi.Neb go rechargeable battery included

Code	WIO2P00	8 018294 009939	To know more
Model	PO905EM		
Master carton	4 pieces		

14513AO © Copyright by FLAEM NUOVA S.p.A., 2020 - Materiale divulgativo riservato alle categorie professionali specializzate. FLAEM NUOVA si riserva di apportare modifiche tecniche e funzionali al prodotto della presente scheda senza alcun preavviso. / FLAEM NUOVA reserves the right to make technical and functional alterations to this product without prior notice.